

TDMS

National Toxicology Support

Data Dictionary: TDMS

NATIONAL TOXICOLOGY SUPPORT

Data Dictionary: TDMS

Prepared for: National Institute of Environmental Health Sciences
Research Triangle Park, NC 2779

Contract no.: 273-95-I-0012

Prepared by: Technology Planning & Management Corporation
Imperial Center
Canterbury Hall, Suite 310
4815 Emperor Boulevard
Durham, NC 27703

Revision History

Date Released	Document Number	Revision Level	Reason for Revision
December 13, 1996	TS-1005	2.2	Regenerated TOC, removed Index and references to Certification
October 7, 1996	TS-1005	2.1	Added missing file description
September 3, 1996	TS-1005	2.0	Updated definitions
January 21, 1992	TS-1005	1.0	New

Table of Contents

1.0 Introduction	1
2.0 Database Descriptions.....	3
Bench.....	3
Production.....	3
Quality Assurance	4
Test.....	4
3.0 File Descriptions.....	7
3.1 TDMS & PCT Tables	7
3.2 Receive Status	15
3.3 Protocol Analysis.....	15
3.4 Receive Transaction.....	26
3.5 Audit History.....	29
3.6 EIS/ PEIS	32
4.0 Field Descriptions.....	37
4.1 Protocol Analysis System Views	38
PAS_ACCT_SITE.....	38
PAS_ANIMAL_ID	41
PAS_CALI_WGT.....	44
PAS_CENSUS.....	47
PAS_CLINICAL_OBS	51
PAS_FAC_CHEM.....	55
PAS_LOAD_WT.....	59
PAS_MASTER.....	63
PAS_MCRO_PTH_TBL.....	64
PAS_MICRO_PROC	67
PAS_MORPHOLOGIES.....	69
PAS_ORGANS.....	73
PAS_PH_P_AST_OPR	75
PAS_PREPARTN_SEG.....	78
PAS_PROCEDUR_SEG.....	82
PAS_PROTO_REQ_TIS.....	86
PAS_RACK_CAGE	89
PAS_REMVL_ACTION.....	93

PAS_RESP_PARTIES	95
PAS_SITES	98
PAS_TREATMNT_SEG	100
PAS_WGT_ITEMS	105
4.2 EIS/ PEIS Views	108
TD-ANIMAL.....	108
TD-ANIMAL-DATA.....	121
TD-BALANCE.....	127
TD-BALANCE-CALIBRATION	128
TD-CAGE	130
TD-CAGE-DATA	133
TD-CORRECTION-REASON	146
TD-MICRO-ANIMAL-NOTE	147
TD-ORGAN	149
TD-ORGAN-DATA	152
TD-TEST	157
TD-TEST-NOTE.....	158
TD-TREATMENT.....	159
TD-TREATMENT-NOTE	161
4.3 Transaction Views	162
ANIMAL-MICRO-NOTES-TRANSACTION.....	162
ANIMAL-NOTE-TRANSACTION	169
ANIMAL-OBSERVATION-TRANSACTION	176
ANIMAL-REMOVAL-TRANSACTION	183
ANIMAL-TRANSFER-TRANSACTION	192
ANIMAL-WEIGHT-TRANSACTION.....	198
ANY-ANIMAL-TRANSACTION	205
ANY-MICRO-OBS-TRANSACTION.....	212
ANY-MICRO-ORGAN-TRANSACTION.....	220
BALANCE-CALIBRATION-TRANSACTION	226
BOTTLE-WEIGHT-TRANSACTION	232
CAGE-ID-TRANSACTION	239
CAGE-NOTE-TRANSACTION	246
FEEDER-WEIGHT-TRANSACTION	252
HISTOLOGY-NUMBER-TRANSACTION	259
MICRO-SITE-STATUS-TRANSACTION.....	265
NOT-EXAMINED-TRANSACTION.....	271
ORGAN-NOTES-TRANSACTION	276
TEST-NOTE-TRANSACTION	283
TREAT-NOTE-TRANSACTION	288
UNCERTAIN-TRANSACTION	294
STATUS_FILE.....	299
4.4 Audit History Views	302
ANIMAL-MICRO-NOTES-CORRECTION	302
ANIMAL-NOTE-CORRECTION	305
ANIMAL-OBSERVATION-CORRECTION.....	307
ANIMAL-REMOVAL-CORRECTION.....	310
ANIMAL-TRANSFER-CORRECTION.....	314
ANIMAL-WEIGHT-CORRECTION	316

ANY-ANIMAL-CORRECTION	319
ANY-MICRO-OBS-CORRECTION	322
ANY-MICRO-ORGAN-CORRECTION	325
BOTTLE-WEIGHT-CORRECTION	327
CAGE-ID-CORRECTION	330
CAGE-NOTE-CORRECTION	333
FEEDER-WEIGHT-CORRECTION	335
HISTOLOGY-NUMBER-CORRECTION	338
MICRO-SITE-STATUS-CORRECTION	340
NOT-EXAMINED-CORRECTION	343
ORGAN-NOTES-CORRECTION	345
TEST-NOTE-CORRECTION	347
TREAT-NOTE-CORRECTION	349
UNCERTAIN-CORRECTION	351
DESTINATION_TBLE	353
4.5 TDMS Table & PCT Views	354
TABLES_DESCRIPTN	354
TBACTONC	355
TBAGENCY	357
TBANCLQF	359
TBANMLSUPPLIER	361
TBCLOBSC	362
TBCLOBSCGRP	365
TBCLOBSI	366
TBCONDTC	368
TBCORREASON	370
TBDOSRTE	371
TBFACTLY	373
TBINLIFETERMS	375
TBORGSTA	376
TBPCT	378
TBQUALGRP	388
TBSTAINC	389
TBSTRAIN	391
TBSUBSTR	393
TBTESTYP	395
TBTRTROLE	397
TBUNITS	398
TBWEIGHTCURVE	400
TBWGTOBJ	401
TBWGTSTATUS	403
TB_OPERATOR_ID	404
4.6 ADABAS Physical Definitions	406
ACCT-SITES-PEIS	406
ACTION-TABLE	408
AGENCY-TABLE	409
ANIMAL-CLASS-QUAL-TABLE	411
ANIMAL-DATA-EIS	412
ANIMAL-EIS	413

ANIMAL-ID-EIS.....	416
ANIMAL-NOTE-PEIS	418
ANIMAL-SUPPLIER-TABLE	419
BALANCE	420
BALANCE-CALIBRATION-EIS	421
CAGE-CENSUS-EIS	422
CAGE-CONDITION-TABLE	424
CAGE-DATA-EIS	426
CAGE-EIS.....	429
CALIB-WGTS-EIS	430
CLIN-OBS-GROUP-TABLE	432
CLIN-OBS-INCOMP-TABLE	433
CLIN-OBSERV-EIS	434
CLINICAL-OBS-TABLE	436
CORRECTION-DATA	438
DATABASE-TRANSACTION	445
DESTINATION-TABLE	468
DOSE-ROUTE-TABLE	469
FACILITY-CHEMICAL-EIS.....	470
FACILITY-TABLE	473
IN-LIFE-TERMS-TABLE	475
LOAD-WEIGHT-EIS	476
MASTER.....	479
MICRO-PROC-PEIS	480
MICRO-TABLE-PEIS	481
MODCOMP-CORRECTION-DATA	483
MORPHOLOGIES-PEIS	484
OPERATOR-ID-TABLE	486
ORGAN-DATA-PEIS.....	487
ORGAN-PEIS.....	488
ORGAN-STATUS-TABLE.....	489
ORGANS-PEIS	491
PATHOLOGY-CODE-TABLE	493
PREPARATION-EIS	503
PROCEDURE-EIS.....	505
PROTO-REQ-TISS-PEIS.....	507
PTH-AST-OPR-PEIS.....	509
RACK-CAGE-EIS.....	511
RECEIVE-STATUS	514
REMOVAL-ACTION-EIS.....	515
RESP-PARTY-EIS	516
SITES-PEIS.....	518
STAIN-TABLE	520
STRAIN-TABLE.....	521
SUBSTRAIN-TABLE.....	523
TABLES_DESCRIPTION	524
TBLDAS	525
TD-NOTES	527
TEST.....	528

TEST-TYPE-TABLE.....	529
TREATMENT-EIS.....	531
TREATMENT-PEIS.....	534
UNITS-TABLE.....	535
WEIGHT-ITEM-EIS.....	536
WEIGHT-OBJECT-TABLE.....	538
Appendix A	540
A-1: Legend for File Listing.....	540
A-2: Legend for Listing of Files with Fields.....	541

1.0 Introduction

This document is used by technical support staff assigned to the Toxicology Data Management System (TDMS). It describes TDMS database parameters, fields, and files. Readers should be technically familiar with TDMS *before* using this manual as a reference. For a general description, refer to the manual, [*TDMS System Overview*](#). The Librarian can help you to obtain a copy.

2.0 Database Descriptions

This chapter describes each of the four databases that compose TDMS.

Bench

The TDMS benchmark database is structured exactly like the production database. It differs only in content. The content comprises copies of production tests and tests created by the testing staff. The purpose of the database is to provide testers with a suitable data environment for acceptance testing of TDMS and LDAS.

DB-ID: TDMS_BENCH

DB-NR: 5

File Parameters		Device Usage		
		Container	Disk	Size
MAX.FILES	250	ASSO	DISK\$ADADSK2	310 Cylinders
CHECKPOINT	1	DATA	DISK\$ADADSK1	550 Cylinders
SECURITY	2	WORK	DISK\$ADADSK6	100 Cylinders
SYSTEM	11			
PREDICT	10			

Production

The production database is a repository for all data generated and used by the TDMS application. Access to, and update of its contents is strictly controlled. Database control is the primary responsibility of the TDMS database administrators.

DB-ID: TDMS_PROD

DB-NR: 7

File Parameters		Device Usage		
		Container	Disk	Size
MAX.FILES	250	ASSO	DSK_ADADSK1,-2,-3 (RA80)	624 cylinders each
CHECKPOINT	1	DATA	DSK_ADADSK4,-5,-6 (RA80)	1246 cylinders each

DB-ID: TDMS_PROD

DB-NR: 7 (continued)

File Parameters		Device Usage		
		Container	Disk	Size
SECURITY	2	WORK	DSK_ADADSK7 (RA80)	190 cylinders
SYSTEM	11			
PREDICT	10			

Quality Assurance

The TDMS QA database is structured exactly like the production database. It differs only in content. The content comprises copies of production tests and tests created by the QA staff. The purpose of the database is to provide testers with a suitable data environment for testing TDMS components.

DB-ID: TDMS_QA DB-NR: 3

File Parameters		Device Usage		
		Container	Disk	Size
MAX.FILES	250	ASSO	DSK_ADADSK3	155 Cylinders
CHECKPOINT	1	DATA	DSK_ADADSK8	600 Cylinders
SECURITY	2	WORK	DSK_ADADSK7	50 Cylinders
SYSTEM	11			
PREDICT	10			

Test

The TDMS test database is structured exactly like the production database. It differs only in content. The content comprises copies of production tests and tests created by the Development staff. The purpose of the database is to provide developers with a suitable data environment for their work.

DB-ID: TDMS_TEST DB-NR: 1

File Parameters		Device Usage		
		Container	Disk	Size
MAX.FILES	250	ASSO	DSK_ADADSK3	855 Cylinders
CHECKPOINT	1	DATA	DSK_ADADSK7	601 Cylinders
SECURITY	2	WORK	DSK_ADADSK1	40 Cylinders
SYSTEM	11			
PREDICT	10			

3.0 File Descriptions

This chapter describes each of the files that compose the databases of TDMS. For a definition of abbreviations used in this chapter, refer to the legend in *Appendix A-1*.

3.1 TDMS & PCT Tables

TB_OPERATOR_ID

Typ Nr U 14

User view for operator table. Corresponds to the LDAS "operator" table. Contains all valid operators of TDMS. Each operator is assigned a unique, sequential number within the facility with which they are associated. All users of the TDMS menu have the first eight characters of their user account entered in the STARUSER field. The menu uses STARUSER to validate users invoking the menu.

OPERATOR-ID-TABLE

Typ Nr A 14

Physical layout for operator table. Supports user view of same name.

TBSTAINC

Typ Nr U 15

Logical layout of stain table. The table contains stains that could be used for tissue section cuts. Stains are a legacy to TDMS, with no use in the current version. PEX software does not require any stain codes to be associated with a test.

STAIN-TABLE

Typ Nr A 15

Physical layout of stain table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBSTAINC. TBGET uses that layout description to form direct calls to retrieve data from the table for PCT and PAS validation reports. This layout corresponds to the view TBSTAINC. (Stain data is not collected or reported by the current version of TDMS.)

TABLES_DESCRIPTOR

Typ Nr U 16

Logical layout of the Tables Description table, which is used by NATURAL.

TABLES_DESCRIPTION

Typ Nr A 16

Physical layout of the Tables Description table. TBGET retrieves data from this file. File contains physical layouts of TDMS tables used by TBGET.

TBSUBSTR

Typ Nr U 17

Logical layout of substrain table. Substrains are a legacy to TDMS, with no use in the current version. PAS software requires a substrain code for inlife test protocol, so the codes for "Not Applicable" or "Not Specified" are entered.

SUBSTRAIN-TABLE

Typ Nr A 17

Physical layout of substrain table. This layout is also stored in the Tables Description table. TBGET uses the layout description to form direct calls and retrieve data from the table for PAS validation reports. This layout corresponds to the view TBSUBSTR. Substrain data is not used in TDMS. Dummy values are entered in PAS.

TBDOSRTE

Typ Nr U 18

Logical layout of Dose Route table. This table contains all routes by which a dose for a test may be administered. Current dose routes are categorized as follows: food, water, gavage, skin paint, and respiration.

DOSE-ROUTE-TABLE

Typ Nr A 18

Physical layout of Dose Route table. This layout is also stored in the Tables Description table. TBGET uses the layout description to form direct calls and retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBDOSRTE.

TBWGTOBJ

Typ Nr U 19

Logical layout of Weight Object table. This table contains objects that may be weighed during the inlife portion of a test and for which there are weight ranges.

WEIGHT-OBJECT-TABLE

Typ Nr A 19

Physical layout of Weight Object table. This layout is also stored in the Tables Description table. TBGET uses the layout description to form direct calls and retrieve data from the table for PAS validation reports. This layout corresponds to the view TBWGTOBJ.

TBUNITS

Typ Nr U 20

Logical layout of Units table. This table lists quantitative units of measure. This table is a legacy to TDMS. The only valid unit of measure in TDMS is decigrams and it is hard-coded into TDMS software.

UNITS-TABLE

Typ Nr A 20

Physical layout of Units table. This layout is described in the Tables Description table. The table has no meaning in TDMS. This file corresponds to the view TBUNITS.

TBSTRAIN

Typ Nr U 21

Logical layout of strain table. Table contains species and species/strains that are used in TDMS. The three current species are mice, rats, and hamsters. Strain information appears on most EIS/ PEIS reports.

STRAIN-TABLE

Typ Nr A 21

Physical layout of strain table. This layout is also stored in the Tables Description table. TBGET uses the layout description to form direct calls and retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBSTRAIN.

TBAGENCY

Typ Nr U 22

Logical layout of agency table. Contains a list of agencies. Table is a legacy to TDMS. Only one entry is valid: the entry for NIEHS/NTP. The code for that entry is '05' and is built into the nine-digit test number entered in PAS. The agency code is hard-coded by the menu system when EIS/PEIS reports are requested.

AGENCY-TABLE

Typ Nr A 22

Physical layout of agency table. This layout is stored in the Tables Description table. TBGET uses the layout description to form direct calls and retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBAGENCY.

TBANCLQF

Typ Nr U 23

Logical layout of animal class qualifier table. Animal class qualifiers are a legacy to TDMS, with no use in its current version. The PEX maintenance software requires an animal class qualifier, so the code for "Not Applicable" is entered.

ANIMAL-CLASS-QUAL-TABLE

Typ Nr A 23

Physical layout of animal class qualifier table. This layout is also stored in the Tables Description table. TBGET uses the layout description to form direct calls and retrieve data from the table for PAS validation reports. This file corresponds to the view TBANCLQF.

TBCONDTC

Typ Nr U 24

Logical layout of cage condition table. This table is a legacy to TDMS and has no use in the current version of TDMS.

CAGE-CONDITION-TABLE

Typ Nr A 24

Physical layout of cage condition table. Layout is stored in the Tables Description table. TBGET uses the information to form direct calls which retrieve data from this table for PAS validation reports. This layout corresponds to the view TBCONDTC.

TBCLOBSC

Typ Nr U 25

Logical layout of the clinical observation table, which contains clinical observations, inlife sites, inlife qualifiers, and removal reasons. Data from this table appears in nearly all EIS/PEIS reports in some manner or fashion.

CLINICAL-OBS-TABLE

Typ Nr A 25

Physical layout of clinical observation table. This layout is stored in the Tables Description table. TBGET uses the information to form direct calls against this table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBCLOBSC.

TBCLOBSI

Typ Nr U 26

Logical layout of clinical incompatibility table. This is a legacy table with no purpose in the current version of TDMS. Clinical observation incompatibility is simply enforced by LDAS software according to the following rules: *unremarkable* is incompatible with any other observation; *hyperactive* and *lethargic* are mutually incompatible.

CLIN-OBS-INCOMP-TABLE

Typ Nr A 26

Physical layout of clinical observation incompatibility table. This layout is stored in the Tables Description table. TBGET uses the layout description to form direct calls and retrieve data from the table for PAS validation reports. This layout corresponds to the view TBCLOBSI.

TBFACLTY

Typ Nr U 27

Logical layout of facility table. Contains all facilities for TDMS, current and former. Used in splash pages of EIS/PEIS reports, PAS validation reports, and almost every other part of the application.

FACILITY-TABLE

Typ Nr A 27

Physical layout of facility table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBFACLTY. TBGET uses the layout description to form direct calls and retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBFACLTY.

TBTESTYP

Typ Nr U 28

Logical layout of test type table. Table contains possible test types by which a test may be classified. Test type information appears on most EIS/PEIS reports, usually on the splash page.

TEST-TYPE-TABLE

Typ Nr A 28

Physical layout of test type table. This layout is also stored in the Tables Description table. TBGET uses the layout description to form direct calls and retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBTESTYP.

TBACTONC

Typ Nr U 29

Logical layout of action table. Table lists actions that may be performed by animal caretakers as specified in protocol procedures. Information is downloaded to LDAS and is used in creating the procedure action schedule.

ACTION-TABLE

Typ Nr A 29

Physical layout of action table. This layout is also stored in the Tables Description table. TBGET uses the layout description to form direct calls and retrieve data from the table for PAS validation reports. This layout corresponds to the view TBACTONC.

TBPCT

Typ Nr U 30

Logical layout of the pathology code table. This layout is the basis on which NATURAL retrieves data for downloads and pathology reports.

PATHOLOGY-CODE-TABLE

Typ Nr A 30

Physical layout of the pathology code table. This layout is the basis on which PCTGET retrieves data for PCT and pathology reports. (PCTGET retrieves all fields in the table for pathology code table reports.) This layout corresponds to the view TBPCT.

TBWGTSTATUS

Typ Nr U 31

Logical view of weight status table. Contains valid statuses for objects whose scheduled weight action is being bypassed. Downloaded to LDAS as part of test protocol. Displayed in all reports that contain weight detail information.

TBTRTROLE

Typ Nr U 31

Logical layout of treatment role table. Contains the possible roles that a treatment group may have on a test. Some older TDMS software still has these roles hard-coded.

TBCORREASON

Typ Nr U 31

Logical view of correction reason table. Contains the valid correction reasons that can be used to classify an entry made in the error correction system.

TBWEIGHTCURVE

Typ Nr U 31

Logical view of animal weight growth table. Contains weight curve upper and lower limits as a function of species, sex, and week-on-test.

TBQUALGRP

Typ Nr U 31

Logical view of the inlife qualifier group table. It contains the possible groups by which current inlife qualifiers are classified.

TBLDAS

Typ Nr A 31

Physical layout of supertable that embraces five tables: weight status, treatment role, correction reason, inlife qualifier group, and weight curve tables. Each of these tables is defined in the Tables Description table. TBGET uses the descriptions to retrieve data via direct calls. Supports views named as TABLE_NAME.

TBANMLSUPPLIER

Typ Nr U 32

Logical view of animal supplier table. Contains valid animal suppliers for tests. (There is a dummy record in the table for tests that preceded entry of animal suppliers in test protocol. This record has a key value of 0 and text of "NOT SPECIFIED".)

ANIMAL-SUPPLIER-TABLE

Typ Nr A 32

Physical layout of animal supplier table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBANMLSU. TBGET uses that layout description to form direct calls to retrieve data from the table for PAS validation reports. This layout corresponds to the view TBANMLSUPPLIER.

TBINLIFETERMS

Typ Nr U 33

Logical view of inlife terms descriptions table. Contains synonyms for clinical observation codes. Used as part of LDAS help system.

IN-LIFE-TERMS-TABLE

Typ Nr A 33

Physical layout of inlife terms table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBINLIFE. TBGET can use that layout description to form direct calls to retrieve data from the table, but does not. This layout corresponds to the view TBINLIFETERMS.

TBORGSTA

Typ Nr U 34

Logical layout of organ (tissue) status table. Contains the possible statuses by which an organ (tissue) that is being microscopically evaluated may be classified.

ORGAN-STATUS-TABLE

Typ Nr A 34

Physical layout of organ (tissue) status table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBORGSTA. TBGET uses that layout description to form direct calls to retrieve data from the table for pathology reports. This layout corresponds to the view TBORGSTA.

DESTINATION_TBLE

Typ Nr U 35

View of valid TDMS output destinations. Used by TDMS menus to validate output destinations on TDMS screens.

DESTINATION-TABLE

Typ Nr A 35

Physical layout of file 35. Supports the DESTINATION_TBLE view. No TDMS software uses physical layout to issue direct calls.

TBCLOBSCGRP

Typ Nr U 36

View of clinical observation groups that classify entries in TBCLOBSC with SUBSETINDICATOR equal to "A".

CLIN-OBS-GROUP-TABLE

Typ Nr A 36

Physical layout of file 36. Supports the TBCLOBSCGRP view. No TDMS software uses this layout to issue direct calls.

3.2 Receive Status**STATUS_FILE**

Typ Nr U 38

View of RECEIVE transaction set status information. Provides a history of set processing.

RECEIVE-STATUS

Typ Nr A 38

Physical layout of file 38. Support the STATUS_FILE view. No TDMS software accesses the file using direct calls.

3.3 Protocol Analysis**PSA_MASTER**

Typ Nr U 39

Supports latent PAS archival feature.

PAS_MASTER

Typ Nr U 40

PAS master table. Indicates whether a test has been registered. Used by report menus to validate a selected test.

PAS_TREATMNT_SEG

Typ Nr U 41

View of treatment, and treatment regimen, protocol. Treatment groups are generally organized according to dose levels and gender.

PAS_PROCEDUR_SEG

Typ Nr U 42

View of in-life procedure action protocol. The procedures are organized into groups of actions. These groups comprise the schedule of activities that apply to cages. Its information is downloaded to LDAS.

PAS_PREPARTN_SEG

Typ Nr U 43

View of in-life, compound preparation protocol. The preparations define the composition of the dose mixtures used on the test. Information in the table is downloaded to LDAS and is used by many of the EIS/PEIS reports.

PAS_CENSUS

Typ Nr U 44

View of in-life, cage census protocol. Used to describe the cages on the test, including weight, treatment, and procedure groups, and number of animals to be allocated. Used by Download to construct what will be the LDAS CAGE table.

PAS_RESP_PARTIES

Typ Nr U 45

View of information on the people who share primary responsibility for the test. Not used by Download or EIS/PEIS reports.

PAS_WGT_ITEMS

Typ Nr U 46

View of weight group in-life protocol. Contains information about the weight groups to which cages are assigned. Specifies absolute minimum and maximum weights for used and full containers.

PAS_LOAD_WT

Typ Nr U 47

View of in-life protocol which describes the characteristics of the animals used on the test. Used throughout the system.

PAS_CLINICAL_OBS

Typ Nr U 48

View of in-life term protocol. Includes the clinical signs, sites, qualifiers, and removal reasons that may be used on the test. (For older MODCOMP tests, it will also include the cage conditions that were used for cage observations.) Processed by Download to generate all in-life term tables.

PAS_ANIMAL_ID

Typ Nr U 49

View of in-life, animal identification protocol. File is not used since implementation of LDAS. There is still the option to maintain the file in

PEX, and validation reports process it, but, otherwise, the system does not reference it. However, the value of CID-FLAG does indicate which of the two methods of constructing TDA-CARCASS-IDENTIFICATION-DE values in TD-ANIMAL was used on the test.

PAS_REMVL_ACTION

Typ Nr U 50

View of in-life, removal action protocol. These actions are those to be performed by the data collector when an animal is removed from the test. Its information is downloaded to LDAS.

PAS_RACK_CAGE

Typ Nr U 51

View of in-life protocol indicating the physical location of cages. Used by download to create the input for the initialization of the CAGE table.

PAS_ORGANS

Typ Nr U 52

View of pathology protocol. Lists the organs that may be used for pathology data collection. This file is initially populated by copying from one of the universal download tests. It is rarely modified after initialization.

PAS_SITES

Typ Nr U 54

View of pathology protocol. Lists the tissue sites that may be used for pathology data collection. This file is initially populated by copying from one of the universal download tests. It is rarely modified after initialization.

PAS_MORPHOLOGIES

Typ Nr U 53

View of pathology protocol. Lists the morphologies that may be used for pathology data collection. This file is initially populated by copying from one of the universal download tests. It is rarely modified after initialization.

PAS_ACCT_SITE

Typ Nr U 55

View of pathology accountable site protocol. Accountable sites pertain only to tests whose pathology data was collected exclusively on the old, MODCOMP based, collection system. They are sites on tissue that the

collectors were obligated to evaluate as if they were organs. (LDAS elevates the tissue sites to the level of tissue.) This view is not used by the current version of the download component.

PAS_MICRO_PROC

Typ Nr U 56

View pertains to pathology protocol for test whose pathology data was collected using the MODCOMP collection system. Described protocol dependent procedures to be used in pathology evaluation.

PAS_PH_P_AST_OPR

Typ Nr U 57

View of operators assigned to a test for pathology data collection. For tests conducted using LDAS, this information is assigned on the local machine and uploaded to the mainframe DBMS; for MODCOMP tests, the information was entered in PEX and downloaded.

PAS_MCRO_PTH_TBL

Typ Nr U 58

View of pathology terms, including qualifiers, stain codes, and cause of death and morbidity codes Only qualifiers are relevant to the current version of TDMS. New records are generally populated by copying from one of the universal download tests. It is generally not modified after creation.

PAS_FAC_CHEM

Typ Nr U 59

Basic information about the test, including the test compound and facility. Used throughout the system. Equivalent to the TEST table on LDAS.

PAS_CALI_WGT

Typ Nr U 60

View of in-life, calibration protocol. It specifies the rules for calibrations of balances used for weight data collection. It is processed by the current version of download. It is always entered for tests on which weights are to be collected.

PSA_TREATMNT_SEG

Typ Nr U 61

Supports latent PAS archival feature.

PSA_PROCEDUR_SEG

Typ Nr U 62

Supports latent PAS archival feature.

PSA_PREPARTN_SEG

Typ Nr U 63

Supports latent PAS archival feature.

PSA_CENSUS

Typ Nr U 64

Supports latent PAS archival feature.

PSA_RESP_PARTIES

Typ Nr U 65

Supports latent PAS archival feature.

PSA_WGT_ITEMS

Typ Nr U 66

Supports latent PAS archival feature.

PSA_LOAD_WT

Typ Nr U 67

Supports latent PAS archival feature.

PSA_CLINICAL_OBS

Typ Nr U 68

Supports latent PAS archival feature.

PSA_ANIMAL_ID

Typ Nr U 69

Supports latent PAS archival feature.

PSA_REMVL_ACTION

Typ Nr U 70

Supports latent PAS archival feature.

PSA_RACK_CAGE

Typ Nr U 71

Supports latent PAS archival feature.

PSA_ORGANS

Typ Nr U 72

Supports latent PAS archival feature.

PSA_MORPHOLOGIES

Typ Nr U 73

Supports latent PAS archival feature.

PSA_SITES

Typ Nr U 74

Supports latent PAS archival feature.

PSA_ACCT_SITE

Typ Nr U 75

Supports latent PAS archiving feature.

PSA_MICRO_PROC

Typ Nr U 76

Supports latent PAS archival feature.

PSA_PH_P_AST_OPR

Typ Nr U 77

Supports latent PAS archival feature.

PSA_MCRO_PTH_TBL

Typ Nr U 78

Supports latent PAS archival feature.

PSA_FAC_CHEM

Typ Nr U 79

Supports latent PAS archival feature.

PSA_CALI_WGT

Typ Nr U 80

Supports latent PAS archival feature.

PAS_PROTO_REQ_TIS

Typ Nr U 84

View of pathology protocol. Defines the criteria by which tissues are required to be evaluated. Criteria are evaluation based on treatment

group, removal reason, and removal date range. Can be created by a liaison and downloaded, or created on the local machine and uploaded. Used by the P14.

MASTER-ARCH

Typ Nr A 39

Physical layout of file 39. Supports the PSA_MASTER view.

MASTER

Typ Nr A 40

Physical layout of file 40. Supports the PAS_MASTER view. This is the layout by which PASGET selects data for PAS validation reports and EIS/ PEIS reports.

TREATMENT-EIS

Typ Nr A 41

Physical layout of file 41. Supports the PAS_TREATMNT_SEG view. This is the layout by which PASGET selects data for PAS validation reports and EIS/PEIS reports.

PROCEDURE-EIS

Typ Nr A 42

Physical layout of file 42. Supports the PAS_PROCEDUR_SEG view. This is the layout by which PASGET selects data for PAS validation reports and EIS/PEIS reports.

PREPARATION-EIS

Typ Nr A 43

Physical layout of file 43. Supports the PAS_PREPARTN_SEG view. This is the layout by which PASGET selects data for PAS validation reports and EIS/PEIS reports.

CAGECENSUSEIS

Typ Nr A 44

Physical layout of file 44. Supports the PAS_CENSUS view. This is the layout by which PASGET selects data for PAS validation reports .

RESP-PARTY-EIS

Typ Nr A 45

Physical layout of file 45. Supports the PAS_RESP_PARTIES view. This is the layout by which PASGET selects data for PAS validation reports.

WEIGHT-ITEM-EIS

Typ Nr A 46

Physical layout of file 46. Supports the PAS_WGT_ITEMS view. This is the layout by which PASGET selects data for PAS validation reports.

LOAD-WEIGHT-EIS

Typ Nr A 47

Physical layout of file 47. Supports the PAS_LOAD_WT view. This is the layout by which PASGET selects data for PAS validation reports.

CLIN-OBSERV-EIS

Typ Nr A 48

Physical layout of file 48. Supports the PAS_CLINICAL_OBS view. This is the layout by which PASGET selects data for PAS validation reports.

ANIMAL-ID-EIS

Typ Nr A 49

Physical layout of file 49. Support the PAS_ANIMAL_ID view. This is the layout used by which PASGET issues direct calls for PAS validation reports.

REMOVAL-ACTION-EIS

Typ Nr A 50

Physical layout of file 50. Supports the PAS_REMVL_ACTION view. PASGET uses this layout to issue direct calls for PAS validation reports.

RACK-CAGE-EIS

Typ Nr A 51

Physical layout of file 51. Supports the PAS_RACK_CAGE view. PASGET uses this layout to issue direct calls for PAS validation reports.

ORGANS-PEIS

Typ Nr A 52

Physical layout of file 52. Supports the PAS_ORGANS view. PASGET uses this layout to issue direct calls for PAS validation reports.

MORPHOLOGIES-PEIS

Typ Nr A 53

Physical layout of file 53. Supports the PAS_MORPHOLOGIES view. PASGET uses this layout to issue direct calls for PAS validation reports.

SITES-PEIS

Typ Nr A 54

Physical layout of file 54. Supports the PAS_SITES view. PASGET uses this layout to issue direct calls for PAS validation reports.

ACCT-SITES-PEIS

Typ Nr A 55

Physical layout of file 55. Support the PAS_ACCT_SITE view. PASGET uses this layout to issue direct calls for EIS/PEIS reports.

MICRO-PROC-PEIS

Typ Nr A 56

Physical layout of file 56. Support the PAS_MICRO_PROC view. No TDMS software uses this layout to issue direct calls.

PTH-AST-OPR-PEIS

Typ Nr A 57

Physical layout of file 57. Support the PAS_PH_P_AST_OPR view. PASGET uses this layout to issue direct calls for PAS validation reports.

MICRO-TABLE-PEIS

Typ Nr A 58

Physical layout of file 58. Support the PAS_MCRO_PTH_TBL view. PASGET uses this layout to issue direct calls for PAS validation reports.

FACILITY-CHEMICAL-EIS

Typ Nr A 59

Physical layout of file 59. Support the PAS_FAC_CHEM view. PASGET uses this layout to issue direct calls for PAS validation and EIS/PEIS reports.

CALIB-WGTS-EIS

Typ Nr A 60

Physical layout of file 60. Support the PAS_CALI_WGT view. PASGET uses this layout to issue direct calls for PAS validation reports.

TREATMENT-ARCH

Typ Nr A 61

Physical layout of file 61. Supports the PSA_TREATMNT_SEG view.

PROCEDURE-ARCH

Typ Nr A 62

Physical layout of file 62. Supports the PSA_PROCEDUR_SEG view.

PREPARATION-ARCH

Typ Nr A 63

Physical layout of file 63. Supports the PSA_PREPARTN_SEG view.

CAGE-CENSUS-ARCH

Typ Nr A 64

Physical layout of file 64. Supports the PSA_CENSUS view.

RESP-PARTY-ARCH

Typ Nr A 65

Physical layout of file 65. Supports the PAS_RESP_PARTIES view.

WEIGHT-ITEM-ARCH

Typ Nr A 66

Physical layout of file 66. Supports the PSA_WGT_ITEMS view.

LOAD-WEIGHT-ARCH

Typ Nr A 67

Physical layout of file 67. Supports the PSA_LOAD_WT view.

CLIN-OBSERV-ARCH

Typ Nr A 68

Physical layout of file 68. Supports the PSA_CLINICAL_OBS view.

ANIMAL-ID-ARCH

Typ Nr A 69

Physical layout of file 69. Support the PSA_ANIMAL_ID view.

REMOVAL-ACTION-ARCH

Typ Nr A 70

Physical layout of file 70. Supports the PSA_REMVL_ACTION view.

RACK-CAGE-ARCH

Typ Nr A 71

Physical layout of file 71. Supports the PSA_RACK_CAGE view.

ORGANS-ARCH

Typ Nr A 72

Physical layout of file 72. Supports the PSA_ORGANS view.

MORPHOLOGIES-ARCH

Typ Nr A 73

Physical layout of file 73. Supports the PSA_MORPHOLOGIES view.

SITES-ARCH

Typ Nr A 74

Physical layout of file 74. Supports the PSA_SITES view.

ACCT-SITES-ARCH

Typ Nr A 75

Physical layout of file 75. Support the PSA_ACCT_SITE view.

MICRO-PROC-ARCH

Typ Nr A 76

Physical layout of file 76. Support the PSA_MICRO_PROC view.

PTH-AST-OPR-ARCH

Typ Nr A 77

Physical layout of file 77. Support the PSA_PH_P_AST_OPR view.

MICRO-TABLE-ARCH

Typ Nr A 78

Physical layout of file 78. Support the PSA_MCRO_PTH_TBL view.

FACILITY-CHEMICAL-ARCH

Typ Nr A 79

Physical layout of file 79. Support the PSA_FAC_CHEM view.

CALIB-WGTS-ARCH

Typ Nr A 80

Physical layout of file 80. Support the PSA_CALI_WGT view.

PROTO-REQ-TISS-PEIS

Typ Nr A 84

Physical layout of file 84. Support the PAS_PROTO_REQ_TIS view.
No direct calls are issued based on this layout.

3.4 Receive Transaction

ANY-ANIMAL-TRANSACTION

Typ Nr U 86

Logical layout for animal allocation transactions. Loaded by RECEIVE and ECS. Processed by Update to add, modify, or delete records to/from the TD-ANIMAL file. For corrections, Update also uses the information to add records to the ANY-ANIMAL-CORRECTION file.

ANIMAL-TRANSFER-TRANSACTION

Typ Nr U 86

Logical layout for animal transfer transactions. Loaded by RECEIVE and ECS. Processed by Update to modify records in the TD-ANIMAL file. For corrections, Update also uses the information to add records to the ANY-ANIMAL-CORRECTION file.

ANIMAL-REMOVAL-TRANSACTION

Typ Nr U 86

Logical layout of animal removal transactions. Loaded by RECEIVE and ECS. Used by Up-date to modify records in TD-ANIMAL. For corrections, the information is also used by Update to add records to ANIMAL-REMOVAL-CORRECTION.

CAGE-ID-TRANSACTION

Typ Nr U 86

Logical layout for cage initialization transactions. Loaded by ECS and RECEIVE. Used by Update to create or modify data in TD-CAGE, and TD-TEST. For corrections, also used to add records to CAGE-ID-CORRECTION.

ANIMAL-OBSERVATION-TRANSACTION

Typ Nr U 86

Logical layout of animal observation transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ANIMAL-DATA file. For corrections, Update also uses the information to add records to ANIMAL-OBSERVATION-CORRECTION.

ANIMAL-WEIGHT-TRANSACTION

Typ Nr U 86

Logical layout of animal weight transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ANIMAL-DATA file. For corrections, Update also uses the information to add records to the ANIMAL-WEIGHT-CORRECTION file.

ANIMAL-NOTE-TRANSACTION

Typ Nr U 86

Logical layout of animal note transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, or delete records to the TD-ANIMAL-DATA file. For corrections, Update also uses the information to add records to ANIMAL-NOTE-CORRECTION.

FEEDER-WEIGHT-TRANSACTION

Typ Nr U 86

Logical layout of feeder weight transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records from the TD-CAGE-DATA file. For corrections, information used to add records to FEEDER-WEIGHT-CORRECTION.

BOTTLE-WEIGHT-TRANSACTION

Typ Nr U 86

Logical layout of bottle weight transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-CAGE-DATA table. For corrections, the information is also used to create records for BOTTLE-WEIGHT-CORRECTION.

CAGE-NOTE-TRANSACTION

Typ Nr U 86

Logical layout of cage note transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-CAGE-DATA file. For corrections, the information is also used to add records to the CAGE-NOTE-TRANSACTION table.

BALANCE-CALIBRATION-TRANSACTION

Typ Nr U 86

Logical layout of balance calibration transactions. Loaded by RECEIVE. Read by Update to add records to the TD-BALANCE-CALIBRATION file.

TEST-NOTE-TRANSACTION

Typ Nr U 86

Logical layout of test note transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-TEST-NOTE file. For corrections, the information is used to add records to TEST-NOTE-CORRECTION.

TREAT-NOTE-TRANSACTION

Typ Nr U 86

Logical layout of treatment group note transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-TREAT-NOTE file. For corrections, the information is used to add records to TREAT-NOTE-CORRECTION.

HISTOLOGY-NUMBER-TRANSACTION

Typ Nr U 86

Logical layout of histology number transactions. Loaded by RECEIVE and ECS. Read by Update to modify records in the TD-ANIMAL file. For corrections, also used to add records to HISTOLOGY-NUMBER CORRECTION.

ANY-MICRO-ORGAN-TRANSACTION

Typ Nr U 86

Logical layout of organ (tissue) status transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ORGAN file. For corrections, used to add records to ANY-MICRO-ORGAN-CORRECTION.

ANY-MICRO-OBS-TRANSACTION

Typ Nr U 86

Logical layout of microscopic observation transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ORGAN-DATA file. For corrections, also used to add records to ANY-MICRO-OBS-CORRECTION.

ORGAN-NOTES-TRANSACTION

Typ Nr U 86

Logical Layout Of Organ Notes transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ORGAN-DATA file. If transactions are audit changes, then Update also uses the information to add records to ORGAN-NOTES CORRECTION.

MICRO-SITE-STATUS-TRANSACTION

Typ Nr U 86

Logical layout for site status transactions. Site status data was a feature of the MODCOMP data collection that was made defunct by LDAS. There is an ECS module that will add records to this table, but it is incorrect procedure to use that module. Correct procedure calls for a test to be converted to LDAS format prior to modification.

ANIMAL-MICRO-NOTES-TRANSACTION

Typ Nr U 86

Logical layout of pathology animal notes transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD- MICRO-ANIMAL-NOTE file. If transactions are corrections or audit changes, then Update uses the information to add records to the ANIMAL-MICRO-NOTES-CORRECTION file.

NOT-EXAMINED-TRANSACTION

Typ Nr U 86

Logical layout of "not examined" transactions. Loaded by RECEIVE and ECS. Read by Update to modify records in the TD-ANIMAL file. If a correction, information also used to add records to NOT-EXAMINED-CORRECTION.

UNCERTAIN-TRANSACTION

Typ Nr U 86

Logical layout for uncertain cause of death transactions. Loaded by RECEIVE and ECS. Read by Update to modify records in the TD-ANIMAL file. For correction, information is used to add records to UNCERTAIN-CORRECTION.

DATABASE-TRANSACTION

Typ Nr A 86

Physical/Logical Layout Of Transaction Super Table, Which embraces 21 user views. Each view corresponds to one of the transaction types that UPDATE will process. This file is updated by both RECEIVE and ECS. The file design is not programmer friendly as it maintains the older VSAM file structure.

3.5 Audit History

ANIMAL-MICRO-NOTES-CORRECTION

Typ Nr U 87

Logical view of pathology animal note audit changes. The master file is TD-MICRO-ANIMAL-NOTES. Data from this view is not retrieved by any current TDMS component.

ANIMAL-NOTE-CORRECTION

Typ Nr U 87

Logical view of animal note corrections. Master file is TD-ANIMAL-DATA. View is used by EIS Report 22.

ANIMAL-OBSERVATION-CORRECTION

Typ Nr U 87

Logical view of animal observation corrections. Master file is TD-ANIMAL-DATA. Used by EIS Report 22.

ANIMAL-WEIGHT-CORRECTION

Typ Nr U 87

Logical layout of animal weight corrections. Master file is TD-ANIMAL-DATA. Used by EIS Report 22.

ANIMAL-REMOVAL-CORRECTION

Typ Nr U 87

Logical view of animal removal corrections. Master file is TD-ANIMAL. Used by EIS Report 22.

ANIMAL-TRANSFER-CORRECTION

Typ Nr U 87

Logical view of animal transfer corrections. Master table is TD-ANIMAL. Used by EIS Report 22.

ANY-ANIMAL-CORRECTION

Typ Nr U 87

Logical view of animal allocation corrections. Master table is TD-ANIMAL. Used by EIS Report 22.

CAGE-ID-CORRECTION

Typ Nr U 87

Logical view of cage ID corrections. Master table is TD-CAGE. Used by EIS Report 22.

FEEDER-WEIGHT-CORRECTION

Typ Nr U 87

Logical view of feeder weight corrections. Master file is TD-CAGE-DATA. Used by EIS Report 22.

BOTTLE-WEIGHT-CORRECTION

Typ Nr U 87

Logical view of bottle weight corrections. Master table is TD-CAGE-DATA. Used by EIS Report 22.

CAGE-NOTE-CORRECTION

Typ Nr U 87

Logical view of cage note corrections. Master file is TD-CAGE-DATA. Used by EIS Report 22.

TEST-NOTE-CORRECTION

Typ Nr U 87

Logical view for test note corrections. Master file is TD-TESTNOTE. Used by EIS Report 22.

TREAT-NOTE-CORRECTION

Typ Nr U 87

Logical view of treatment note corrections. Master file is TD-TREATMENTNOTE. Used by EIS Report 22.

HISTOLOGY-NUMBER-CORRECTION

Typ Nr U 87

Logical view of histology number audit changes. Master file is TD-ANIMAL. Not used for retrieval by any TDMS component.

ANY-MICRO-ORGAN-CORRECTION

Typ Nr U 87

Logical view for tissue status observation audit changes. Master file is TD-ORGAN. Not used for retrieval by any TDMS component.

ANY-MICRO-OBS-CORRECTION

Typ Nr U 87

Logical view of pathology microscopic observation audit changes. Master file is TD-ORGAN-DATA. Not referenced by any TDMS retrieval component.

ORGAN-NOTES-CORRECTION

Typ Nr U 87

Logical view of organ note audit changes. Master file is TD-ORGAN-DATA. Not used for retrieval by any TDMS component.

MICRO-SITE-STATUS-CORRECTION

Typ Nr U 87

Logical view of accountable site status observation corrections. Master file is TD-ORGAN. Not used for retrieval by any TDMS component.

NOT-EXAMINED-CORRECTION

Typ Nr U 87

Logical view of not examined audit changes. Master file is TD-ANIMAL. Not used for retrieval by any TDMS component.

UNCERTAIN-CORRECTION

Typ Nr U 87

Logical view of uncertain cause of death audit changes. Master table is TD-ANIMAL. Not used for retrieval by any TDMS component.

CORRECTION-DATA

Typ Nr A 87

This is the physical layout of correction data. The physical file supports the views for each of the 20 correction segments. In-life correction segments are reported in EIS report 22; pathology correction segments are not reported at all. File contains only correction data since the introduction of LDAS. Prior correction data is found in TD-CORRECTION-REASON.

3.6 EIS/ PEIS**TD-TREATMENT**

Typ Nr U 100

View of pre LDAS protocol required tissues. Only used by the P14 for those tests that were not converted to LDAS. Replaced by the PAS protocol required tissues segment.

TD-OTHER-TERMS

Typ Nr U 101

Contains pathology terms assigned by the lab during collection using MODCOMP collection system. No longer used.

TD-BALANCE-CALIBRATION

Typ Nr U 102

View of balance calibration data. Although this file is updated daily, no data in it is retrieved by any TDMS software.

TD-ANIMAL

Typ Nr U 103

View embracing animal allocation, animal transfer, animal removal, histology number, "not examined", and uncertain cause of death segments. These segments are unique observations of animals.

TD-BALANCE

Typ Nr U 104

View Of Balance Information. Some reports access this view.

TD-CAGE

Typ Nr U 105

View embracing the cage segment and the defunct "cage aborted" segment.

TD-ANIMAL-DATA

Typ Nr U 106

View embracing animal weight, animal observation, and in-life animal note segments.

TD-CAGE-DATA

Typ Nr U 107

View embracing feeder weight, bottle weight, cage note, feed consumption, water consumption, cage observation, and misidentified animals. (The last four are inactive.)

TD-ORGAN-DATA

Typ Nr U 108

View of microscopic observation, organ note, and organ stain segments, the last of which is inactive.

TD-CORRECTION-REASON

Typ Nr U 109

File contains MODCOMP corrections. In-life information is reported in E13; pathology information is not reported. Since LDAS' introduction, no new records are added to this file.

TD-MICRO-ANIMAL-NOTE

Typ Nr U 110

File contains pathology animal notes. Only used in P14.

TD-ORGAN

Typ Nr Typ Nr U 111

View of organ (tissue) status and accountable site status segments. (The last only applies to data not converted for LDAS.)

TD-TEST

Typ Nr U 112

View of test information.

TD-TRANSACTION

Typ Nr U 113

Formerly used by MODCOMP ECS and UPDATE, but no longer used.

TD-TEST-NOTE

Typ Nr U 114

View of test notes. Used by E10 only.

TD-TREATMENT-NOTE

Typ Nr U 114

View of treatment notes. Used by E10 only.

TREATMENT-PEIS

Typ Nr A 100

Physical layout of file 100. Supports TD-TREATMENT view. No current TDMS software accesses this file using direct calls.

OTHER-TERMS-PEIS

Typ Nr A 101

Supports TD-OTHER-TERMS view, which is no longer used.

BALANCE-CALIBRATION-EIS

Typ Nr A 102

Physical layout of file 102. Supports the TD-BALANCE-CALIBRATION view. No current TDMS software accesses this file using direct calls.

ANIMAL-EIS

Typ Nr A 103

Physical layout of file 103. Supports the TD-ANIMAL view. No current TDMS software accesses this file using direct calls.

BALANCE

Typ Nr A 104

Physical layout of file 104. Supports the TD-BALANCE view. No current TDMS software accesses this file using direct calls.

CAGE-EIS

Typ Nr A 105

Physical layout of file 105. Supports the TD-CAGE view. No current TDMS software accesses this file using direct calls.

ANIMAL-DATA-EIS

Typ Nr A 106

Physical layout of file 106. Supports the TD-ANIMAL-DATA view. No TDMS software accesses this file using direct calls.

CAGE-DATA-EIS

Typ Nr A 107

Physical layout for file 107. Supports the TD-CAGE-DATA view. No current TDMS software accesses this file using direct calls.

ORGAN-DATA-PEIS

Typ Nr A 108

Physical layout of file 108. Supports the TD-ORGAN-DATA view. No current TDMS software accesses this file using direct calls.

MODCOMP-CORRECTION-DATA

Typ Nr A 109

Physical layout of file 109. Supports the TD-CORRECTION-REASON view. No TDMS software accesses this file using direct calls.

ANIMAL-NOTE-PEIS

Typ Nr A 110

Physical layout of file 110. Supports TDMICROANIMALNOTE view. No current TDMS software accesses this file using direct calls.

ORGAN-PEIS

Typ Nr A 111

Physical layout of file 111. Supports the TD-ORGAN view. No TDMS software accesses this file using direct calls.

TEST

Typ Nr A 112

Physical layout for file 112. Supports the TD-TEST view. No TDMS software accesses this file using direct calls.

MODCOMP-TRANSACTION

Typ Nr A 113

Supports TD-TRANSACTION view, which is no longer used.

TD-NOTES

Typ Nr A 114

Physical layout of file 114. Supports TD-TESTNOTE and TD-TREATMENTNOTE views. No TDMS software accesses this file using direct calls.

4.0 Field Descriptions

This chapter describes each of the fields of TDMS. There are approximately 1,355 fields. For a definition of abbreviations used in this chapter, refer to the legend in *Appendix A-2*.

4.1 Protocol Analysis System Views

PAS_ACCT_SITE

View of pathology accountable-site protocol. Accountable sites pertain only to tests whose pathology data was collected exclusively on the old, MODCOMP-based, collection system. They are sites on tissue that the collectors were obligated to evaluate as if they were organs. (LDAS elevates the tissue-sites to the level of tissue.) This view is not used by the current version of the download component.

FILE...: PAS_ACCT_SITE

TYPE...: USER VIEW

FILE-NR: 55

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER.				
1	TE		DATA-DATE	A	8.0		N
			The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF		DATA-TIME	A	6.0		N
			The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG		DATA-FACILITY-NUMBER	N	5.0		N
			The facility with which the entry operator is associated. Foreign key to TBFACILITY.				
1	TH		OPERATOR-ID	N	5.0		N
			The entry or modification operator. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				

FILE...: PAS_ACCT_SITE

TYPE...: USER VIEW

FILE-NR: 55

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TL	ORG-OR-COR	Flag relating to a latent PEX/ DOWNLOAD feature. Always set to "0".	A	1.0	F	
1	TM	REASON	Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	78.0	N	
1	TN	COR-DATE	Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	8.0	N	
1	TO	COR-TIME	Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	6.0	N	
1	BB	ORGAN-CODE	The organ to which the ac- countable sites pertain. For- eign key to TBPCT.	N	5.0		
1	CC	SITE-GRP-CT	The number of entries in the multiple-entry field SITE- CODE. Should be equal to C*SITE-CODE.	N	4.0	N	
M 1	DD	SITE-CODE	The tissue-sites being desig- nated as accountable. Foreign key to TBPCT.	N	5.0	N	
T	L	DB	NAME	F	LENG	S	DE
-	-	-	- - - -				
1	SD	SUPER-KEY	SOURCE FIELD(S) --- -START- --END- DATA-BASE-KEY ORGAN-CODE	A	11.0		SP
					1	9	
					1	2	
			Used to perform discrete sel- ections based on DATA-BASE- AGENCY and ORGAN-CODE. To build the search value, ORGAN- CODE must be converted to a 2- byte binary, and the appended				

FILE...: PAS_ACCT_SITE

TYPE...: USER VIEW

FILE-NR: 55

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME

F LENG S DE

to DATA-BASE-KEY. Primary
key for the table.

PAS_ANIMAL_ID

View of in-life, animal identification protocol. File is not used since implementation of LDAS. There is still the option to maintain the file in PEX, and validation reports process it, but, otherwise, the system does not reference it. However, the value of CID-FLAG does indicate which of the two methods of constructing TDA-CARCASS-IDENTIFICATION-DE values in TD-ANIMAL was used on the test.

FILE...: PAS_ANIMAL_ID

TYPE...: USER VIEW

FILE-NR: 49

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MASTER.	A	9.0		DE
1	TE	DATA-DATE	The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)	A	8.0		N
1	TF	DATA-TIME	The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)	A	6.0		N
1	TG	DATA-FACILITY-NUMBER	The facility with which the entry operator is associated. Foreign key to TBFACILITY.	N	5.0		N
1	TH	OPERATOR-ID	The entry or modification operator. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.	N	5.0		N
1	TL	ORG-OR-COR	Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".	A	1.0		F

FILE...: PAS_ANIMAL_ID

TYPE...: USER VIEW

FILE-NR: 49

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
1 TM REASON Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	78.0 N
1 TN COR-DATE Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	8.0 N
1 TO COR-TIME Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	6.0 N
1 BB ID-FLAG Flag indicating the method of identifying an animal for in- life data collection. Always "A", for animal number.	A	1.0 F
1 CC CID-FLAG Flag indicating the method of identifying an animal for pathology data collection, (carcass ID). C indicates that CIDs were assigned based on the 4-digit cage number, and a 1-digit sequential iden- tifier within the cage; E in- dicates that the animal num- ber was used. C was not used after the fall of 1986.	A	1.0 F
T L DB NAME	F	LENG S DE
- - - - -		
1 DD ID-CT Number of entries used in the ANY-ID-INFO. Same value as C*ANY-ID-INFO. Since ANY-ID- INFO was not used in TDMS, this value is always 0.	N	2.0 N
P 1 EE ANY-ID-INFO Repeating group contain animal identification information. Never used in TDMS.		

FILE...: PAS_ANIMAL_ID

TYPE...: USER VIEW

FILE-NR: 49

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

2 FF ANIMAL-ID	A	1.0 F
----------------	---	-------

ID indicating the alphabetic
characters used to different-
iate among animals in a cage.
Never used in TDMS.

2 GG DESCRIPTION	A	16.0 N
------------------	---	--------

Description of each ID. Never
used in TDMS.

PAS_CALI_WGT

View of in-life, calibration protocol. It specifies the rules for calibrations of balances used for weight data collection. It is processed by the current version of download. It is always entered for tests on which weights are to be collected.

FILE...: PAS_CALI_WGT

TYPE...: USER VIEW

FILE-NR: 60

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number.	A	9.0		DE
1	TE	DATA-DATE	The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)	A	8.0		N
1	TF	DATA-TIME	The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)	A	6.0		N
1	TG	DATA-FACILITY-NUMBER	The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS.	N	5.0		N
1	TH	OPERATOR-ID	ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.)	N	5.0		N
1	TL	ORG-OR-COR	Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".	A	1.0		F

FILE...: PAS_CALI_WGT

TYPE...: USER VIEW

FILE-NR: 60

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
1 TM REASON Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	78.0 N
1 TN COR-DATE Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	8.0 N
1 TO COR-TIME Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	6.0 N
1 BB NO-OF-WEIGHINGS Number of weighings that may be taken in a data collection session before recalibration is required. (Calibration is always required at the be- ginning of a data collection session.) Loaded by LDAS into the CALI_NUM field on the TEST table.	N	4.0 N
1 CC CALIBRATION-CT Number of entries in the CALI- BRATION-WT repeating group. Should be equal to C*CALIBRA- TION-WT.	N	4.0 N
P 1 DD CALIBRATION-WT Repeating group that describes each weight to be used for balance calibration.		
T L DB NAME	F	LENG S DE
2 EE CAL-WT-DES Free text describing calibra- weight, and used by LDAS as a prompt. Generally, is the weight expressed in grams. Loaded by LDAS into the WEIGHT_TEXT field on the CALI_WGT table.	A	32.0 N

FILE...: PAS_CALI_WGT

TYPE...: USER VIEW

FILE-NR: 60

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

2 FF MAX-WT	N	8.0 N
-------------	---	-------

Weight, in decigrams, that is the maximum acceptable value for the associated weight. Loaded by LDAS into the MAX_-WEIGHT field on the CALI_WGT table.

2 GG MIN-WT	N	8.0 N
-------------	---	-------

Weight, in decigrams, that is the minimum acceptable value for the associated weight. Loaded by LDAS into the MIN_-WEIGHT field on the CALI_WGT table.

PAS_CENSUS

View of in-life, cage census protocol. Used to describe the cages on the test, including weight, treatment, and procedure groups, and number of animals to be allocated. Used by DOWNLOAD to construct what will be the LDAS CAGE table.

FILE...: PAS_CENSUS

TYPE...: USER VIEW

FILE-NR: 44

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER.	A	9.0		DE
1	TE	DATA-DATE	The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)	A	8.0		N
1	TF	DATA-TIME	The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)	A	6.0		N
1	TG	DATA-FACILITY-NUMBER	The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS. Foreign key to TBFACLT.	N	5.0		F
1	TH	OPERATOR-ID	ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.) When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.	N	5.0		F

FILE...: PAS_CENSUS

TYPE...: USER VIEW

FILE-NR: 44

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
1 TL ORG-OR-COR Flag relating to a latent PEX/ DOWNLOAD feature. Always set to "O".	A	1.0 F
1 TM REASON Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	78.0 N
1 TN COR-DATE Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	8.0 N
1 TO COR-TIME Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	6.0 N
1 BB STARTING-CAGE The starting cage number for the range of cages being description- cribed in the record.	N	4.0 F
1 CC ENDING-CAGE The ending cage number for the range of cages being description- cribed in the record.	N	4.0 F
1 DD TREATMENT-ID The ID of the treatment group to which the cage range be- longs. Loaded by LDAS into the TREAT_NUM field on the CAGE table. When appended to DATA-BASE-KEY, a foreign key to PAS_TREATMNT_SEG.	A	3.0 N
1 EE PROCEDURE-ID The procedure group to which the cages in the range belong. Loaded by LDAS into the PROC_- NUM field in the CAGE table. When appended to DATA-BASE-KEY a foreign key to PAS_PROCEDUR- _SEG.	A	3.0 N

FILE...: PAS_CENSUS

TYPE...: USER VIEW

FILE-NR: 44

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	FF	WGT-ITEM	N	3.0	N	
<p>The weight group to which the cages in the range belong.</p> <p>Loaded by LDAS into the WEIGHT_GROUP field in the TREATMENT table. (MODCOMP supported multiple weight groups per treatment group;</p>						
1	FF	WGT-ITEM	N	3.0	N	
<p>LDAS does not.) When appended to DATA-BASE-KEY, a foreign key to PAS_WGT_ITEMS.</p>						
1	GG	RELATIVE-START	N	4.0	F	
<p>The date, relative to the cage start date, on which dosing in the cage is initiated. Overridden in TDMS by the first value of REL-START in the related PAS_TREATMNT_SEG record.</p>						
P 1	HH	ANY-ANIMAL-IDS				
<p>Defined as a repeating group, but TDMS, by procedure, has only one entry per record.</p>						
2	II	SSSC-ID	N	3.0	F	
<p>ID of the Strain/Substrain/Sex/Class group to which the cages in the range belong.</p> <p>When appended to DATA-BASE-KEY, a foreign key to PAS_LOAD_WT.</p>						
2	JJ	NO-ANIMALS	N	2.0	N	
<p>The number of animals that are to be allocated to cages within the cage ranges. Loaded by LDAS into the NUM_ANIMALS field on the CAGE table.</p>						
1	SD	SUPER-KEY	A	11.0	SP	
<p>SOURCE FIELD(S) --- -START- --END-</p>						
		DATA-BASE-KEY		1	9	
		STARTING-CAGE		1	2	

FILE...: PAS_CENSUS

TYPE...: USER VIEW

FILE-NR: 44

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

Used in PEX for indexed reads based on test number and starting cage. To build a search value, STARTING-CAGE must be converted to a 2-byte binary, and then appended to DATA-BASE-KEY.

1	SL	AGENCY-EXP-TEST-ENDING-CAGE	A	11.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY		1		9
		ENDING-CAGE		1		2

Used by RECEIVE for unique record selection.

PAS_CLINICAL_OBS

View of in-life term protocol. Includes the clinical signs, sites, qualifiers, and removal reasons that may be used on the test. (For older MODCOMP tests, it will also include the cage conditons that were used for cage observations.) Processed by Download to generate all in-life term tables.

FILE...: PAS_CLINICAL_OBS

TYPE...: USER VIEW

FILE-NR: 48

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	A	9.0	DE	
		Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER.				
1	TE	DATA-DATE	A	8.0	N	
		The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF	DATA-TIME	A	6.0	N	
		The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG	DATA-FACILITY-NUMBER	N	5.0	N	
		Facility with which the entry operator is associated. Foreign key to TBFACLT.				
1	TH	OPERATOR-ID	N	5.0	N	
		Entry or modification operator. When appended to DATA-FACILITY-NUMBER, foreign key to TB_OPERATOR_ID.				
1	TL	ORG-OR-COR	A	1.0	F	
		Flag relating to a latent PEX/DOWNLOAD feature. Always set to "0".				
1	TM	REASON	A	78.0	N	
		Text relating to a latent PEX/				

FILE...: PAS_CLINICAL_OBS

TYPE...: USER VIEW

FILE-NR: 48

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
			DOWNLOAD feature. No values assigned.				
1	TN		COR-DATE	A	8.0	N	
			Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO		COR-TIME	A	6.0	N	
			Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB		PAS-OBS-CODE-CT	N	3.0		
			Number of entries in the ANY-OBSERVATION repeating group. Should be equal to C*ANY-OBSERVATION.				
1	CC		PAS-SELECT-SITE-CT	N	3.0	N	
			Number of entries in the multiple-entry field, SITE-CODE. Should be equal to C*SITE-CODE.				
1	DD		PAS-SELECT-SIZE-CT	N	3.0	N	
			Number of entries in the multiple-entry field, SIZE-CODE. Should be equal to C*SIZE-CODE.				
1	EE		PAS-REM-REASON-CT	N	3.0	N	
			Number of entries in the multiple-entry field, REMOVAL-CODE. Should be equal to C*REMOVAL-CODE.				
1	FF		PAS-CAGE-COND-CT	N	3.0	N	
			Number of entries in the multiple-entry field, CAGE-COND-CODE. Should be equal to C*CAGE-COND-CODE.				
P 1	GG		ANY-OBSERVATION				
			Repeating group for in-life observation codes.				
2	HH		OBSERVATION-CODE	N	5.0	N	
			Valid clinical signs for animal and removal observations. Loaded into the CLNOBS_NUM				

FILE...: PAS_CLINICAL_OBS

TYPE...: USER VIEW

FILE-NR: 48

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

field in the LDAS tables CLNOBS, CLNOBS_QUAL, and CLNOBS-SITE, and into TERM_NUM in INLIFE_TERMS. Also used to build the CLNOBS_GROUP and INLIFE_QUAL tables. Foreign key to TBCLOBSC.

2 II SITE-REQ	A	1.0
---------------	---	-----

Flag indicating whether the corresponding observation requires that a site be associated with it. Possible values are R for required, and O for optional, but only R is currently used. The only current clinical signs for which this flag is set are 20 for mass, and 22 for ulcer. DOWNLOAD uses this flag to build what LDAS will load into CLNOBS_SITE table.

2 JJ SIZE-REQ	A	1.0
---------------	---	-----

Flag indicating whether the corresponding observation requires that a size be associated with it. Possible values are R for required, and O for optional, but only R is currently used. The only current clinical sign for which this flag is set is 20 for mass. However, it is not used by DOWNLOAD.

M 1 LL SITE-CODE	N	5.0
------------------	---	-----

The valid in-life sites which may be used for animal and removal observations. Loaded by LDAS into SITE_NUM in the INLIFE_SITE table, and for each clinical site for which the SITE-REQ flag is set, into SITE_NUM in the CLNOBS-

FILE...: PAS_CLINICAL_OBS

TYPE...: USER VIEW

FILE-NR: 48

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
			SITE table. Foreign key to TBCLOBSC.				
M	1	MM	SIZE-CODE	N	5.0		
			The valid size codes that may be used for animal and removal observations. Loaded by LDAS into the QUAL_NUM field of the INLIFE_QUAL table. Foreign key to TBCLOBSC.				
M	1	NN	REMOVAL-CODE	N	5.0		
			Valid removal reasons that may be used to define an animal removal. Loaded by LDAS into the REASON_CODE field of the REMOVL_REASON table. Foreign key to TBCLOBSC.				
M	1	OO	CAGE-COND-CODE	N	5.0		
			Valid cage condition codes. Only used for older MODCOMP tests. Foreign key to TBCONDTC.				

PAS_FAC_CHEM

Basic information about the test, including the test compound and facility. Used throughout the system. Equivalent to the TEST table on LDAS.

FILE...: PAS_FAC_CHEM

TYPE...: USER VIEW

FILE-NR: 59

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number.				
1	TE		DATA-DATE	A	8.0		N
			The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF		DATA-TIME	A	6.0		N
			The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG		DATA-FACILITY-NUMBER	N	5.0		N
			The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS.				
1	TH		OPERATOR-ID	N	5.0		N
			ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.)				
1	TL		ORG-OR-COR	A	1.0		F
			Flag relating to a latent PEX/DOWNLOAD feature. Always set to "0".				
1	TM		REASON	A	78.0		N
			Text relating to a latent PEX/DOWNLOAD feature. No values assigned.				
1	TN		COR-DATE	A	8.0		N

FILE...: PAS_FAC_CHEM

TYPE...: USER VIEW

FILE-NR: 59

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.		
1 TO COR-TIME	A	6.0 N
Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.		
1 BB TEST-TYPE	N	2.0 N DE
Code for the test type by which the test has been class- ified. Loaded by LDAS into TYPE_NUM on the TEST table.		
1 CC DOSE-RT	N	2.0 N
Code for the dose route as- signed to the test. Used by DOWNLOAD to retrieve the as- sociated dose route text.		
1 DD FACILITY-ID	N	5.0 N DE
The ID for the facility con- ducting the text. Loaded by LDAS into the FACLTY_NUM field on the TEST table.		
1 EE FACILITY-NAME	A	54.0 N
Free text describing the facility. Should correspond to the LONG-TEXT value in TBFACLTY. Loaded by LDAS into the FACILITY field in the TEST table.		
1 FF PRIN-INV	A	32.0 N
The name of the principal in- vestigator entered as free text. Not used by DOWNLOAD.		
1 GG CONTRACT-1	A	12.0 N
The number of the contract under which the test is fund- ed. Loaded by LDAS into the CONTRACT_NUM field on the TEST table.		

FILE...: PAS_FAC_CHEM

TYPE...: USER VIEW

FILE-NR: 59

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
1 HH CONTRACT-2 Secondary contract number. Not used in TDMS.	A	12.0 N
M 1 II HAZARD-TXT Free text narrative concerning safety issues relative to ex- posure to the test chemicals. Up to 400 characters are load- ed by LDAS into the HAZARD_- TEXT field on the TEST table.	A	64.0 N
1 JJ TEST-COMPOUND-CT The number of compounds description- cribed in the ANY-TEST-COM- POUND repeating group. Should be equal to C*ANY-TEST-COM- POUND.	N	3.0 N
P 1 KK ANY-TEST-COMPOUND Repeating group containing descriptions of each of the chemicals being tested.		
2 LL NTP The NTP chemical number. Also referred to as C number. The first one in the repeating group appears on the header or splash page for all in-life reports, and some pathology reports.	A	6.0 N
2 MM CAS The Chemical Abstract Services (CAS) number for the test com- pound. The first one in the repeating group appears on the header or splash page for all in-life and pathology reports.	A	12.0 N
2 NN COMPOUND Name of compound entered as free text. First one in the repeating group appears on in- life and pathology report headers and splash pages.	A	54.0 N DE

FILE...: PAS_FAC_CHEM

TYPE...: USER VIEW

FILE-NR: 59

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	OO	LOCK-DATE	P	5.0	N	
The date when the test is submitted to NTP. Stored in TDMS TOX format. Loaded by LDAS into the LOCK_DATE field in the TEST table.						
1	PP	LOCK-ENTRY-DATE	P	5.0	N	
The system date when the lock date is set, in TDMS TOX format. Processed by DOWNLOAD, but not used by LDAS.						
1	QQ	LOCK-OPER-ID	P	5.0	N	
The ID of the NTP operator who set the lock date.						
1	RR	ANML-SUPPLR	P	5.0	N	
The code for the animal supplier for the test. Used by DOWNLOAD to extract the corresponding supplier text.						
1	SS	PROTO-REQ-TIS-FLAG	A	1.0	N	
Flag indicating whether protocol-required tissues may be defined on the local machine. L-may be defined on a local machine; M-must be defined on VAX. Loaded by LDAS into the PRT_FLAG field in the TEST table.						

PAS_LOAD_WT

View of in-life protocol which describes the characteristics of the animals used on the test. Used throughout the system.

FILE...: PAS_LOAD_WT

TYPE...: USER VIEW

FILE-NR: 47

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER	A	9.0		DE
1	TE	DATA-DATE	The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)	A	8.0		N
1	TF	DATA-TIME	The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)	A	6.0		N
1	TG	DATA-FACILITY-NUMBER	The facility with which the entry operator is associated. Foreign key to TBFACLT.	N	5.0		N
1	TH	OPERATOR-ID	Entry or modification operator. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.	N	5.0		N
1	TL	ORG-OR-COR	Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".	A	1.0		F
1	TM	REASON	Text relating to a latent PEX/DOWNLOAD feature. No values assigned.	A	78.0		N
1	TN	COR-DATE	Date relating to a latent PEX/	A	8.0		N

FILE...: PAS_LOAD_WT

TYPE...: USER VIEW

FILE-NR: 47

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
DOWNLOAD feature. No values assigned.		
1 TO COR-TIME Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	6.0 N
1 BB SSSC-ID Sequentially assigned ID which serves to distinguish between records for the same test. Since there is never more than 2 records per test, this value is either 1 or 2.	N	3.0 N
1 CC SPECIES-CODE The species of the animal profile. Foreign key to TBSTRAIN. Loaded by LDAS into the ANML_SPECIES field in the TEST table. Also used by DOWNLOAD to determine species-dependant tissues and lesions.	N	5.0 N
1 DD STRAIN-CODE The strain of the animal profile. Foreign key to TBSTRAIN. Loaded by LDAS into the STRAIN_NUM field in the TREATMENT table.	N	5.0 N
1 EE SUBSTRAIN-CODE Sub-strain of the animal profile. (Sub-strains are a legacy to TDMS, with no current use. New tests are assigned "Not Applicable.") Foreign key to TBSUBSTR.	N	5.0 N
1 FF CLASS-CODE Animal class qualifier associated with the animal profile. (Animal class qualifiers are a legacy to TDMS, with no current use. New tests are assigned "Not Applicable".) Foreign key to TBANCLQF.	N	5.0 N

1 GG	SEX-OF-ANIMAL	A	1.0 F
	<p>The sex of the animals in the SSSC group. M for Male; F for Female. Loaded by LDAS into the TREAT_SEX field of the TREATMENT table. (This is the only field that should have a different values for records of the same test.)</p>		
1 HH	NO-ANIMALS	N	4.0 N
	<p>The number of animals in the SSSC group. Essentially, the number of animals of the sex identified in SEX-OF-ANIMAL.</p>		
1 II	DESIRED-AGE	N	3.0 N
	<p>Desired age of animals at test start. Current liaison procedures dictate no value for this field, although it is loaded by LDAS into the TREAT_AGE field of the TREATMENT table.</p>		
1 JJ	HI-ALLOCATION-WEIGHT	N	6.0 N
	<p>Maximum weight, in decigrams, of an animal at allocation. Pertains only to MODCOMP tests. Current liaison procedures indicate that no value should be entered for new tests.</p>		
1 KK	LO-ALLOCATION-WEIGHT	N	6.0 N
	<p>Minimum weight, in decigrams, of an animal at allocation. Pertains only to MODCOMP tests. Current liaison procedures indicate that no value should be entered for new tests.</p>		
1 LL	WT-DEVIATION	N	2.0 N
	<p>The maximum acceptable difference between two consecutive weights of the same animal, expressed as a percentage. Pertains only to MODCOMP tests. Function replaced in LDAS by body weight growth curve. Current liaison procedures indicate that no value should be entered.</p>		

FILE...: PAS_LOAD_WT

TYPE...: USER VIEW

FILE-NR: 47

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

1	MM	MAX-TEST-WT	N	6.0	N	
---	----	-------------	---	-----	---	--

Maximum weight of animal on test, in decigrams. Pertains only to MODCOMP tests. Current liaison procedures indicate that no value should be entered.

1	NN	MIN-TEST-WT	N	6.0	N	
---	----	-------------	---	-----	---	--

Minimum weight of animal on test, in decigrams. Pertains only to MODCOMP tests. Current liaison procedures indicate that no value should be entered.

1	SD	SUPER-KEY	A	11.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-

DATA-BASE-KEY	1	9
SSSC-ID	1	2

Used for single-record select and indexed reads in PEX, DOWNLOAD, and some reports. Primary key for table.

PAS_MASTER

PAS master table. Indicates whether a test has been registered. Used by report menus to validate a selected test.

FILE...: PAS_MASTER

TYPE...: USER VIEW

FILE-NR: 40

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
Composite field whose elements are: the 2-digit agency number which is always set to "05", and the 7-digit experiment-test number. This field is assigned when an experiment-test is registered in PEX and used throughout the menu system to validate experiment-test values entered on the screens.							
P 1	TD		DOWNLOAD-INFO				
Periodic group relating to a latent PEX/DOWNLOAD feature.							
2	TE		DOWNLOAD-DATE	A	8.0		N
Date relating to a latent PEX/DOWNLOAD feature. No values assigned.							
2	TF		DOWNLOAD-TIME	A	6.0		N
Time relating to a latent PEX/DOWNLOAD feature. No values assigned.							

PAS_MCRO_PTH_TBL

View of pathology terms, including qualifiers, stain codes, and cause-of-death and morbidity codes. Only qualifiers are relevant to the current version of TDMS. New records are generally populated by copying from one of the universal download tests. It is generally not modified after creation.

FILE...: PAS_MCRO_PTH_TBL

TYPE...: USER VIEW

FILE-NR: 58

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number.	A	9.0		DE
1	TE	DATA-DATE	The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)	A	8.0		N
1	TF	DATA-TIME	The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)	A	6.0		N
1	TG	DATA-FACILITY-NUMBER	The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS.	N	5.0		N
1	TH	OPERATOR-ID	ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.)	N	5.0		N
1	TL	ORG-OR-COR	Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".	A	1.0		F

FILE...: PAS_MCRO_PTH_TBL

TYPE...: USER VIEW

FILE-NR: 58

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	TM	REASON	A	78.0	N	
		Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TN	COR-DATE	A	8.0	N	
		Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO	COR-TIME	A	6.0	N	
		Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB	P-STAIN-CT	N	4.0	N	
		The number of entries in the multiple-entry field STAIN- CODE. Should be equal to C*STAIN-CODE.				
1	CC	P-COD-COM-CT	N	4.0	N	
		The number of entries in the multiple-entry field COD-COM- CODE. Should be equal to C*COD-COM-CODE.				
1	DD	P-QUALIFIER-CT	N	4.0	N	
		The number of entries in the multiple-entry field QUALIFI- ER-CODE. Should be equal to C*QUALIFIER-CODE.				
M 1	EE	STAIN-CODE	N	5.0	N	
		Stain codes that could be used to create stain transactions. Pertains only to older MODCOMP tests.				
M 1	FF	COD-COM-CODE	N	5.0	N	
		Codes for cause-of-death/mor- bidity. No use in TDMS.				
M 1	GG	QUALIFIER-CODE	N	5.0	N	
		The qualifier codes that may be used to build histopatho- logy observations. Loaded by LDAS into the QUAL_NUM field of the PATH_QUAL table, and				

FILE...: PAS_MCRO_PTH_TBL
TYPE...: USER VIEW
FILE-NR: 58
PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

is also used by DOWNLOAD to
constuct what LDAS will load
into the P_QUAL_GROUP table.

PAS_MICRO_PROC

View pertains to pathology protocol for test whose pathology data was collected using the MODCOMP collection system. Described protocol-dependant procedures to be used in pathology evaluation.

FILE...: PAS_MICRO_PROC

TYPE...: USER VIEW

FILE-NR: 56

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	A	9.0	DE	
		Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER.				
1	TE	DATA-DATE	A	8.0	N	
		The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF	DATA-TIME	A	6.0	N	
		The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG	DATA-FACILITY-NUMBER	N	5.0	N	
		The facility which the entry operator is associated. Foreign key to TBFACILITY.				
1	TH	OPERATOR-ID	N	5.0	N	
		The entry or modification operator. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
1	TL	ORG-OR-COR	A	1.0	F	
		Flag relating to a latent PEX/DOWNLOAD feature. Always set to "0".				
1	TM	REASON	A	78.0	N	
		Text relating to a latent PEX/DOWNLOAD feature. No values assigned.				

FILE...: PAS_MICRO_PROC

TYPE...: USER VIEW

FILE-NR: 56

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	TN	COR-DATE	A	8.0	N	
		Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO	COR-TIME	A	6.0	N	
		Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB	MICRO-PROC-CT	N	4.0	N	
		The number of entries in the multiple-entry field MICRO- PROC-CHAR. Should be equal to C*MICRO-PROC-CHAR.				
M 1	CC	MICRO-PROC-CHAR	A	65.0	N	
		Free text describing proce- dures used during histopatho- logy evaluation. Pertains only to MODCOMP tests.				

PAS_MORPHOLOGIES

View of pathology protocol. Lists the morphologies that may be used for pathology data collection. This file is initially populated by copying from one of the universal download tests. It is rarely modified after initialization.

FILE...: PAS_MORPHOLOGIES

TYPE...: USER VIEW

FILE-NR: 53

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MASTER.	A	9.0		DE
1	TE	DATA-DATE	The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)	A	8.0		N
1	TF	DATA-TIME	The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)	A	6.0		N
1	TG	DATA-FACILITY-NUMBER	The facility with which the entry operator is associated. Foreign key to TBFACILITY.	N	5.0		N
1	TH	OPERATOR-ID	The entry or modification operator. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.	N	5.0		N
1	TL	ORG-OR-COR	Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".	A	1.0		F
1	TM	REASON	Text relating to a latent PEX/	A	78.0		N

FILE...: PAS_MORPHOLOGIES

TYPE...: USER VIEW

FILE-NR: 53

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
DOWNLOAD feature. No values assigned.		
1 TN COR-DATE Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	8.0 N
1 TO COR-TIME Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	6.0 N
1 BB ORGAN-CODE The organ to which the morpho- logies apply. Foreign key to TBPCT. Loaded by LDAS into the TISSUE_NUM field of the TISSUE_MORPH table.	N	5.0
1 CC NON-NEO-MORPH-CT The number of entries in the repeating group ANY-NON-NEO- PLASTIC. Should be equal to C*ANY-NON-NEOPLASTIC.	N	4.0 N
1 DD NEO-MORPH-CT The number of entries in the repeating group ANY-NEOPLAS- TIC-GROUP. Should be equal to C*ANY-NEOPLASTIC-GROUP.	N	4.0 N
P 1 EE ANY-NON-NEOPLASTIC Group embracing non-neoplastic morphologies and their assoc- iated qualifier flags.		
2 FF MORPHOLOGY-CODE-NON The non-neoplastic morpholo- gies that may be used to build pathology observations. Load- ed by LDAS into the MORPH_NUM fields in the MORPHOLOGY, TISS_MORPH, and MORPH_QUAL tables. Foreign key to TBPCT.	N	5.0 N
2 GG QUALIFIER-FLAG-NON Flag indicating whether the associated morphology can or must be qualified. The value	A	1.0 F

FILE...: PAS_MORPHOLOGIES

TYPE...: USER VIEW

FILE-NR: 53

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
must be moved into an integer to be evaluated since that is the way the value is stored. The possible values are: 1 - Must be qualified 32 - May be qualified 78 - May not be qualified 81 - May be qualified Loaded by LDAS into the FORCED_GROUP field in the MORPHOLOGY table.		
P 1 HH ANY-NEOPLASTIC-GRP Repeating group embracing neoplastic morphology codes and their associated qualifier flags.		
2 II MORPHOLOGY-CODE The neoplastic morphologies that may be used to build pathology observations. Load- ed by LDAS into the MORPH_NUM fields in the MORPHOLOGY, TISS_MORPH, and MORPH_QUAL tables. Foreign key to TBPCT.	N	5.0 N
2 JJ QUALIFIER-FLAG Flag indicating whether the associated morphology can or must be qualified. The value must be moved into an integer to be evaluated since that is the way the value is stored. The possible values are: 1 - Must be qualified 32 - May be qualified 78 - May not be qualified 81 - May be qualified Loaded by LDAS into the FORCED_GROUP field in the MORPHOLOGY table.	A	1.0 F
1 SD SUPER-KEY SOURCE FIELD(S) --- -START- --END- DATA-BASE-KEY	A	11.0 SP 1 9

FILE...: PAS_MORPHOLOGIES

TYPE...: USER VIEW

FILE-NR: 53

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
ORGAN-CODE	1	2
Used in PEX for indexed reads based on test and organ val- ues.		

PAS_ORGANS

View of pathology protocol. Lists the organs that may be used for pathology data collection. This file is initially populated by copying from one of the universal download tests. It is rarely modified after initialization.

FILE...: PAS_ORGANS

TYPE...: USER VIEW

FILE-NR: 52

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MASTER.	A	9.0		DE
1	TE	DATA-DATE	The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)	A	8.0		N
1	TF	DATA-TIME	The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)	A	6.0		N
1	TG	DATA-FACILITY-NUMBER	The facility with which the entry operator is associated. Foreign key to TBFACILITY.	N	5.0		N
1	TH	OPERATOR-ID	Entry or modification operator. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.	N	5.0		N
1	TL	ORG-OR-COR	Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".	A	1.0		F
1	TM	REASON	Text relating to a latent PEX/	A	78.0		N

FILE...: PAS_ORGANS

TYPE...: USER VIEW

FILE-NR: 52

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
		DOWNLOAD feature. No values assigned.				
1	TN	COR-DATE	A	8.0	N	
		Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO	COR-TIME	A	6.0	N	
		Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB	SSSC-ID-CT	N	4.0	N	
		The number of entries in the multiple-field SSSC-ID. Always 1.				
1	CC	ORGAN-GROUP-CT	N	4.0	N	
		The number of the entries in the multiple field ORGAN-CODE. Should be equal to C*ORGAN-CODE.				
M 1	DD	SSSC-ID	N	3.0		
		The SSSC ID to which the list of organ codes pertains. Although a multiple-entry field, always only one entry. When appended to DATA-BASE-KEY, a foreign key to PAS_LOAD_WT.				
M 1	EE	ORGAN-CODE	N	5.0		
		The PCT codes for the organs that may be used to build pathology observations loaded by LDAS into the TISSUE_NUM fields in TISSUE, TISS_SITE, and TISS_MORPH tables. Foreign key to TBPCT.				
1	FF	STUDY-ORGAN	A	11.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY	1		9	
		ORGAN-CODE	1		2	
		Used for indexed finds based on test number and organ code values.				

PAS_PH_P_AST_OPR

View of operators assigned to a test for pathology data collection. For tests conducted using LDAS, this information is assigned on the local machine and uploaded to the mainframe DBMS; for MODCOMP tests, the information was entered in PEX and downloaded.

FILE...: PAS_PH_P_AST_OPR

TYPE...: USER VIEW

FILE-NR: 57

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MASTER.				
1	TE		DATA-DATE	A	8.0		N
			The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF		DATA-TIME	A	6.0		N
			The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG		DATA-FACILITY-NUMBER	N	5.0		N
			The facility with which the entry operator is associated. A foreign key to TBFACILITY.				
1	TH		OPERATOR-ID	N	5.0		N
			The entry or modification operator. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
1	TL		ORG-OR-COR	A	1.0		F
			Flag relating to a latent PEX/DOWNLOAD feature. Always set to "0".				

FILE...: PAS_PH_P_AST_OPR

TYPE...: USER VIEW

FILE-NR: 57

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
1 TM REASON	A	78.0 N
Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.		
1 TN COR-DATE	A	8.0 N
Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.		
1 TO COR-TIME	A	6.0 N
Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.		
1 BB PERSONNEL-ID-CT	N	4.0 N
Number of entries in the re- peating group ANY-ID-GRP. Should be equal to C*ANY-ID- GRP.		
P 1 CC ANY-ID-GRP		
Repeating group embracing description- criptions of each operator assigned to a test for patho- logy data collection.		
2 DD NAME	A	32.0 N
The name of the operator, ex- tracted from the operator tab- le at the time of assignment by LDAS. (For MODCOMP tests, name was entered as free text on PEX.)		
2 EE PERSONNEL-CODE	A	1.0 F
Two different sets of codes, depending on whether code was assigned for MODCOMP or LDAS. Both code systems may exist for the same record, if both MODCOMP and LDAS were used for the pathology portion of the test. For MODCOMP tests the codes are: A-Pathologist,B-		

FILE...: PAS_PH_P_AST_OPR

TYPE...: USER VIEW

FILE-NR: 57

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

Path. Assistant, C-Entry Operator, D-Audit Path.; for LDAS tests the codes are: P-Pathologist, T-Technician, D-Data Administrator, S-Security Administrator.

2	FF	PATH-CODE	A	12.0	N	
---	----	-----------	---	------	---	--

Operator assigned to the test. Constructed in one of two ways, depending on the subsystem where it was assigned. For MODCOMP, the first 2 digits are 09, followed by the 5-digit facility number, the 4-digit operator number, and a 1-digit check byte; for LDAS, the first 8 digits are 0, followed by the 4-digit operator number, which must be appended to FACILITY-ID from PAS_FAC_-CHEM. Foreign key to TB_OPERATOR_ID.

1	GG	SUPER-KEY	A	22.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-

DATA-BASE-KEY	1	9
PERSONNEL-CODE	1	1
PATH-CODE	1	12

Used by RECEIVE for discrete record selection.

PAS_PREPARTN_SEG

View of in-life, compound preparation protocol. The preparations define the composition of the dose mixtures used on the test.

Information in the table is downloaded to LDAS and is used by many of the EIS/PEIS reports.

FILE...: PAS_PREPARTN_SEG

TYPE...: USER VIEW

FILE-NR: 43

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MASTER.	A	9.0		DE
1	TE		DATA-DATE The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)	A	8.0		N
1	TF		DATA-TIME The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)	A	6.0		N
1	TG		DATA-FACILITY-NUMBER The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS. Foreign key to TBFACILITY.	N	5.0		N
1	TH		OPERATOR-ID ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.) Appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.	N	5.0		N

FILE...: PAS_PREPARTN_SEG

TYPE...: USER VIEW

FILE-NR: 43

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	TL	ORG-OR-COR	A	1.0	F	
		Flag relating to a latent PEX/ DOWNLOAD feature. Always set to "O".				
1	TM	REASON	A	78.0	N	
		Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TN	COR-DATE	A	8.0	N	
		Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO	COR-TIME	A	6.0	N	
		Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB	ID-NO	N	3.0	N	
		Unique, sequentially-assigned, identification number for a preparation. Loaded by LDAS into the PREP_NUM field in the PREPARATION and TREAT_PREP tables.				
1	CC	COMPOUND-CT	N	1.0	N	
		The number of entries in the repeating group ANY-COMPOUND- INFO. Should be equal to C*ANY-COMPOUND-INFO.				
P 1	DD	ANY-COMPOUND-INFO				
		Repeating group defining each of the constituent elements of a dose mixture.				
2	EE	NTP-NO	A	6.0	N	
		The NTP number associated with the chemical. Not entered for vehicles.				
2	FF	COMPOUND-NAME	A	54.0	N	
		The name of the compound en-				

FILE...: PAS_PREPARTN_SEG

TYPE...: USER VIEW

FILE-NR: 43

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

tered as free text. Loaded by
LDAS into the COMP_TEXT field
of the PREPARATION table.

2	GG	COMP-TYPE	A	1.0	F	
---	----	-----------	---	-----	---	--

Flag indicating the type of
compound being defined. The
possible flag values, and
their meanings, are:

- A - Vehicle
- B - Test
- C - Initiator
- D - Promotor

The associated text, not the
flag itself, is loaded by
LDAS into the ROLE_TEXT field
in the PREPARATION table.

2	HH	COMP-PERCENT	N	3.6	N	
---	----	--------------	---	-----	---	--

The percentage of the compound
relative to the whole prepara-
. Loaded by LDAS into the
COMP_PERCENT of the PREPARA-
TION table. Also used in many
EIS/PEIS reports.

2	II	MEAS-TYPE	A	2.0	N	
---	----	-----------	---	-----	---	--

Flag indicating the basis for
the percentage calculation.
The possible values, and their
meanings, are:

- WW - Weight/Weight
- VV - Volume/Volume
- WV - Weight/Volume
- VW - Volume/Weight

The flag itself is loaded by
LDAS into the COMP_MEAS field
in the PREPARATION table.

1	SD	SUPER-KEY	A	11.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-

DATA-BASE-KEY	1	9
---------------	---	---

ID-NO	1	2
-------	---	---

Used for indexed reads to the

FILE...: PAS_PREPARTN_SEG
TYPE...: USER VIEW
FILE-NR: 43
PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

file in test and preparation
number order. Primary key for
the table.

PAS_PROCEDUR_SEG

View of in-life, procedure-action protocol. The procedures are organized into groups of actions. These groups comprise the schedule of activities that apply to cages. Its information is downloaded to LDAS.

FILE...: PAS_PROCEDUR_SEG

TYPE...: USER VIEW

FILE-NR: 42

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER.				
1	TE		DATA-DATE	A	8.0		N
			The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF		DATA-TIME	A	6.0		N
			The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG		DATA-FACILITY-NUMBER	N	5.0		N
			The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS. Foreign key to TBFACILITY.				
1	TH		OPERATOR-ID	N	5.0		N
			ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.) Appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				

FILE...: PAS_PROCEDUR_SEG

TYPE...: USER VIEW

FILE-NR: 42

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
1 TL ORG-OR-COR Flag relating to a latent PEX/ DOWNLOAD feature. Always set to "O".	A	1.0 F
1 TM REASON Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	78.0 N
1 TN COR-DATE Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	8.0 N
1 TO COR-TIME Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.	A	6.0 N
1 BB ID-NO A sequentially-assigned, un- ique number that identifies a procedure set on a test. Loaded by LDAS into the PROC_- NUM field in the PROCEDURE table.	N	3.0 N
1 CC PROCEDURE-GRP-CNT The number of actions defined in the repeating group ANY- PROCEDURE-GROUP. Should be equal to C*ANY-PROCEDURE-GROUP	N	3.0 N
P 1 DD ANY-PROCEDURE-GROUP Repeating group that defines each actions in a procedure set.		
2 EE ACTION-CODE The action being scheduled. Loaded by LDAS into the ACTION_NUM field in the PROCEDURE table. Foreign key to TBACTONC.	N	4.0 N DE

FILE...: PAS_PROCEDUR_SEG

TYPE...: USER VIEW

FILE-NR: 42

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
2	FF	REL-START	I	2.0	F	
The day-on-dose on which the schedule for the action begins. Loaded by LDAS into the START_DAY field in the PROCEDURE table.						
2	GG	REL-STOP	I	2.0	F	
The day-on-dose on which the schedule for the action ends. Loaded by LDAS into the END_DAY field in the PROCEDURE table.						
2	HH	WEEKLY-INTERVALS	N	2.0	N	
The number of weeks between occurrences of the action. Loaded by LDAS into the WEEK_INTRVL field in the PROCEDURE table.						
2	II	REL-OR-SPEC	A	1.0		
Flag indicating the basis for the DAY-OF-WEEK field. R - test week, S - calendar week. Loaded by LDAS into the REL_OR_SPEC field in the PROCEDURE table.						
2	JJ	DAY-OF-WEEK	A	7.0	N	
Each character is a flag for the day corresponding to the position. If the character has a value of 'X', then the action is scheduled for the corresponding day; if blank, it isn't scheduled. Loaded by LDAS into the DAY_OF_WEEK field in the PROCEDURE table.						
1	SD	SUPER-KEY	A	11.0	SP	
SOURCE FIELD(S) --- -START- --END-						
		DATA-BASE-KEY		1	9	
		ID-NO		1	2	

FILE...: PAS_PROCEDUR_SEG

TYPE...: USER VIEW

FILE-NR: 42

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME

F LENG S DE

Used for indexed reads by
test and procedure set number.
Primary key for table.

PAS_PROTO_REQ_TIS

View of pathology protocol. Defines the criteria by which tissues are required to be evaluated. Criteria are evaluation based on treatment group, removal reason, and removal date range. Can be created by a liaison and downloaded, or created on the local machine and uploaded. Used by the P14.

FILE...: PAS_PROTO_REQ_TIS

TYPE...: USER VIEW

FILE-NR: 84

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number.				
1	TE		DATA-DATE	A	8.0		N
			The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF		DATA-TIME	A	6.0		N
			The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG		DATA-FACILITY-NUMBER	N	5.0		N
			The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS.				
1	TH		OPERATOR-ID	N	5.0		N
			ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.)				
1	TL		ORG-OR-COR	A	1.0		N
			Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".				
1	TM		REASON	A	78.0		N
			Text relating to a latent PEX/DOWNLOAD feature. No values assigned.				

FILE...: PAS_PROTO_REQ_TIS

TYPE...: USER VIEW

FILE-NR: 84

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TN		COR-DATE	A	8.0	N	
			Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO		COR-TIME	A	6.0	N	
			Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB		TREAT-NUM	N	3.0		
			The treatment group to which the PRT set pertains. Loaded by LDAS into the TREATMENTS field of the PRO_REQ_TISS tab- le.				
1	CC		BEG-REMOVAL-DATE	N	5.0		
			The beginning removal date, in TDMS TOX format, for the PRT set. Loaded by LDAS, as a gregorian date, into REM_- BEG_DATE in the PRO_REQ_TISS table.				
1	DD		END-REMOVAL-DATE	N	5.0		
			The ending removal date, in TDMS TOX format, for the PRT set. Loaded by LDAS, as a gregorian date, into REM_- END_DATE in the PRO_REQ_TISS table.				
1	EE		PROTO-REQ-TIS-PATHOLOGIST	N	5.0	N	
			The operator ID of the patho- logist assigned to the test. Loaded by LDAS into the PATH_- NUM field in the PRO_REQ_TISS table.				
M 1	FF		REASON-NUM	P	5.0	N	
			Codes for the removal reasons that apply to the PRT set. Loaded by LDAS into the REM_- REASONS field in the PRO_REQ_- TISS table.				
M 1	GG		TISSUE-NUM	P	5.0	N	
			Codes for the organs that ap- ply to the PRT set. Loaded by				

FILE...: PAS_PROTO_REQ_TIS

TYPE...: USER VIEW

FILE-NR: 84

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
			LDAS into the TISSUE_NUM field of the PRO_REQ_TISS table.				
1	HH		ALL-TREAT-FLAG	A	1.0	F	
			Flag indicating whether the PRT set applies to all treat- ment groups on the test.				
1	II		ALL-REM-REASON-FLAG	A	1.0	F	
			Flag indicating whether the PRT set applies to all removal reasons on the test.				
1	AB		SYSTEM-DATE-TIME	A	13.0	DE	
			The system date and time when the set was created. In the form MMDDYYHHMMSSH.				
1	SD		PROTO-REQ-TIS-SUPERKEY	A	22.0	SP	
			----- SOURCE FIELD(S) -----				
			DATA-BASE-KEY		1	9	
			BEG-REMOVAL-DATE		1	5	
			END-REMOVAL-DATE		1	5	
			TREAT-NUM		1	3	
			Used for indexed reads of the file.				
1	SE		AGENCY-EXP-TEST-TREAT-REASON	A	17.0	SP	
			----- SOURCE FIELD(S) -----				
			DATA-BASE-KEY		1	9	
			TREAT-NUM		1	3	
			REASON-NUM		1	5	
			Used by E18 and P14 for dis- crete record selection.				

PAS_RACK_CAGE

View of in-life protocol indicating the physical location of cages. Used by download to create the input for the initialization of the CAGE table.

FILE...: PAS_RACK_CAGE

TYPE...: USER VIEW

FILE-NR: 51

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER.	A	9.0	DE	
1	TE	DATA-DATE	The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)	A	8.0	N	
1	TF	DATA-TIME	The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)	A	6.0	N	
1	TG	DATA-FACILITY-NUMBER	The facility with which the entry operator is associated. Foreign key to TBFACILITY.	N	5.0	N	
1	TH	OPERATOR-ID	The entry or modification operator. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.	N	5.0	N	
1	TL	ORG-OR-COR	Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".	A	1.0	F	
1	TM	REASON	Text relating to a latent PEX/DOWNLOAD feature. No values assigned.	A	78.0	N	

FILE...: PAS_RACK_CAGE

TYPE...: USER VIEW

FILE-NR: 51

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	TN	COR-DATE	A	8.0	N	
		Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO	COR-TIME	A	6.0	N	
		Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB	C-AREA	N	5.0		
		The location at the facility of the rack (e.g. room number) Loaded by LDAS into the AREA_- NUM field of the CAGE table.				
1	CC	RACK-NO	N	2.0		
		The number of the rack being defined, which number must be unique on the test. Loaded by LDAS into the RACK_NUM field of the CAGE table.				
1	DD	C-SIDE	A	1.0	F	
		The side of the rack which is being described. Either A or B. Loaded by LDAS into the SIDE_FLAG field of the CAGE table.				
1	EE	RACK-ROW	N	2.0	N	
		The number of rows in the rack-side.				
1	FF	RACK-COL	N	2.0	N	
		The number of columns in the rack side.				
P 1	GG	CAGES				
		Repeating group defining phy- sical location of cages on the rack-side. Each group entry represents a row.				

FILE...: PAS_RACK_CAGE

TYPE...: USER VIEW

FILE-NR: 51

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
2	HA	COLUMN-1		N	4.0	N	
<p>The cage numbers in column 1. Loaded by LDAS into the CAGE_- NUM field in the CAGE table. When appended to right-most 7 characters in DATA-BASE-KEY, a foreign key to TD-CAGE.</p>							
2	HB	COLUMN-2		N	4.0	N	
<p>The cage numbers in column 2. Loaded by LDAS into the CAGE_- NUM field in the CAGE table. When appended to right-most 7 characters in DATA-BASE-KEY, a foreign key to TD-CAGE.</p>							
2	HC	COLUMN-3		N	4.0	N	
<p>The cage numbers in column 3. Loaded by LDAS into the CAGE_- NUM field in the CAGE table. When appended to right-most 7 characters in DATA-BASE-KEY, a foreign key to TD-CAGE.</p>							
2	HD	COLUMN-4		N	4.0	N	
<p>The cage numbers in column 4. Loaded by LDAS into the CAGE_- NUM field in the CAGE table. When appended to right-most 7 characters in DATA-BASE-KEY, a foreign key to TD-CAGE.</p>							
2	HE	COLUMN-5		N	4.0	N	
<p>The cage numbers in column 5. Loaded by LDAS into the CAGE_- NUM field in the CAGE table. When appended to right-most 7 characters in DATA-BASE-KEY, a foreign key to TD-CAGE.</p>							
2	HF	COLUMN-6		N	4.0	N	
<p>The cage numbers in column 6. Loaded by LDAS into the CAGE_- NUM field in the CAGE table. When appended to right-most 7 characters in DATA-BASE-KEY, a foreign key to TD-CAGE.</p>							

FILE...: PAS_RACK_CAGE

TYPE...: USER VIEW

FILE-NR: 51

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

2	HG	COLUMN-7	N	4.0	N	
---	----	----------	---	-----	---	--

The cage numbers in column 7.
Loaded by LDAS into the CAGE_-
NUM field in the CAGE table.
When appended to right-most
7 characters in DATA-BASE-KEY,
a foreign key to TD-CAGE.

2	HH	COLUMN-8	N	4.0	N	
---	----	----------	---	-----	---	--

The cage numbers in column 8.
Loaded by LDAS into the CAGE_-
NUM field in the CAGE table.
When appended to right-most
7 characters in DATA-BASE-KEY,
a foreign key to TD-CAGE.

2	HI	COLUMN-9	N	4.0	N	
---	----	----------	---	-----	---	--

The cage numbers in column 9.
Loaded by LDAS into the CAGE_-
NUM field in the CAGE table.
When appended to right-most
7 characters in DATA-BASE-KEY,
a foreign key to TD-CAGE.

2	HJ	COLUMN-10	N	4.0	N	
---	----	-----------	---	-----	---	--

The cage numbers in column 10.
Loaded by LDAS into the CAGE_-
NUM field in the CAGE table.
When appended to right-most
7 characters in DATA-BASE-KEY,
a foreign key to TD-CAGE.

1	SD	SUPER-KEY	A	17.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S)	---	-START-	--END-
DATA-BASE-KEY	1		9
C-AREA	1		5
RACK-NO	1		2
C-SIDE	1		1

Used for indexed reads of the
file. Primary key for table.

PAS_REMVL_ACTION

View of in-life, removal-action protocol. These actions are those to be performed by the data collector when an animal is removed from the test. Its information is downloaded to LDAS.

FILE...: PAS_REMVL_ACTION

TYPE...: USER VIEW

FILE-NR: 50

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MASTER.				
1	TE		DATA-DATE	A	8.0		N
			The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF		DATA-TIME	A	6.0		N
			The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG		DATA-FACILITY-NUMBER	N	5.0		N
			The facility with which the entry operator is associated. Foreign key to TBFACILITY.				
1	TH		OPERATOR-ID	N	5.0		N
			The entry or modification operator. When appended to DATA-FACILITY-NUMBER, foreign key to TB_OPERATOR_ID.				
1	TL		ORG-OR-COR	A	1.0		F
			Flag relating to a latent PEX/DOWNLOAD feature. Always set to "0".				

FILE...: PAS_REMVL_ACTION

TYPE...: USER VIEW

FILE-NR: 50

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	TM	REASON	A	78.0	N	
		Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TN	COR-DATE	A	8.0	N	
		Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO	COR-TIME	A	6.0	N	
		Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB	REMOVAL-ACT-CT	N	2.0	N	
		The number of entries in the multiple-entry field REMOVAL- ACT-CODE. Should be equal to C*REMOVAL-ACT-CODE.				
M 1	CC	REMOVAL-ACT-CODE	N	5.0	N	
		Action to performed when an animal is removed from the test. Loaded by LDAS into the ACTION_NUM field of the REMVL_ACTION table. Foreign key to TBACTONC.				

PAS_RESP_PARTIES

View of information on the people who share primary responsibility for the test. Not used by download or EIS/PEIS reports.

FILE...: PAS_RESP_PARTIES

TYPE...: USER VIEW

FILE-NR: 45

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER.				
1	TE		DATA-DATE	A	8.0		N
			The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF		DATA-TIME	A	6.0		N
			The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG		DATA-FACILITY-NUMBER	N	5.0		N
			The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS. Foreign key to TBFACLT.				
1	TH		OPERATOR-ID	N	5.0		N
			ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.) When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
1	TL		ORG-OR-COR	A	1.0		F
			Flag relating to a latent PEX/DOWNLOAD feature. Always set				

FILE...: PAS_RESP_PARTIES

TYPE...: USER VIEW

FILE-NR: 45

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
		to "O".				
1	TM	REASON	A	78.0	N	
		Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TN	COR-DATE	A	8.0	N	
		Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO	COR-TIME	A	6.0	N	
		Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB	RESP-PARTY	A	32.0	N	
		Name of one of the people with primary responsibility for the test, entered as free text.				
1	CC	JOB-POSITION	A	2.0	N	
		Two-letter identification for the RP's role on the test. The possible values are: PI - Principal Investigator PO - Project Manager CM - Chemical Manager (Defined in code, not data model.)				
1	DD	FACILITY-CODE	N	5.0	N	
		Code for the facility with which the RP is associated. Foreign key to TBFACTY.				
1	EE	EFFECTIVE-DATE	N	6.0	N	
		Date, entered in the form of MMDDYY, on which the RP's role on the test becomes effective.				
1	FF	PHONE-NUMBER	N	10.0	N	
		Phone number of the RP. En- tered as free numerical text.				

FILE...: PAS_RESP_PARTIES

TYPE...: USER VIEW

FILE-NR: 45

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
1 GG EXTENSION	N	4.0 N
The extension for the RP, if there is one. Entered as free numeric text.		
1 HH FACILITY-ADDRESS	A	128.0 N
The address for the facility with which the RP is associated. For entry purposes, the field is broken into 4 32-character subfields, into which the address is entered as free text.		

PAS_SITES

View of pathology protocol. Lists the tissue-sites that may be used for pathology data collection. This file is initially populated by copying from one of the universal download tests. It is rarely modified after initialization.

FILE...: PAS_SITES

TYPE...: USER VIEW

FILE-NR: 54

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	TC	DATA-BASE-KEY	A	9.0		DE
		Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER.				
1	TE	DATA-DATE	A	8.0		N
		The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF	DATA-TIME	A	6.0		N
		The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG	DATA-FACILITY-NUMBER	N	5.0		N
		The facility with which the entry operator is associated. Foreign key for TBFACILITY.				
1	TH	OPERATOR-ID	N	5.0		N
		Entry or modification operator. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
1	TL	ORG-OR-COR	A	1.0		F
		Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".				
1	TM	REASON	A	78.0		N
		Text relating to a latent PEX/				

FILE...: PAS_SITES

TYPE...: USER VIEW

FILE-NR: 54

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

DOWNLOAD feature. No values assigned.

1	TN	COR-DATE	A	8.0	N	
---	----	----------	---	-----	---	--

Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.

1	TO	COR-TIME	A	6.0	N	
---	----	----------	---	-----	---	--

Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.

1	BB	ORGAN-CODE	N	5.0		
---	----	------------	---	-----	--	--

The organ (tissue) to which the sites apply. Foreign key to TBPCT. Loaded by LDAS into the TISSUE_NUM field in the TISS_SITE table.

1	CC	SITE-GRP-CT	N	4.0	N	
---	----	-------------	---	-----	---	--

The number of tissue-sites in the multiple-entry field SITE-CODE. Should be equal to C*SITE-CODE.

M 1	DD	SITE-CODE	N	5.0	N	
-----	----	-----------	---	-----	---	--

The tissue-sites which may be used to further define observations of the associated organ. Foreign key to TBPCT. Loaded by LDAS into the SITE_NUM field of the TISS_SITE and PATH_SITE tables.

1	SD	SUPER-KEY	A	11.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-
DATA-BASE-KEY 1 9
ORGAN-CODE 1 2

Used for indexed reads of the table based on test number and organ code. Primary key for table.

PAS_TREATMNT_SEG

View of in-life, treatment, and treatment-regimen, protocol. Treatment groups are generally organized according to dose levels and gender.

FILE...: PAS_TREATMNT_SEG

TYPE...: USER VIEW

FILE-NR: 41

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MAS-TER.				
1	TE		DATA-DATE	A	8.0		N
			The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF		DATA-TIME	A	6.0		N
			The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG		DATA-FACILITY-NUMBER	N	5.0		N
			The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS. Foreign key to TBFACILITY.				
1	TH		OPERATOR-ID	N	5.0		N
			ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.) Foreign key to TB_OPERATOR_ID.				
1	TL		ORG-OR-COR	A	1.0		F
			Flag relating to a latent PEX/DOWNLOAD feature. Always set to "O".				

FILE...: PAS_TREATMNT_SEG

TYPE...: USER VIEW

FILE-NR: 41

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	TM	REASON	A	78.0	N	
		Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TN	COR-DATE	A	8.0	N	
		Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO	COR-TIME	A	6.0	N	
		Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB	ID-NO	N	3.0	N	
		Sequentially assigned number for the treatment group. Uni- que on the test. Values range from 1 to 32. Loaded by LDAS into the TREAT_NUM field in the TREATMENT and TREAT_PREP tables.				
1	CC	CONTROL-DES	A	1.0	F	
		The treatment role. Foreign key to TBTRTROLE.				
1	DD	TEXT-ID	A	16.0	N	
		Free text describing the treatment groups. Entered as two 8-character sub-fields due to design of some report headers. Loaded by LDAS into the TREAT_TEXT field of the TREATMENT table.				
1	EE	REGIMEN-CT	N	3.0	N	
		Number of entries in the re- peating group ANY-REGIMEN. Should be equal to C*ANY-REGI- MEN.				
P 1	FF	ANY-REGIMEN				
		Repeating group embracing each				

FILE...: PAS_TREATMNT_SEG

TYPE...: USER VIEW

FILE-NR: 41

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
		regiment associated with a treatment group.				
2	GG	PREP-ID	N	3.0	N	
		The preparation associated with the regimen. Loaded by LDAS into the PREP_NUM field in the TREAT_PREP table. Appended to DATA-BASE-KEY, a foreign key to PAS_PREPARTN_SEG.				
2	HH	ROUTE	N	5.0	N	
		Dose route associated with regimen. Foreign key to TBDOSRTE. Loaded by LDAS into the ROUTE_NUM field in the PREPARATION table.				
2	II	PREP-VOL	N	5.0	N	
		Ratio of preparation to animal weight, expressed as .000cc/g. Only for gavage and injection tests. Loaded by LDAS into the DOSE_VOLUME field in the TREAT_PREP table.				
2	JJ	DOSE-CALC	A	1.0		
		Legacy field since implementation of LDAS. Formerly, a flag indicating to the animal caretaker how to determine the method of determining the preparation to be given to each animal. Values: M - Manual Entry R - Real-Time Computation C - Computer-Generated Table				
2	KK	REL-START	N	4.0	N	
		The day-on-dose on which the schedule for the regimen begins. Loaded by LDAS into the REL_START field in the TREAT_PREP table.				

FILE...: PAS_TREATMNT_SEG

TYPE...: USER VIEW

FILE-NR: 41

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
2	LL	REL-STOP	N	4.0	N	
<p>The day-on-dose on which the schedule for the regimen ends. Loaded by LDAS into the REL_STOP field in the TREAT_PREP table.</p>						
2	MM	INTERVALS	N	2.0	N	
<p>The number of weeks between occurrences of the regimen. Loaded by LDAS into the WEEK_INTRVL field in the TREAT_PREP table.</p>						
2	NN	REL-OR-SPEC	A	1.0		
<p>Flag indicating the basis for the DAY-OF-WEEK field. R - test week, S - calendar week. Loaded by LDAS into the REL_OR_SPEC field in the TREAT_PREP table.</p>						
2	OO	DAY-OF-WEEK	A	7.0	N	
<p>Each character is a flag for the day corresponding to the position. If the character has a value of 'X', then the action is scheduled for the corresponding day; if blank, it isn't scheduled. Loaded by LDAS into the DAY_OF_WEEK field in the TREAT_PREP table.</p>						
1	PP	MIN-WATER	P	5.0	N	
<p>Minimum water consumption value, expressed in decigrams. Loaded by LDAS into the MIN_WATER field in the the TREATMENT table.</p>						
1	QQ	MAX-WATER	P	5.0	N	
<p>Maximum water consumption value, expressed in decigrams. Loaded by LDAS into the MAX_WATER field in the the TREATMENT table.</p>						

FILE...: PAS_TREATMNT_SEG

TYPE...: USER VIEW

FILE-NR: 41

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

1	RR	MIN-FOOD	P	5.0	N	
---	----	----------	---	-----	---	--

Minimum feed consumption value, expressed in decigrams. Loaded by LDAS into the MIN_FEED field in the the TREATMENT table.

1	SS	MAX-FOOD	P	5.0	N	
---	----	----------	---	-----	---	--

Maximum feed consumption value, expressed in decigrams. Loaded by LDAS into the MAX_FEED field in the the TREATMENT table.

1	TT	CURVE-NO	P	5.0	N	
---	----	----------	---	-----	---	--

The body weight growth curve that applies to the treatment group. Loaded by LDAS into the WEIGHT_CURVE field in the TREATMENT table. Foreign key to TBWEIGHTCURVE.

1	SD	SUPER-KEY	A	11.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-

DATA-BASE-KEY	1	9
---------------	---	---

ID-NO	1	2
-------	---	---

Used for indexed reads and discrete selections based on test number and treatment group ID. Primary key for table.

PAS_WGT_ITEMS

View of weight group in-life protocol. Contains information about the weight groups to which cages are assigned. Specifies absolute minimum and maximum weights for used and full containers.

FILE...: PAS_WGT_ITEMS

TYPE...: USER VIEW

FILE-NR: 46

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Composite field whose elements are the 2-digit agency number, which is always "05", and the 7-digit experiment-test number. Foreign key to PAS_MASTER.				
1	TE		DATA-DATE	A	8.0		N
			The date, in TDMS TOX format, on which the record was added, or last modified. (The TOX date is in the last 5 digits.)				
1	TF		DATA-TIME	A	6.0		N
			The time, in TDMS TOX format, at which the record was added, or last modified. (The TOX time is in the last 5 digits.)				
1	TG		DATA-FACILITY-NUMBER	N	5.0		N
			The ID number of the facility to which the entering operator is associated. Historically, has always been 20 for NIEHS. Foreign key to TBFACILITY.				
1	TH		OPERATOR-ID	N	5.0		N
			ID of operator who entered, or last modified, the record. (By procedure, the ID should always belong to a liaison.) When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
1	TL		ORG-OR-COR	A	1.0		F
			Flag relating to a latent PEX/				

FILE...: PAS_WGT_ITEMS

TYPE...: USER VIEW

FILE-NR: 46

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
		DOWNLOAD feature. Always set to "O".				
1	TM	REASON	A	78.0	N	
		Text relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TN	COR-DATE	A	8.0	N	
		Date relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	TO	COR-TIME	A	6.0	N	
		Time relating to a latent PEX/ DOWNLOAD feature. No values assigned.				
1	BB	WGT-GRP-ID	N	5.0	N	
		Sequentially assigned ID, which is unique on a test, for a group of container weight items, and to which a cage is assigned in the cage census. Loaded by LDAS into the GROUP_NUM field in the WEIGHT_GROUP table.				
1	CC	WGT-CT	N	4.0	N	
		Number of items detailed in the repeating group ANY-WGT-ITEM. Should be equal to C*ANY-WGT-ITEM.				
P 1	DD	ANY-WGT-ITEM				
		Repeating group embracing each weight object in the weight group.				
2	EE	ITEM-ID	N	2.0	N	
		The weight object to which a weight range is assigned. Loaded by LDAS into the ITEM_NUM field in the WEIGHT_GROUP table. Foreign key to TBWGTOBJ.				
2	FF	LOW-WEIGHT	N	6.0	N	
		The absolute minimum weight				

FILE...: PAS_WGT_ITEMS

TYPE...: USER VIEW

FILE-NR: 46

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

for the weight object, expressed in decigrams. Loaded by LDAS into the MIN_WEIGHT field in the WEIGHT_GROUP table.

2	GG	HI-WEIGHT		N	6.0	N	
---	----	-----------	--	---	-----	---	--

The absolute maximum weight for the weight object, expressed in decigrams. Loaded by LDAS into the MAX_WEIGHT field in the WEIGHT_GROUP table.

1	SD	SUPER-KEY		A	11.0	SP	
---	----	-----------	--	---	------	----	--

SOURCE FIELD(S) --- -START- --END-							
		DATA-BASE-KEY		1		9	
		WGT-GRP-ID		1		2	

Used for indexed reads of the file, based on test number and weight group ID. Primary key for the table.

4.2 EIS/ PEIS Views

TD-ANIMAL

View embracing animal allocation, animal transfer, animal removal, histology number, "not examined", and uncertain cause-of-death segments. These segments are unique observations of animals.

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AA	TDA-EXP-TEST-ANIMAL-UQ	Composite field whose elements are the 7-digit test number, the 4-digit number of the original cage to which the animal was allocated, and the 5-digit animal number assigned at allocation. The test and cage portions form a foreign key to TD-CAGE.	A	16.0	UQ	
1	AB	TDA-ORIGINAL-CAGE-NUMBER-DE	The number of the cage to which the animal was originally allocated. Duplicate of cage number portion of TDA-EXP-TEST-ANIMAL-UQ. Used by report queries for cage range selection.	N	4.0		
1	CA	TDA-CARCASS-IDENTIFICATION-DE	Alternative primary key assigned at removal and used for pathology record selection. Made depending on CID-FLAG in PAS_ANIMAL_ID. In all cases, the first 9-digits are 05 and the 7-digit test number. When CID-FLAG is 'C', the last 5 digits are the 4-digit cage number and a 1-digit sequen-	A	14.0	N	DE

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T	L	DB NAME	F	LENG	S	DE
		tial identifier for each animal- mal in the cage; when CID-FLAG is 'E', it is merely the animal- mal number. (Animal numbers were not necessarily unique on older tests.)				
1	CB	TDA-BALANCE-NUMBER-DE	N	5.0		
		Legacy field. Formerly, the unique number of the balance used to determine the removal weight. Still used in ECS and Update, although Update merely assigns the default value, which is 0. Foreign key to TD-BALANCE.				
1	CC	TDA-REASON-FOR-REMOVAL-DE	N	5.0	N	
		The code for the removal rea- son that was assigned by the operator recording the remov- al. Foreign key to TBCLOBSC.				
1	CD	TDA-HOSPITAL-CAGE-DE	N	4.0	N	
		If an animal has been trans- ferred from its original cage, this will contain the number of the cage to which it was transferred. Used by EIS/PEIS reports for cage range selec-				
1	AC	TDA-ANIMAL-IDENTIFICATION	A	1.0	F	
		Legacy field. No values exist for field in database. Refer- enced by Update and various report queries.				
1	AD	TDA-ANIMAL-CLASS	N	5.0		
		Animal class code associated with the test in protocol. (Animal class qualifiers are a legacy to TDMS, with no cur- rent use. New tests are as- signed the code for "Not Ap- plicable.") Referenced, but				

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T	L	DB NAME	F	LENG	S	DE
		not used, in report software. Foreign key to TBANCLQF.				
1	AE	TDA-SEX	A	1.0	F	
		Gender identifier. F for females; M for males.				
1	AF	TDA-SPECIES	N	5.0		
		Code for species as assigned in test protocol. Foreign key to TBSTRAIN.				
1	AG	TDA-STRAIN	N	5.0		
		Code for species strain, as defined in test protocol. Foreign key to TBSTRAIN.				
1	AH	TDA-SUBSTRAIN	N	5.0		
		Substrain code associated with the test in protocol. (Substrains are a legacy to TDMS, with no current use. New tests are assigned the code for "Not Applicable.") Referenced, but not used, in report software. Foreign key to TBSUBSTR.				
1	AI	TDA-ANIMAL-DATE	N	5.0		
		The date, in TDMS TOX format, on which the transaction was created. Not used in report software, except for correction history report, where changes to it are audited.				
1	AJ	TDA-ANIMAL-TIME	N	5.0		
		The time, in TDMS TOX format, at which the transaction was created. Not used in report software, except for correction history report, where changes to it are audited.				
1	AK	TDA-ANIMAL-OPERATOR-ID	N	5.0		
		The ID of the operator who performed the allocation. Not used in report software.				

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	AL	TDA-ANIMAL-ERROR-FLAG	A	1.0	F	
---	----	-----------------------	---	-----	---	--

Flag indicating whether animal allocation segment was flagged as being in error upon update. N - No; Y - Yes. Applies only to early MODCOMP data.

1	AM	TDA-ACTION-OVERRIDE-FLAG-CID	A	1.0	F	
---	----	------------------------------	---	-----	---	--

Flag indicating whether the scheduled CID assignment was overridden. 0 - Not overridden; 1 - Overridden. Pertains only to MODCOMP tests.

1	AN	TDA-MANUAL-ENTRY-FLAG-CID	A	1.0	F	
---	----	---------------------------	---	-----	---	--

Flag indicating whether CID assignment was made as part of a manual entry session. 0 - No; 1 - Yes. Applies only to MODCOMP tests.

1	AO	TDA-REMOVAL-ACTION-FLAG-CID	A	1.0	F	
---	----	-----------------------------	---	-----	---	--

Flag indicating whether the CID assignment was a removal action. 0 - No; 1 - Yes. CID assignment, by definition, is a removal action. Hence, the value is always 1 when a CID has been assigned. Pertains only to MODCOMP tests.

1	AP	TDA-BACKWARD-FORWARD-FLAG-CID	A	1.0	F	
---	----	-------------------------------	---	-----	---	--

Flag indicating whether CID assignment was performed in scheduled order. 0 - Yes; 1 - No. Pertains only to MODCOMP tests.

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AQ		TDA-EXPECTED-SELECTED-FLAG-CID	A	1.0	F	
			Flag indicating whether CID assignment was a scheduled ac-				
			. 0 - Yes; 1 - No. Per-				
			tains only to MODCOMP tests.				
1	AS		TDA-ACTION-OVERRIDE-FLAG-REM	A	1.0	F	
			Flag indicating whether the scheduled removal action was				
			overridden. 0 - Not overrid-				
			den; 1 - Overridden. Pertains				
			only to MODCOMP tests.				
1	AT		TDA-MANUAL-ENTRY-FLAG-REM	A	1.0	F	
			Flag indicating whether the removal was recorded as part				
			of a manual data entry ses-				
			sion. 0 - No; 1 - Yes. Not				
			used in report software.				
1	AU		TDA-REMOVAL-ACTION-FLAG-REM	A	1.0	F	
			Flag indicating whether the removal was a removal action.				
			0 - No; 1 - Yes. An animal				
			removal, by definition, is a				
			removal action. Hence, the				
			value is always 1 when a re-				
			moval is recorded. Pertains				
			only to MODCOMP tests.				
1	AV		TDA-BACKWARD-FORWARD-FLAG-REM	A	1.0	F	
			Flag indicating whether animal removal was performed in the				
			scheduled order. 0 - Yes;				
			1 - No. Pertains only to				
			MODCOMP tests.				
1	AW		TDA-EXPECTED-SELECTED-FLAG-REM	A	1.0	F	
			Flag indicating whether animal removal was a scheduled ac-				
			. 0 - Yes; 1 - No. Per-				
			tains only to MODCOMP tests.				
1	AX		TDA-REMOVAL-DATE	N	5.0	N	
			Date, in TDMS TOX format, on				
			which the removal was per-				
			formed.				

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T L DB NAME	F	LENG S DE
1 CF TDA-REMOVAL-TIME Time, in TDMS TOX format, at which the removal was performed.	N	5.0 N
1 CE TDA-REMOVAL-OPERATOR-ID The ID of the operator who performed the removal. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.	N	5.0 N
1 CG TDA-REMOVAL-ERROR-FLAG Flag indicating whether animal removal segment was flagged as being in error upon update. N - No; Y - Yes. Applies only to early MODCOMP data.	A	1.0 F
1 CM TDA-REMOVED-DATE Date, in TDMS TOX format, on which the removal occurred. Duplicate of TDA-REMOVAL-DATE. Any historical distinction between the two fields is no longer known.	N	5.0 N
1 AY TDA-DAYS-ON-EXPERIMENT Number of days on experiment. The difference between the removal and allocation dates.	N	5.0
1 AZ TDA-ACTION-OVERRIDE-FLAG-WGT Flag indicating whether the scheduled removal weight was overridden. 0 - Not overridden; 1 - Overridden. Pertains only to MODCOMP tests.	A	1.0 F
1 BA TDA-MANUAL-ENTRY-FLAG-WGT Flag indicating whether removal weight was recorded as part of a manual entry session. 0 - No; 1 - Yes. Not used in report software, except Correct History, where any changes to it are audited.	A	1.0 F

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T L	DB NAME	F	LENG	S DE
1 BB	TDA-REMOVAL-ACTION-FLAG-WGT	A	1.0	F
	Flag indicating whether the removal weight was a removal action. 0 - No; 1 - Yes. A removal weight, by definition, is a removal action. Hence, the value is always when a removal weight is recorded. Pertains only to MODCOMP tests.			
1 BC	TDA-BACKWARD-FORWARD-FLAG-WGT	A	1.0	F
	Flag indicating whether removal weight was performed in scheduled order. 0 - Yes; 1 - No. Pertains only to MODCOMP tests.			
1 BD	TDA-EXPECTED-SELECTED-FLAG-WGT	A	1.0	F
	Flag indicating whether removal weight was a scheduled action. 0 - Yes; 1 - No. Pertains only to MODCOMP tests.			
1 BE	TDA-ANIMAL-REM-WEIGHT	N	10.0	N
	Weight, in decigrams, of the animal at removal. Pre-LDAS removals where a weight was not recorded will show a value of 99999, which value is flagged by reports and displayed as "OVERRIDE"; post-LDAS data where a weight was not recorded will show a value of 0, and have a value assigned to TDA-REM-WEIGHT-STATUS.			
1 BF	TDA-WEIGHT-OVERRIDE	A	1.0	F
	Flag indicating weight range validation. A = Within range; + = Above range; - = Below			

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T L DB NAME	F	LENG S DE
range. Applies only to MOD- COMP tests.		
1 BI TDA-ACTION-OVERRIDE-FLAG-OBS	A	1.0 F
Flag indicating whether the scheduled removal observation action was overridden. 0 - Not overridden; 1 - Overrid- den. Pertains only to MODCOMP tests.		
1 BJ TDA-MANUAL-ENTRY-FLAG-OBS	A	1.0 F
Flag indicating whether remov- al observations were recorded as part of a manual entry ses- sion. 0 - No; 1 - Yes. Not used in reports, except Corresponds- rection History, where any changes to it are audited.		
1 BK TDA-REMOVAL-ACTION-FLAG-OBS	A	1.0 F
Flag indicating whether the removal observations were re- moval actions. 0 - No; 1 - Yes. Removal observations, by definition, are removal ac- tions. Hence, the value is always 1 when removal obser- vations are recorded. Per- tains only to MODCOMP tests.		
1 BL TDA-BACKWARD-FORWARD-FLAG-OBS	A	1.0 F
Flag indicating whether remov- al observations were performed in scheduled order. 0 - Yes; 1 - No. Pertains only to MODCOMP tests.		
1 BM TDA-EXPECTED-SELECTED-FLAG-OBS	A	1.0 F
Flag indicating whether remov- al observations were scheduled actions. 0 - Yes; 1 - No. Pertains only to MODCOMP tests.		
M 1 BN TDA-DISPOSITION-CODE	N	5.0 N
Legacy field with no meaning		

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

in TDMS. Intended to be codes indicating the disposition of animals after removal. No such codes exists in the TDMS tables.

1	BP	TDA-TRANSFER-DATE	N	5.0	N		
---	----	-------------------	---	-----	---	--	--

Date, in TDMS TOX format, on which animal was transferred to a hospital cage.

1	BQ	TDA-TRANSFER-TIME	N	5.0	N		
---	----	-------------------	---	-----	---	--	--

Time, in TDMS TOX format, at which animal was transferred to a hospital cage.

1	BR	TDA-TRANSFER-OPERATOR-ID	N	5.0	N		
---	----	--------------------------	---	-----	---	--	--

ID of the operator who recorded the transfer. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	BS	TDA-TRANSFER-ERROR-FLAG	A	1.0	F		
---	----	-------------------------	---	-----	---	--	--

Flag indicating whether animal transfer segment was flagged as being in error upon update. N - No; Y - Yes. Applies only to early MODCOMP data.

1	BT	TDA-PATHOLOGIST-IDENTIFICATION	N	5.0			
---	----	--------------------------------	---	-----	--	--	--

Operator ID of the pathologist assigned to the test when the histology number was assigned. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	BU	TDA-HISTOLOGY-NUMBER-DATE	N	5.0	N		
---	----	---------------------------	---	-----	---	--	--

Date, in TDMS TOX format, on which the histology number was assigned. Not used in report

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T L DB NAME	F	LENG S DE
software.		
1 BV TDA-HISTOLOGY-NUMBER-TIME	N	5.0 N
Time, in TDMS TOX format, at which the histology number was assigned. Not used in report software.		
1 BW TDA-HISTOLOGY-NUM-OPERATOR-ID	N	5.0 N
ID of the operator who recorded the histology number assignment. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.		
1 BY TDA-HISTOLOGY-NUM-ERROR-FLAG	A	1.0 F
Flag indicating whether transaction was flagged as being in error by Update. N - No; Y - Yes. Applies only to early MODCOMP data.		
1 BX TDA-HISTOLOGY-NUMBER	A	16.0 N
Lab-specific identification number for carcasses.		
1 BZ TDA-PATHOLOGIST-TYPE	A	1.0 F
Code identifying operator type of pathologist assigned to the test. If TDMS were to have assigned values, they would have been 'A' (MODCOMP) or 'P' (LDAS), but no values have ever been assigned.		
1 EA TDA-REM-WEIGHT-STATUS	P	5.0 F
Code for weight status that describes the reason why a removal weight wasn't recorded. (TDA-ANIMAL-REM-WEIGHT should have a value of 0 if this field has a non-0 value.) A foreign key to TBWGTSTATUS.		
P 1 FA TDA-OBSERVATION-DATA		
Repeating group embracing each removal observation.		

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T L DB NAME	F	LENG S DE
2 FB TDA-OBSERVATION Code for clinical sign object- served at removal. A foreign key to TBCLOBSC.	P	5.0 N
2 FC TDA-SITE Code for in-life site where clinical sign was observed. (Only applicable for certain clinical signs, as defined in the TDMS tables.) A foreign key to TBCLOBSC.	P	5.0 N
2 FD TDA-EIQUAL1 Code for in-life size description- cribing clinical sign. (Only applicable for certain clini- cal signs, as defined in the TDMS tables.) A foreign key to TBCLOBSC.	P	5.0 N
2 FE TDA-EIQUAL2 Code for in-life count description- cribing clinical sign. (Only applicable for certain clini- cal signs, as defined in the TDMS tables.) A foreign key to TBCLOBSC.	P	5.0 N
1 AR TDA-NO-EXAM-PATHOLOGIST Operator ID of pathologist assigned to the test when the "not examined" status is re- corded. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR- _ID.	P	5.0 N
1 A0 TDA-NO-EXAM-OPERATOR-ID ID of operator who recorded the "not examined" status. Not used in report software. When appended to the appropri- ate facility code, a foreign key to TB_OPERATOR_ID.	P	5.0 N

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T	L	DB	NAME	F	LENG	S	DE
1	A1		TDA-NO-EXAM-DATE	P	5.0	N	DE
			Date, in TDMS TOX format, on which the "not examined" status was recorded. Used by reports to flag whether such a status was assigned.				
1	A2		TDA-NO-EXAM-TIME	P	5.0	N	
			Time, in TDMS TOX format, at which the "not examined" status was recorded. Not used by reports.				
1	A3		TDA-UNCERTAIN-PATHOLOGIST	P	5.0	N	
			Operator ID of pathologist assigned to the test when the uncertain cause-of-death is recorded. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	A4		TDA-UNCERTAIN-OPERATOR-ID	P	5.0	N	
			ID of operator who recorded the uncertain cause-of-death. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	A5		TDA-UNCERTAIN-DATE	P	5.0	N	
			Date, in TDMS TOX format, on which the uncertain cause-of-death was recorded. Used by reports to flag whether an uncertain cause-of-death was recorded.				
1	A6		TDA-UNCERTAIN-TIME	P	5.0	N	
			Time, in TDMS TOX format, at which the uncertain cause-of-death was recorded. Not used in report software.				

FILE...: TD-ANIMAL

TYPE...: USER VIEW

FILE-NR: 103

PRIMARY SEQUENCE FIELD: TDA-EXP-TEST-ANIMAL-UQ

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

1	DP	TDA-EXP-TEST-ANIMAL-SP		A	12.0		SP
		SOURCE FIELD(S) --- -START- --END-					
		EXP-TEST-ANIMAL		1		7	
		EXP-TEST-ANIMAL		12		16	

Superdescriptor based on
test number and animal number
portions of the primary key.
Refer to ANIMAL-EIS field
descriptions for source field
details.

1	CN	TDA-EXP-TEST-HOSPITAL-CAGE-SP		A	11.0		SP
		SOURCE FIELD(S) --- -START- --END-					
		EXP-TEST-ANIMAL		1		7	
		HOSPITAL-CAGE		1		4	

Superkey used to select trans-
ferred animals. A foreign
key to TD-CAGE. See ANIMAL-
EIS field descriptions for
source field information.

TD-ANIMAL-DATA

View embracing animal weight, animal observation, and in-life animal note segments.

FILE...: TD-ANIMAL-DATA

TYPE...: USER VIEW

FILE-NR: 106

PRIMARY SEQUENCE FIELD: TDD-EXP-TEST-ANIMAL-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
1	AA	TDD-EXP-TEST-ANIMAL-DATE-DE	A	21.0	DE		
Composite field whose elements are the 7-digit test number, the 4-digit number of the original cage to which the animal was allocated, the 5-digit animal number assigned at allocation, and the 5-digit date, in TDMS TOX format, on which the observation was recorded. The test, cage, and animal portions form a foreign key to TD-ANIMAL.							
1	AL	TDD-BALANCE-NUMBER-DE	N	5.0	N		
Legacy field. Formerly, the unique balance number on which the weight was recorded; currently, Update assigns a 0. A foreign key to TD-BALANCE.							
1	AB	TDD-WEIGHT-TIME	N	5.0	N		
The time, in TDMS TOX format, at which the transaction was created.							
1	AC	TDD-ANIMAL-WEIGHT	N	10.0	N		
The weight of the animal, in decigrams. Pre-LDAS records where a weight was not recorded will show a value of 99999, which value is flagged by reports and displayed as "OVERRIDE"; post-LDAS data where a weight was not recorded will show a value of 0, and have a value assigned to TDD-							

FILE...: TD-ANIMAL-DATA

TYPE...: USER VIEW

FILE-NR: 106

PRIMARY SEQUENCE FIELD: TDD-EXP-TEST-ANIMAL-DATE-DE

T	L	DB NAME	F	LENG	S	DE
		WEIGHT-STATUS.				
1	AD	TDD-WEIGHT-OPERATOR-ID	N	5.0	N	
		The ID of the operator who performed the allocation. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	AE	TDD-WEIGHT-ERROR-FLAG	A	1.0	F	
		Flag indicating whether weight segment record was flagged as being in error at time of Update. N - No; Y - Yes. Applies only to early MODCOMP data.				
1	AF	TDD-ACTION-OVERRIDE-FLAG-WGT	A	1.0	F	
		Flag indicating whether the scheduled weight action was overridden. 0 - No; 1 - Yes. Pertains only to MODCOMP tests.				
1	AG	TDD-MANUAL-ENTRY-FLAG-WGT	A	1.0	F	
		Flag indicating whether the weight was recorded as part of a manual data entry session. 0 - No; 1 - Yes.				
1	AH	TDD-REMOVAL-ACTION-FLAG-WGT	A	1.0	F	
		Flag indicating whether weight was a removal action. 0 - No; 1 - Yes. Pertains only to MODCOMP tests.				
1	AI	TDD-BACKWARD-FORWARD-FLAG-WGT	A	1.0	F	
		Flag indicating whether weight was recorded in scheduled order. 0 - Yes; 1 - No. Pertains only to MODCOMP tests.				
1	AJ	TDD-EXPECTED-SELECTED-FLAG-WGT	A	1.0	F	
		Flag indicating whether weight				

FILE...: TD-ANIMAL-DATA

TYPE...: USER VIEW

FILE-NR: 106

PRIMARY SEQUENCE FIELD: TDD-EXP-TEST-ANIMAL-DATE-DE

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

		was a scheduled action. 0 - Yes; 1 - No. Pertains only to MODCOMP tests.				
1	AK	TDD-WEIGHT-OVERRIDE	A	1.0	F	
		Flag indicating whether weight was within acceptable range. A = Within range; + = Above range; - = Below range. Applies only to MODCOMP tests.				
1	BB	TDD-OBSERVATION-TIME	N	5.0	N	
		Time, in TDMS TOX format, at which the observations were recorded.				
1	AM	TDD-OBSERVATION-OPERATOR-ID	N	5.0	N	
		ID of the operator who recorded the observations. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	AN	TDD-OBSERVATION-ERROR-FLAG	A	1.0	F	
		Flag indicating whether observation record was flagged as being in error at time of Update. N - No; Y - Yes. Applies only to early MODCOMP data.				
1	AO	TDD-ACTION-OVERRIDE-FLAG-OBS	A	1.0	F	
		Flag indicating whether the scheduled observation action was overridden. 0 - No; 1 - Yes. Pertains only to MODCOMP tests.				
1	AP	TDD-MANUAL-ENTRY-FLAG-OBS	A	1.0	F	
		Flag indicating whether the observation was recorded as part of a manual data entry session. 0 - No; 1 - Yes.				
1	AQ	TDD-REMOVAL-ACTION-FLAG-OBS	A	1.0	F	
		Flag indicating whether observation was a removal action. 0 - No; 1 - Yes. Pertains on-				

FILE...: TD-ANIMAL-DATA

TYPE...: USER VIEW

FILE-NR: 106

PRIMARY SEQUENCE FIELD: TDD-EXP-TEST-ANIMAL-DATE-DE

T	L	DB NAME	F	LENG	S	DE
		ly to MODCOMP tests.				
1	AR	TDD-BACKWARD-FORWARD-FLAG-OBS	A	1.0	F	
		Flag indicating whether obser- vation was recorded in sched- uled order. 0 - Yes; 1 - No. tains only to MODCOMP tests.				
1	AS	TDD-EXPECTED-SELECTED-FLAG-OBS	A	1.0	F	
		Flag indicating whether obser- vation was a scheduled action. 0 - Yes; 1 - No. Pertains on- ly to MODCOMP tests.				
1	AT	TDD-NOTE-OPERATOR-ID	N	5.0	N	
		ID of the operator who entered the in-life animal note. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	BC	TDD-NOTE-TIME	N	5.0	N	
		Time, in TDMS TOX format, at which the note was entered.				
1	AU	TDD-NOTE-ERROR-FLAG	A	1.0	F	
		Flag indicating whether note segment record was flagged as being in error at time of Up- date. N - No; Y - Yes. Ap- plies only to early MODCOMP tests.				
1	AV	TDD-ACTION-NUMBER	N	5.0	N	
		The code for the action assoc- iated with the observation. Only relevant for pre-removal observation version of MOD- COMP. A foreign key to TBACTONC.				
1	AW	TDD-ANIMAL-NOTE	A	80.0	N	
		Free text describing, explain- ing, or commenting on the				

FILE...: TD-ANIMAL-DATA

TYPE...: USER VIEW

FILE-NR: 106

PRIMARY SEQUENCE FIELD: TDD-EXP-TEST-ANIMAL-DATE-DE

T L DB NAME	F	LENG S DE
animal, or observations of the animal.		
P 1 DA TDD-OBSERVATION-DATA Repeating group contain clinical observations.		
2 DB TDD-OBSERVATION Codes for clinical signs observed on the animal. A foreign key to TBCLOBSC.	N	5.0 N
2 DC TDD-SITE Codes for in-life sites where clinical sign was observed. (Only applicable for certain clinical signs, as defined in the TDMS tables.) A foreign key to TBCLOBSC.	N	5.0 N
2 DD TDD-EIQUAL1 Codes for in-life sizes describing clinical signs. (Only applicable for certain clinical signs, as defined in the TDMS tables.) A foreign key to TBCLOBSC.	N	5.0 N
2 DE TDD-EIQUAL2 Codes for in-life counts describing clinical signs. (Only applicable for certain clinical signs, as defined in the TDMS tables.) A foreign key to TBCLOBSC.	N	5.0 N
1 DF TDD-MANUAL-ENTRY-FLAG-NOTE Flag indicating whether note was recorded as part of a manual entry session. 0-No; 1-Yes.	A	1.0 F
1 DG TDD-WEIGHT-STATUS Code for weight status that describes the reason why a	P	5.0 N

FILE...: TD-ANIMAL-DATA

TYPE...: USER VIEW

FILE-NR: 106

PRIMARY SEQUENCE FIELD: TDD-EXP-TEST-ANIMAL-DATE-DE

T L DB NAME

F LENG S DE

weight wasn't recorded. (TDD-
ANIMAL-WEIGHT should be 0 if
this field has a non-0 value.)
Foreign key to TBWGTSTATUS.

TD-BALANCE

View of balance information. Some reports access this view.

FILE...: TD-BALANCE

TYPE...: USER VIEW

FILE-NR: 104

PRIMARY SEQUENCE FIELD: TDB-BALANCE-NUMBER-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AA	TDB-BALANCE-NUMBER-UQ		N	5.0	F	DE
			The unique ID number for a balance used to weigh animals and containers.				
1	AB	TDB-AGENCY-NUMBER		N	2.0	F	
			The agency. Foreign key to TBAGENCY.				
1	AC	TDB-WEIGHT-UNIT-CODE		N	5.0	F	
			The unit of measure. Foreign key to TBUNITS.				

TD-BALANCE-CALIBRATION

View of balance calibration data. Although this file is updated daily, no data in it is retrieved by any TDMS software.

FILE...: TD-BALANCE-CALIBRATION

TYPE...: USER VIEW

FILE-NR: 102

PRIMARY SEQUENCE FIELD: TDF-EXP-TEST-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
1	AA		TDF-EXP-TEST-DATE-DE	N	12.0		DE
			Composite field whose elements are the 7-digit test number and the 5-digit date, in TDMS TOX format, on which the calibration data was recorded. The test portion is a foreign key to TD-TEST.				
1	AB		TDF-BALANCE-NUMBER-DE	N	5.0		
			The balance on which the calibration data was recorded. Foreign key to TD-BALANCE.				
P 1	AC		TDF-CALIBRATION-WEIGHT-DATA				
			Repeating group embracing calibration weight details.				
2	AD		TDF-CALIBRATION-WEIGHT	N	5.0		N
			The weight, in grams, of the calibration object.				
2	AE		TDF-WEIGHT-OVERRIDE	A	1.0		N
			Flag indicating whether the weight was within the specified range. A = Within range; + = Above range; - = Below range. Applies only to MOD-COMP tests.				
1	AF		TDF-TIME	N	5.0		
			Time, in TDMS TOX format, at which the calibration data was recorded.				
1	AG		TDF-OPERATOR-IDENTIFICATION	N	5.0		
			The operator who recorded the calibration data. When appended to the appropriate facility code, a foreign key to				

FILE...: TD-BALANCE-CALIBRATION

TYPE...: USER VIEW

FILE-NR: 102

PRIMARY SEQUENCE FIELD: TDF-EXP-TEST-DATE-DE

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

TB_OPERATOR_ID.

1	AH	TDF-ERROR-FLAG	A	1.0	F	
---	----	----------------	---	-----	---	--

Flag indicating whether balance calibration transaction was flagged by Update as being in error. N - No; Y - Yes. Applies only to early MODCOMP data.

1	AI	TDF-CALI-MANUAL-FLAG	A	1.0	F	
---	----	----------------------	---	-----	---	--

Flag indicating whether calibration data was recorded as part of a manual entry session. 0 - No; 1 - Yes.

TD-CAGE

View embracing the cage segment and the defunct "cage aborted" segment.

FILE...: TD-CAGE

TYPE...: USER VIEW

FILE-NR: 105

PRIMARY SEQUENCE FIELD: TDC-EXP-TEST-CAGE-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AA		TDC-EXP-TEST-CAGE-UQ	N	11.0	F	DE
			Composite field whose elements are the 7-digit test number and the 4-digit cage-number. The test number portion is a foreign key to TD-TEST.				
1	AB		TDC-CAGE-OPERATOR-ID	N	5.0	N	
			ID of the operator who performed the cage initialization. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	AC		TDC-CAGE-ERROR-FLAG	A	1.0	F	
			Flag indicating whether transaction was flagged by Update as being in error. N - No; Y - Yes. Pertains only to early MODCOMP data.				
1	AD		TDC-CAGE-DATE	N	5.0	N	
			Date, in TDMS TOX format, on which cage initialization was performed. Not used in report software, except in the correction history, where changes to it are audited.				
1	AE		TDC-CAGE-TIME	N	5.0	N	
			Time, in TDMS TOX format, at which cage initialization was performed. Not used in report software, except in the correction history, where changes to it are audited.				

FILE...: TD-CAGE

TYPE...: USER VIEW

FILE-NR: 105

PRIMARY SEQUENCE FIELD: TDC-EXP-TEST-CAGE-UQ

T	L	DB NAME	F	LENG	S	DE
1	AF	TDC-TREATMENT-NUMBER	N	3.0	F	
		ID of treatment group to which the cage belongs, as defined in test protocol. When appended to the appropriate agency code and test number, a foreign key to PAS_TREATMNT-_SEG.				
1	AG	TDC-PROCEDURE-ACTION-SET-NUM	N	5.0	N	
		Number of the procedure action set to which the cage was assigned in test protocol. Not used in report software. When appended to the appropriate agency number and test number, a foreign key to PAS_PROCEDUR-_SEG.				
1	AH	TDC-START-DATE	N	5.0	N	
		Date, in TDMS TOX format, on which the cage started. For pre-LDAS studies, this should be the same as TDC-CAGE-DATE; for post-LDAS studies, this is the first day-on-dose.				
1	AI	TDC-AREA-NUMBER	N	5.0	N	
		Number of the area to which the cage was assigned in test protocol. Not used in report software.				
1	AJ	TDC-ABORT-DATE	N	5.0	N	
		Date, in TDMS TOX format, when cage was aborted. (Cage abort data has never been part of TDMS.)				
1	AK	TDC-ABORT-TIME	N	5.0	N	
		Time, in TDMS TOX format, when cage was aborted. (Cage Abort				

FILE...: TD-CAGE

TYPE...: USER VIEW

FILE-NR: 105

PRIMARY SEQUENCE FIELD: TDC-EXP-TEST-CAGE-UQ

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

data has never been part of
TDMS.)

1	AL	TDC-ABORT-OPERATOR-ID	N	5.0	N		
---	----	-----------------------	---	-----	---	--	--

ID of operator performing cage
abort. (Cage abort data has
never been part of TDMS.)
When appended to the appropri-
ate facility code, a foreign
key to TB_OPERATOR_ID.

1	AM	TDC-ABORT-ERROR-FLAG	A	1.0	F		
---	----	----------------------	---	-----	---	--	--

Flag indicating whether cage
abort transaction was flagged
by Update as being in error.
(Cage abort data has never
been part of TDMS.)

1	AN	TDC-ABORTED-DATE	N	5.0	N		
---	----	------------------	---	-----	---	--	--

Date, in TDMS TOX format, when
cage was aborted. (Cage abort
data has never been part of
TDMS.)

1	AO	TDC-EXP-TEST-TREATMENT-SP	N	10.0	SP		
---	----	---------------------------	---	------	----	--	--

SOURCE FIELD(S) --- -START- --END-
>>> SOURCE FIELD(S) NOT SPECIFIED <<<

Super-key used to select all
cages in a treatment group.

1	AP	TDC-EXP-TEST-TREATMENT-CAGE-SP	N	14.0	SP		
---	----	--------------------------------	---	------	----	--	--

SOURCE FIELD(S) --- -START- --END-
>>> SOURCE FIELD(S) NOT SPECIFIED <<<

Used to read cage records
logically in treatment-cage
sequence.

1	AQ	TDC-EXP-TEST-SP	N	7.0	SP		
---	----	-----------------	---	-----	----	--	--

SOURCE FIELD(S) --- -START- --END-
>>> SOURCE FIELD(S) NOT SPECIFIED <<<

Used to select all cage re-
cords for a specific test
number.

TD-CAGE-DATA

View embracing feeder weight, bottle weight, cage note, feed consumption, water consumption, cage observation, and misidentified animals. (The last four are inactive.)

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
1	AA		TDG-EXP-TEST-CAGE-DATE-DE	N	16.0		DE
			Composite field whose elements are the 7-digit test number, the 4-digit number of the cage, and the 5-digit date, in TDMS TOX format, on which the weights or note was recorded. The test and cage number portion form a foreign key to TD-CAGE.				
1	AY		TDG-FEEDER-BALANCE-NUMBER-DE	N	5.0		N
			Legacy field. Formerly, the balance number on which the weight was recorded; currently, Update assigns a 0. A foreign key TD-BALANCE.				
1	BV		TDG-BOTTLE-BALANCE-NUMBER-DE	N	5.0		N
			Legacy field. Formerly, the balance number on which the weight was recorded; currently, Update assigns a 0. A foreign key to TD-BALANCE.				
1	AB		TDG-FEEDER-TIME	N	5.0		N
			Time, in TDMS TOX format, at which the feeder weight was recorded.				
1	AC		TDG-FEEDER-OPERATOR-ID	N	5.0		N
			The operator who recorded the feeder weight. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

1	AD	TDG-FEEDER-ERROR-FLAG		A	1.0	F	
---	----	-----------------------	--	---	-----	---	--

Flag indicating whether the feeder weight transaction was flagged in error by Update.
N - No; Y - Yes. Applies only to early MODCOMP data.

1	AE	TDG-FEED-CONSUMED-BY-CAGE		N	10.0	N	
---	----	---------------------------	--	---	------	---	--

The weight of food consumed by all animals in the cage during monitoring period, measured in decigrams. Calculated by subtracting used feeder weight from full feeder weight. Pertains only to MODCOMP tests. (Since 1989, consumption values have been dynamically calculated by reports themselves, due to design-related problems maintaining integrity of stored consumption values whose calculation was based on other stored values.)

1	AF	TDG-HOURS-FEEDER-IN-CAGE		N	5.0	N	
---	----	--------------------------	--	---	-----	---	--

Hours feeder was in the cage. Calculated by subtracting used feeder weight date and time from full feeder weight date and time. Pertains only to MODCOMP data. (Since 1989, hours-in-cage values have been dynamically calculated by reports themselves, due to design-related problems maintaining integrity of stored calculated values whose calculations were based on other stored values.)

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T L DB NAME	F	LENG S DE
1 AG TDG-AVG-FEED-PER-ANIMAL-HOUR	N	5.5 N
<p>The average amount of feed consumed by an animal, measured as decigrams/hour. Calculated by dividing feed consumption by animal hours-in-cage. Pertains only to MODCOMP data. (Since 1989, average feed per animal hour values have been dynamically calculated by reports themselves due to design-related problems maintaining integrity of stored calculated values whose calculations were based on other stored values.)</p>		
1 AH TDG-FEEDER-ANIMAL-HOURS-USED	N	5.0 N
<p>The sum of the hours each animal was in the cage during the monitoring period. Calculated by multiplying the hours the feeder was in the cage by the number of animals in the cage. Pertains only to MODCOMP data. (Since 1989, animal hour values have been dynamically calculated by the reports themselves due to design-related problems maintaining integrity of stored, calculated values whose calculations were based on other stored values.)</p>		
1 AI TDG-ACT-OVERRIDE-FLAG-OLD-FEED	A	1.0 F
<p>Flag indicating whether the scheduled used feeder weight action was overridden. 0 - No; 1 - Yes. Applies only to MODCOMP data.</p>		
1 AJ TDG-MANUAL-ENTRY-FLAG-OLD-FEED	A	1.0 F
<p>Flag indicating whether used</p>		

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
			feeder weight was recorded as part of a manual entry session. 0 - No; 1 - Yes.				
1	AK		TDG-REMOVAL-ACT-FLAG-OLD-FEED	A	1.0	F	
			Flag indicating whether the used feeder weight is a removal action. 0 - No; 1 - Yes. Pertains only to MODCOMP data.				
1	AL		TDG-BACK-FORWARD-FLAG-OLD-FEED	A	1.0	F	
			Flag indicating whether used feeder weight was recorded in scheduled order. 0 - Yes; 1 - No. Pertains only to MODCOMP data.				
1	AM		TDG-EXP-SELECTED-FLAG-OLD-FEED	A	1.0	F	
			Flag indicating whether used feeder weight was a scheduled action. 0 - Yes; 1 - No. Pertains only to MODCOMP data.				
1	AN		TDG-OLD-FEEDER-MIX-NUMBER	N	10.0	N	
			Legacy field with no meaning in TDMS.				
1	AO		TDG-OLD-FEEDER-WEIGHT	N	10.0	N	
			The weight of the used feeder, in decigrams. Pre-LDAS records where a weight was not recorded will show a value of 99999, which value is flagged by reports and displayed as "OVERRIDE"; post-LDAS data where a weight was not recorded will show a value of 0, and have a value assigned to TDG-OLD-FEEDER-STATUS.				
1	AP		TDG-OLD-FEEDER-WEIGHT-OVERRIDE	A	1.0	F	
			Flag indicating whether used feeder weight was with acceptable range. A = Within range; + = Above range; - = Below range. Pertains only to MODCOMP data.				

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
1	AQ	TDG-ACT-OVERRIDE-FLAG-NEW-FEED	A	1.0	F		
		Flag indicating whether the scheduled full feeder weight action was overridden. 0 - No; 1 - Yes. Applies only to MODCOMP data.					
1	AR	TDG-MANUAL-ENTRY-FLAG-NEW-FEED	A	1.0	F		
		Flag indicating whether the full feeder weight was recorded as part of manual entry session. 0 - No; 1 - Yes.					
1	AS	TDG-REMOVAL-ACT-FLAG-NEW-FEED	A	1.0	F		
		Flag indicating whether the full feeder weight is a removal action. 0 - No; 1 - Yes. Pertains only to MODCOMP data.					
1	AT	TDG-BACK-FORWARD-FLAG-NEW-FEED	A	1.0	F		
		Flag indicating whether full feeder weight was recorded in scheduled order. 0 - Yes; 1 - No. Pertains only to MODCOMP data.					
1	AU	TDG-EXP-SELECTED-FLAG-NEW-FEED	A	1.0	F		
		Flag indicating whether full feeder weight was a scheduled action. 0 - Yes; 1 - No. Pertains only to MODCOMP data.					
1	AV	TDG-NEW-FEEDER-MIX-NUMBER	N	10.0	N		
		Legacy field with no meaning in TDMS.					
1	AW	TDG-NEW-FEEDER-WEIGHT	N	10.0	N		
		Weight, in decigrams, of the full feeder.					
1	AX	TDG-NEW-FEEDER-WEIGHT-OVERRIDE	A	1.0	F		
		Flag indicating whether full feeder weight was with acceptable range. A = Within range; + = Above range; - = Below range. Pertains only to MODCOMP data.					
1	CL	TDG-WATER-TIME	N	5.0	N		
		Time, in TDMS TOX format, at					

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
			which the bottle weight was recorded.				
1	AZ	TDG-WATER-OPERATOR-ID		N	5.0	N	
			The operator who recorded the bottle weight. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	BA	TDG-WATER-ERROR-FLAG		A	1.0	F	
			Flag indicating whether bottle weight transaction was flagged as being in error by Update. N - No; Y - Yes. Applies only to early MODCOMP data.				
1	BB	TDG-WATER-CONSUMED-BY-CAGE		N	10.0	N	
			The weight of water consumed by all animals in the cage during monitoring period, measured in decigrams. Calculated by subtracting used bottle weight from full bottle weight. Pertains only to MODCOMP tests. (Since 1989, consumption values have been dynamically calculated by reports themselves, due to design-related problems maintaining integrity of stored consumption values whose calculation was based on other stored values.)				
1	BC	TDG-AVG-WATER-PER-ANIMAL-HOUR		N	5.5	N	
			The average amount of water consumed by an animal, measured as decigrams/hour. Calculated by dividing water consumption by animal hours-in-cage. Pertains only to MODCOMP tests. (Since 1989, average water per animal hour values have been dynamically calcul-				

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

ated by reports themselves due to design-related problems maintaining integrity of stored calculated values whose calculations were based on other stored values.)

1	BD	TDG-HOURS-BOTTLE-IN-CAGE	N	5.0	N	
---	----	--------------------------	---	-----	---	--

Hours bottle was in the cage. Calculated by subtracting used bottle weight date and time from full bottle weight date and time. Pertains only to MODCOMP data. (Since 1989, hours-in-cage values have been dynamically calculated by reports themselves, due to design-related problems maintaining stored calculated values whose calculations were based on other stored values.)

1	BE	TDG-BOTTLE-ANIMAL-HOURS-USED	N	5.0	N	
---	----	------------------------------	---	-----	---	--

The sum of the hours each animal was in the cage during the monitoring period. Calculated by multiply the hours the bottle was in the cage by the number of animals in the cage. Pertains only to MODCOMP data. (Since 1989, animal hour values have been dynamically calculated by the reports themselves due to design-related problems maintaining integrity of stored, calculated values whose calculations were based on other stored values.)

1	BF	TDG-ACT-OVERRIDE-FLAG-OLD-BOTT	A	1.0	F	
---	----	--------------------------------	---	-----	---	--

Flag indicating whether the scheduled used bottle weight action was overridden. 0 - No;

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB NAME	F	LENG	S	DE
		1 - Yes. Applies only to MOD- COMP data.				
1	BG	TDG-MANUAL-ENTRY-FLAG-OLD-BOTT	A	1.0	F	
		Flag indicating whether the used bottle weight was record- ed as part of a manual entry session. 0 - No; 1 - Yes.				
1	BH	TDG-REMOVAL-ACT-FLAG-OLD-BOTT	A	1.0	F	
		Flag indicating whether the used bottle weight is a remov- al action. 0 - No; 1 - Yes. Pertains only to MODCOMP data.				
1	BI	TDG-BACK-FORWARD-FLAG-OLD-BOTT	A	1.0	F	
		Flag indicating whether used bottle weight was recorded in scheduled order. 0 - Yes; 1 - No. Pertains only to MODCOMP data.				
1	BJ	TDG-EXP-SELECTED-FLAG-OLD-BOTT	A	1.0	F	
		Flag indicating whether used bottle weight was a scheduled action. 0 - Yes; 1 - No. Pertains only to MODCOMP data.				
1	BK	TDG-OLD-BOTTLE-MIX-NUMBER	N	10.0	N	
		Legacy field with no meaning in TDMS.				
1	BL	TDG-OLD-BOTTLE-WEIGHT	N	10.0	N	
		The weight of the used bottle, in decigrams. Pre-LDAS re- cords where a weight was not recorded will show a value of 99999, which value is flagged by reports and displayed as "OVERRIDE"; post-LDAS data where a weight was not record- ed will show a value of 0, and have a value assigned to TDG- OLD-BOTTLE-STATUS.				
1	BM	TDG-OLD-BOTTLE-WEIGHT-OVERRIDE	A	1.0	F	
		Flag indicating whether used bottle weight was with accep-				

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
			table range. A = Within range; + = Above range; - = Below range. Pertains only to MODCOMP data.				
1	BN		TDG-ACT-OVERRIDE-FLAG-NEW-BOTT	A	1.0	F	
			Flag indicating whether the scheduled full bottle weight action was overridden. 0 - No; 1 - Yes. Applies only to MOD- COMP data.				
1	BO		TDG-MANUAL-ENTRY-FLAG-NEW-BOTT	A	1.0	F	
			Flag indicating whether the full bottle weight was record- ed as part of manual entry session. 0 - No; 1 - Yes.				
1	BP		TDG-REMOVAL-ACT-FLAG-NEW-BOTT	A	1.0	F	
			Flag indicating whether the full bottle weight is a remov- al action. 0 - No; 1 - Yes. Pertains only to MODCOMP data.				
1	BQ		TDG-BACK-FORWARD-FLAG-NEW-BOTT	A	1.0	F	
			Flag indicating whether full bottle weight was recorded in scheduled order. 0 - Yes; 1 - No. Pertains only to MODCOMP data.				
1	BR		TDG-EXP-SELECTED-FLAG-NEW-BOTT	A	1.0	F	
			Flag indicating whether full bottle weight was a scheduled action. 0 - Yes; 1 - No. Pertains only to MODCOMP data.				
1	BS		TDG-NEW-BOTTLE-MIX-NUMBER	N	10.0	N	
			Legacy field with no meaning in TDMS.				
1	BT		TDG-NEW-BOTTLE-WEIGHT	N	10.0	N	
			Weight, in decigrams, of the full bottle.				

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB NAME	F	LENG	S	DE
1	BU	TDG-NEW-BOTTLE-WEIGHT-OVERRIDE	A	1.0	F	
		Flag indicating whether full bottle weight was with acceptable range. A = Within range; + = Above range; - = Below range. Pertains only to MODCOMP data.				
1	CM	TDG-OBSERVATION-TIME	N	5.0	N	
		Time, in TDMS TOX format, on which cage condition observations were made. Pertains only to MODCOMP tests.				
1	BW	TDG-OBSERVATION-OPERATOR-ID	N	5.0	N	
		The operator making cage condition observations. Pertains only to MODCOMP tests. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	BX	TDG-OBSERVATION-ERROR-FLAG	A	1.0	F	
		Flag indicating whether cage observation was flagged as being in error by Update. N - No; Y - Yes. Pertains only to MODCOMP data.				
1	BY	TDG-ACTION-OVERRIDE-FLAG-OBS	A	1.0	F	
		Flag indicating whether the scheduled cage observation action was overridden. 0 - No; 1 - Yes. Applies only to MODCOMP data.				
1	BZ	TDG-MANUAL-ENTRY-FLAG-OBS	A	1.0	F	
		Flag indicating whether cage observation was recorded as part of a manual entry session. 0 - No; 1 - Yes.				
1	CA	TDG-REMOVAL-ACTION-FLAG-OBS	A	1.0	F	
		Flag indicating whether the cage observation is a removal action. 0 - No; 1 - Yes. Pertains only to MODCOMP data.				

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB NAME	F	LENG	S	DE
1	CB	TDG-BACKWARD-FORWARD-FLAG-OBS	A	1.0	F	
		Flag indicating whether cage observation was recorded in scheduled order. 0 - Yes; 1 - No. Pertains only to MODCOMP data.				
1	CC	TDG-EXPECTED-SELECTED-FLAG-OBS	A	1.0	F	
		Flag indicating whether cage observation was a scheduled action. 0 - Yes; 1 - No. Pertains only to MODCOMP data.				
1	CE	TDG-MIS-ANIMAL-ID-TYPE	A	1.0	F	
		Indicates if animals in cage are identified by numeric or alphabetic characters. N - Numeric, A - Alphabetic. Only numeric is valid for TDMS. (Mis-identified animals was an early feature of MODCOMP that was abandoned due to design-related problems.)				
1	CN	TDG-NOTE-TIME	N	5.0	N	
		Time, in TDMS TOX format, at which the note was entered.				
1	CF	TDG-NOTE-OPERATOR-ID	N	5.0	N	
		The operator who recorded the note. Whe appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	CG	TDG-NOTE-ERROR-FLAG	A	1.0	F	
		Flag indicating whether the cage note transaction was flagged in error by Update. N - No; Y - Yes. Applies only to MODCOMP data.				
1	CH	TDG-NOTE-ACTION-NUMBER	N	5.0	N	
		Action with which the note is associated. Applies only to MODCOMP data. Foreign key to TBACTONC.				

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
1	CI	TDG-CAGE-NOTE		A	80.0	N	
		Free text describing, explaining, or commenting on the cage, or observations of the cage.					
M 1	CJ	TDG-CAGE-OBSERVATION		N	5.0	N	
		Cage condition that are being observed. Pertains only to MODCOMP tests. Foreign key to TBCONDTC.					
M 1	CK	TDG-MIS-ANIMAL-IDENTIFICATION		N	5.0	N	
		Mis-identified animal that has been incorrectly removed from the census. (Mis-identified animals were an early feature of MODCOMP that was abandoned due to design-related problems.)					
1	CO	TDG-MISID-ANIMAL-ERROR-FLAG		A	1.0	F	
		Flag indicating whether mis-identified animal transaction was flagged by Update as being in error. (Mis-identified animals were an early feature of MODCOMP that was abandoned due to design-related problems.)					
1	CP	TDG-MISID-ANIMAL-TIME		N	5.0	N	
		Time, in TDMS TOX format, at which mis-identified animal was recorded. (Mis-identified animals were an early feature of MODCOMP that was abandoned due to design-related problems.)					
1	CQ	TDG-MISID-ANIMAL-OPERATOR-ID		N	5.0	N	
		ID of operator who recorded the mis-identified animal					

FILE...: TD-CAGE-DATA

TYPE...: USER VIEW

FILE-NR: 107

PRIMARY SEQUENCE FIELD: TDG-EXP-TEST-CAGE-DATE-DE

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

transaction. (Mis-identified animals were an early feature of MODCOMP that was abandoned due to design-related problems.) When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	EA	TDG-OLD-BOTTLE-STATUS	P	5.0	N		
---	----	-----------------------	---	-----	---	--	--

Code for weight status that describes why the used weight was not recorded. Foreign key to TBWGTSTATUS.

1	EB	TDG-OLD-FEEDER-STATUS	P	5.0	N		
---	----	-----------------------	---	-----	---	--	--

Weight status that describes the used bottle weight was not recorded. Foreign key to TBWGTSTATUS.

1	EC	TDG-CAGE-NOTE-MANUAL-FLAG	A	1.0	N		
---	----	---------------------------	---	-----	---	--	--

Flag indicating whether the cage note was recorded as part of a manual entry session. 0 - No; 1 - Yes.

1	DA	TDG-FEEDER-SUPERKEY	N	21.0	SP		
---	----	---------------------	---	------	----	--	--

SOURCE FIELD(S) --- -START- --END-
EXP-TEST-CAGE-DATE 1 16
FEEDER-TIME 1 5

Super-descriptor used to select feeder weight records.

1	DB	TDG-WATER-SUPERKEY	N	21.0	SP		
---	----	--------------------	---	------	----	--	--

SOURCE FIELD(S) --- -START- --END-
EXP-TEST-CAGE-DATE 1 16
WATER-TIME 1 5

Super-descriptor used to select bottle weights.

1	DC	TDG-NOTES-SUPERKEY	A	21.0	SP		
---	----	--------------------	---	------	----	--	--

SOURCE FIELD(S) --- -START- --END-
EXP-TEST-CAGE-DATE 1 16
NOTE-TIME 1 5

TD-CORRECTION-REASON

File contains MODCOMP corrections. In-life information is reported in E13; pathology information is not reported. Since LDAS's introduction, no new records are added to this file.

FILE...: TD-CORRECTION-REASON

TYPE...: USER VIEW

FILE-NR: 109

PRIMARY SEQUENCE FIELD: TDJ-KEY-VALUE-DE

T	L	DB	NAME	F	LENG	S	DE
1	AA		TDJ-KEY-VALUE-DE	A	29.0		DE
			Composite field whose elements are the 2-byte segment ID and the primary key of the record that was modified, which can be from 11 to 27 characters. The used characters from positions 3 to 27 form a foreign key to the view associated with the segment ID.				
1	AC		TDJ-CORRECTION-SEQUENCE-NUMBER	N	2.0		
			Number that identifies the corrections of a particular key.				
M 1	AD		TDJ-CORRECTION-REASON	A	80.0	N	
			Free text describing the corresponds-rection.				
1	AE		TDJ-CORRECTION-DATE	N	5.0		
			Date, in TDMS TOX format, on which the correction was made.				
1	AF		TDJ-CORRECTION-TIME	N	5.0		
			Time, in TDMS TOX format, at which the correction was made.				
1	AG		TDJ-CORRECTION-REASON-OP-ID	N	5.0		
			The operator who made the corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	AH		TDJ-EXP-TEST-CAGE-DE	A	11.0		SB
			SOURCE FIELD(S) --- -START- --END-				
			CORRECTION-KEY		3		13
			Key used to select data from the file by E13. Foreign key to TD-CAGE.				

TD-MICRO-ANIMAL-NOTE

File contains pathology animal notes. Only used in P14.

FILE...: TD-MICRO-ANIMAL-NOTE**TYPE...: USER VIEW****FILE-NR: 110****PRIMARY SEQUENCE FIELD: TDM-EXP-TST-ANUM-PATHT-DATE-DE**

T	L	DB	NAME	F	LENG	S	DE
1	AA	TDM-EXP-TST-ANUM-PATHT-DATE-DE	A	22.0	DE		
		Composite field whose elements are the 7-digit test number, the 4-digit number of the original cage to which the animal was allocated, the 5-digit animal number assigned at allocation, the letter 'A' for pathologist type, and the 5-digit date, in TDMS TOX format, on which the note was recorded. The first 16 characters form a foreign key to TD-ANIMAL.					
1	AB	TDM-MICRO-ANIMAL-NOTE-SEQUENCE	N	5.0			
		Time, in TDMS TOX format, at which the note was recorded. Not used in report software.					
1	AC	TDM-OPERATOR-IDENTIFICATION	N	5.0			
		Operator who recorded the note. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.					
1	AD	TDM-ERROR-FLAG	A	1.0	F		
		Flag indicating whether micro-animal note transaction was flagged by Update as being in error. N - No; Y - Yes. Applies only to early MODCOMP data.					
1	AF	TDM-PATHOLOGIST-IDENTIFICATION	N	5.0			
		The pathologist who was assigned to the test when the					

FILE...: TD-MICRO-ANIMAL-NOTE

TYPE...: USER VIEW

FILE-NR: 110

PRIMARY SEQUENCE FIELD: TDM-EXP-TST-ANUM-PATHT-DATE-DE

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

note was recorded. Not used
in report software. When ap-
pended to the appropriate faci-
lity code, a foreign key to
TB_OPERATOR_ID.

M	1	AE	TDM-MICRO-ANIMAL-NOTE	A	80.0	N
---	---	----	-----------------------	---	------	---

Free text describing, comment-
ing on, or explaining the
animal, or any observations of
the animal.

1	AG	TDM-COD-FLAG		A	1.0	F
---	----	--------------	--	---	-----	---

Introduced with LDAS, but nev-
er implemented.

TD-ORGAN

View of organ (tissue) status and accountable site status segments. (The last only applies to data not converted for LDAS.)

FILE...: TD-ORGAN

TYPE...: USER VIEW

FILE-NR: 111

PRIMARY SEQUENCE FIELD: TDO-EXP-TST-ANUM-ORG-PATHT-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AA	TDO-EXP-TST-ANUM-ORG-PATHT-UQ	A	22.0	UQ		
<p>Composite field whose elements are the 7-digit test number, the 4-digit number of the original cage to which the animal was allocated, the 5-digit animal number assigned at allocation, the organ being object-served, and the pathologist type, which is always 'A'. The test, cage, and animal number portions form a foreign key to TD-ANIMAL; the organ is a foreign key to TBPCT.</p>							
1	AB	TDO-ORGAN-STATUS-DATE	N	5.0	N		
<p>Date, in TDMS TOX format, on which the status was recorded. Not used in report software.</p>							
1	AC	TDO-ORGAN-STATUS-TIME	N	5.0	N		
<p>Time, in TDMS TOX format, at which the status was recorded.</p>							
1	AD	TDO-ORGAN-STATUS-OPERATOR-ID	N	5.0	N		
<p>The operator who recorded the status. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.</p>							
1	AE	TDO-ORGAN-STATUS-ERROR-FLAG	A	1.0	F		
<p>Flag indicating whether organ status transaction was flagged by Update as being in error. N - No; Y - Yes. Applies only to early MODCOMP data.</p>							

FILE...: TD-ORGAN

TYPE...: USER VIEW

FILE-NR: 111

PRIMARY SEQUENCE FIELD: TDO-EXP-TST-ANUM-ORG-PATHT-UQ

T	L	DB NAME	F	LENG	S	DE
1	AF	TDO-ORGAN-STATUS	N	5.0	N	
		The status code of the organ. Foreign key to TBORGSTA.				
1	AG	TDO-ORGAN-STAT-PATHOLOGIST-ID	N	5.0	N	
		The pathologist assigned to the test when the status was recorded. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ ID.				
P 1	AH	TDO-SITE-STATUS-INFORMATION				
		Repeating group for account- able site status information. This group only applies for tests not converted for LDAS data collection. Also, only applies to the following or- gans: adrenal gland, large and small intestines, stomach, lymph node, and skin.				
2	AK	TDO-SITE-STATUS-DATE	N	5.0	N	
		Date, in TDMS TOX format, on which the accountable site status was recorded. Not used in report software.				
2	AL	TDO-SITE-STATUS-TIME	N	5.0	N	
		Time, in TDMS TOX format, at which the accountable site status was recorded.				
2	AM	TDO-SITE-STATUS-PATHOLOGIST-ID	N	5.0	N	
		The pathologist assigned to the test when the accountable site status was recorded. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				

FILE...: TD-ORGAN

TYPE...: USER VIEW

FILE-NR: 111

PRIMARY SEQUENCE FIELD: TDO-EXP-TST-ANUM-ORG-PATHT-UQ

T	L	DB NAME	F	LENG	S	DE
2	AN	TDO-SITE-STATUS-OPERATOR-ID	N	5.0	N	
		The operator who recorded the accountable site status. Not used in report software. When appended to appropriate facility code, a foreign key to TB_OPERATOR_ID.				
2	AO	TDO-SITE-STATUS-ERROR-FLAG	A	1.0	N	
		Flag indicating whether site status transaction was flagged by Update as being in error. N - No; Y - Yes.				
2	AI	TDO-SITE-CODE	N	5.0	N	
		The accountable site to which the status applies. A foreign key to TBPCT.				
2	AJ	TDO-SITE-STATUS	N	5.0	N	
		Organ status assigned to the accountable site. Foreign key to TBORGSTA.				

TD-ORGAN-DATA

View of microscopic observation, organ note, and organ stain segments,
the last of which is inactive.

FILE...: TD-ORGAN-DATA

TYPE...: USER VIEW

FILE-NR: 108

PRIMARY SEQUENCE FIELD: TDH-ORGAN-OBSERVATION-ID-DE

T	L	DB	NAME	F	LENG	S	DE
1	AA	TDH-ORGAN-OBSERVATION-ID-DE	A	27.0	DE		
<p>Composite field whose elements are the 7-digit test number, the 4-digit number of the original cage to which the animal was allocated, the 5-digit animal number assigned at allocation, the 5-digit code for the organ being observed, the letter 'A' for pathologist type, and the 5-digit date, in TDMS TOX format, on which the observation was recorded. The test, cage, and animal number portions form a foreign key to TD-ANIMAL; the organ code is a foreign key to TBPCT.</p>							
1	AB	TDH-MICRO-OBSERVATION-SEQUENCE	N	5.0	N		
<p>The time, in TDMS TOX format, at which the observation was made. Not used in report software.</p>							
1	AC	TDH-MICRO-OBS-OPERATOR-ID	N	5.0	N		
<p>The ID of the operator who recorded the observation. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.</p>							
1	AO	TDH-MICRO-OBS-PATHOLOGIST-ID	N	5.0	N		
<p>The operator ID of the pathologist assigned to the test when the observation was recorded. Not used in report</p>							

FILE...: TD-ORGAN-DATA

TYPE...: USER VIEW

FILE-NR: 108

PRIMARY SEQUENCE FIELD: TDH-ORGAN-OBSERVATION-ID-DE

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	AF	TDH-OBSERVATION-MORPHOLOGY-DE	N	5.0	N		
---	----	-------------------------------	---	-----	---	--	--

The code for the morphology that was observed. A foreign key to TBPCT.

1	AE	TDH-PRIMARY-SITE	N	5.0	N		
---	----	------------------	---	-----	---	--	--

The primary site for a morphology that has metastasized to the organ being observed. (If blank there is no metastasis.) Foreign key to TBPCT.

1	AW	TDH-NEOPLASTIC-FLAG	A	1.0	F		
---	----	---------------------	---	-----	---	--	--

Flag indicating whether the observed morphology is neoplastic. Y - Yes; N - No.

M	1	AG	TDH-OBSERVATION-QUALIFIER	N	5.0	N	
---	---	----	---------------------------	---	-----	---	--

The qualifiers that further describe the observed morphology. Up to 4 may be designated. Foreign key to TBPCT.

M	1	AH	TDH-OBSERVATION-SITE	N	5.0	N	
---	---	----	----------------------	---	-----	---	--

The sites on the organ being observed where the morphology detected. Up to 3 may be designated. Foreign key to TBPCT

1	BA	TDH-COD-FLAG	A	1.0	F		
---	----	--------------	---	-----	---	--	--

Flag indicating the observed morphology's role in the animal's death. The field will only have a value if the animal was removed due to natural death or moribund sacrifice.
P - Primary cause-of-death;
C - Contributory cause-of-

FILE...: TD-ORGAN-DATA

TYPE...: USER VIEW

FILE-NR: 108

PRIMARY SEQUENCE FIELD: TDH-ORGAN-OBSERVATION-ID-DE

T	L	DB	NAME	F	LENG	S	DE
			death.				
1	AV	TDH-MORPH-GENERIC-TYPE		N	5.0	N	
			The code for the generic morphology to which the object-served morphology is associated. Applies only to MODCOMP tests. Foreign key to TBPCT.				
1	AD	TDH-MICRO-OBS-ERROR-FLAG		A	1.0	F	
			Flag indicating whether micro observation transaction was flagged by Update as being in error. N - No; Y - Yes. Applies only to MODCOMP tests.				
1	AX	TDH-TRACE-LESION-NUMBER		N	5.0	N	
			A sequentially assigned gross lesion number which is used to associate a gross lesion with a microscopic observation. Never used in TDMS.				
1	AI	TDH-ORGAN-NOTE-OPERATOR-ID		N	5.0	N	
			ID of operator recording the note. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	AP	TDH-ORGAN-NOTE-PATHOLOGIST-ID		N	5.0	N	
			Operator ID of the pathologist assigned to the test when the note was recorded. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	AR	TDH-ORGAN-NOTE-SEQUENCE		N	5.0	N	
			The time, in TDMS TOX format, at which the note was recorded. Not used in report software.				
M 1	AK	TDH-ORGAN-NOTE		A	80.0	N	
			Free text describing, comment-				

FILE...: TD-ORGAN-DATA

TYPE...: USER VIEW

FILE-NR: 108

PRIMARY SEQUENCE FIELD: TDH-ORGAN-OBSERVATION-ID-DE

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

ing or explaining, the organ
or observations of it. Also,
may contain TGL information.
TGL information can be identi-
fied by character positions
22-27 having the value
"TGLs =".

1	AJ	TDH-ORGAN-NOTE-ERROR-FLAG	A	1.0	F		
---	----	---------------------------	---	-----	---	--	--

Flag indicating whether the
organ note transaction was
flagged by Update as being in
error. N - No; Y - Yes. Ap-
plies only to early MODCOMP
data.

1	AL	TDH-ORGAN-STAIN-OPERATOR-ID	N	5.0	N		
---	----	-----------------------------	---	-----	---	--	--

The ID of the operator who
created the organ stain trans-
action. Pertains only to MOD-
COMP tests. When appended to
the appropriate facility code,
a foreign key to TB_OPERATOR_ID.

1	AQ	TDH-ORGAN-STAIN-PATHOLOGIST-ID	N	5.0	N		
---	----	--------------------------------	---	-----	---	--	--

Operator ID of the pathologist
assigned to test. Pertains
only to MODCOMP tests. When
appended to the appropriate
facility code, a foreign key
to TB_OPERATOR_ID.

1	AS	TDH-ORGAN-STAIN-SEQUENCE	N	5.0	N		
---	----	--------------------------	---	-----	---	--	--

Time, in TDMS TOX format, at
which the organ stain transac-
tion was created. Pertains
only to MODCOMP tests.

1	AM	TDH-ORGAN-STAIN-ERROR-FLAG	A	1.0	F		
---	----	----------------------------	---	-----	---	--	--

Flag indicating whether the
organ stain transaction was
flagged by Update as being in
error. N - No; Y - Yes. Per-
tains only to early MODCOMP
data.

P	1	AT	TDH-ORGAN-STAIN-INFORMATION				
---	---	----	-----------------------------	--	--	--	--

FILE...: TD-ORGAN-DATA

TYPE...: USER VIEW

FILE-NR: 108

PRIMARY SEQUENCE FIELD: TDH-ORGAN-OBSERVATION-ID-DE

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

Repeating group embracing
each stain being observed.
Recording of stain observation
applies only to MODCOMP tests.
Recording of such data was
discouraged by NTP. No re-
ports access stain informa-

2	AN	TDH-ORGAN-STAIN	N	5.0	N	
---	----	-----------------	---	-----	---	--

The code for the stain
being observed. Pertains on-
ly to MODCOMP tests. Foreign
key to TBSTAIN.

2	AU	TDH-ORGAN-STAIN-CONTRIBUT-FLAG	A	1.0	N	
---	----	--------------------------------	---	-----	---	--

Flag indicating the results of
the stain. + = Positive; - =
Negative; ? = Non-contributory
Applies only to MODCOMP tests.

1	UB	TDH-EXP-TEST-ORGAN-SP	A	12.0	SP	
---	----	-----------------------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-
>>> SOURCE FIELD(S) NOT SPECIFIED <<<

Used to select all observa-
tions for an organ on a test

TD-TEST

View of test information.

FILE...: TD-TEST**TYPE...: USER VIEW****FILE-NR: 112****PRIMARY SEQUENCE FIELD: TDT-EXPERIMENT-TEST-DE**

T	L	DB	NAME	F	LENG	S	DE
1	AA	TDT-EXPERIMENT-TEST-DE		N	7.0	UQ	
		Composite field whose elements are the 5-digit experiment (or study) number, and the 2-digit test number. When appended to TDT-AGENCY-NUMBER, a foreign key to PAS_MASTER.					
1	AC	TDT-AGENCY-NUMBER		N	2.0		
		Then agency. A foreign key to TBAGENCY. Not used in report software.					
1	AD	TDT-TEST-TYPE		N	5.0		
		Test type code, as assigned in test protocol. Foreign key to TBTESTYP.					

TD-TEST-NOTE

View of test notes. Used by E10 only.

FILE...: TD-TEST-NOTE**TYPE...: USER VIEW****FILE-NR: 114****PRIMARY SEQUENCE FIELD: TDS-EXP-TEST-DATE-TIME-UQ**

T L DB NAME	F	LENG S DE
1 AB TDS-EXP-TEST-DATE-TIME-UQ	A	17.0 N UQ
Composite field whose elements are the 7-digit test number, the 5-digit date, in TDMS TOX format, on which the note was recorded, and the 5-digit time, in TDMS TOX format, at which the note was recorded. The test number portions forms a foreign key to TD-TEST.		
1 AC TDS-OPER-NUM	P	5.0 F
The ID of the operator who recorded the test note. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.		
1 AD TDS-NOTE-TEXT	A	160.0 N
Free text describing, commenting on, or explaining the test or some observation of the test.		
1 BB TDS-TEST	A	7.0 SP
SOURCE FIELD(S) --- -START- --END-		
TDS-EXP-TEST-DATE-TIME-UQ	1	7

TD-TREATMENT

View of pre-LDAS protocol-required tissues. Only used by the P14 for those tests that were not converted to LDAS. Replaced by the PAS protocol-required tissues segment.

FILE...: TD-TREATMENT

TYPE...: USER VIEW

FILE-NR: 100

PRIMARY SEQUENCE FIELD: TDI-EXP-TEST-TREATMENT-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AA		TDI-EXP-TEST-TREATMENT-UQ	N	10.0		UQ
Composite field whose elements are the 7-digit test number and the 3-digit treatment group ID number. The test number portion forms a foreign key to TD-TEST. When the test number and the treatment number (converted to a 2-byte binary) are appended to the appropriate agency number, a foreign key to PAS_TREATMENT_SEG.							
1	AB		TDI-TREATMENT-CODE	A	1.0		F
Code for treatment role assigned to the treatment group in test protocol. Not used in report software. Foreign key to TBTRTROLE.							
M 1	AD		TDI-PROTOCOL-REQ-TISSUE-CODE	N	5.0		N
Codes for the organs that are required to be assigned an organ status. Used only in P14. Foreign key to TBPCT.							
1	AE		TDI-PATHOLOGIST-IDENTIFICATION	N	5.0		
Operator ID of the pathologist assigned to the test when the PRTs for the treatment group are assigned. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.							

FILE...: TD-TREATMENT

TYPE...: USER VIEW

FILE-NR: 100

PRIMARY SEQUENCE FIELD: TDI-EXP-TEST-TREATMENT-UQ

T	L	DB NAME	F	LENG	S	DE
1	AF	TDI-OPERATOR-IDENTIFICATION	N	5.0		
		ID of the operator who assigns the PRTs to the treatment group. Not used in report software. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	AG	TDI-DATE	N	5.0		
		Date, in TDMS TOX format, on which the PRTs were assigned. Not used in report software.				
1	AH	TDI-TREATMENT-SEQUENCE	N	5.0		
		Time, in TDMS TOX format, at which the PRTs were assigned. Not used in report software.				
1	AI	TDI-ERROR-FLAG	A	1.0	F	
		Flag indicating whether treatment transaction was flagged by Update as being in error. N - No; Y - Yes. Applies only to early MODCOMP data.				

TD-TREATMENT-NOTE

View of treatment notes. Used by E10 only.

FILE...: TD-TREATMENT-NOTE**TYPE...: USER VIEW****FILE-NR: 114****PRIMARY SEQUENCE FIELD: TDT-EXP-TEST-TRT-DATE-TIME-UQ**

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

1	AA		TDT-EXP-TEST-TRT-DATE-TIME-UQ	A	20.0	N	UQ
---	----	--	-------------------------------	---	------	---	----

Composite field whose elements are the 7-digit test number, the 3-digit treatment group ID, the 5-digit date, in TDMS TOX format, on which the note was recorded, and the 5-digit time, in TDMS TOX format, at which the note was recorded.

The test number portions forms a foreign key to TD-TEST.

When the the treatment group is converted to a 2-byte binary and append to the agency number and test number, a foreign key to PAS_TREAT-MNT_SEG.

1	AC		TDT-OPER-NUM	P	5.0	F	
---	----	--	--------------	---	-----	---	--

The ID of the operator who recorded the treatment group note. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	AD		TDT-NOTE-TEXT	A	160.0	N	
---	----	--	---------------	---	-------	---	--

Free text describing, commenting on, or explaining the test or some observation of the treatment group.

4.3 Transaction Views

ANIMAL-MICRO-NOTES-TRANSACTION

Logical layout of pathology animal notes transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-MICRO-ANIMAL-NOTE file. If transactions are corrections or audit changes, then Update uses the information to add records to the ANIMAL-MICRO-NOTES-CORRECTION file.

FILE...: ANIMAL-MICRO-NOTES-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
G	1	AC HEADER				
		Group embracing fields that are common to most segments.				
2	AA	SEGMENT-ID	A	2.0	F	
		Two-letter segment identifier. Not referenced by Update; informational only. Set to MA. Used to form RECEIVE-UPDATE-KEY.				
2	AB	MICRO-OR-ECS-FLAG	A	1.0	F	
		Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE	TOX-MICRO-DATE	B	2.0	F	
		The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF	TOX-MICRO-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				

FILE...: ANIMAL-MICRO-NOTES-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	AG	MICRO-TRANS-SET	B	2.0	F	
		The number of the transac- tion's source set. Not refer- enced by Update.				
2	AI	TOX-MAINFRAME-DATE	B	2.0	F	
		The date, in TDMS TOX date format, on which the source set was received. Not refer- enced by Update. Used to form RECEIVE-UPDATE-KEY.				
2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set was received. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY.				
2	AL	AGENCY-NUMBER	B	2.0	F	
		Not referenced by Update. Foreign key to TBAGENCY.				
2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
		When the transaction is a corresponds- rection or an audit change, this indicates the responsible operator. Loaded by Update into the MAC-CORRECTION- OPERATOR-ID field of the ANI- MAL-MICRO-ANIMAL-NOTES-CORREC- TION file. When appended to DATA-FACILITY-NUMBER, a for- eign key to TB_OPERATOR_ID.				
2	AS	SOFTWARE-VERSION	B	2.0	F	
		The version of LDAS under which the source set was generated. Not referenced by Update.				
2	AT	TOX-DATA-DATE	B	2.0	F	
		The date, in TDMS TOX date format, on which the trans- action was created. Used by Update in the formation of the field TDM-EXP-TST-ANUM-PATHT-				

FILE...: ANIMAL-MICRO-NOTES-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		DATE-DE of the TD-MICRO-ANI-MAL-NOTES file and the MAC-MICRO-ANI-NOTE-CORR-KEY-UQ field of the ANIMAL-MICRO-NOTES-CORRECTION file. Also used to form RECEIVE-UPDATE-KEY.				
2	AU	TOX-DATA-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, at which the transaction was created. Loaded by Update into the TDM-MICRO-ANIMAL-NOTE-SEQUENCE field of the TD-MICRO-ANIMAL-NOTES file. For corrections or audit changes, used to form MAC-MICRO-ANI-NOTE-CORR-KEY-UQ of the ANIMAL-MICRO-NOTES-CORRECTION file. Also used to form RECEIVE-UPDATE-KEY.				
2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
		The facility associated with the machine on which the source set was created. Not referenced by Update. A foreign key to TBFACLTy.				
2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
		The operator ID for the pathologist who is assigned to the test. Loaded by Update into the TDM-PATHOLOGIST-IDENTIFICATION field of the TD-MICRO-ANIMAL-NOTES file. If an add with audit, the value is loaded into the MAC-PATHOLOGIST-ID field of the ANIMAL-MICRO-NOTES-CORRECTION file. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				

FILE...: ANIMAL-MICRO-NOTES-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	AX	PROCESS-FLAG	A	2.0	F	
Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.						
2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 14.						
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 3.						
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
Indicates sub-system. Set to 3 for pathology.						
2	A3	DATA-BASE-KEY	A	14.0	F	
The carcass ID (CID) for the animal to which the note applies. Foreign key to TD-ANIMAL. If a correction or audit transaction, then Update uses values from the corresponding TD-ANIMAL master record to form MICRO-ANI-NOTE-CORR-KEY-UQ in ANIMAL-MICRO-NOTES-CORRECTION.						
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
This field is parsed into 10 2-byte binary sub-fields. Only the first three are used. The first is set to 80; the second is for the pathologist type and is set to "A"; the third is set to 30. Update						

FILE...: ANIMAL-MICRO-NOTES-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

						uses the pathologist type value in forming the TDM-EXP-TST-ANUM-PATHT-DATE-DE field of the TD-MICRO-ANIMAL-NOTES file. If a correction or audit transaction, Update uses the value to form MAC-MICRO-ANI-NOTE-CORR-KEY-UQ in ANIMAL-MICRO-NOTES-CORRECTION.
G 1	D5	MICRO-ANIMAL-NOTES-DATA				Group embracing segment-specific fields of pathology animal note transactions.
2	D6	MICRO-ANIMAL-NOTES-OPERATOR	B	2.0	N	ID of operator recording the note. Loaded by Update into the TDM-OPERATOR-IDENTIFICATION field of the TD-MICRO-ANIMAL-NOTE file. If an add with audit, loaded by Update into MAC-OPERATOR-ID of ANIMAL-MICRO-NOTES-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.
M 2	D8	MICRO-ANIMAL-NOTE	A	80.0	N	Free text describing a microscopically evaluated animal. Loaded by Update into the TDM-MICRO-ANIMAL-NOTE field of the TD-MICRO-ANIMAL-NOTE file. If add with audit, loaded by Update into MAC-MICRO-ANIMAL-NOTE of ANIMAL-MICRO-NOTES-CORRECTION.
2	D4	CAUSE-OF-DEATH-FLAG	A	1.0	N	
G 1	EP	CORRECTION-PORTION				Group embracing segment-common correction fields.
2	EQ	CORRECTION-DATE	P	5.0	N	Date, in TDMS TOX date format,

FILE...: ANIMAL-MICRO-NOTES-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
on which correction was made. Used by Update in forming MAC-MICRO-ANI-NOTE-CORR-KEY-UQ for ANIMAL-MICRO-NOTE-CORRECTION.		
2 ER CORRECTION-TIME	P	5.0 N
Time of day, in TDMS TOX time format, at which correction was made. Used by Update in forming MAC-MICRO-ANI-NOTE-CORR-KEY-UQ of ANIMAL-MICRO-NOTE-CORRECTION.		
2 ES CORRECTION-REASON	A	2.0 N
Correction reason. Loaded by Update into MAC-CORRECTION-REASON of ANIMAL-MICRO-NOTE-CORRECTION. Foreign key to TBCORREASON.		
2 ET CORRECTION-NOTE	A	79.0 N
Free text describing the corresponds-rection. Loaded by Update into MAC-CORRECTION-NOTE in ANIMAL-MICRO-NOTE-CORRECTION.		
2 EU CORRECTION-TYPE	A	1.0 N
Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into MAC-CORRECTION-TYPE in ANIMAL-MICRO-NOTE-CORRECTION.		
2 FD CORRECTION-KEY	A	37.0 N DE
If either the date portion of TDM-EXP-TST-ANUM-PATHT-DATE-DE or TDM-MICRO-ANIMAL-NOTE-SEQUENCE is being changed, then this field contains the original value of both, along with CORRECTION-DATE and CORRECTION-TIME. Update loads it into MAC-MICRO-ANI-NTE-COR-KEY-LNK-UQ of ANIMAL-MICRO-NOTE-CORRECTION. First 22 characters are a foreign key		

FILE...: ANIMAL-MICRO-NOTES-TRANSACTION
TYPE...: USER VIEW
FILE-NR: 86
PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

to TD-MICRO-ANIMAL-NOTE.

1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
---	----	--------------------	---	------	---	----

Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDM-EXP-TST-ANUM-PATH-DATE-DE, 5 spaces, TOX-DATA-DATE, TOX-DATA-TIME, TOX-MAINFRAME-DATE, TOX-MAINFRAME-TIME, and a 4-digit sequence number for the transaction in its segment within the source set.

ANIMAL-NOTE-TRANSACTION

Logical layout of animal note transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, or delete records to the TD-ANIMAL-DATA file. For corrections, Update also uses the information to add records to ANIMAL-NOTE-CORRECTION.

FILE...: ANIMAL-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
G 1 AC HEADER		
Group embracing fields that are common to most segments.		
2 AA SEGMENT-ID	A	2.0 F
Two-letter segment identifier. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY. Set to AN.		
2 AB MICRO-OR-ECS-FLAG	A	1.0 F
Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.		
2 AE TOX-MICRO-DATE	B	2.0 F
The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.		
2 AF TOX-MICRO-TIME	B	2.0 F
The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.		
2 AG MICRO-TRANS-SET	B	2.0 F
The number of the transaction's source set. Not referenced by Update.		
2 AI TOX-MAINFRAME-DATE	B	2.0 F
The date, in TDMS TOX date		

FILE...: ANIMAL-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

format, on which the source set was received. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
---	----	--------------------	---	-----	---	--

The time of day, in TDMS TOX time format, on which the source set was received. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY.

2	AL	AGENCY-NUMBER	B	2.0	F	
---	----	---------------	---	-----	---	--

Foreign key to TBAGNECY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
---	----	-----------------------	---	-----	---	--

When the transaction is a correspondence, this indicates the operator who made the correction and is loaded by Update into ANC-CORRECTION-OPERATOR of ANIMAL-NOTE-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.

2	AS	SOFTWARE-VERSION	B	2.0	F	
---	----	------------------	---	-----	---	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F	
---	----	---------------	---	-----	---	--

The date, in TDMS TOX date format, on which the transaction was created. Used by Update to form the TDD-EXP-TEST-ANIMAL-DATE-DE value that will be added to, or selected for modification from, TD-ANIMAL-DATA. If a correction,

FILE...: ANIMAL-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

also used to form ANC-ANIMAL-NOTE-CORR-KEY-UQ of ANIMAL-NOTE-CORRECTION.

2	AU	TOX-DATA-TIME	B	2.0	F	
---	----	---------------	---	-----	---	--

The time of day, in TDMS TOX time format, at which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Loaded by Update into TDD-NOTE-TIME of TD-ANIMAL-DATA. If a corresponds-rection, also used to form ANC-ANIMAL-NOTE-CORR-KEY-UQ of ANIMAL-NOTE-CORRECTION.

2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
---	----	----------------------	---	-----	---	--

The facility associated with the machine on which the source set was created. Not referenced by Update. A foreign key to TBFACLT.

2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
---	----	-------------------------	---	-----	---	--

For an original addition, the the operator who created the transaction. Loaded by Update TDD-NOTE-OPERATOR-ID, and, if an add with audit, also into ANC-NOTE-OPERATOR-ID of ANIMAL-NOTE-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.

2	AX	PROCESS-FLAG	A	2.0	F	
---	----	--------------	---	-----	---	--

Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indica-

FILE...: ANIMAL-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		ting an original transaction, or 'A', indicating an audit transaction.				
2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
		The number of digits being used in the DATA-BASE-KEY field. Not referenced by Up- date. Set to 13.				
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
		The number of two-byte binary segments of the SEGMENT-TREE- KEY field being used. Set to 3.				
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
		Indicates sub-system. Set to 2 for in-life.				
2	A3	DATA-BASE-KEY	A	14.0	F	
		Indicates the number of the agency-test-cage where the animal being observed is lo- cated. (AATTTTTTCCCC). Used by Update to form the TDD-EXP- TEST-ANIMAL-DATE-DE value that will be added to, or selected for modification from, TD-ANI- MAL-DATA. Also, if a correc- , used to form ANC-ANIMAL- NOTE-CORR-KEY-UQ of ANIMAL- NOTE-CORRECTION. The test- cage portion is a foreign key to TD-CAGE.				
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
		This field is parsed into 10 2-byte binary sub-fields. Only 3 are used for animal note transactions. The first is set to 10; the second is set to the number of the animal- mal being observed; the third is set to 60. Update uses the animal number to form				

FILE...: ANIMAL-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

TDD-EXP-TEST-ANIMAL-DATE-DE
for TD-ANIMAL-DATA, and, if a
correction, to form ANC-ANI-
MAL-NOTE-CORR-KEY-UQ of ANI-
MAL-NOTE-CORRECTION.

G 1 CT ANIMAL-NOTE-DATA

Group embracing segment-speci-
fic fields of in-life animal
notes.

2	CV	ANIMAL-NOTE	A	80.0	N
---	----	-------------	---	------	---

Free text entered by operator
modifying or explaining animal
data. Loaded by Update into
the TDD-ANIMAL-NOTE of the
TD-ANIMAL-DATA file. If an
add-with-audit transaction,
loaded by Update into ANC-ANI-
MAL-NOTE of ANIMAL-NOTE-COR-
RECTION.

2	CW	ANIMAL-NOTE-MANUAL-FLAG	A	1.0	N
---	----	-------------------------	---	-----	---

Flag indicating whether note
was entered as part of a man-
ual entry session. 0 - No;
1 - Yes. Loaded by Update
into the TDD-MANUAL-ENTRY-
FLAG-NOTE field of the TD-
ANIMAL-DATA file. If an add-
with-audit transaction, load-
ed by Update into ANC-NOTE-
MANUAL-FLAG of ANIMAL-NOTE-
CORRECTION.

G 1 EP CORRECTION-PORITION

Group embracing segment-common
correction fields.

2	EQ	CORRECTION-DATE	P	5.0	N
---	----	-----------------	---	-----	---

Date, in TDMS TOX date format,
on which correction was made.
Used by Update in forming
ANC-ANIMAL-NOTE-CORR-KEY-UQ of
ANIMAL-NOTE-CORRECTION.

FILE...: ANIMAL-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	ER	CORRECTION-TIME	P	5.0	N	
Time of day, in TDMS TOX time format, at which correction was made. Used by Update in forming ANC-ANIMAL-NOTE-CORR-KEY-UQ of ANIMAL-NOTE-CORRECTION.						
2	ES	CORRECTION-REASON	A	2.0	N	
Correction reason. Foreign key to TBCORREASON. Loaded into ANC-CORRECTION-REASON in ANIMAL-NOTE-CORRECTION.						
2	ET	CORRECTION-NOTE	A	79.0	N	
Free text describing the corresponds-rection. Loaded by Update into ANC-CORRECTION-NOTE of ANIMAL-NOTE-CORRECTION.						
2	EU	CORRECTION-TYPE	A	1.0	N	
Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into ANC-CORRECTION-TYPE of ANIMAL-NOTE-CORRECTION.						
2	FC	CORRECTION-KEY	A	36.0	N	DE
If either the date portion of TDD-EXP-TEST-ANIMAL-DATE-DE or TDD-NOTE-TIME is being changed then RECEIVE or ECS will concatenate the original TDD-EXP-TEST-ANIMAL-DATE-DE, the original TDD-NOTE-TIME, CORRECTION-DATE, and CORRECTION-TIME to form a string that it loads into this field. Loaded by Update into AWC-ANIMAL-NOTE-CORR-KEY-LINK-UQ in ANIMAL-NOTE-CORRECTION. The first 21 characters are a foreign key to TD-ANIMAL-DATA.						
1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
Unique key for database trans-						

FILE...: ANIMAL-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME

F LENG S DE

actions. Elements of key are:
SEGMENT-ID, MICRO-OR-ECS-FLAG,
TDD-EXP-TEST-ANIMAL-DATE-DE
for/from TD-ANIMAL-DATA left-
justified into 27-byte string,
TOX-DATA-DATE, TOX-DATA-TIME,
TOX-MAINFRAME-DATE, TOX-MAIN-
FRAME-TIME, and a 4-digit se-
quence number for the trans-
action in its segment within
the source set.

ANIMAL-OBSERVATION-TRANSACTION

Logical layout of animal observation transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ANIMAL-DATA file. For corrections, Update also uses the information to add records to ANIMAL-OBSERVATION-CORRECTION.

FILE...: ANIMAL-OBSERVATION-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY. Set to AO.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				

FILE...: ANIMAL-OBSERVATION-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	AI	TOX-MAINFRAME-DATE	B	2.0	F	
		The date, in TDMS TOX date format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AL	AGENCY-NUMBER	B	2.0	F	
		Foreign key to TBAGENCY. Not referenced by Update.				
2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
		When the transaction is a correspondence, this indicates the operator who made the correction and is loaded by Update into AOC-CORRECTION-OPERATOR of ANIMAL-OBSERVATION-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
2	AS	SOFTWARE-VERSION	B	2.0	F	
		The version of LDAS under which the source set was generated. Not referenced by Update.				
2	AT	TOX-DATA-DATE	B	2.0	F	
		The date, in TDMS TOX date format, on which the transaction was created. Used by RECEIVE and ECS in the forming of RECEIVE-UPDATE-KEY. Used by Update to form the EXP-TEST-ANIMAL-DATE-DE value that will be added, or select-				

FILE...: ANIMAL-OBSERVATION-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
ed for modification, in TD-ANIMAL-DATA. If a correction, also used to form AOC-ANIMAL-OBS-CORR-KEY-UQ in ANIMAL-OBSERVATION-CORRECTION.		
2 AU TOX-DATA-TIME	B	2.0 F
The time of day, in TDMS TOX time format, at which the transaction was created. Used by RECEIVE and ECS in forming RECEIVE-UPDATE-KEY. Loaded by Update into TDD-OBSERVATION-TIME of TD-ANIMAL-DATA. If a correction, also used to form AOC-ANIMAL-OBS-CORR-KEY-UQ in ANIMAL-OBSERVATION-CORRECTION.		
2 AV DATA-FACILITY-NUMBER	B	2.0 F
The facility associated with the machine on which the source set was created. Not referenced by Update. A foreign key to TBFACLT.		
2 AW OPERATOR-OR-PATHOLOGIST	B	2.0 F
The operator who created the transaction. Loaded by Update into TDD-OBSERVATION-OPERATOR-ID, and, if an add-with-audit, into AOC-OBSERVATION-OPERATOR-ID in ANIMAL-OBSERVATION-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.		
2 AX PROCESS-FLAG	A	2.0 F
Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.		

FILE...: ANIMAL-OBSERVATION-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
2 AZ DATA-BASE-KEY-LENGTH The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 13.	B	2.0 F
2 A0 NUMBER-OF-TREE-KEY-LEVELS The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 3.	B	2.0 F
2 A2 DATA-BASE-INDICATOR Indicates sub-system. Set to 2 for in-life.	A	1.0 F
2 A3 DATA-BASE-KEY Indicates the cage in which the animal being observed is located. (AATTTTTTCCCC). Used by Update to form the TDD-EXP-TEST-ANIMAL-DATE-DE value that will be added to, or selected for modification from, TD-ANIMAL-DATA. If a correction, also used to form AOC-ANIMAL-OBS-CORR-KEY-UQ of ANIMAL-OBSERVATION-CORRECTION. The test-cage portion is a foreign key to TD-CAGE.	A	14.0 F
2 A4 SEGMENT-TREE-KEY This field is parsed into 10 2-byte binary sub-fields. The first is set to 10; the second is set to the number of the animal being observed; the third is set to 40; the rest are unused. Update uses the animal number to form TDD-EXP-TEST-ANIMAL-DATE-DE of TD-ANIMAL-DATA and, if a correction, to form AOC-ANIMAL-OBS-CORR-KEY-UQ of ANIMAL-OBSERVATION-CORRECTION.	B	20.0 F

FILE...: ANIMAL-OBSERVATION-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME F LENG S DE

G 1 BT OBSERVATION-DATA
 Group embracing segment-specific, non-periodic group fields of animal observation transactions. (There is only one field in this group.)

2 BV MANUAL-ANIMAL-OBSERVATION-FLAG A 1.0 N
 Flag indicating whether the observation was entered during a manual entry session. 0 - No; 1 - Yes. Loaded by Update into the TDD-MANUAL-ENTRY-FLAG-OBS field of the TD-ANIMAL-DATA file, and, if an add-with-audit, also into AOC-MANUAL-ENTRY-FLAG of ANIMAL-OBSERVATION-CORRECTION.

P 1 BY ALL-ANIMAL-OBSERVATION-SETS
 Repeating group embracing segment-specific transactions of animal observations.

2 BZ CLINICAL-ANIMAL-OBSERVATION-NUM B 2.0 N
 Clinical sign observed on the animal. Loaded by Update into the TDD-OBSERVATION fields of the TD-ANIMAL-DATA file, and, if an add-with-audit, into AOC-OBSERVATION-DE of ANIMAL-OBSERVATION-CORRECTION. A foreign key to TBCLOBSC.

2 B0 CLINICAL-ANIMAL-OBSERVATION-SITE B 2.0 N
 Site where clinical sign was observed. Foreign key to TBCLOBSC. Loaded by Update into TDD-SITE, and, if an add-with-audit, into AOC-SITE of ANIMAL-OBSERVATION-CORRECTION.

2 B1 ANIMAL-OBSERVATION-QUALIFIER-1 B 2.0 N
 Size of the clinical sign. Foreign key to TBCLOBSC. Loaded by Update into TDD-

FILE...: ANIMAL-OBSERVATION-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		<p>EIQUAL1, and, if an add-with-audit, into AOC-EIQUAL1 of ANIMAL-OBSERVATION-CORRECTION.</p>				
2	B2	ANIMAL-OBSERVATION-QUALIFIER-2	B	2.0	N	<p>The count for the clinical sign. A foreign key for TBCLOBSC. Loaded by Update by Update into TDD-EIQUAL2, and, if an add-with-audit, into AOC-EIQUAL2 of ANIMAL-OBSERVATION-CORRECTION.</p>
G 1	EP	CORRECTION-PORITION				<p>Group embracing segment-common correction fields.</p>
2	EQ	CORRECTION-DATE	P	5.0	N	<p>Date, in TDMS TOX date format, on which correction was made. Used by Update in forming AOC-ANIMAL-OBS-CORR-KEY-UQ of ANIMAL-OBSERVATION-CORRECTION.</p>
2	ER	CORRECTION-TIME	P	5.0	N	<p>Time of day, in TDMS TOX time format, at which correction was made. Used by Update in forming AOC-ANIMAL-OBS-CORR-KEY-UQ of ANIMAL-OBSERVATION-CORRECTION.</p>
2	ES	CORRECTION-REASON	A	2.0	N	<p>The correction reason. Foreign key to TBCORREASON. Loaded by Update into AOC-CORRECTION-REASON in ANIMAL-OBSERVATION-CORRECTION.</p>
2	ET	CORRECTION-NOTE	A	79.0	N	<p>Free text describing the corresponds-rection. Loaded by Update into AOC-CORRECTION-NOTE of ANIMAL-OBSERVATION-CORRECTION.</p>
2	EU	CORRECTION-TYPE	A	1.0	N	<p>Single letter indicating type of correction. A-Add, C-Change</p>

FILE...: ANIMAL-OBSERVATION-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

D-Delete. Loaded by Update
 into AOC-CORRECTION-TYPE of
 ANIMAL-OBSERVATION-CORRECTION.

2	FC	CORRECTION-KEY	A	36.0	N	DE
---	----	----------------	---	------	---	----

If either the date portion of
 TDD-EXP-TEST-ANIMAL-DATE-DE or
 TDD-OBSERVATION-TIME is being
 changed then RECEIVE or ECS
 will concatenate the original
 TDD-EXP-TEST-ANIMAL-DATE-DE,
 the original TDD-OBSERVATION-
 TIME, CORRECTION-DATE, and
 CORRECTION-TIME to form a
 string that it loads into this
 field. Loaded by Update into
 AOC-ANIMAL-OBS-CORR-KEY-LINK-
 UQ in ANIMAL-OBSERVATION-COR-
 RECTION. The first 21 char-
 acters are a foreign key to
 TD-ANIMAL-DATA.

1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
---	----	--------------------	---	------	---	----

Unique key for database trans-
 actions. Elements of key are:
 SEGMENT-ID, MICRO-OR-ECS-FLAG,
 TDD-EXP-TEST-ANIMAL-DATE-DE
 for/from TD-ANIMAL-DATA left-
 justified into 27-byte string,
 TOX-DATA-DATE, TOX-DATA-TIME,
 TOX-MAINFRAME-DATE, TOX-MAIN-
 FRAME-TIME, and a 4-digit se-
 quence number for the trans-
 action in its segment within
 the source set.

ANIMAL-REMOVAL-TRANSACTION

Logical layout of animal removal transactions. Loaded by RECEIVE and ECS. Used by Update to modify records in TD-ANIMAL. For corrections, the information is also used by Update to add records to ANIMAL-REMOVAL-CORRECTION.

FILE...: ANIMAL-REMOVAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY. Set to AR.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date				

FILE...: ANIMAL-REMOVAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
---	----	--------------------	---	-----	---	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AL	AGENCY-NUMBER	B	2.0	F	
---	----	---------------	---	-----	---	--

Foreign key to TBAGENCY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
---	----	-----------------------	---	-----	---	--

When the transaction is a correspondence, this indicates the operator making the correction and is loaded by Update into ARC-CORRECTION-OPERATOR-ID of ANIMAL-REMOVAL-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.

2	AS	SOFTWARE-VERSION	B	2.0	F	
---	----	------------------	---	-----	---	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F	
---	----	---------------	---	-----	---	--

The date, in TDMS TOX date format, on which the transaction was created. Loaded by Update into TDA-REMOVAL-DATE and TDA-REMOVED-DATE of TD-ANIMAL. If an add-with-

FILE...: ANIMAL-REMOVAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

audit, also loaded into ARC-REMOVAL-DATE of ANIMAL-REMOVAL-CORRECTION.

2	AU	TOX-DATA-TIME	B	2.0	F	
---	----	---------------	---	-----	---	--

The time of day, in TDMS TOX time format, at which the transaction was created. Loaded by Update into TDA-REMOVAL-TIME of TD-ANIMAL. If an add-with-audit, also loaded into ARC-REMOVAL-TIME of ANIMAL-REMOVAL-CORRECTION.

2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
---	----	----------------------	---	-----	---	--

The facility associated with the machine on which the source set was created. Not referenced by Update. A foreign key to TBFACILITY.

2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
---	----	-------------------------	---	-----	---	--

The operator who created the transaction. Loaded by Update into TDA-REMOVAL-OPERATOR-ID. If an add-with-audit, also loaded into ARC-REMOVAL-OPERATOR-ID of ANIMAL-REMOVAL-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.

2	AX	PROCESS-FLAG	A	2.0	F	
---	----	--------------	---	-----	---	--

Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.

2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
---	----	----------------------	---	-----	---	--

The number of digits being used in the DATA-BASE-KEY

FILE...: ANIMAL-REMOVAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

field. Not referenced by Update. Set to 13.

2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
---	----	---------------------------	---	-----	---	--

The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 3.

2	A2	DATA-BASE-INDICATOR	A	1.0	F	
---	----	---------------------	---	-----	---	--

Indicates sub-system. Set to 2 for in-life.

2	A3	DATA-BASE-KEY	A	14.0	F	
---	----	---------------	---	------	---	--

Indicates the cage of the animal being removed, in the form of agency, test, and cage number. (AATTTTTTCCCC) Used by Update to form the TDA-EXP-TEST-ANIMAL-UQ value that it uses to select the TD-ANIMAL record to which it will apply the removal transaction. If a correction, it is used to form ARC-ANIMAL-REM-CORR-KEY-UQ of ANIMAL-REMOVAL-CORRECTION. The test-cage portion is a foreign key to TD-CAGE.

2	A4	SEGMENT-TREE-KEY	B	20.0	F	
---	----	------------------	---	------	---	--

This field is parsed into 10 2-byte binary sub-fields. Only 3 are used for animal removals. The first is set to 10; the second is the number of the animal being removed; the third is set to 20. Update uses the animal number in forming the TDA-EXP-TEST-ANIMAL-UQ value that is used to select the TD-ANIMAL record to which it applies the transaction. For corrections, also used to form ARC-ANIMAL-REM-CORR-KEY-LINK-UQ for ANIMAL-

FILE...: ANIMAL-REMOVAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

REMOVAL-CORRECTION.

G	1	C1	ANIMAL-REMOVAL-DATA				
---	---	----	---------------------	--	--	--	--

Group embracing segment-specific, non-repeating group, fields in animal removal transactions.

2	C3	CARCASS-IDENTIFICATION		A	14.0	N	
---	----	------------------------	--	---	------	---	--

Unique numeric code that identifies the animal for pathology reports. Loaded by Update into TDA-CARCASS-IDENTIFICATION-DE of TD-ANIMAL. If an add-with-audit, also loaded into ARC-CARCASS-IDENTIFICATION-DE of ANIMAL-REMOVAL-CORRECTION.

2	C5	MANUAL-REMOVAL-FLAG		A	1.0	N	
---	----	---------------------	--	---	-----	---	--

Flag indicating whether removal was recorded as part of a manual entry session. 0 - No; 1 - Yes. Loaded by Update into TDA-MANUAL-ENTRY-FLAG-REM of TD-ANIMAL. If an add-with-audit, also loaded into ARC-MANUAL-FLAG of ANIMAL-REMOVAL-CORRECTION.

2	C8	DAYS-ON-EXPERIMENT		B	2.0	N	
---	----	--------------------	--	---	-----	---	--

Dynamically calculated value equal to the removal date minus the start date of the animal's cage. For LDAS tests, this value is the same as days-on-dose. Loaded by Update into TDA-DAYS-ON-EXPERIMENT of TD-ANIMAL. If an add-with audit, also loaded into ARC-DAYS-ON-EXPERIMENT of ANIMAL-REMOVAL-CORRECTION.

2	DA	MANUAL-REMOVAL-WEIGHT-FLAG		A	1.0	N	
---	----	----------------------------	--	---	-----	---	--

Flag indicating whether removal was recorded as part of a manual entry session.

FILE...: ANIMAL-REMOVAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

al weight was recorded as part of a manual entry session. 0-No; 1 - Yes. Loaded by Update into TDA-MANUAL-ENTRY-FLAG-WGT of TD-ANIMAL. If an add-with-audit, also loaded into ARC-MANUAL-ENTRY-FLAG-WGT of ANIMAL-REMOVAL-CORRECTION.

2	DC	REMOVAL-WEIGHT-STATUS	P	5.0	N	
---	----	-----------------------	---	-----	---	--

Weight status indicating the reason for which a removal weight was not taken. Loaded by Update into TDA-REM-WEIGHT-STATUS of TD-ANIMAL. If an add-with-audit, also loaded into ARC-REM-WEIGHT-STATUS of ANIMAL-REMOVAL-CORRECTION. Foreign key to TBWGTSTATUS.

2	DD	REMOVAL-WEIGHT	B	4.0	N	
---	----	----------------	---	-----	---	--

Weight, in decigrams, of the animal at removal. Loaded by Update into TDA-ANIMAL-REM-WEIGHT of TD-ANIMAL. If an add-with-audit, also loaded into ARC-REMOVAL-WEIGHT of ANIMAL-REMOVAL-CORRECTION.

2	DF	REMOVAL-REASON	B	2.0	N	
---	----	----------------	---	-----	---	--

Reason for removal. Loaded by Update into TDA-REASON-FOR-REMOVAL-DE. f an add-with-audit, also loaded into ARC-REASON-FOR-REMOVAL-DE of ANIMAL-REMOVAL-CORRECTION. A foreign key to TBCLOBSC.

2	DH	MANUAL-REMOVAL-OBSERVATION-FLAG	A	1.0	N	
---	----	---------------------------------	---	-----	---	--

Flag indicating whether removal observations were recorded as part of a manual entry session. 0 - No; 1 - Yes. Loaded by Update into TDA-MANUAL-ENTRY-FLAG-OBS of TD-ANIMAL. If

FILE...: ANIMAL-REMOVAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

an add-with-audit, also loaded
into ARC-MANUAL-ENTRY-FLAG-OBS
of ANIMAL-REMOVAL-CORRECTION.

P 1 BY ALL-REMOVAL-OBSERVATION-SETS

Repeating group embracing an-
imal removal observations.

2 BZ REMOVAL-OBSERVATION-NUMBER	B	2.0 N
---------------------------------	---	-------

Code for clinical sign object-
served on animal being re-
moved. Loaded by Update into
TDA-OBSERVATION of TD-ANIMAL.
If an add-with-audit, also
loaded into ARC-OBSERVATION-DE
of ANIMAL-REMOVAL-CORRECTION.
Foreign key to TBCLOBSC.

2 B0 REMOVAL-OBSERVATION-SITE	B	2.0 N
-------------------------------	---	-------

Code for site of clinical sign
when that sign requires that
a site be specified. Loaded
by Update into TDA-SITE of TD-
ANIMAL. If an add-with-audit,
also loaded into ARC-SITE of
ANIMAL-REMOVAL-CORRECTION.
Foreign key to TBCLOBSC.

2 B1 REMOVAL-QUALIFIER-1	B	2.0 N
--------------------------	---	-------

Code for the size of the clin-
ical sign when that sign re-
quires that a size be speci-
fied. Loaded by Update into
TDA-EIQUAL1 of TD-ANIMAL. If
an add-with-audit, also loaded
into ARC-EIQUAL1 of ANIMAL-RE-
MOVAL-CORRECTION. Foreign key
to TBCLOBSC.

2 B2 REMOVAL-QUALIFIER-2	B	2.0 N
--------------------------	---	-------

Code for the count of the
clinical sign when that sign
requires that count be speci-
fied. Loaded by Update into
TDA-EIQUAL2 of TD-ANIMAL. If
an add-with-audit, also loaded

FILE...: ANIMAL-REMOVAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

into ARC-EIQUAL2 of ANIMAL-REMOVAL-CORRECTION. Foreign key to TBCLOBSC.

G 1 EP CORRECTION-PORTION

Group embracing segment-common correction fields.

2	EQ	CORRECTION-DATE	P	5.0	N		
---	----	-----------------	---	-----	---	--	--

Date, in TDMS TOX date format, on which correction was made. Used by Update to form ARC-ANIMAL-REM-CORR-KEY-UQ of ANIMAL-REMOVAL-CORRECTION.

2	ER	CORRECTION-TIME	P	5.0	N		
---	----	-----------------	---	-----	---	--	--

Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form ARC-ANIMAL-REM-CORR-KEY-UQ of ANIMAL-REMOVAL-CORRECTION.

2	ES	CORRECTION-REASON	A	2.0	N		
---	----	-------------------	---	-----	---	--	--

Correction reason. Foreign key to TBCORREASON. Loaded by Update into ARC-CORRECTION-REASON of ANIMAL-REMOVAL-CORRECTION.

2	ET	CORRECTION-NOTE	A	79.0	N		
---	----	-----------------	---	------	---	--	--

Free text describing the corresponds-rection. Loaded by Update into ARC-CORRECTION-NOTE of ANIMAL-REMOVAL-CORRECTION.

2	EU	CORRECTION-TYPE	A	1.0	N		
---	----	-----------------	---	-----	---	--	--

Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into ARC-CORRECTION-TYPE of ANIMAL-REMOVAL-CORRECTION.

2	EX	CORRECTION-KEY	A	26.0	N	DE	
---	----	----------------	---	------	---	----	--

No value assigned by RECEIVE or ECS. Update loads the

FILE...: ANIMAL-REMOVAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

field into ANIMAL-REM-CORR-
KEY-LINK-UQ in ANIMAL-REMOVAL-
CORRECTION.

1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
---	----	--------------------	---	------	---	----

Unique key for database trans-
actions. Elements of key are:
SEGMENT-ID, MICRO-OR-ECS-FLAG,
TDA-EXP-TEST-ANIMAL-UQ from
TD-ANIMAL left-justified into
a 27-byte string, TOX-DATA-
DATE, TOX-DATA-TIME, TOX-MAIN-
FRAME-DATE, TOX-MAINFRAME-
TIME, and a 4-digit sequence
number of the transaction for
its segment within the source
set.

ANIMAL-TRANSFER-TRANSACTION

Logical layout for animal transfer transactions. Loaded by RECEIVE and ECS. Processed by Update to modify records in the TD-ANIMAL file. For corrections, Update also uses the information to add records to the ANY-ANIMAL-CORRECTION file.

FILE...: ANIMAL-TRANSFER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY. Set to AT.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date				

FILE...: ANIMAL-TRANSFER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
			format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AJ		TOX-MAINFRAME-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AL		AGENCY-NUMBER	B	2.0	F	
			Foreign key to TBAGENCY. Not referenced by Update.				
2	AP		ENTRY-OPERATOR-NUMBER	B	2.0	F	
			When the transaction is a corresponds-rection, this indicates the operator ID for the person making the correction and is loaded by Update into ATC-CORRECTION-OPERATOR-ID of ANIMAL-TRANSFER-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
2	AS		SOFTWARE-VERSION	B	2.0	F	
			The version of LDAS under which the source set was generated. Not referenced by Update.				
2	AT		TOX-DATA-DATE	B	2.0	F	
			The date, in TDMS TOX date format, on which the transaction was created. Duplicate of TOX-ALLOCATION-DATE. If an add-with-audit, loaded by Update into ATC-TRANSFER-DATE of ANIMAL-TRANSFER-CORRECTION.				
2	AU		TOX-DATA-TIME	B	2.0	F	
			The time of day, in TDMS TOX				

FILE...: ANIMAL-TRANSFER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

time format, at which the transaction was created. Duplicate of TOX-ALLOCATION-TIME. If an add-with-audit, loaded by Update into ATC-TRANSFER-TIME of ANIMAL-TRANSFER-CORRECTION.

2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
---	----	----------------------	---	-----	---	--

The facility associated with the machine on which the source set was created. Not referenced by Update. A foreign key to TBFACILITY.

2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
---	----	-------------------------	---	-----	---	--

The ID of the operator who created the transaction. Loaded by Update into TDA-TRANSFER-OPERATOR-ID of TD-ANIMAL, and, if a correction, into ATC-TRANSFER-OPERATOR-ID of ANIMAL-TRANSFER-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.

2	AX	PROCESS-FLAG	A	2.0	F	
---	----	--------------	---	-----	---	--

Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.

2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
---	----	----------------------	---	-----	---	--

The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 13.

2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
---	----	---------------------------	---	-----	---	--

The number of two-byte binary

FILE...: ANIMAL-TRANSFER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

segments of the SEGMENT-TREE-KEY field being used. Set to 3.

2	A2	DATA-BASE-INDICATOR	A	1.0	F	
---	----	---------------------	---	-----	---	--

Indicates sub-system. Set to 2 for in-life.

2	A3	DATA-BASE-KEY	A	14.0	F	
---	----	---------------	---	------	---	--

Indicates the agency-test-cage number to which the transaction applies. (AATTTTTTCCCC)
Used by Update to form the value of TDA-EXP-TEST-ANIMAL-UQ of TD-ANIMAL that identifies the record being transferred (updated). If also a correction, used by Update to form ATC-ANIMAL-TRAN-CORR-KEY-UQ of ANIMAL-TRANSFER-CORRECTION. The test-cage portion is a foreign key to TD-CAGE.

2	A4	SEGMENT-TREE-KEY	B	20.0	F	
---	----	------------------	---	------	---	--

This field is parsed into 10 2-byte binary sub-fields. Only 3 are used. The first is set to 10; the second is set to the animal number being transferred; the third is set to 15. The animal number is used to form a TDA-EXP-TEST-ANIMAL-UQ value to select the record in TD-ANIMAL to update with transfer information. If a correction, used by Update to form ATC-ANIMAL-TRAN-CORR-KEY-UQ of ANIMAL-TRANSFER-CORRECTION.

G	1	CX	ANIMAL-TRANSFER-DATA			
---	---	----	----------------------	--	--	--

Group embracing segment-specific fields in animal transfer transactions.

2	CY	NEW-CAGE-NUMBER	P	5.0	N	
---	----	-----------------	---	-----	---	--

FILE...: ANIMAL-TRANSFER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

Number of cage to which the animal is being transferred. Loaded by Update into the TDA-HOSPITAL-CAGE-DE field of the TD-ANIMAL file. If an add-with-audit, loaded by Update into ATC-HOSPITAL-CAGE of ANIMAL-TRANSFER-CORRECTION.

2	CZ	TOX-ALLOCATION-DATE	B	2.0	N	
---	----	---------------------	---	-----	---	--

The date, in TDMS TOX date format, on which the animal was transferred to the new cage. Duplicate of TOX-DATA-DATE. Loaded by Update into TDA-TRANSFER-DATE of TD-ANIMAL

2	C0	TOX-ALLOCATION-TIME	B	2.0	N	
---	----	---------------------	---	-----	---	--

The time of day, in TDMS TOX time format, at which the animal was transferred. Duplicate of TOX-DATA-TIME. Loaded by Update into TDA-TRANSFER-TIME of TD-ANIMAL.

G	1	EP	CORRECTION-PORTION			
---	---	----	--------------------	--	--	--

Group embracing segment-common correction fields.

2	EQ	CORRECTION-DATE	P	5.0	N	
---	----	-----------------	---	-----	---	--

Date, in TDMS TOX date format, on which correction was made. Used by Update to form ATC-ANIMAL-TRAN-CORR-KEY-UQ of ANIMAL-TRANSFER-CORRECTION.

2	ER	CORRECTION-TIME	P	5.0	N	
---	----	-----------------	---	-----	---	--

Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form ATC-ANIMAL-TRAN-CORR-KEY-UQ of ANIMAL-TRANSFER-CORRECTION.

2	ES	CORRECTION-REASON	A	2.0	N	
---	----	-------------------	---	-----	---	--

Correction reason. Loaded by Update into ATC-CORRECTION-

FILE...: ANIMAL-TRANSFER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		REASON of ANIMAL-TRANSFER-COR-				
		RECTION. Foreign key to				
		TBCORREASON.				
2	ET	CORRECTION-NOTE	A	79.0	N	
		Free text describing the corresponds-				
		rection. Loaded by Update in-				
		to ATC-CORRECTION-NOTE of ANI-				
		MAL-TRANSFER-CORRECTION.				
2	EU	CORRECTION-TYPE	A	1.0	N	
		Single letter indicating type				
		of correction. A-Add, C-Change				
		D-Delete. Loaded by Update				
		into ATC-CORRECTION-TYPE of				
		ANIMAL-TRANSFER-CORRECTION.				
2	EX	CORRECTION-KEY	A	26.0	N	DE
		Not assigned a value by ECS or				
		RECEIVE. Loaded by Update in-				
		to ATC-ANIMAL-TRAN-CORR-KEY-				
		LINK-UQ in ANIMAL-TRANSFER-				
		CORRECTION.				
1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
		Unique key for database trans-				
		actions. Elements of key are:				
		SEGMENT-ID, MICRO-OR-ECS-FLAG,				
		TDA-EXP-TEST-ANIMAL-UQ from				
		TD-ANIMAL left-justified into				
		a 27-byte string, TOX-DATA-				
		DATE, TOX-DATA-TIME, TOX-MAIN-				
		FRAME-DATE, TOX-MAINFRAME-				
		TIME, and a 4-digit sequence				
		number of the transaction for				
		its segment within the source				
		set.				

ANIMAL-WEIGHT-TRANSACTION

Logical layout of animal weight transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ANIMAL-DATA file. For corrections, Update also uses the information to add records to the ANIMAL-WEIGHT-CORRECTION file.

FILE...: ANIMAL-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY. Set to AW.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	

FILE...: ANIMAL-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
			The date, in TDMS TOX date format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AJ		TOX-MAINFRAME-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AL		AGENCY-NUMBER	B	2.0	F	
			A foreign key to TBAGENCY. Not referenced by Update.				
2	AP		ENTRY-OPERATOR-NUMBER	B	2.0	F	
			When the transaction is a corresponds-rection, this indicates the operator ID for the person making the correction, and is loaded by Update into AWC-CORRECTION-OPERATOR-ID of ANIMAL-WEIGHT-CORRECTION. When appended to DATA-FACILITY-NUMBER a foreign key to TB_OPERATOR_ID.				
2	AS		SOFTWARE-VERSION	B	2.0	F	
			The version of LDAS under which the source set was generated. Not referenced by Update.				
2	AT		TOX-DATA-DATE	B	2.0	F	
			The date, in TDMS TOX date format, on which the transaction was created. Used by RECEIVE and ECS in forming RE-				

FILE...: ANIMAL-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

CEIVE-UPDATE-KEY. Used by Update to form the TDD-EXP-TEST-ANIMAL-DATE-DE value that will be added or modified to TD-ANIMAL-DATA. If a correction, also used to form AWC-ANIMAL-WGT-CORR-KEY-UQ of ANIMAL-WEIGHT-CORRECTION.

2	AU	TOX-DATA-TIME	B	2.0	F	
---	----	---------------	---	-----	---	--

The time of day, in TDMS TOX time format, at which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Loaded by Update into TDD-WEIGHT-TIME of TD-ANIMAL. If a correction, also used to form AWC-ANIMAL-WGT-CORR-KEY-UQ of ANIMAL-WEIGHT-CORRECTION.

2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
---	----	----------------------	---	-----	---	--

The facility associated with the machine on which the source set was created. Not referenced by Update. A foreign key to TBFACILITY.

2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
---	----	-------------------------	---	-----	---	--

The ID of the operator who created the transaction. Loaded by Update into TDD-WEIGHT-OPERATOR-ID of TD-ANIMAL-DATA, and, if an add-with-audit, into AWC-WEIGHT-OPERATOR-ID of ANIMAL-WEIGHT-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.

2	AX	PROCESS-FLAG	A	2.0	F	
---	----	--------------	---	-----	---	--

Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or del-

FILE...: ANIMAL-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

ete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.

2	AZ		DATA-BASE-KEY-LENGTH	B	2.0	F	
---	----	--	----------------------	---	-----	---	--

The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 13.

2	A0		NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
---	----	--	---------------------------	---	-----	---	--

The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 2.

2	A2		DATA-BASE-INDICATOR	A	1.0	F	
---	----	--	---------------------	---	-----	---	--

Indicates sub-system. Set to 2 for in-life.

2	A3		DATA-BASE-KEY	A	14.0	F	
---	----	--	---------------	---	------	---	--

Indicates the number of the agency-test-cage number where the animal being observed is located. (AATTTTTTCCCC)
Used by Update to form the EXP-TEST-ANIMAL-DATE-DE value that will be added to, or selected from, TD-ANIMAL-DATE. If a correction, also used to form AWC-ANIMAL-WGT-CORR-KEY-UQ of ANIMAL-WEIGHT-CORRECTION. The test-cage portion is a foreign key to TD-CAGE.

2	A4		SEGMENT-TREE-KEY	B	20.0	F	
---	----	--	------------------	---	------	---	--

This field is parsed into 10 2-byte binary sub-fields. Only 3 are used for animal weight transactions. The first is set to 10; the second is set to the number of the animal being observed; the third is set to 30. Update

FILE...: ANIMAL-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME F LENG S DE

used the animal number to form
TDD-EXP-TEST-ANIMAL-DATE-DE FO
for TD-ANIMAL-DATA, and, if a
correction, to form AWC-ANI-
MAL-WGT-CORR-KEY-UQ of ANIMAL-
WEIGHT-CORRECTION.

G 1 B3 ANIMAL-WEIGHT-DATA

Group embracing segment-speci-
fic fields of animal weight
transactions.

2 B5 ANIMAL-MANUAL-WEIGHT-FLAG A 1.0 N

Flag indicating whether animal
weight was manually entered.
0 - Automatically recorded
from attached scale 1 - Man-
ually entered. Loaded by
Update into the TDD-MANUAL-
ENTRY-FLAG-WGT field of the
TD-ANIMAL-DATA file. If an
add-with-audit, also loaded
into AWC-MANUAL-ENTRY-FLAG-WGT
of ANIMAL-WEIGHT-CORRECTION.

2 B7 ANIMAL-WEIGHT-STATUS P 5.0 N

When no weight is recorded,
status code that explains why
weight was not recorded.
Loaded by Update into the
field TDD-WEIGHT-STATUS of the
TD-ANIMAL-DATA file. If an
add-with-audit, also loaded
into AWC-WEIGHT-STATUS of ANI-
MAL-WEIGHT-CORRECTION. For-
eign key to TBWGSTATUS.

2 B8 ANIMAL-WEIGHT B 4.0 N

Weight of animal, expressed in
decigrams. Loaded by Update
into the TDD-ANIMAL-WEIGHT of
the TD-ANIMAL-DATA file. If
an add-with-audit, also loaded
into AWC-ANIMAL-WEIGHT of ANI-
MAL-WEIGHT-CORRECTION.

FILE...: ANIMAL-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
G 1 EP CORRECTION-PORTION		
Group embracing segment-common correction fields.		
2 EQ CORRECTION-DATE	P	5.0 N
Date, in TDMS TOX date format, on which correction was made. Used by Update in forming AWC-ANIMAL-WGT-CORR-KEY-UQ of ANIMAL-WEIGHT-CORRECTION.		
2 ER CORRECTION-TIME	P	5.0 N
Time of day, in TDMS TOX time format, at which correction was made. Used by Update in forming AWC-ANIMAL-WGT-CORR-KEY-UQ of ANIMAL-WEIGHT-CORRECTION.		
2 ES CORRECTION-REASON	A	2.0 N
Correction reason. Loaded by Update into AWC-CORRECTION-REASON of ANIMAL-WEIGHT-CORRECTION. Foreign key to TBCORRREASON.		
2 ET CORRECTION-NOTE	A	79.0 N
Free text describing the corresponds-rection. Loaded by Update into AWC-CORRECTION-NOTE of ANIMAL-WEIGHT-CORRECTION.		
2 EU CORRECTION-TYPE	A	1.0 N
Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into AWC-CORRECTION-TYPE of ANIMAL-WEIGHT-CORRECTION.		
2 FC CORRECTION-KEY	A	36.0 N DE
If either the date portion of TDD-EXP-TEST-ANIMAL-DATE-DE or TDD-WEIGHT-TIME is being changed then RECEIVE or ECS will concatenate the original TDD-EXP-TEST-ANIMAL-DATE-DE, the original TDD-WEIGHT-TIME,		

FILE...: ANIMAL-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F LENG S DE
-------------	-------------

CORRECTION-DATE, and CORREC-
TION-TIME to form a string
that it loads into this field.
Loaded by Update into AWC-ANI-
MAL-WGT-CORR-KEY-LINK-UQ in
ANIMAL-WEIGHT-CORRECTION.
First 21 characters form a
foreign key to TD-ANIMAL-
DATA.

1 EW RECEIVE-UPDATE-KEY	A 54.0 N DE
-------------------------	-------------

Unique key for database trans-
actions. Elements of key are:
SEGMENT-ID, MICRO-OR-ECS-FLAG,
TDD-EXP-TEST-ANIMAL-DATE-DE
for/from TD-ANIMAL-DATA left-
justified into 27-byte string,
TOX-DATA-DATE, TOX-DATA-TIME,
TOX-MAINFRAME-DATE, TOX-MAIN-
FRAME-TIME, and a 4-digit se-
quence number for the trans-
action in its segment within
the source set.

ANY-ANIMAL-TRANSACTION

Logical layout for animal allocation transactions. Loaded by RECEIVE and ECS. Processed by Update to add, modify, or delete records to/from the TD-ANIMAL file. For corrections, Update also uses the information to add records to the ANY-ANIMAL-CORRECTION file.

FILE...: ANY-ANIMAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Set to AA. Used to form RECEIVE-UPDATE-KEY.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date format, on which the source				

FILE...: ANY-ANIMAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
---	----	--------------------	---	-----	---	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AL	AGENCY-NUMBER	B	2.0	F	
---	----	---------------	---	-----	---	--

Foreign key to TBAGENCY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
---	----	-----------------------	---	-----	---	--

When the transaction is a corresponds-rection, this indicates the operator ID for the person making the correction and is loaded by Update into AAC-CORRECTION-OPERATOR-ID of ANY-ANIMAL-CORRECTION. When appended to DATA-FACILITY-NUMBER it is a foreign key to TB_OPERATOR_ID.

2	AS	SOFTWARE-VERSION	B	2.0	F	
---	----	------------------	---	-----	---	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F	
---	----	---------------	---	-----	---	--

The date, in TDMS TOX date format, on which the transaction was created. Loaded by Update into TDA-ANIMAL-DATE of TD-ANIMAL. If an add-with-audit, also loaded into AAC-

FILE...: ANY-ANIMAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		ANIMAL-DATE of ANY-ANIMAL-COR- RECTION.				
2	AU	TOX-DATA-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, at which the transaction was created. Loaded by Update into TDA-ANI- MAL-TIME. If an add-with-aud- it, also loaded into AAC-ANI- MAL-TIME of ANY-ANIMAL-CORREC- TION.				
2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
		The facility associated with the machine on which the source set was created. Not referenced by Update. A for- eign key to TBFACLTy.				
2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
		The ID of the operator who created the transaction. Loaded by Update into TDA- ANIMAL-OPERATOR-ID. If an add add-with-audit, loaded into AAC-ANIMAL-OPERATOR-ID of ANY-ANIMAL-CORRECTION. When appended to DATA-FACILITY- NUMBER, a foreign key to TB_OPERATOR_ID.				
2	AX	PROCESS-FLAG	A	2.0	F	
		Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or del- ete, respectively. The second digit is either blank, indica- ting an original transaction, or 'A', indicating an audit transaction.				
2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
		The number of digits being used in the DATA-BASE-KEY				

FILE...: ANY-ANIMAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		field. Set to 13. Not referenced by Update.				
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
		The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 2.				
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
		Indicates sub-system. Set to 1 for in-life.				
2	A3	DATA-BASE-KEY	A	14.0	F	
		Indicates the agency, test, and cage to which the transaction applies. Used by Update to form TDA-EXP-TEST-ANIMAL-UQ of TD-ANIMAL, and, if a correction, AAC-ANY-ANIMAL-CORR-KEY-UQ of ANY-ANIMAL-CORRECTION. The test and cage portion forms a foreign key to TD-CAGE.				
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
		This field is parsed into 10 2-byte binary sub-fields. Only the first two are used. The first is set to 10; the second is set to the animal number being allocated. Update uses the animal number to form TDA-EXP-TEST-ANIMAL-UQ of TD-ANIMAL, and, if a correction, AAC-ANY-ANIMAL-CORR-KEY-UQ of ANY-ANIMAL-CORRECTION.				
G	1	BC ANY-ANIMAL-DATA				
		Group embracing all fields specific to the animal allocation segment.				
2	BE	ANIMAL-NUMBER	B	2.0	N	
		The number of the animal being allocated. Same as the second item in SEGMENT-TREE-KEY. Update				

FILE...: ANY-ANIMAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

date does not reference this field, using SEGMENT-TREE-KEY instead.

2	BF	ANIMAL-ID		A	1.0	N	
---	----	-----------	--	---	-----	---	--

Loaded by ECS only for legacy purposes. Loaded by Update to TDA-ANIMAL-IDENTIFICATION of TD-ANIMAL. If an add-with-audit, also loaded to AAC-ANIMAL-IDENTIFICATION of ANY-ANIMAL-CORRECTION.

2	BG	SEX		A	1.0	N	
---	----	-----	--	---	-----	---	--

Sex of animal as defined in test protocol. Loaded by Update into TDA-SEX of TD-ANIMAL. If an add-with-audit, also loaded into AAC-SEX of ANY-ANIMAL-CORRECTION.

2	BH	SPECIES		B	2.0	N	
---	----	---------	--	---	-----	---	--

The species code of the animal as defined in test protocol. Loaded by Update into TDA-SPECIES of TD-ANIMAL. If an add-with-audit, also loaded into AAC-SPECIES of ANY-ANIMAL-CORRECTION. Foreign key to TBSTRAIN.

2	BI	STRAIN		B	2.0	N	
---	----	--------	--	---	-----	---	--

The strain code of the animal as defined in test protocol. Loaded by Update into TDA-STRAIN of TD-ANIMAL. If an add-with-audit, also loaded into AAC-STRAIN of ANY-ANIMAL-CORRECTION. Foreign key to TBSTRAIN.

2	BJ	SUBSTRAIN		B	2.0	N	
---	----	-----------	--	---	-----	---	--

The substrain code of the animal as defined in test protocol. Loaded by Update into TDA-SUBSTRAIN of TD-ANIMAL.

FILE...: ANY-ANIMAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
If an add-with-audit, also loaded into AAC-SUBSTRAIN of ANY-ANIMAL-CORRECTION. Foreign key to TBSUBSTR.		
2 BK ANIMAL-CLASS	B	2.0 N
The animal class code as defined in test protocol. Loaded by Update into TDA-ANIMAL-CLASS of TD-ANIMAL. If an add-with-audit, also loaded into AAC-ANIMAL-CLASS of ANY-ANIMAL-CORRECTION. Foreign key to TBANCLQF.		
G 1 EP CORRECTION-PORTRION		
Group embracing segment-common correction fields.		
2 EQ CORRECTION-DATE	P	5.0 N
Date, in TDMS TOX date format, on which correction was made. Used by Update to form AAC-ANY-ANIMAL-CORR-KEY-UQ of ANY-ANIMAL-CORRECTION.		
2 ER CORRECTION-TIME	P	5.0 N
Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form AAC-ANY-ANIMAL-CORR-KEY-UQ of ANY-ANIMAL-CORRECTION.		
2 ES CORRECTION-REASON	A	2.0 N
Code for correction reason. Loaded by Update into AAC-CORRECTION-REASON of ANY-ANIMAL-CORRECTION. Foreign key to TBCORREASON.		
2 ET CORRECTION-NOTE	A	79.0 N
Free text describing the corresponds-rection. Loaded by Update into AAC-CORRECTION-NOTE of ANY-ANIMAL-CORRECTION.		
2 EU CORRECTION-TYPE	A	1.0 N
Single letter indicating type		

FILE...: ANY-ANIMAL-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
of correction. A-Add, C-Change		
D-Delete. Loaded by Update		
into AAC-CORRECTION-TYPE of		
ANY-ANIMAL-CORRECTION.		
2 EX CORRECTION-KEY	A	26.0 N DE
If the animal number is being		
changed by ECS, then the orig-		
inal value of TDA-EXP-TEST-		
ANIMAL-UQ from TD-ANIMAL, COR-		
RECTION-DATE, and CORRECTION-		
TIME are concatenated into		
this field. Loaded by Update		
into AAC-ANY-ANIMAL-CORR-KEY-		
LINK-UQ of ANY-ANIMAL-CORREC-		
TION. The first 16 characters		
form a foreign key to TD-ANI-		
MAL.		
1 EW RECEIVE-UPDATE-KEY	A	54.0 N DE
Unique key for database trans-		
actions. Elements of key are:		
SEGMENT-ID, MICRO-OR-ECS-FLAG,		
TDA-EXP-TEST-ANIMAL-UQ for/		
from TD-ANIMAL left-justified		
into 27-byte string, TOX-DATA-		
DATE, TOX-DATA-TIME, TOX-MAIN-		
FRAME-DATE, TOX-MAINFRAME-		
TIME, and 4-digit sequence		
number of the transaction for		
its segment in the source set.		
1 GA ANY-ANIMAL-KEY	A	9.0 N SP
SOURCE FIELD(S) --- -START- --END-		
DATA-BASE-KEY	3	9
ANIMAL-NUMBER	1	2
Super-key whose elements are		
the test number (7) and the		
animal (B2).		

ANY-MICRO-OBS-TRANSACTION

Logical layout of microscopic observation transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ORGAN-DATA file. For corrections, also used to add records to ANY-MICRO-OBS-CORRECTION.

FILE...: ANY-MICRO-OBS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
G	1	AC	HEADER				
			Group embracing fields that				
			are common to most segments.				
2	AA	SEGMENT-ID		A	2.0	F	
			Two-letter segment identifier.				
			Not referenced by Update. Set				
			to MO. Used to form RECEIVE-				
			UPDATE-KEY.				
2	AB	MICRO-OR-ECS-FLAG		A	1.0	F	
			Identifier for source of				
			transaction. A - LDAS,				
			E - Mainframe ECS. Used by				
			Update to generate Update				
			statistics report. Used to				
			form RECEIVE-UPDATE-KEY.				
2	AE	TOX-MICRO-DATE		B	2.0	F	
			The date, in TDMS TOX format,				
			on which the source set for				
			the transaction was created.				
			Not referenced by Update.				
2	AF	TOX-MICRO-TIME		B	2.0	F	
			The time of day, in TDMS TOX				
			time format, on which the				
			source set for the transac-				
			was created. Not refer-				
			enced by Update.				
2	AG	MICRO-TRANS-SET		B	2.0	F	
			The number of the transac-				
			tion's source set. Not refer-				
			enced by Update.				
2	AI	TOX-MAINFRAME-DATE		B	2.0	F	
			The date, in TDMS TOX date				

FILE...: ANY-MICRO-OBS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F		
---	----	--------------------	---	-----	---	--	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AL	AGENCY-NUMBER	B	2.0	F		
---	----	---------------	---	-----	---	--	--

Foreign key to TBAGENCY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F		
---	----	-----------------------	---	-----	---	--	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction and is loaded by Update into MOC-CORRECTION-OPERATOR-ID in ANY-MICRO-OBS-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.

2	AS	SOFTWARE-VERSION	B	2.0	F		
---	----	------------------	---	-----	---	--	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F		
---	----	---------------	---	-----	---	--	--

The date, in TDMS TOX date format, on which the transaction was created. Used by Update to form the value of TDH-ORGAN-OBSERVATION-ID-DE that will be added to, or selected for modification from, TD-ORGAN-DATA. For a correction, used to form MOC-MICRO-

FILE...: ANY-MICRO-OBS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		OBS-CORR-KEY-UQ for ANY-MICRO-OBS-CORRECTION.				
2	AU	TOX-DATA-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, at which the transaction was created. Loaded by Update into TDH-MICRO-OBSERVATION-SEQUENCE in TD-ORGAN-DATA. For a correction, used to form MOC-MICRO-OBS-CORR-KEY-UQ for ANY-MICRO-OBS-CORRECTION.				
2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
		The facility associated with the machine on which the source set was created. Not referenced by Update. Foreign key to TBFACILITY.				
2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
		The operator number of the pathologist who is assigned to the test. Loaded by Update into TDH-MICRO-OBS-PATHOLOGIST-ID in TD-ORGAN-DATA. If an add-with-audit, also loaded into MOC-MICRO-OBS-PATHOLOGIST-ID in ANY-MICRO-OBS-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
2	AX	PROCESS-FLAG	A	2.0	F	
		Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.				

FILE...: ANY-MICRO-OBS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
		The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 14.				
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
		The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 6.				
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
		Indicates sub-system. Set to 3 for pathology.				
2	A3	DATA-BASE-KEY	A	14.0	F	
		The carcass ID (CID) for the animal from which the organ being observed was taken. Used by Update to select the master TD-ANIMAL record for the transaction.				
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
		This field is parsed into 10 2-byte binary sub-fields. The second is for pathologist type and is set to "A"; the fourth is the organ; the sixth is the lesion. The organ code and pathologist type are used by Update to form TDH-MICRO-OBSERVATION-ID-DE, and MOC-MICRO-OBS-CORR-KEY-UQ for ANY-MICRO-OBS-CORRECTION. The organ and lesion codes are foreign keys to TBPCT. Update does not use the lesion code, using the MORPHOLOGY-NUMBER value, which is the same.				
G 1	DV	MICRO-OBSERVATION-DATA				
		Group embracing segment-specific fields of microscopic observation transactions.				

FILE...: ANY-MICRO-OBS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
2 DW MICRO-OBSERVATION-OPERATOR ID of operator who recorded the microscopic observation. Loaded by Update into TDH-MIC- RO-OBS-OPERATOR-ID in TD-ORG- AN-DATA. If an add-with-aud- it, also loaded into MOC-MIC- RO-OBS-OPERATOR-ID in ANY-MIC- RO-OBS-CORRECTION. When ap- pended to DATA-FACILITY-NUMBER it is a foreign key to TB_OPERATOR_ID.	B	2.0 N
2 DY PRIMARY-SITE Organ from which the micro- scopic observation metasta- sized. (Blank indicates that organ is original site.) Loaded by Update into TDH-PRI- MARY-SITE in TD-ORGAN-DATA. If an add-with-audit, also loaded into MOC-PRIMARY-SITE in ANY-MICRO-OBS-CORRECTION. Foreign key to TBPCT.	B	2.0 N
2 FR TRACE-LESION-NUMBER No meaning in TDMS. However, it is referenced by ECS and Update.	B	2.0 N
2 D1 MORPHOLOGY-NUMBER Morphology being microscopic- ally observed. Loaded by Update into TDH-OBSERVATION- MORPHOLOGY-DE in TD-ORGAN- DATA. If an add-with-audit, also loaded into MOC-OBSERVA- TION-MORPHOLOGY-DE in ANY- MICRO-ORGAN-CORRECTION. For- eign key to TBPCT.	B	2.0 N
M 2 D2 QUALIFIER-NUMBER Multiple-entry field contain- ing qualifiers that describe morphology being observed.	B	2.0 N

FILE...: ANY-MICRO-OBS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME F LENG S DE

Loaded by Update into TDH-OBSERVATION-QUALIFIER. If an add-with-audit, also loaded into MOC-OBSERVATION-QUALIFIER in ANY-MICRO-OBS-CORRECTION. (Qualifiers, by design, are limited to 4 entries.) Foreign key to TBPCT.

M 2 D3 SITE-NUMBER B 2.0 N

Multiple-entry field containing sites on the organ being observed where the morphology was detected. Loaded by Update into TDH-OBSERVATION-SITE in TD-ORGAN-DATA. For an add-with-audit, also loaded into MOC-OBSERVATION-SITE in ANY-MICRO-OBS-CORRECTION. (Limited by design to 3 entries.) Foreign key to TBPCT.

2 D4 CAUSE-OF-DEATH-FLAG A 1.0 N

Flag indicating whether the observation caused the death of the animal. C for contributory; P for primary. Loaded by Update the TDH-COD-FLAG in TD-ORGAN-DATA. For an add-with-audit, also loaded into MOC-COD-FLAG in ANY-MICRO-OBS-CORRECTION.

G 1 EP CORRECTION-PORTION

Group embracing segment-common correction fields.

2 EQ CORRECTION-DATE P 5.0 N

Date, in TDMS TOX date format, on which correction was made. If the observation date or time are being changed, then used by RECEIVE and ECS to

FILE...: ANY-MICRO-OBS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		form CORRECTION-KEY. Used by Update to form MOC-MICRO-OBS- CORR-KEY-UQ for ANY-MICRO-OBS- CORRECTION.				
2	ER	CORRECTION-TIME	P	5.0	N	
		Time of day, in TDMS TOX time format, at which correction was made. If the observation date or time are being changed used by RECEIVE and ECS to form CORRECTION-KEY. Used by Update to form MOC-MICRO-OBS- CORR-KEY-UQ for ANY-MICRO-OBS- CORRECTION.				
2	ES	CORRECTION-REASON	A	2.0	N	
		Code for correction reason. Loaded by Update into MOC-COR- RECTION-REASON in ANY-MICRO- OBS-CORRECTION. Foreign key to TBCORREASON.				
2	ET	CORRECTION-NOTE	A	79.0	N	
		Free text describing the corresponds- rection. Loaded by Update in- to MOC-CORRECTION-NOTE in ANY- MICRO-OBS-CORRECTION.				
2	EU	CORRECTION-TYPE	A	1.0	N	
		Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into MOC-CORRECTION-TYPE in ANY-MICRO-OBS-CORRECTION.				
2	FE	CORRECTION-KEY	A	42.0	N	DE
		If the observation date or time is being changed, ECS and RECEIVE combines into it TDO-EXP-TST-ANUM-ORG-PATHT-UQ from TD-ORGAN, the original observation date and time, both in TDMS TOX format, COR- RECTION-DATE, and CORRECTION- TIME. Loaded by Update into				

FILE...: ANY-MICRO-OBS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

MOC-MICRO-OBS-CORR-KEY-LINK-
UQ in ANY-MICRO-OBS-CORREC-
TION. The first 27 charac-
ters form a foreign key to
TD-ORGAN-DATA.

1 EW RECEIVE-UPDATE-KEY	A	54.0 N DE
-------------------------	---	-----------

Unique key for database trans-
actions. Elements of key are:
SEGMENT-ID, MICRO-OR-ECS-FLAG,
TDH-ORGAN-OBSERVATION-ID-DE,
TOX-DATA-DATE, TOX-DATA-TIME,
TOX-MAINFRAME-DATE, TOX-MAIN-
FRAME-TIME, and a 4-digit se-
quence number for the trans-
action in its segment within
the source set.

ANY-MICRO-ORGAN-TRANSACTION

Logical layout of organ (tissue) status transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ORGAN file. For corrections, used to add records to ANY-MICRO-ORGAN-CORRECTION.

FILE...: ANY-MICRO-ORGAN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Set to AM. Used to form RECEIVE-UPDATE-KEY.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date				

FILE...: ANY-MICRO-ORGAN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
---	----	--------------------	---	-----	---	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AL	AGENCY-NUMBER	B	2.0	F	
---	----	---------------	---	-----	---	--

Foreign key to TBAGENCY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
---	----	-----------------------	---	-----	---	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction and is loaded into AMC-CORRECTION-OPERATOR-ID in ANY-MICRO-ORGAN-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID>

2	AS	SOFTWARE-VERSION	B	2.0	F	
---	----	------------------	---	-----	---	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F	
---	----	---------------	---	-----	---	--

The date, in TDMS TOX date format, on which the transaction was created. Loaded by Update into TDO-ORGAN-STATUS-DATE in TD-ORGAN. If an add-with-audit, also loaded into AMC-ORGAN-STATUS-DATE in ANY-MICRO-ORGAN-CORRECTION.

2	AU	TOX-DATA-TIME	B	2.0	F	
---	----	---------------	---	-----	---	--

The time of day, in TDMS TOX

FILE...: ANY-MICRO-ORGAN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

time format, at which the transaction was created. Loaded by Update into TDO-ORGAN-STATUS-TIME in TD-ORGAN. If an add-with-audit, also loaded into AMC-ORGAN-STATUS-TIME in ANY-MICRO-ORGAN-CORRECTION.

2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
---	----	----------------------	---	-----	---	--

The facility associated with the machine on which the source set was created. Not referenced by Update. Foreign key to TBFACLT.

2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
---	----	-------------------------	---	-----	---	--

The operator ID of the pathologist who is assigned to the test. Loaded by Update into TDO-ORGAN-STAT-PATHOLOGIST-ID in TD-ORGAN. If an add-with-audit, also loaded into AMC-ORGAN-STAT-PATHOLOGIST-ID in ANY-MICRO-ORGAN-CORRECTION. When appended to DATA-FACILITY-CODE, a foreign key to TB_OPERATOR_ID.

2	AX	PROCESS-FLAG	A	2.0	F	
---	----	--------------	---	-----	---	--

Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.

2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
---	----	----------------------	---	-----	---	--

The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 14.

FILE...: ANY-MICRO-ORGAN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
2 A0 NUMBER-OF-TREE-KEY-LEVELS The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 4.	B	2.0 F
2 A2 DATA-BASE-INDICATOR Indicates sub-system. Set to 3 for pathology.	A	1.0 F
2 A3 DATA-BASE-KEY The carcass ID (CID) for the animal from which the organ being observed was removed. Used by Update to select the master TD-ANIMAL record for the transaction.	A	14.0 F
2 A4 SEGMENT-TREE-KEY This field is parsed into 10 2-byte binary sub-fields. Only the first four are used for organ status transactions. The first is set to 80; the second is for pathologist type and is set to "A"; the third is set to 10; the fourth is the organ code (foreign key to TBPCT). The organ code and pathologist type values are used by Update to form TDO-EXP-TST-ANUM-ORG-PATHT-UQ for TD-ANIMAL, and AMC-MICRO-ORG-CORR-KEY-UQ in ANY-MICRO-ORGAN-CORRECTION.	B	20.0 F
G 1 DR ANY-MICRO-ORGAN-DATA Group embracing segment-specific fields in organ status transactions.		

FILE...: ANY-MICRO-ORGAN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
2 DS MICRO-ORGAN-OPERATOR	B	2.0 N
ID of operator assigning organ status. Loaded by Update into TDO-ORGAN-STATUS-OPERATOR-ID in TD-ORGAN. If an add-with-audit, also loaded into AMC-ORGAN-STATUS-OPERATOR-ID in ANY-MICRO-ORGAN-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.		
2 DU STATUS-NUMBER	B	2.0 N
Status of organ being microscopically evaluated. Loaded by Update into TDO-ORGAN-STATUS of TD-ORGAN. If an add-with audit, also loaded into AMC-ORGAN-STATUS in ANY-MICRO-ORGAN-CORRECTION. Foreign key to TBORGSTA.		
G 1 EP CORRECTION-PORION		
Group embracing segment-common correction fields.		
2 EQ CORRECTION-DATE	P	5.0 N
Date, in TDMS TOX date format, on which correction was made. Used by Update to form AMC-MICRO-ORGAN-CORR-KEY-UQ for ANY-MICRO-ORGAN-CORRECTION.		
2 ER CORRECTION-TIME	P	5.0 N
Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form AMC-MICRO-ORG-CORR-KEY-UQ for ANY-MICRO-ORGAN-CORRECTION.		
2 ES CORRECTION-REASON	A	2.0 N
Code for correction reason. Loaded by Update into AMC-CORRECTION-REASON in ANY-MICRO-ORGAN-CORRECTION. Foreign key to TBCORREASON.		

FILE...: ANY-MICRO-ORGAN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	ET	CORRECTION-NOTE	A	79.0	N	
Free text describing the corresponds-rection. Loaded by Update into AMC-CORRECTION-NOTE in ANY-MICRO-ORGAN-CORRECTION.						
2	EU	CORRECTION-TYPE	A	1.0	N	
Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into AMC-CORRECTION-TYPE in ANY-MICRO-ORGAN-CORRECTION.						
2	FB	CORRECTION-KEY	A	32.0	N	DE
Not assigned a value by ECS or RECEIVE. Loaded by Update into AMC-MICRO-ORG-CORR-KEY-LINK-UQ in ANY-MICRO-ORGAN-CORRECTION.						
1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDO-EXP-TST-ANUM-ORG-PATHT-UQ, 5 spaces, TOX-DATA-DATE, TOX-DATA-TIME, TOX-MAINFRAME-DATE, TOX-MAINFRAME-TIME, and a 4-digit sequence number for the transaction in its segment within the source set.						

BALANCE-CALIBRATION-TRANSACTION

Logical layout of balance calibration transactions. Loaded by RECEIVE.
Read by Update to add records to the TD-BALANCE-CALIBRATION
file.

FILE...: BALANCE-CALIBRATION-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Set to BC. Used to form RECEIVE- UPDATE-KEY.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transac- tion was created. Not refer- enced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transac- tion's source set. Not refer- enced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date format, on which the source				

FILE...: BALANCE-CALIBRATION-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		set was received. Used by RE- CEIVE and ECS to form RECEIVE- UPDATE-KEY. Not referenced by Update.				
2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RE- CEIVE-UPDATE-KEY. Not refer- enced by Update.				
2	AL	AGENCY-NUMBER	B	2.0	F	
		Foreign key to TBAGENCY. Not referenced by Update.				
2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
		The ID of the operator who created the transaction. Same value as in OPERATOR-OR-PATHO- LOGIST. Update references, but does not use, this field, using OPERATOR-OR-PATHOLOGIST instead. When appended to DATA-FACILITY-NUMBER, a for- eign key to TB_OPERATOR_ID.				
2	AS	SOFTWARE-VERSION	B	2.0	F	
		The version of LDAS under which the source set was generated. Not referenced by Update.				
2	AT	TOX-DATA-DATE	B	2.0	F	
		The date, in TDMS TOX date format, on which the trans- action was created. Used by RECEIVE and ECS to form RE- CEIVE-UPDATE-KEY. Used by Up- date to from the value of TDF- EXP-TEST-DATE-DE that will be				

FILE...: BALANCE-CALIBRATION-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
added to, or selected for modification from, TD-BALANCE-CALIBRATION.		
2 AU TOX-DATA-TIME	B	2.0 F
The time of day, in TDMS TOX time format, at which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Loaded by Update into TDF-TIME of TD-BALANCE-CALIBRATION.		
2 AV DATA-FACILITY-NUMBER	B	2.0 F
The facility associated with the machine on which the source set was created. Not referenced by Update. Foreign key to TBFACLT.		
2 AW OPERATOR-OR-PATHOLOGIST	B	2.0 F
The ID of the operator who created the transaction. Contains same value as in ENTRY-OPERATOR-NUMBER. Loaded by Update into TDF-OPERATOR-IDENTIFICATION in TD-BALANCE-CALIBRATION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.		
2 AX PROCESS-FLAG	A	2.0 F
Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.		
2 AZ DATA-BASE-KEY-LENGTH	B	2.0 F
The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 9.		

FILE...: BALANCE-CALIBRATION-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
		The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 4.				
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
		Indicates sub-system. Set to 1 for balance calibration.				
2	A3	DATA-BASE-KEY	A	14.0	F	
		Indicates the agency-test to which the transaction applies. (AATTTTTT). Used by Update to form the value of TDF-EXP-TEST-DATE-DE that will be added to, or selected for modification from, TD-BALANCE-CALIBRATION.				
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
		This field is parsed into 10 2-byte binary sub-fields. Only the first four are used. The first is set to 100; the second is set to 4; the third and fourth are also set to 100				
G 1	BL	BALANCE-CALIBRATION-DATA				
		Group embracing segment-specific fields of balance calibration transactions.				
2	BN	BALANCE-NUMBER	B	2.0	N	
		Number of the balance whose calibration weights are recorded in the transaction. Loaded by Update into the TDF-BALANCE-NUMBER-DE field of the TD-BALANCE-CALIBRATION file and, if it is a new balance, into the TDB-BALANCE-				

FILE...: BALANCE-CALIBRATION-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
NUMBER-UQ field of the TD-BALANCE file. Foreign key to TD-BALANCE.		
2 BP CALIBRATION-WEIGHT-MANUAL	A	1.0 N
Flag indicating whether the calibration weights were entered manually. Loaded by Update into the TDF-CALIBRATION-WEIGHT-MANUAL-FLAG field of the TD-BALANCE-CALIBRATION file.		
P 1 BQ BALANCE-CALIBRATION-WEIGHT-DATA		
Repeating group, which is not a group at all, containing only calibration weights.		
2 BR CALIBRATION-WEIGHT	B	2.0 N
The calibration weights recorded automatically or manually on LDAS. Loaded into the TDF-CALIBRATION-WEIGHT field of the TD-BALANCE-CALIBRATION file.		
G 1 EP CORRECTION-PORZION		
Group embracing segment-common correction fields. (There is no correction capability for calibration records either on LDAS or mainframe ECS.)		
2 EQ CORRECTION-DATE	P	5.0 N
2 ER CORRECTION-TIME	P	5.0 N
2 ES CORRECTION-REASON	A	2.0 N
2 ET CORRECTION-NOTE	A	79.0 N
2 EU CORRECTION-TYPE	A	1.0 N
2 FC CORRECTION-KEY	A	36.0 N DE
1 EW RECEIVE-UPDATE-KEY	A	54.0 N DE
Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDF-EXP-TEST-DATE-DE for/from TD-BALANCE-CALIBRATION, TOX-DATA-DATE, TOX-DATA-TIME, TOX-MAINFRAME-DATE, TOX-MAINFRAME-TIME, and a 4-digit sequence		

FILE...: BALANCE-CALIBRATION-TRANSACTION
 TYPE...: USER VIEW
 FILE-NR: 86
 PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

number for the transaction in
 its segment within the source
 set.

1	GB	BAL-CAL-SUPERKEY	A	37.0	N	SP
SOURCE FIELD(S) --- -START- --END-						
		DATA-BASE-KEY		3		9
		BALANCE-NUMBER		1		2
		TOX-DATA-DATE		1		2
		TOX-DATA-TIME		1		2
		RECEIVE-UPDATE-KEY		31		54

Super-key used to identify
 whether there is an existing
 record to which the calibra-
 should be appended.

BOTTLE-WEIGHT-TRANSACTION

Logical layout of bottle weight transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-CAGE-DATA table. For corrections, the information is also used to create records for BOTTLE-WEIGHT-CORRECTION.

FILE...: BOTTLE-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Set to BW. Used to form RECEIVE-UPDATE-KEY.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date				

FILE...: BOTTLE-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
---	----	--------------------	---	-----	---	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AL	AGENCY-NUMBER	B	2.0	F	
---	----	---------------	---	-----	---	--

Foreign key to TBAGENCY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
---	----	-----------------------	---	-----	---	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction and is loaded by Update into BWC-CORRECTION-OPERATOR-ID in BOTTLE-WEIGHT-CORRECTION. When appended to DATA-FACILITY-NUMBER a foreign key to TB_OPERATOR_ID.

2	AS	SOFTWARE-VERSION	B	2.0	F	
---	----	------------------	---	-----	---	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F	
---	----	---------------	---	-----	---	--

The date, in TDMS TOX date format, on which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Used by

FILE...: BOTTLE-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		Update to from the TDG-EXP-TEST-CAGE-DATE-DE value that will be added to, or selected for modification from, TD-CAGE-DATA. If a correction, also used to form BWC-BOTTLE-WGT-CORR-KEY-UQ of BOTTLE-WEIGHT-CORRECTION.				
2	AU	TOX-DATA-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, at which the transaction was created. Loaded by Update into TDG-WATER-TIME of TD-CAGE-DATA. If a correction, used to form BWC-BOTTLE-WGT-CORR-KEY-UQ of BOTTLE-WEIGHT-CORRECTION.				
2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
		The facility associated with the machine on which the source set was created. Not referenced by Update. Foreign key to TBFACILITY.				
2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
		The ID of the operator that created the transaction. Loaded by Update into TDG-WATER-OPERATOR-ID of TD-CAGE-DATA. If an add-with-audit, also loaded into BWC-BOTTLE-OPERATOR-ID in BOTTLE-WEIGHT-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
2	AX	PROCESS-FLAG	A	2.0	F	
		Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction,				

FILE...: BOTTLE-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		or 'A', indicating an audit transaction.				
2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
		The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 13.				
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
		The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 3.				
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
		Indicates sub-system. Set to 2 for in-life.				
2	A3	DATA-BASE-KEY	A	14.0	F	
		The number of the agency-test-cage to which the transaction applies. (AATTTTTTTCCCC) Used by Update to form the TDG-EXP-TEST-CAGE-DATE-DE value that will be added to, or selected for modification from, TD-CAGE-DATA. If a correspondence, also used to form BWC-BOTTLE-WGT-CORR-KEY-UQ for BOTTLE-WEIGHT-CORRECTION. The test-cage portions forms a foreign key to TD-CAGE.				
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
		This field is parsed into 10 2-byte binary sub-fields. Only the first three are used for bottle weight transactions. The first is set to 1000; the second is set to 10; the third is set to 15.				

FILE...: BOTTLE-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
G	1	B9 BOTTLE-WEIGHT-DATA				
		Group embracing segment-specific fields of bottle weight transactions.				
2	CB	OLD-BOTTLE-MANUAL-WEIGHT-FLAG	A	1.0	N	
		Flag indicating whether used bottle weight was recorded manually. 0 - No; 1 - Yes. Loaded by Update into TDG-MANUAL-ENTRY-FLAG-OLD-BOTT of TD-CAGE-DATA. If an add-with-audit, also loaded into BWC-MANUAL-ENTRY-FLAG-OLD-BOTT in BOTTLE-WEIGHT-CORRECTION.				
2	CD	OLD-BOTTLE-WEIGHT	B	4.0	N	
		Weight, in decigrams, of used bottle. Loaded by Update into TDG-OLD-BOTTLE-WEIGHT in TD-CAGE-DATA. If an add-with-audit, also loaded into BWC-OLD-BOTTLE-WEIGHT in BOTTLE-WEIGHT-CORRECTION.				
2	CF	NEW-BOTTLE-MANUAL-WEIGHT-FLAG	A	1.0	N	
		Flag indicating whether full bottle weight was recorded manually. 0 - No; 1 - Yes. Loaded by Update into TDG-MANUAL-ENTRY-FLAG-NEW-BOTT of TD-CAGE-DATA. If an add-with-audit, also loaded into BWC-MANUAL-ENTRY-FLAG-NEW-BOTT in BOTTLE-WEIGHT-CORRECTION.				
2	CH	NEW-BOTTLE-WEIGHT	B	4.0	N	
		Weight, in decigrams, of the full bottle. Loaded by Update into TDG-NEW-BOTTLE-WEIGHT in TD-CAGE-DATA. If an add-with-audit, also loaded into BWC-NEW-BOTTLE-WEIGHT in BOTTLE-WEIGHT-CORRECTION.				

FILE...: BOTTLE-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
2 CJ OLD-BOTTLE-STATUS	P	5.0 N
Weight status code indicating reason why used weight was not recorded. Loaded by Update into TDG-OLD-BOTTLE-STATUS in TD-CAGE-DATA. If an add-with audit, also loaded into BWC-OLD-BOTTLE-STATUS in BOTTLE-WEIGHT-CORRECTION. Foreign key to TBWGSTATUS.		
G 1 EP CORRECTION-PORION		
Group embracing segment-common correction fields.		
2 EQ CORRECTION-DATE	P	5.0 N
Date, in TDMS TOX date format, on which correction was made. Used by Update in forming BWC-BOTTLE-WGT-CORR-KEY-UQ for BOTTLE-WEIGHT-CORRECTION.		
2 ER CORRECTION-TIME	P	5.0 N
Time of day, in TDMS TOX time format, at which correction was made. Used by Update in forming BWC-BOTTLE-WGT-CORR-KEY-UQ for BOTTLE-WEIGHT-CORRECTION.		
2 ES CORRECTION-REASON	A	2.0 N
Code for correction reason. Loaded by Update into BWC-CORRECTION-REASON in BOTTLE-WEIGHT-CORRECTION. Foreign key to TBCORREASON.		
2 ET CORRECTION-NOTE	A	79.0 N
Free text describing the corresponds-rection. Loaded by Update into BWC-CORRECTION-NOTE in BOTTLE-WEIGHT-CORRECTION.		
2 EU CORRECTION-TYPE	A	1.0 N
Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update		

FILE...: BOTTLE-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

into BWC-CORRECTION-TYPE in
BOTTLE-WEIGHT-CORRECTION.

2	FA	CORRECTION-KEY	A	31.0	N	DE
---	----	----------------	---	------	---	----

If either the date portion of
TDG-EXP-TEST-CAGE-DATE-DE or
TDG-WATER-TIME is being
changed then RECEIVE or ECS
will concatenate the original
TDG-EXP-TEST-CAGE-DATE-DE,
the original TDD-WATER-TIME,
CORRECTION-DATE, and CORREC-
TION-TIME to form a string
that it loads into this field.
Loaded by Update into BWC-BOT-
TLE-WGT-CORR-KEY-LINK-UQ in
BOTTLE-WEIGHT-CORRECTION. The
first 16 characters form a
foreign key to TD-CAGE-DATA.

1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
---	----	--------------------	---	------	---	----

Unique key for database trans-
actions. Elements of key are:
SEGMENT-ID, MICRO-OR-ECS-FLAG,
TDG-EXP-TEST-CAGE-DATE-DE for/
from TD-CAGE-DATA left-justi-
fied into a 27-byte string,
TOX-DATA-DATE, TOX-DATA-TIME,
TOX-MAINFRAME-DATE, TOX-MAIN-
FRAME-TIME, and a 4-digit se-
quence number for the trans-
action in its segment within
the source set.

CAGE-ID-TRANSACTION

Logical layout for cage initialization transactions. Loaded by ECS and RECEIVE. Used by Update to create or modify data in TD-CAGE, and TD-TEST. For corrections, also used to add records to CAGE-ID-CORRECTION.

FILE...: CAGE-ID-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Set to 'CI'. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date				

FILE...: CAGE-ID-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
			format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AJ		TOX-MAINFRAME-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AL		AGENCY-NUMBER	B	2.0	F	
			Foreign key to TBAGENCY. Not referenced by Update.				
2	AP		ENTRY-OPERATOR-NUMBER	B	2.0	F	
			When the transaction is a corresponds-rection, this indicates the operator ID for the person making the correction and is loaded by Update into CIC-CORRECTION-OPERATOR-ID of CAGE-ID-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
2	AS		SOFTWARE-VERSION	B	2.0	F	
			The version of LDAS under which the source set was generated. Not referenced by Update.				
2	AT		TOX-DATA-DATE	B	2.0	F	
			The date, in TDMS TOX date format, on which the transaction was created. Loaded by Update into TDC-CAGE-DATE of TD-CAGE. If an add-with-audit, also loaded into CIC-CAGE-				

FILE...: CAGE-ID-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
DATE of CAGE-ID-CORRECTION.		
2 AU TOX-DATA-TIME	B	2.0 F
The time of day, in TDMS TOX time format, at which the transaction was created.		
Loaded by Update into TDC-CAGE-TIME. If an add-with-audit, also loaded into CIC-CAGE-TIME of CAGE-ID-CORRECTION.		
2 AV DATA-FACILITY-NUMBER	B	2.0 F
The facility associated with the machine on which the source set was created. Not referenced by Update. Foreign key to TBFACILITY.		
2 AW OPERATOR-OR-PATHOLOGIST	B	2.0 F
The ID of the operator who created the transaction.		
Loaded by Update into TDC-CAGE-OPERATOR-ID of TD-CAGE.		
If an add-with-audit, also loaded into CIC-CAGE-OPERATOR-ID of CAGE-ID-CORRECTION.		
When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.		
2 AX PROCESS-FLAG	A	2.0 F
Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.		
2 AZ DATA-BASE-KEY-LENGTH	B	2.0 F
The number of digits being used in the DATA-BASE-KEY field. Set to 13. Not referenced by Update.		

FILE...: CAGE-ID-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
		The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 2.				
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
		Indicates sub-system. Set to 2 for in-life.				
2	A3	DATA-BASE-KEY	A	14.0	F	
		The agency-test-cage number in the form AATTTTTTTTCCCC. The last 11 digits represent the TDC-EXP-TEST-CAGE-UQ value that will be added to, or selected for modification from, TD-CAGE. If a correction, also used to form CIC-CAGE-ID-CORR-KEY-UQ of CAGE-ID-CORRECTION. The test-cage portion is a foreign key to TD-CAGE.				
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
		This field is parsed into 10 2-byte binary sub-fields. Only 2 are used for cage ID transaction. The first is set to 10; the second is set to 1000. Neither is used by Update.				
G 1	A5	CAGE-ID-DATA				
		Group embracing fields specific to the cage ID segment.				
2	A6	TREATMENT-NUMBER	B	2.0	N	
		The number of the treatment group to which the cage is assigned. Loaded by Update into TDC-TREATMENT-NUMBER of TD-CAGE. If an add-with-aud-				

FILE...: CAGE-ID-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		it, also loaded into CIC-TREATMENT-NUMBER of CAGE-ID-CORRECTION. When appended to the agency-test portion of DATA-BASE-KEY, a foreign key to PAS_TREATMNT_SEG.				
2	A7	CONTROL-TREATMENT-FLAG	A	1.0	N	
		Alphabetic identifier of the role of the cage's treatment group. Not processed by Update. Foreign key to TBTRT-ROLE.				
2	A8	PROCEDURE-ACTION-SET-NUMBER	B	2.0	N	
		The procedure action set that applies to the cage, as defined in test protocol. Loaded by Update into TDC-PROCEDURE-ACTION-SET-NUM in TD-CAGE. If an add-with-audit, also loaded into CIC-PROCEDURE-ACTION-SET-NUM of CAGE-ID-CORRECTION. When appended to the agency-test portion of DATA-BASE-KEY, a foreign key to PAS_PROCEDUR_SEG.				
2	A9	TEST-TYPE	B	2.0	N	
		The test type for the test. Loaded by Update into TDT-TEST-TYPE of TD-TEST. A foreign key to TBTESTYP.				
2	BA	DOSE-DATE	B	2.0	N	
		The first dosing date for the cage, in TDMS TOX date format. Loaded by Update into TDC-START-DATE of TD-CAGE. If an add-with-audit, also loaded into CIC-START-DATE of CAGE-ID-CORRECTION.				
2	BB	AREA-NUMBER	B	2.0	N	
		The lab room where the cage is located, as defined in test protocol. Loaded by Update				

FILE...: CAGE-ID-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

into TDC-AREA-NUMBER of TD-CAGE. If an add-with-audit, also loaded into CIC-AREA-NUMBER of CAGE-ID-CORRECTION.

G 1 EP CORRECTION-PORION

Group embracing segment-common correction fields.

2 EQ CORRECTION-DATE P 5.0 N

Date, in TDMS TOX date format, on which correction was made. Used by Update to form CIC-CAGE-ID-CORR-KEY-UQ of CAGE-ID-CORRECTION.

2 ER CORRECTION-TIME P 5.0 N

Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form CIC-CAGE-ID-CORR-KEY-UQ in CAGE-ID-CORRECTION.

2 ES CORRECTION-REASON A 2.0 N

Code for correction reason. Loaded by Update into CIC-CORRECTION-REASON of CAGE-ID-CORRECTION. Foreign key to TBCORREASON.

2 ET CORRECTION-NOTE A 79.0 N

Free text describing the corresponds-rection. Loaded by Update into CIC-CORRECTION-NOTE of CAGE-ID-CORRECTION.

2 EU CORRECTION-TYPE A 1.0 N

Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into CIC-CORRECTION-TYPE of CAGE-ID-CORRECTION.

FILE...: CAGE-ID-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L	DB NAME	F	LENG	S	DE
2	EV CORRECTION-KEY	A	21.0	N	DE
	Not assigned a value by ECS or RECEIVE. Loaded by Update into CIC-CAGE-ID-CORR-LINK-UQ in CAGE-ID-CORRECTION.				
1	EW RECEIVE-UPDATE-KEY	A	54.0	N	DE
	Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDC-EXP-TEST-CAGE-UQ for/from TD-CAGE left-justified into a 27-byte string, TOX-DATA-DATE, TOX-DATA-TIME, TOX-MAINFRAME-DATE, TOX-MAINFRAME-TIME, and a 4-digit sequence number for the transactions in its segment-ment within the source set.				

CAGE-NOTE-TRANSACTION

Logical layout of cage note transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-CAGE-DATA file For corrections, the information is also used to add records to the CAGE-NOTE-TRANSACTION table.

FILE...: CAGE-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Set to CN. Used to form RECEIVE-UPDATE-KEY.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date				

FILE...: CAGE-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F		
---	----	--------------------	---	-----	---	--	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AL	AGENCY-NUMBER	B	2.0	F		
---	----	---------------	---	-----	---	--	--

Foreign key to TBAGENCY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F		
---	----	-----------------------	---	-----	---	--	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction and is loaded by Update into CNC-CORRECTION-OPERATOR-ID in CAGE-NOTE-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.

2	AS	SOFTWARE-VERSION	B	2.0	F		
---	----	------------------	---	-----	---	--	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F		
---	----	---------------	---	-----	---	--	--

The date, in TDMS TOX date format, on which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Used by Update to form the TDG-EXP-TEST-

FILE...: CAGE-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
CAGE-DATE-DE value that will be added to, or selected for modification from, TD-CAGE-DATA. If a correction, it is also used to form CNC-CAGE-NOTE-CORR-KEY-UQ for CAGE-NOTE-CORRECTION.		
2 AU TOX-DATA-TIME	B	2.0 F
The time of day, in TDMS TOX time format, at which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Loaded by Update into TDG-NOTE-TIME in TD-CAGE-DATA. If a correction, used by Update to form CNC-CAGE-NOTE-CORR-KEY-UQ for CAGE-NOTE-CORRECTION.		
2 AV DATA-FACILITY-NUMBER	B	2.0 F
The facility associated with the machine on which the source set was created. Not referenced by Update. A foreign key to TBFACILITY.		
2 AW OPERATOR-OR-PATHOLOGIST	B	2.0 F
The ID of the operator who created the transaction. Loaded by Update into TDG-NOTE-OPERATOR-ID. If an add-with-audit, also loaded into CNC-CAGE-NOTE-OPERATOR-ID. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.		
2 AX PROCESS-FLAG	A	2.0 F
Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction,		

FILE...: CAGE-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

or 'A', indicating an audit transaction.

2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
---	----	----------------------	---	-----	---	--

The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 13.

2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
---	----	---------------------------	---	-----	---	--

The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 2.

2	A2	DATA-BASE-INDICATOR	A	1.0	F	
---	----	---------------------	---	-----	---	--

Indicates sub-system. Set to 2 for in-life.

2	A3	DATA-BASE-KEY	A	14.0	F	
---	----	---------------	---	------	---	--

The number of the agency-test-cage to which the transaction applies. (AATTTTTTCCCC)
Used by Update to form the TDG-EXP-TEST-CAGE-DATE-DE value that will be added to, or selected for modification from, TD-CAGE-DATA. If a correspondence, also used to form CNC-CAGE-NOTE-CORR-KEY-UQ for CAGE-NOTE-CORRECTION. The test-cage portion forms a foreign key to TD-CAGE.

2	A4	SEGMENT-TREE-KEY	B	20.0	F	
---	----	------------------	---	------	---	--

This field is parsed into 10 2-byte binary sub-fields. Only the first two are used for cage note transactions. The first is set to 1000; the second is set to 80.

FILE...: CAGE-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
G 1 CT CAGE-NOTE-DATA		
Group embracing segment-specific fields of cage note transactions.		
2 CV CAGE-NOTE	A	80.0 N
Free text entered by operator modifying or explaining cage data. Loaded by Update into TDG-CAGE-NOTE of TD-CAGE-DATA. If an add-with-audit, also loaded into CNC-CAGE-NOTE in CAGE-NOTE-CORRECTION.		
2 CW CAGE-NOTE-MANUAL-FLAG	A	1.0 N
Flag indicating whether note was entered as part of a manual entry session. 0 - No; 1 - Yes. Loaded by Update into TDG-CAGE-NOTE-MANUAL-FLAG of TD-CAGE-DATA. If an add-with-audit, also loaded into CNC-CAGE-NOTE-MANUAL-FLAG in CAGE-NOTE-CORRECTION.		
G 1 EP CORRECTION-PORITION		
Group embracing segment-common correction fields.		
2 EQ CORRECTION-DATE	P	5.0 N
Date, in TDMS TOX date format, on which correction was made. Used by Update to form CNC-CAGE-NOTE-CORR-KEY-UQ for CAGE-NOTE-CORRECTION.		
2 ER CORRECTION-TIME	P	5.0 N
Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form CNC-CAGE-NOTE-CORR-KEY-UQ for CAGE-NOTE-CORRECTION.		
2 ES CORRECTION-REASON	A	2.0 N
Code for correction reason. Loaded by Update into CNC-CORRECTION-REASON in CAGE-NOTE-CORRECTION. A foreign key to		

FILE...: CAGE-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		TBCORREASON.				
2	ET	CORRECTION-NOTE	A	79.0	N	
		Free text describing the corresponds- rection. Loaded by Update in- to CNC-CORRECTION-NOTE in CAGE-NOTE-CORRECTION.				
2	EU	CORRECTION-TYPE	A	1.0	N	
		Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into CNC-CORRECTION-TYPE in CAGE-NOTE-CORRECTION.				
2	FA	CORRECTION-KEY	A	31.0	N	DE
		If either the date portion of TDG-EXP-TEST-CAGE-DATE-DE or TDG-NOTE-TIME is being changed then RECEIVE or ECS will concatenate the original TDG-EXP-TEST-CAGE-DATE-DE, the original TDD-NOTE-TIME, CORRECTION-DATE, and CORREC- TION-TIME to form a string that it loads into this field. Loaded by Update into CNC- CAGE-NOTE-CORR-KEY-LINK-UQ in CAGE-NOTE-CORRECTION.				
1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
		Unique key for database trans- actions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDG-EXP-TEST-CAGE-DATE-DE for/ from TD-CAGE-DATA, TOX-DATA- DATE, TOX-DATA-TIME, TOX-MAIN- FRAME-DATE, TOX-MAINFRAME- TIME, and a 4-digit sequence number for the transaction in its segment within the source set.				

FEEDER-WEIGHT-TRANSACTION

Logical layout of feeder weight transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records from the TD-CAGE-DATA file. For corrections, information used to add records to FEEDER-WEIGHT-CORRECTION.

FILE...: FEEDER-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Set to FW. Used to form RECEIVE-UPDATE-KEY.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date				

FILE...: FEEDER-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
---	----	--------------------	---	-----	---	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AL	AGENCY-NUMBER	B	2.0	F	
---	----	---------------	---	-----	---	--

Foreign key to TBAGENCY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
---	----	-----------------------	---	-----	---	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction and is loaded by Update into FWC-CORRECTION-OPERATOR-ID in FEEDER-WEIGHT-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.

2	AS	SOFTWARE-VERSION	B	2.0	F	
---	----	------------------	---	-----	---	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F	
---	----	---------------	---	-----	---	--

The date, in TDMS TOX date format, on which the transaction was created. Used by Update to form the TDG-EXP-TEST-CAGE-DATE-DE value that

FILE...: FEEDER-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
will be added to, or selected for modification from, TD-CAGE-DATA. If a correction, it is also used to form FWC-FEEDER-WGT-CORR-KEY-UQ for FEEDER-WEIGHT-CORRECTION.		
2 AU TOX-DATA-TIME	B	2.0 F
The time of day, in TDMS TOX time format, at which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Loaded by Update into TDG-FEEDER-TIME in TD-CAGE-DATA. If an add-with-audit, also used to form FWC-FEEDER-WGT-CORR-KEY-UQ of FEEDER-WEIGHT-CORRECTION.		
2 AV DATA-FACILITY-NUMBER	B	2.0 F
The facility associated with the machine on which the source set was created. Not referenced by Update. Foreign key to TBFACILITY.		
2 AW OPERATOR-OR-PATHOLOGIST	B	2.0 F
The ID of the operator who created the transaction. Loaded by Update into TDG-FEEDER-OPERATOR-ID of TD-CAGE-DATA. If an add-with-audit, also loaded into FWC-FEEDER-OPERATOR-ID in FEEDER-WEIGHT-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.		
2 AX PROCESS-FLAG	A	2.0 F
Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction,		

FILE...: FEEDER-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		or 'A', indicating an audit transaction.				
2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
		The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 13.				
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
		The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 3.				
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
		Indicates sub-system. Set to 2 for in-life.				
2	A3	DATA-BASE-KEY	A	14.0	F	
		The agency-test-cage number to which the transaction applies. (AATTTTTTTTCCCC) Used by Update to form the TDG-EXP-TEST-CAGE-DATE-DE value that will be added to, or selected for modification from, TD-CAGE-DATA. If a correction, also used to form FWC-FEEDER-WGT-CORR-KEY-UQ for FEEDER-WEIGHT-CORRECTION. The test-cage portion forms a foreign key to TD-CAGE.				
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
		This field is parsed into 10 2-byte binary sub-fields. Only 3 are used for feeder weight transactions. The first is set to 1000; the second is set to 10; the third is set to 15.				

FILE...: FEEDER-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
G	1	B9 FEEDER-WEIGHT-DATA				
		Group embracing segment-specific fields of feeder weight transactions.				
2	CB	OLD-FEEDER-MANUAL-WEIGHT-FLAG	A	1.0	N	
		Flag indicating whether used feeder weight was recorded manually. 0 - Yes; 1 - No. Loaded by Update into TDG-MANUAL-ENTRY-FLAG-OLD-FEED in TD-CAGE-DATA. If an add-with-audit, also loaded into FWC-MANUAL-ENTRY-FLAG-OLD-FEED in FEEDER-WEIGHT-CORRECTION.				
2	CD	OLD-FEEDER-WEIGHT	B	4.0	N	
		Weight, in decigrams, of used feeder. Loaded by Update into TDG-OLD-FEEDER-WEIGHT in TD-CAGE-DATA. If an add-with-audit, also loaded into FWC-OLD-FEEDER-WEIGHT in FEEDER-WEIGHT-CORRECTION.				
2	CF	NEW-FEEDER-MANUAL-WEIGHT-FLAG	A	1.0	N	
		Flag indicating whether full feeder weight was recorded manually. 0 - No; 1 - Yes. Loaded by Update into TDG-MANUAL-ENTRY-FLAG-NEW-FEED in TD-CAGE-DATA. If an add-with-audit, also loaded into FWC-MANUAL-ENTRY-FLAG-NEW-FEED in FEEDER-WEIGHT-CORRECTION.				
2	CH	NEW-FEEDER-WEIGHT	B	4.0	N	
		Weight, in decigrams, of the full feeder. Loaded by Update into TDG-NEW-FEEDER-WEIGHT in TD-CAGE-DATA. If an add-with-audit, also loaded into FWC-NEW-FEEDER-WEIGHT in FEEDER-WEIGHT-CORRECTION.				
2	CJ	OLD-FEEDER-STATUS	P	5.0	N	
		Weight status code indicating				

FILE...: FEEDER-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		reason why used weight was not recorded. Loaded by Update into TDG-OLD-FEEDER-STATUS in TD-CAGE-DATA. If an add-with-audit, also loaded into FWC-OLD-FEEDER-STATUS in FEEDER-WEIGHT-CORRECTION. Foreign key to TBWGTSTATUS.				
G 1	EP	CORRECTION-PORZION				
		Group embracing segment-common correction fields.				
2	EQ	CORRECTION-DATE	P	5.0	N	
		Date, in TDMS TOX date format, on which correction was made. Used by Update to form FWC-FEEDER-WGT-CORR-KEY-UQ for FEEDER-WEIGHT-CORRECTION.				
2	ER	CORRECTION-TIME	P	5.0	N	
		Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form FWC-FEEDER-WGT-CORR-KEY-UQ for FEEDER-WEIGHT-CORRECTION.				
2	ES	CORRECTION-REASON	A	2.0	N	
		Code for correction reason. Loaded by Update into FWC-CORRECTION-REASON in FEEDER-WEIGHT-CORRECTION. Foreign key to TBCORREASON.				
2	ET	CORRECTION-NOTE	A	79.0	N	
		Free text describing the corresponds-rection. Loaded by Update into FWC-CORRECTION-NOTE in FEEDER-WEIGHT-CORRECTION.				
2	EU	CORRECTION-TYPE	A	1.0	N	
		Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into FWC-CORRECTION-TYPE in FEEDER-WEIGHT-CORRECTION.				
2	FA	CORRECTION-KEY	A	31.0	N	DE

FILE...: FEEDER-WEIGHT-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F LENG S DE
-------------	-------------

If either the date portion of TDG-EXP-TEST-CAGE-DATE-DE or TDG-FEEDER-TIME is being changed then RECEIVE or ECS will concatenate the original TDG-EXP-TEST-CAGE-DATE-DE, the original TDD-FEEDER-TIME, CORRECTION-DATE, and CORRECTION-TIME to form a string that it loads into this field. Loaded by Update into FWC-FEEDER-WGT-CORR-KEY-LINK-UQ in FEEDER-WEIGHT-CORRECTION.

1 EW RECEIVE-UPDATE-KEY	A 54.0 N DE
-------------------------	-------------

Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDG-EXP-TEST-CAGE-DATE-DE for/ from TD-CAGE-DATA, TOX-DATA-DATE, TOX-DATA-TIME, TOX-MAIN-FRAME-DATE, TOX-MAINFRAME-TIME, and a 4-digit sequence number for the transaction in its segment within the source set.

HISTOLOGY-NUMBER-TRANSACTION

Logical layout of histology number transactions. Loaded by RECEIVE and ECS. Read by Update to modify records in the TD-ANIMAL file. For corrections, also used to add records to HISTOLOGY-NUMBER-CORRECTION.

FILE...: HISTOLOGY-NUMBER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Set to HN. Used to form RECEIVE-UPDATE-KEY.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date				

FILE...: HISTOLOGY-NUMBER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
			format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AJ		TOX-MAINFRAME-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AL		AGENCY-NUMBER	B	2.0	F	
			Foreign key to TBAGENCY. Not referenced by Update.				
2	AP		ENTRY-OPERATOR-NUMBER	B	2.0	F	
			When the transaction is a corresponds-rection, this indicates the operator ID for the person making the correction and is loaded by Update into HIC-CORRECTION-OPERATOR-ID in HISTOLOGY-NUMBER-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
2	AS		SOFTWARE-VERSION	B	2.0	F	
			The version of LDAS under which the source set was generated. Not referenced by Update.				
2	AT		TOX-DATA-DATE	B	2.0	F	
			The date, in TDMS TOX date format, on which the transaction was created. Loaded by Update into TDA-HISTOLOGY-NUMBER-DATE in TD-ANIMAL. If an				

FILE...: HISTOLOGY-NUMBER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
add-with-audit, also loaded into HIC-HISTOLOGY-NUMBER-DATE in HISTOLOGY-NUMBER-CORREC- TION.		
2 AU TOX-DATA-TIME	B	2.0 F
The time of day, in TDMS TOX time format, at which the transaction was created. Loaded by Update into TDA-HIS- TOLOGY-NUMBER-TIME in TD-ANI- MAL. If a correction, also loaded into HIC-HISTOLOGY-NUM- BER-TIME.		
2 AV DATA-FACILITY-NUMBER	B	2.0 F
The facility associated with the machine on which the source set was created. Not referenced by Update. A for- eign key to TBFACLTy.		
2 AW OPERATOR-OR-PATHOLOGIST	B	2.0 F
The ID of the pathologist as- signed to the test. Loaded by Update into TDA-PATHOLOGIST- IDENTIFICATION in TD-ANIMAL. If an add-with-audit, also loaded into HIC-PATHOLOGIST- IDENTIFICATION in HISTOLOGY- NUMBER-CORRECTION. When ap- pended to DATA-FACILITY-NUM- BER, a foreign key to TB_OPERATOR_ID.		
2 AX PROCESS-FLAG	A	2.0 F
Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or del- ete, respectively. The second digit is either blank, indica- ting an original transaction, or 'A', indicating an audit transaction.		
2 AZ DATA-BASE-KEY-LENGTH	B	2.0 F

FILE...: HISTOLOGY-NUMBER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update. Set to 14.				
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
		The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 1.				
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
		Indicates sub-system. Set to 3 for pathology.				
2	A3	DATA-BASE-KEY	A	14.0	F	
		The carcass ID (CID) for the animal being assigned an histology number. Used by Update to select the TD-ANIMAL record to which the transaction applies. For corrections, also used to form HIC-HISTOLOGY-CORR-KEY-UQ for HISTOLOGY-NUMBER-CORRECTION.				
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
		This field is parsed into 10 2-byte binary sub-fields. Only the first one is used for histology number transactions, and it is assigned a value of 2000.				
G	1	DJ HISTOLOGY-NUMBER-DATA				
		Group embracing segment-specific fields of histology number transactions.				
2	DK	HISTOLOGY-OPERATOR	B	2.0	N	
		ID for operator assigning histology number. Loaded by Update into TDA-HISTOLOGY-NUM-				

FILE...: HISTOLOGY-NUMBER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
<p>OPERATOR-ID in TD-ANIMAL. If an add-with-audit, also loaded into HIC-HISTOLOGY-NUM-OPERATOR-ID of HISTOLOGY-NUMBER-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.</p>		
2 DL HISTOLOGY-NUMBER	A	16.0 N
<p>Lab-specific identification number for carcasses. Loaded by Update into TDA-HISTOLOGY-NUMBER in TD-ANIMAL. If an add-with-audit, also loaded into HIC-HISTOLOGY-NUMBER of HISTOLOGY-NUMBER-CORRECTION.</p>		
G 1 EP CORRECTION-PORION		
<p>Group embracing segment-common correction fields.</p>		
2 EQ CORRECTION-DATE	P	5.0 N
<p>Date, in TDMS TOX date format, on which correction was made. Used by Update to form HIC-HISTOLOGY-CORR-KEY-UQ for HISTOLOGY-NUMBER-CORRECTION.</p>		
2 ER CORRECTION-TIME	P	5.0 N
<p>Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form HIC-HISTOLOGY-CORR-KEY-UQ for HISTOLOGY-NUMBER-CORRECTION.</p>		
2 ES CORRECTION-REASON	A	2.0 N
<p>Code for correction reason. Loaded by Update into HIC-CORRECTION-REASON in HISTOLOGY-NUMBER-CORRECTION. A foreign key to TBCORREASON.</p>		
2 ET CORRECTION-NOTE	A	79.0 N
<p>Free text describing the corresponds-rection. Loaded by Update into HIC-CORRECTION-NOTE in HISTOLOGY-NUMBER-CORRECTION.</p>		

FILE...: HISTOLOGY-NUMBER-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	EU	CORRECTION-TYPE	A	1.0	N	
		Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into HIC-CORRECTION-TYPE in HISTOLOGY-NUMBER-CORRECTION.				
2	EX	CORRECTION-KEY	A	26.0	N	DE
		Not loaded by RECEIVE or ECS. Loaded by Update into HIC-HISTOLOGY-CORR-KEY-LINK-UQ in HISTOLOGY-NUMBER-CORRECTION.				
1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
		Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDA-EXP-TEST-ANIMAL-UQ, 11 spaces, TOX-DATA-DATE, TOX-DATA-TIME, TOX-MAINFRAME-DATE, TOX-MAINFRAME-TIME, and a 4-digit sequence number for the transaction in its segment within the source set.				

MICRO-SITE-STATUS-TRANSACTION

Logical layout for site status transactions. Site status data was a feature of the MODCOMP data collection that was made defunct by LDAS. There is an ECS module that will add records to this table, but it is incorrect procedure to use that module. Correct procedure calls for a test to be converted to LDAS format prior to modification.

FILE...: MICRO-SITE-STATUS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
- - - - -						
G	1	AC HEADER				
		Group embracing fields that are common to most segments.				
2	AA	SEGMENT-ID	A	2.0	F	
		Two-letter segment identifier. Not referenced by Update. Used by ECS to form RECEIVE-UPDATE-KEY. Set to SS.				
2	AB	MICRO-OR-ECS-FLAG	A	1.0	F	
		Identifier for source of transaction. Set to E for Mainframe ECS. Used by Update to generate Update statistics report. Used by ECS to form RECEIVE-UPDATE-KEY.				
2	AE	TOX-MICRO-DATE	B	2.0	F	
		The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF	TOX-MICRO-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG	MICRO-TRANS-SET	B	2.0	F	
		The number of the transaction's source set. Not referenced by Update.				

FILE...: MICRO-SITE-STATUS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	AI	TOX-MAINFRAME-DATE	B	2.0	F	
		The date, in TDMS TOX date format, on which the source set was received. Used by ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.				
2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set was received. Used by ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update				
2	AL	AGENCY-NUMBER	B	2.0	F	
		Foreign key to TBAGENCY. Not referenced by Update.				
2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
		The operator ID for the person making the correction. It is loaded by Update into SSC-CORRECTION-OEPRATOR-ID in MICRO-SITE-STATUS-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
2	AS	SOFTWARE-VERSION	B	2.0	F	
		The version of LDAS under which the source set was generated. Not referenced by Update.				
2	AT	TOX-DATA-DATE	B	2.0	F	
		The date, in TDMS TOX date format, on which the transaction was created. Loaded by Update into TDO-SITE-STATUS-DATE in TD-ORGAN. Used to form SSC-MICRO-SITE-CORR-KEY-UQ for MICRO-SITE-STATUS-CORRECTION. Used by ECS to form				

FILE...: MICRO-SITE-STATUS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
RECEIVE-UPDATE-KEY.		
2 AU TOX-DATA-TIME	B	2.0 F
<p>The time of day, in TDMS TOX time format, at which the transaction was created.</p> <p>Loaded by Update into TDO-SITE-STATUS-TIME in TD-ORGAN.</p> <p>Used to form SCC-MICRO-SITE-CORR-KEY-UQ for MICRO-SITE-CORRECTION. Used by RECEIVE to form RECEIVE-UDPATE-KEY.</p>		
2 AV DATA-FACILITY-NUMBER	B	2.0 F
<p>The facility associated with the machine on which the source set was created. Not referenced by Update. A foreign key to TBFACLT.</p>		
2 AW OPERATOR-OR-PATHOLOGIST	B	2.0 F
<p>The operator ID for the pathologist assigned to the test.</p> <p>Loaded by Update into TDO-SITE-STATUS-PATHOLOGIST-ID in TD-ORGAN. If an add-with-audit, also loaded into SSC-SITE-STATUS-PATHOLOGIST-ID in MICRO-SITE-STATUS-CORRECTION.</p> <p>When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.</p>		
2 AX PROCESS-FLAG	A	2.0 F
<p>Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.</p>		
2 AZ DATA-BASE-KEY-LENGTH	B	2.0 F
<p>The number of digits being used in the DATA-BASE-KEY</p>		

FILE...: MICRO-SITE-STATUS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

field. Not referenced by Update. Set to 14.

2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
---	----	---------------------------	---	-----	---	--

The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used. Set to 5.

2	A2	DATA-BASE-INDICATOR	A	1.0	F	
---	----	---------------------	---	-----	---	--

Indicates sub-system. Set to 3 for pathology.

2	A3	DATA-BASE-KEY	A	14.0	F	
---	----	---------------	---	------	---	--

The carcass ID (CID) for the animal from which the organ being observed was taken. Used by Update to select the master TD-ANIMAL record for the transaction.

2	A4	SEGMENT-TREE-KEY	B	20.0	F	
---	----	------------------	---	------	---	--

This field is parsed into 10 2-byte binary sub-fields. The second is for pathologist type and is set to "A"; the fourth is the organ (a foreign key to TBPCT) where the site is. Both values are used to form TDO-EXP-TST-ANUM-ORG-PATHT-UQ that will be used by Update to select the record to which the site status transaction will be applied. Also used to from the SSC-MICRO-SITE-CORR-KEY-UQ for MICRO-SITE-STATUS-CORRECTION.

G	1	EL	MICRO-SITE-STATUS-DATA			
---	---	----	------------------------	--	--	--

Legacy group embracing all segment-specific fields in site status transactions.

FILE...: MICRO-SITE-STATUS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	EM	MICRO-SITE-STATUS-OPERATOR-ID	B	2.0	N	
Id for operator assigning site status. Loaded by Update into TDO-SITE-STATUS-OPERATOR-ID of TD-ORGAN. If an add-with-audit, also loaded into SSC-SITE-STATUS-OPERATOR-ID in MICRO-SITE-STATUS-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.						
2	EN	MICRO-SITE	B	2.0	N	
Organ-Site for which a status is being assigned. Loaded by Update into TDO-SITE-CODE in TD-ORGAN. If an add-with-audit, also loaded into SSC-SITE-CODE in MICRO-SITE-STATUS-CORRECTION. A foreign key to TBPCT.						
2	EO	MICRO-SITE-STATUS	B	2.0	N	
Status of the organ-site being evaluated. Loaded by Update into TDO-SITE-STATUS in TD-ORGAN. If an add-with-audit, also loaded into SSC-SITE-STATUS in MICRO-SITE-STATUS-CORRECTION. A foreign key to TBORGSTA.						
G 1	EP	CORRECTION-PORTION				
Group embracing segment-common correction fields.						
2	EQ	CORRECTION-DATE	P	5.0	N	
Date, in TDMS TOX date format, on which correction was made. Used by Update to form SSC-MICRO-SITE-CORR-KEY-UQ for MICRO-SITE-STATUS-CORRECTION.						
2	ER	CORRECTION-TIME	P	5.0	N	
Time of day, in TDMS TOX time format, at which correction was made. Used by Update to						

FILE...: MICRO-SITE-STATUS-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		form SSC-MICRO-SITE-CORR-KEY- UQ for MICRO-SITE-STATUS-COR- RECTION.				
2	ES	CORRECTION-REASON	A	2.0	N	
		Code for correction reason. Loaded by Update into SSC-COR- RECTION-REASON of MICRO-SITE- STATUS-CORRECTION. Foreign key to TBCORREASON.				
2	ET	CORRECTION-NOTE	A	79.0	N	
		Free text describing the corresponds- rection. Loaded by Update in- to SSC-CORRECTION-NOTE of MIC- RO-SITE-STATUS-CORRECTION.				
2	EU	CORRECTION-TYPE	A	1.0	N	
		Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into SSC-CORRECTION-TYPE of MICRO-SITE-STATUS-CORRECTION.				
2	FB	CORRECTION-KEY	A	32.0	N	DE
		Not assigned a value by ECS or RECEIVE. Loaded by Update in- to AMC-MICRO-ORG-CORR-KEY- LINK-UQ in ANY-MICRO-ORGAN- CORRECTION.				
1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
		Unique key for database trans- actions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDO-EXP-TST-ANUM-ORG-PATHT-UQ, 5 spaces, TOX-DATA-DATE, TOX- DATA-TIME, TOX-MAINFRAME-DATE, TOX-MAINFRAME-TIME, and a 4- digit sequence number for the transaction in its segment within the source set.				

NOT-EXAMINED-TRANSACTION

Logical layout of "not examined" transactions. Loaded by RECEIVE and ECS. Read by Update to modify records in the TD-ANIMAL file. If a correction, information also used to add records to NOT-EXAMINED-CORRECTION.

FILE...: NOT-EXAMINED-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-

G 1 AC HEADER

Group embracing fields that are common to most segments.

2	AA	SEGMENT-ID	A	2.0	F
---	----	------------	---	-----	---

Two-letter segment identifier. Not referenced by Update. Set to NE. Used to form RECEIVE-UPDATE-KEY.

2	AB	MICRO-OR-ECS-FLAG	A	1.0	F
---	----	-------------------	---	-----	---

Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.

2	AE	TOX-MICRO-DATE	B	2.0	F
---	----	----------------	---	-----	---

The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.

2	AF	TOX-MICRO-TIME	B	2.0	F
---	----	----------------	---	-----	---

The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.

2	AG	MICRO-TRANS-SET	B	2.0	F
---	----	-----------------	---	-----	---

The number of the transaction's source set. Not referenced by Update.

2	AI	TOX-MAINFRAME-DATE	B	2.0	F
---	----	--------------------	---	-----	---

The date, in TDMS TOX date

FILE...: NOT-EXAMINED-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F		
---	----	--------------------	---	-----	---	--	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F		
---	----	-----------------------	---	-----	---	--	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction and is loaded by Update into NEC-CORRECTION-OPERATOR-ID in NOT-EXAMINED-CORRECTION. When appended to the appropriate facility code, a foreign code to TB_OPERATOR_ID.

2	AS	SOFTWARE-VERSION	B	2.0	F		
---	----	------------------	---	-----	---	--	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F		
---	----	---------------	---	-----	---	--	--

The date, in TDMS TOX date format, on which the transaction was created. Loaded by Update into TDA-NO-EXAM-DATE in TD-ANIMAL. If an add-with-audit, also loaded into NEC-NO-EXAM-DATE in NOT-EXAMINED-CORRECTION.

FILE...: NOT-EXAMINED-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	AU	TOX-DATA-TIME	B	2.0	F	
<p>The time of day, in TDMS TOX time format, at which the transaction was created.</p> <p>Loaded by Update into TDA-NO-EXAM-TIME in TD-ANIMAL. If an add-with-audit, also loaded into NEC-NO-EXAM-TIME in NOT-EXAMINED-TRANSACTION.</p>						
2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
<p>The ID for the operator who created the transaction.</p> <p>Loaded by Update into TDA-NO-EXAM-OPERATOR-ID in TD-ANIMAL.</p> <p>If an add-with-audit, also loaded into NEC-NO-EXAM-OPERATOR-ID in NOT-EXAMINED-CORRECTION. When appended to the appropriate facility code, a foreign key to TB_OPERATOR-ID.</p>						
2	AX	PROCESS-FLAG	A	2.0	F	
<p>Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.</p>						
G 1	FF	NO-EXAM-DATA				
<p>Group embracing segment-specific fields in "not examined" transactions.</p>						
2	FG	NO-EXAM-STUDY	A	7.0	N	
<p>TDMS test number to which the observation applies. Used by Update to select the TD-ANIMAL records to which the transaction applies. A foreign key to TD-TEST.</p>						

FILE...: NOT-EXAMINED-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
2 FH NO-EXAM-ANIMAL-NUM Number of animal being classified as not examined. Not referenced by Update.	P	5.0 N
2 FI NO-EXAM-PATHOLOGIST Operator ID of the pathologist assigned to the test. Loaded by Update into TDA-NO-EXAM-PATHOLOGIST in TD-ANIMAL. If an add-with-audit, also loaded into NEC-NO-EXAM-PATH-ID in NOT-EXAMINED-CORRECTION. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.	P	5.0 N
G 1 EP CORRECTION-PORION Group embracing segment-common correction fields.		
2 EQ CORRECTION-DATE Date, in TDMS TOX date format, on which correction was made. Used by Update in forming NEC-NEC-NO-EXAM-CORR-KEY-UQ for NOT-EXAMINED-CORRECTION.	P	5.0 N
2 ER CORRECTION-TIME Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form NEC-NO-EXAM-CORR-KEY-UQ for NOT-EXAMINED-CORRECTION.	P	5.0 N
2 ES CORRECTION-REASON Code for correction reason. Loaded by Update into NEC-CORRECTION-REASON in NOT-EXAMINED-CORRECTION. A foreign key to TBCORREASON.	A	2.0 N

FILE...: NOT-EXAMINED-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	ET	CORRECTION-NOTE	A	79.0	N	
Free text describing the corresponds- rection. Loaded by Update in- to NEC-CORRECTION-NOTE in NOT- EXAMINED-CORRECTION.						
2	EU	CORRECTION-TYPE	A	1.0	N	
Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into NEC-CORRECTION-TYPE in NOT-EXAMINED-CORRECTION.						
2	EX	CORRECTION-KEY	A	26.0	N	DE
Not loaded by RECEIVE and ECS. Loaded by Update into NEC-NO- EXAM-CORR-KEY-LINK-UQ in NOT- EXAMINED-CORRECTION.						
1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
Unique key for database trans- actions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDA-EXP-TEST-ANIMAL-UQ, 11 spaces, TOX-DATA-DATE, TOX- DATA-TIME, TOX-MAINFRAME-DATE, TOX-MAINFRAME-TIME, and a 4- digit sequence number for the transaction in its segment within the source set. Char- acters 4-19 form a foreign key to TD-ANIMAL.						

ORGAN-NOTES-TRANSACTION

Logical layout of organ notes transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-ORGAN-DATA file. If transactions are audit changes, then Update also uses the information to add records to ORGAN-NOTES-CORRECTION.

FILE...: ORGAN-NOTES-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
- - - - -						
G	1	AC HEADER				
		Group embracing fields that are common to most segments.				
2	AA	SEGMENT-ID	A	2.0	F	
		Two-letter segment identifier. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY. Set to ON.				
2	AB	MICRO-OR-ECS-FLAG	A	1.0	F	
		Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE	TOX-MICRO-DATE	B	2.0	F	
		The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF	TOX-MICRO-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG	MICRO-TRANS-SET	B	2.0	F	
		The number of the transaction's source set. Not referenced by Update.				

FILE...: ORGAN-NOTES-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	AI	TOX-MAINFRAME-DATE	B	2.0	F	
		The date, in TDMS TOX date format, on which the source set was received. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY.				
2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set was received. Not referenced by Update. Used to form RECEIVE-UPDATE-KEY.				
2	AL	AGENCY-NUMBER	B	2.0	F	
		Always 5, for NIH/NIEHS/NTP Not referenced by Update.				
2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
		When the transaction is a corresponds- rection, this indicates the operator ID for the person making the correction and is loaded by Update into OMC- CORRECTION-OPERATOR-ID in ORGAN-NOTE-CORRECTION. When appended to DATA-FACILITY-NUM- BER, a foreign key to TB_OPER- ATOR_ID.				
2	AS	SOFTWARE-VERSION	B	2.0	F	
		The version of LDAS under which the source set was generated. Not referenced by Update.				
2	AT	TOX-DATA-DATE	B	2.0	F	
		The date, in TDMS TOX date format, on which the trans- action was created. Used to form RECEIVE-UPDATE-KEY. Used by Update to form the value of TDH-ORGAN-OBSERVA-				

FILE...: ORGAN-NOTES-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		TION-ID-DE that will be added to, or selected for modification, TD-ORGAN-DATA. For corrections, this value will be used to form ONC-ORGAN-NOTE-CORR-KEY-UQ for ORGAN-NOTES-CORRECTION.				
2	AU	TOX-DATA-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, at which the transaction was created. Used to form RECEIVE-UPDATE-KEY. Loaded by Update into TDH-ORGAN-NOTE-SEQUENCE. For correspondences, it is also used to form the value of ONC-ORGAN-NOTE-CORR-KEY-UQ for ORGAN-NOTES-CORRECTION.				
2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
		The facility associated with the machine on which the source set was created. Not referenced by Update. A foreign key to TBFACILITY.				
2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
		The operator ID for the pathologist assigned to the test. Loaded by Update into TDH-ORGAN-NOTE-PATHOLOGIST-ID in TD-ORGAN-DATA. For an add-with-audit, Update loads the field into ONC-ORGAN-NOTE-PATHOLOGIST-ID in ORGAN-NOTES-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR_ID.				
2	AX	PROCESS-FLAG	A	2.0	F	
		Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second				

FILE...: ORGAN-NOTES-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.				
2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F	
		The number of digits being used in the DATA-BASE-KEY field. Not referenced by Update.				
2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F	
		The number of two-byte binary segments of the SEGMENT-TREE-KEY field being used.				
2	A2	DATA-BASE-INDICATOR	A	1.0	F	
		Indicates sub-system. Set to 3 for pathology.				
2	A3	DATA-BASE-KEY	A	14.0	F	
		The carcass ID for the animal to which the observation applies. Used by Update to retrieve the master record in TD-ANIMAL.				
2	A4	SEGMENT-TREE-KEY	B	20.0	F	
		This field is parsed into 10 2-byte binary sub-fields. The second is set to 'A' for pathologist type; the fourth is the organ (foreign key to TBPCT). Both are used to form the value of TDH-ORGAN-OBSERVATION-ID-DE that will be added or selected for modification, TD-ORGAN-DATA. For corrections, also used to form ONC-ORGAN-NOTE-CORR-KEY-UQ for ORGAN-NOTES-CORRECTION.				

FILE...: ORGAN-NOTES-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
G	1	D5 ORGAN-NOTES-DATA				
		Group embracing fields specific to organ note transactions				
	2	D6 ORGAN-NOTES-OPERATOR	B	2.0	N	
		The ID for the operator who recorded the transaction. Loaded by Update into TDH-ORGAN-NOTE-OPERATOR-ID in TD-ORGAN-DATA. If an add-with-audit, also loaded into ONC-ORGAN-NOTE-OPERATOR-ID in ORGAN-NOTES-CORRECTION. When appended to DATA-FACILITY-NUMBER, a foreign key to TB_OPERATOR-ID.				
M	2	D8 ORGAN-NOTE	A	80.0	N	
		Free text describing a microscopically evaluated animal or organ. Loaded by Update into TDH-ORGAN-NOTE in TD-ORGAN-DATA. If an add-with-audit, also loaded into ONC-ORGAN-NOTE in ORGAN-NOTES-CORRECTION.				
G	1	EP CORRECTION-PORION				
		Group embracing segment-common correction fields.				
	2	EQ CORRECTION-DATE	P	5.0	N	
		Date, in TDMS TOX date format, on which correction was made. Used by Update to form ONC-ORGAN-NOTE-CORR-KEY-UQ for ORGAN-NOTES-CORRECTION.				
	2	ER CORRECTION-TIME	P	5.0	N	
		Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form ONC-ORGAN-NOTE-CORR-KEY-UQ for ORGAN-NOTES-CORRECTION.				
	2	ES CORRECTION-REASON	A	2.0	N	
		Code for correction reason. Loaded by Update into ONC-COR-				

FILE...: ORGAN-NOTES-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
		RECTION-REASON in ORGAN-NOTES-CORRECTION. A foreign key to TBCORREASON.				
2	ET	CORRECTION-NOTE	A	79.0	N	
		Free text describing the corresponds-rection. Loaded by Update into ONC-CORRECTION-NOTE in ORGAN-NOTES-CORRECTION.				
2	EU	CORRECTION-TYPE	A	1.0	N	
		Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into ONC-CORRECTION-TYPE in ORGAN-NOTES-CORRECTION.				
2	FE	CORRECTION-KEY	A	42.0	N	DE
		If the date portion of TDH-ORGAN-OBSERVATION-ID-DE or TDH-ORGAN-NOTE-SEQUENCE are being changed, then this field is build by concatenating the original value of both TDH-ORGAN-OBSERVATION-ID-DE and TDH-ORGAN-NOTE-SEQUENCE, with CORRECTION-DATE and CORRECTION-TIME. It is loaded by Update into ONC-ORGAN-NOTE-CORR-KEY-LINK-UQ in ORGAN-NOTES-CORRECTION.				
1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
		Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDH-ORGAN-OBSERVATION-ID-DE, TOX-DATA-DATE, TOX-DATE-TIME, TOX-MAINFRAME-DATE, TOX-MAINFRAME-TIME, and a 4-digit sequence number for the transaction in its segment within				

FILE...: ORGAN-NOTES-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
the source set.		

TEST-NOTE-TRANSACTION

Logical layout of test note transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-TEST-NOTE file. For corrections, the information is used to add records to TEST-NOTE-CORRECTION.

FILE...: TEST-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
- - - - -						
G	1	AC HEADER				
		Group embracing fields that are common to most segments.				
2	AA	SEGMENT-ID	A	2.0	F	
		Two-letter segment identifier. Not referenced by Update. Set to TN. Used to form RECEIVE-UPDATE-KEY.				
2	AB	MICRO-OR-ECS-FLAG	A	1.0	F	
		Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE	TOX-MICRO-DATE	B	2.0	F	
		The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF	TOX-MICRO-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG	MICRO-TRANS-SET	B	2.0	F	
		The number of the transaction's source set. Not referenced by Update.				
2	AI	TOX-MAINFRAME-DATE	B	2.0	F	
		The date, in TDMS TOX date				

FILE...: TEST-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
---	----	--------------------	---	-----	---	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
---	----	-----------------------	---	-----	---	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction and is loaded by Update into TNC-CORRECTION-OPERATOR-ID in TEST-NOTE-CORRECTION.

2	AS	SOFTWARE-VERSION	B	2.0	F	
---	----	------------------	---	-----	---	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F	
---	----	---------------	---	-----	---	--

The date, in TDMS TOX date format, on which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

date to form the value of TDS-EXP-TEST-DATE-TIME-UQ that will be added to, or selected for modification from, TD-TEST-NOTE. If a correction, used to form TNC-TEST-NOTE-CORR-KEY-UQ for TEST-NOTE-COR-

FILE...: TEST-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
RECTION.		
2 AU TOX-DATA-TIME	B	2.0 F
The time of day, in TDMS TOX time format, at which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Used by Update to form the value of TDS-EXP-TEST-DATE-TIME-UQ that will be added to, or selected for modification from, TD-TEST-NOTE. If a correction, also used to form TNC-TEST-NOTE-CORR-KEY-UQ of TEST-NOTE-CORRECTION.		
2 AW OPERATOR-OR-PATHOLOGIST	B	2.0 F
The ID of the operator who create the transaction. Loaded by Update into TDS-OPER-NUM in TD-TEST-NOTE. If an add-with-audit, also loaded into TNC-TEST-OPERATOR-ID in TEST-NOTE-CORRECTION.		
2 AX PROCESS-FLAG	A	2.0 F
Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.		
G 1 FK TEST-NOTE-DATA		
Group embracing segment-specific fields in test note transactions.		
2 FL TEST-STUDY	A	7.0 N
The TDMS test number to which the note applies. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Used by Up-		

FILE...: TEST-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

date to form the value of TDS-
EXP-TEST-DATE-TIME-UQ that
will be added to, or selected
for modification from, TD-
TEST-NOTE. If a correction,
also used to form TNC-TEST-
NOTE-CORR-KEY-UQ for TEST-
NOTE-CORRECTION.

2 FM TEST-NOTE	A	160.0 N
----------------	---	---------

Free text commenting on some
aspect of a test. Loaded by
Update into TDS-NOTE-TEXT in
TD-TEST-NOTE. If an add-with-
audit, also loaded into TNC-
TEST-NOTE in TEST-NOTE-CORREC-
TION.

G 1 EP CORRECTION-PORTION

Group embracing segment-common
correction fields.

2 EQ CORRECTION-DATE	P	5.0 N
----------------------	---	-------

Date, in TDMS TOX date format,
on which correction was made.
Used by Update to form TNC-
TEST-NOTE-CORR-KEY-UQ for
TEST-NOTE-CORRECTION.

2 ER CORRECTION-TIME	P	5.0 N
----------------------	---	-------

Time of day, in TDMS TOX time
format, at which correction
was made. Used by Update to
form TNC-TEST-NOTE-CORR-KEY-UQ
for TEST-NOTE-CORRECTION.

2 ES CORRECTION-REASON	A	2.0 N
------------------------	---	-------

Code for correction reason.
Loaded by Update into TNC-COR-
RECTION-REASON in TEST-NOTE-
CORRECTION.

2 ET CORRECTION-NOTE	A	79.0 N
----------------------	---	--------

FILE...: TEST-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
Free text describing the corresponds- rection. Loaded by Update in- to TNC-CORRECTION-NOTE in TEST-NOTE-CORRECTION.		
2 EU CORRECTION-TYPE	A	1.0 N
Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into TNC-CORRECTION-TYPE in TEST-NOTE-CORRECTION.		
2 EY CORRECTION-KEY	A	27.0 N DE
If either the date or time portion of TDS-EXP-TEST-DATE- TIME-UQ is being changed then then RECEIVE or ECS will con- catenate the original TDS-EXP- -TEST-DATE-TIME-UQ, CORREC- TION-DATE, and CORRECTION- TIME to form a string that it loads into this field. Loaded by Update into TNC-TEST-NOTE- CORR-KEY-LINK-UQ in TEST-NOTE- CORRECTION.		
1 EW RECEIVE-UPDATE-KEY	A	54.0 N DE
Unique key for database trans- actions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDS-EXP-TEST-DATE-TIME-UQ for/ from TD-TEST-NOTE, TOX-DATA- DATE, TOX-DATA-TIME, TOX-MAIN- FRAME-DATE, TOX-MAINFRAME- TIME, and a 4-digit sequence number for the transaction in its segment within the source set.		

TREAT-NOTE-TRANSACTION

Logical layout of treatment group note transactions. Loaded by RECEIVE and ECS. Read by Update to add, modify, and delete records in the TD-TREAT-NOTE file. For corrections, the information is used to add records to TREAT-NOTE-CORRECTION.

FILE...: TREAT-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
G	1	AC	HEADER				
			Group embracing fields that				
			are common to most segments.				
2	AA	SEGMENT-ID		A	2.0	F	
			Two-letter segment identifier.				
			Not referenced by Update. Set				
			to TG. Used to form RECEIVE-				
			UPDATE-KEY.				
2	AB	MICRO-OR-ECS-FLAG		A	1.0	F	
			Identifier for source of				
			transaction. A - LDAS,				
			E - Mainframe ECS. Used by				
			Update to generate Update				
			statistics report. Used to				
			form RECEIVE-UPDATE-KEY.				
2	AE	TOX-MICRO-DATE		B	2.0	F	
			The date, in TDMS TOX format,				
			on which the source set for				
			the transaction was created.				
			Not referenced by Update.				
2	AF	TOX-MICRO-TIME		B	2.0	F	
			The time of day, in TDMS TOX				
			time format, on which the				
			source set for the transac-				
			was created. Not refer-				
			enced by Update.				
2	AG	MICRO-TRANS-SET		B	2.0	F	
			The number of the transac-				
			tion's source set. Not refer-				
			enced by Update.				
2	AI	TOX-MAINFRAME-DATE		B	2.0	F	
			The date, in TDMS TOX date				

FILE...: TREAT-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
---	----	--------------------	---	-----	---	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
---	----	-----------------------	---	-----	---	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction and is loaded by Update into TGC-CORRECTION-OPERATOR-ID of TREAT-NOTE-CORRECTION.

2	AS	SOFTWARE-VERSION	B	2.0	F	
---	----	------------------	---	-----	---	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F	
---	----	---------------	---	-----	---	--

The date, in TDMS TOX date format, on which the transaction was created. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Used by Update to form the value of TDT-EXP-TEST-TRT-DATE-TIME-UQ that will be added to, or selected for modification from, TD-TREATMENT-NOTE. If a correction, also used to form TGC-

FILE...: TREAT-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
TREAT-NOTE-CORR-KEY-UQ for TREAT-NOTE-CORRECTION.		
2 AU TOX-DATA-TIME The time of day, in TDMS TOX time format, at which the transaction was created. Used by RECEIVE and ECS to form RE- CEIVE-UPDATE-KEY. Used by Up- date to form the value of TDG- EXP-TEST-TRT-DATE-TIME-UQ that will be added to, or selected for modification from, TD- TREATMENT-NOTE. If a correc- , also used to form TGC- TREAT-NOTE-CORR-KEY-UQ for TREAT-NOTE-CORRECTION.	B	2.0 F
2 AW OPERATOR-OR-PATHOLOGIST The ID of the operator who created the transaction. Loaded by Update into TDT- OPER-NUM of TD-TREATMENT-NOTE. If an add-with-audit, also loaded into TGC-TREAT-OPERA- TOR-ID in TREAT-NOTE-CORREC- TION.	B	2.0 F
2 AX PROCESS-FLAG Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or del- ete, respectively. The second digit is either blank, indica- ting an original transaction, or 'A', indicating an audit transaction.	A	2.0 F
G 1 FN TREAT-NOTE-DATA Group embracing segment-speci- fic fields of treatment group note transactions.		
2 FO TREAT-STUDY The study to which the note applies. Used by RECEIVE and	A	7.0 N

FILE...: TREAT-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

ECS to form RECEIVE-UPDATE-KEY. Used by Update to form the value of TDT-EXP-TEST-TRT-DATE-TIME-UQ for TD-TREATMENT-NOTE. If a correction, also used to form TGC-TREAT-NOTE-CORR-KEY-UQ of TREAT-NOTE-CORRECTION.

2 FP TREAT-GROUP	A	3.0 N
------------------	---	-------

The ID of the treatment group to which the note applies. Used by Update to form the value of TDT-EXP-TEST-TRT-DATE-TIME-UQ of TD-TREATMENT-NOTE. If a correction, also used to form TGC-TREAT-NOTE-CORR-KEY-UQ for TREAT-NOTE-CORRECTION.

2 FQ TREAT-NOTE	A	160.0 N
-----------------	---	---------

Free text commenting on some aspect of a treatment group. Loaded by Update into TDT-NOTE-TEXT of TD-TREATMENT-NOTE. If an add-with-audit, also loaded into TGC-TREAT-NOTE in TREAT-NOTE-CORRECTION.

G 1 EP CORRECTION-PORION		
--------------------------	--	--

Group embracing segment-common correction fields.

2 EQ CORRECTION-DATE	P	5.0 N
----------------------	---	-------

Date, in TDMS TOX date format, on which correction was made. Used by Update to form TGC-TREAT-NOTE-CORR-KEY-UQ for TREAT-NOTE-CORRECTION.

2 ER CORRECTION-TIME	P	5.0 N
----------------------	---	-------

Time of day, in TDMS TOX time

FILE...: TREAT-NOTE-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
format, at which correction was made. Used by Update to form TGC-TREAT-NOTE-CORR-KEY-UQ of TREAT-NOTE-CORRECTION.		
2 ES CORRECTION-REASON	A	2.0 N
Code for correction reason. Loaded by Update into TGC-CORRECTION-REASON in TREAT-NOTE-CORRECTION.		
2 ET CORRECTION-NOTE	A	79.0 N
Free text describing the corresponds-rection. Loaded by Update into TGC-CORRECTION-NOTE of TREAT-NOTE-CORRECTION.		
2 EU CORRECTION-TYPE	A	1.0 N
Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into TGC-CORRECTION-TYPE in TREAT-NOTE-CORRECTION.		
2 EZ CORRECTION-KEY	A	30.0 N DE
If either the date, time, or treatment group number portion of TDT-EXP-TEST-TRT-DATE-TIME-UQ is changing then RECEIVE or ECS concatenates the original value of TDT-EXP-TEST-TRT-DATE-TIME-UQ, CORRECTION-DATE, and CORRECTION-TIME into a string that it loads into this field. Loaded by Update into TGC-TREAT-NOTE-CORR-KEY-LINK-UQ of TREAT-NOTE-CORRECTION.		
1 EW RECEIVE-UPDATE-KEY	A	54.0 N DE
Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDT-EXP-TEST-TRT-DATE-TIME-UQ for/from TD-TREATMENT-NOTE left-justified into a 27-byte string, TOX-DATA-DATE, TOX-DATA-TIME, TOX-MAINFRAME-DATE,		

FILE...: TREAT-NOTE-TRANSACTION
TYPE...: USER VIEW
FILE-NR: 86
PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

TOX-MAINFRAME-TIME, and a 4-
digit sequence number for the
transaction in its segment
within the source set.

UNCERTAIN-TRANSACTION

Logical layout for uncertain cause-of-death transactions. Loaded by RECEIVE and ECS. Read by Update to modify records in the TD-ANIMAL file. For correction, information is used to add records to UNCERTAIN-CORRECTION.

FILE...: UNCERTAIN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	AC	HEADER				
			Group embracing fields that are common to most segments.				
2	AA		SEGMENT-ID	A	2.0	F	
			Two-letter segment identifier. Not referenced by Update. Set to UC. Used to form RECEIVE-UPDATE-KEY.				
2	AB		MICRO-OR-ECS-FLAG	A	1.0	F	
			Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report. Used to form RECEIVE-UPDATE-KEY.				
2	AE		TOX-MICRO-DATE	B	2.0	F	
			The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF		TOX-MICRO-TIME	B	2.0	F	
			The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG		MICRO-TRANS-SET	B	2.0	F	
			The number of the transaction's source set. Not referenced by Update.				
2	AI		TOX-MAINFRAME-DATE	B	2.0	F	
			The date, in TDMS TOX date				

FILE...: UNCERTAIN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F		
---	----	--------------------	---	-----	---	--	--

The time of day, in TDMS TOX time format, on which the source set was received. Used by RECEIVE and ECS to form RECEIVE-UPDATE-KEY. Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F		
---	----	-----------------------	---	-----	---	--	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction and is loaded by Update into UCC-CORRECTION-OPERATOR-ID.

2	AS	SOFTWARE-VERSION	B	2.0	F		
---	----	------------------	---	-----	---	--	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F		
---	----	---------------	---	-----	---	--	--

The date, in TDMS TOX date format, on which the transaction was created. Loaded by Update into TDA-UNCERTAIN-OPERATOR-ID in TD-ANIMAL. If an add-with-audit, also loaded into UCC-UNCERTAIN-DATE in UNCERTAIN-CORRECTION.

2	AU	TOX-DATA-TIME	B	2.0	F		
---	----	---------------	---	-----	---	--	--

The time of day, in TDMS TOX time format, at which the transaction was created.

FILE...: UNCERTAIN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
Loaded by Update into TDA-UNCERTAIN-TIME in TD-ANIMAL. If an add-with-audit, also loaded into UCC-UNCERTAIN-TIME in UNCERTAIN-CORRECTION.		
2 AW OPERATOR-OR-PATHOLOGIST	B	2.0 F
The ID of the operator who created the transaction. Loaded by Update into TDA-UNCERTAIN-OPERATOR-ID in TD-ANIMAL. If an add-with-audit, also loaded into UCC-UNCERTAIN-OPERATOR-ID in UNCERTAIN-CORRECTION.		
2 AX PROCESS-FLAG	A	2.0 F
Two-letter identifier for the type of transaction. The first digit can be A, C, or D, which signifies add, change, or delete, respectively. The second digit is either blank, indicating an original transaction, or 'A', indicating an audit transaction.		
G 1 FF UNCERTAIN-DATA		
Group embracing segment-specific fields in "uncertain cause-of-death" transactions.		
2 FG UNCERTAIN-STUDY	A	7.0 N
TDMS test number to which the observation applies. Used by Update to select the TD-ANIMAL record to which the transaction applies.		
2 FH UNCERTAIN-ANIMAL-NUM	P	5.0 N
Number of animal being classified as having uncertain cause-of-death. Not referenced by Update.		
2 FI UNCERTAIN-PATHOLOGIST	P	5.0 N
Operator ID of the pathologist assigned to the test. Loaded		

FILE...: UNCERTAIN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

by Update into TDA-UNCERTAIN-PATHOLOGIST in TD-ANIMAL. If an add-with-audit, also loaded into UCC-UNCERTAIN-PATH-ID in UNCERTAIN-CORRECTION.

G 1 EP CORRECTION-PORION

Group embracing segment-common correction fields.

2 EQ CORRECTION-DATE P 5.0 N

Date, in TDMS TOX date format, on which correction was made. Used by Update to form UCC-UNCERTAIN-CORR-KEY-UQ for UNCERTAIN-CORRECTION.

2 ER CORRECTION-TIME P 5.0 N

Time of day, in TDMS TOX time format, at which correction was made. Used by Update to form UCC-UNCERTAIN-CORR-KEY-UQ for UNCERTAIN-CORRECTION.

2 ES CORRECTION-REASON A 2.0 N

Code for correction reason. Loaded by Update into UCC-CORRECTION-REASON in UNCERTAIN-CORRECTION.

2 ET CORRECTION-NOTE A 79.0 N

Free text describing the corresponds-rection. Loaded by Update into UNC-CORRECTION-NOTE in UNCERTAIN-CORRECTION.

2 EU CORRECTION-TYPE A 1.0 N

Single letter indicating type of correction. A-Add, C-Change D-Delete. Loaded by Update into UNC-CORRECTION-TYPE in UNCERTAIN-CORRECTION.

FILE...: UNCERTAIN-TRANSACTION

TYPE...: USER VIEW

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
2 EX CORRECTION-KEY	A	26.0 N DE
Not assigned a value by RECEIVE or ECS. Loaded by Update into UCC-UNCERTAIN-CORR-KEY-LINK-UQ in UNCERTAIN-CORRECTION.		
1 EW RECEIVE-UPDATE-KEY	A	54.0 N DE
Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, TDA-EXP-TEST-ANIMAL-UQ, 11 spaces, TOX-DATA-DATE, TOX-DATA-TIME, TOX-MAINFRAME-DATE, TOX-MAINFRAME-TIME, and a 4-digit sequence number for the transaction in its segment within the source set.		

STATUS_FILE

View of RECEIVE transaction set status information. Provides a history of set processing.

FILE...: STATUS_FILE

TYPE...: USER VIEW

FILE-NR: 38

PRIMARY SEQUENCE FIELD: STATUS_EXP_AND_TEST_NUMBER

T	L	DB	NAME	F	LENG	S	DE
1	CD		STATUS_AGENCY_FACILT_DISK_SUBSYS	A	11	0	
			Composite field whose elements are: the 2-digit agency number, which is a foreign key to TBAGENCY, the facility code, which is a foreign key to TBFACLT and is right-justified and padded into a 5-digit string, the diskette ID as a 2-byte binary string (not assigned to LDAS sets); and a 2-digit subsystem ID, EI for in-life, and HI for pathology.				
1	IJ		STATUS_EXP_AND_TEST_NUMBER	N	7	0	DE
			The experiment and test number to which the set applies. When appended to the agency number, a foreign key to PAS_MASTER.				
1	KK		STATUS_TRANSACTION_SET_NUM	N	5	0	
			The sequentially assigned set number. In-life and pathology sets on a test both start from 1, and are unique for that test.				
1	LL		STATUS_HEADER_TOX_DATE	A	8	0	
			The date on which the set was created, in the form YYYYMMDD.				
1	MM		STATUS_HEADER_TOX_TIME	A	6	0	
			The time at which the set was created, in the form HHMMSS.				
1	OO		STATUS_RECEIVE_DATE	A	8	0	N
			The date on which the set was successfully processed by RE-				

FILE...: STATUS_FILE

TYPE...: USER VIEW

FILE-NR: 38

PRIMARY SEQUENCE FIELD: STATUS_EXP_AND_TEST_NUMBER

T	L	DB	NAME	F	LENG	S	DE
			CEIVE, in the form YYYYMMDD.				
1	PP		STATUS_RECEIVE_TIME	A	6.0	N	
			The time at which the set was successfully processed by RE-CEIVE, in the form HHMMSS.				
1	QQ		STATUS_RECEIVE_INDICATOR	A	1.0	F	
			Flag used to track interim status of MODCOMP RECEIVE processing. Assigned a value by current RECEIVE for compatibility with STATTRAN reports, but has no meaning since there is no longer interim status.				
1	RR		STATUS_VERIFY_INDICATOR	A	1.0	F	
			Flag used to track interim status of MODCOMP RECEIVE processing. Assigned a value by current RECEIVE for compatibility with STATTRAN reports, but has no meaning since there is no longer interim status.				
1	SS		STATUS_CONVERT_INDICATOR	A	1.0	F	
			Flag used to track interim status of MODCOMP RECEIVE processing. Assigned a value by current RECEIVE for compatibility with STATTRAN reports, but has no meaning since there is no longer interim status.				
1	TT		STATUS_NOTIFY_INDICATOR	A	1.0	N	
			Flag used to track interim status of MODCOMP RECEIVE processing. Assigned a value by current RECEIVE for compatibility with STATTRAN reports, but has no meaning since there is no longer interim status.				

FILE...: STATUS_FILE

TYPE...: USER VIEW

FILE-NR: 38

PRIMARY SEQUENCE FIELD: STATUS_EXP_AND_TEST_NUMBER

T	L	DB NAME	F	LENG	S	DE
1	UU	STATUS_TRANSACTION_COUNT	N	5.0	N	

The number of database transactions RECEIVE created. This total is different than the number of transaction sets themselves, since RECEIVE consolidates and splits out transactions in certain circumstances.

1	KE	STATUS_KEY	A	34.0	SP	
SOURCE FIELD(S) --- -START- --END-						
		STATUS_AGENCY_FACILT_DISK_SUBSYS	1		11	
		STATUS_EXP_AND_TEST_NUMBER	1		7	
		STATUS_TRANSACTION_SET_NUM	1		2	
		STATUS_HEADER_TOX_DATE	1		8	
		STATUS_HEADER_TOX_TIME	1		6	

Used to select data for STAT-TRAN and RECEIVE. The primary key for the table.

1	VV	STATUS_EXP_TEST_TRANS_SET_NUM	B	9.0	SP	
SOURCE FIELD(S) --- -START- --END-						
		STATUS_EXP_AND_TEST_NUMBER	1		7	
		STATUS_TRANSACTION_SET_NUM	1		2	

Used to select data based on experiment-test and transaction set number

4.4 Audit History Views

ANIMAL-MICRO-NOTES-CORRECTION

Logical view of pathology animal note audit changes. The master file is TD-MICRO-ANIMAL-NOTES. Data from this view is not retrieved by any current TDMS component.

FILE...: ANIMAL-MICRO-NOTES-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: MAC-MICRO-ANI-NOTE-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AY		MAC-MICRO-ANI-NOTE-CORR-KEY-UQ	A	37.0	N	UQ
			Composite field formed at Up-date by the concatenation of: the values of TDM-EXP-TST-ANUM-PATHT-DATE-DE and TDM-MICRO-ANIMAL-NOTE-SEQUENCE after the record change, and CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 22 characters constitute a foreign key to TD-MICRO-ANIMAL-NOTE.				
1	ED		MAC-OPERATOR-ID	P	5.0	N	
			Historic value of TDM-OPERATOR-IDENTIFICATION. If MAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. Appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	EE		MAC-PATHOLOGIST-ID	P	5.0	N	
			Historic value of TDM-PATHOLOGIST-IDENTIFICATION. If MAC-CORRECTION-TYPE is 'A', it is the original value; otherwise, it is the value prior to modification. Appended to				

FILE...: ANIMAL-MICRO-NOTES-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: MAC-MICRO-ANI-NOTE-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
		the appropriate facility code, a foreign key to TB_OPERATOR_ ID.				
1	DA	MAC-COD-FLAG	A	1.0	F	Historic value of TDM-COD- FLAG. If MAC-CORRECTION-TYPE is 'A', then it is the origi- nal value; otherwise, it is the value prior to modifica- tion.
M 1	DI	MAC-MICRO-ANIMAL-NOTE	A	80.0	N	Historic value of TDM-MICRO- ANIMAL-NOTE. If MAC-CORREC- TION-TYPE is 'A', then this is the original value; otherwise, it is the value prior to modi- fication.
1	CA	MAC-CORRECTION-OPERATOR-ID	P	5.0	F	The ID of the operator who made the correction. Append- ed to the appropriate facility code, a foreign key to TB_OPERATOR_ID.
1	CB	MAC-CORRECTION-REASON	A	2.0	F	The reason for the correction. A foreign key to TBCORREASON.
1	CC	MAC-CORRECTION-NOTE	A	79.0	N	Free text entered by correc- operator describing the correction (optional).
1	CD	MAC-CORRECTION-TYPE	A	1.0	F	Type of correction, A for add, C for change, D for delete.
1	AZ	MAC-MICRO-ANI-NTE-COR-KEY-LNK-UQ	A	37.0	N	UQ If MAC-CORRECTION-TYPE is 'C' and either TDM-EXP-TST-ANUM- PATHT-DATE-DE or TDM-MICRO- ANIMAL-NOTE-SEQUENCE was changed, then the value of this is field built by concat- enating the prior values of

FILE...: ANIMAL-MICRO-NOTES-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: MAC-MICRO-ANI-NOTE-CORR-KEY-UQ

T L DB NAME	F L E N G S D E
TDM-EXP-TST-ANUM-PATHT-DATE-DE	
and TDM-MICRO-ANIMAL-NOTE-SE-	
QUENCE, and CORRECTION-DATE	
and CORRECTION-TIME. A for-	
eign key to ANIMAL-MICRO-	
NOTES-CORRECTION itself.	

ANIMAL-NOTE-CORRECTION

Logical view of animal note corrections. Master file is TD-ANIMAL-DATA. View is used by EIS Report 22.

FILE...: ANIMAL-NOTE-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ANC-ANIMAL-NOTE-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AA		ANC-ANIMAL-NOTE-CORR-KEY-UQ	A	36.0	N	UQ
			A composite field formed at Update by concatenating: the values of TDD-EXP-TEST-ANIMAL-DATA-DE and TDD-NOTE-TIME after the record change, and CORRECTION-DATE and CORRECTION-TIME from the transaction record. The first 21 characters are a foreign key to TD-ANIMAL-DATA.				
1	ED		ANC-NOTE-OPERATOR-ID	P	5.0	N	
			Historic value of TDD-NOTE-OPERATOR-ID. If ANC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	DA		ANC-NOTE-MANUAL-FLAG	A	1.0	F	
			Historic value of TDD-MANUAL-ENTRY-FLAG-NOTE. If ANC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
M 1	DI		ANC-ANIMAL-NOTE	A	80.0	N	
			Historic value of TDD-ANIMAL-NOTE. If ANC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modifica-				

FILE...: ANIMAL-NOTE-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ANC-ANIMAL-NOTE-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
		tion.				
1	CA	ANC-CORRECTION-OPERATOR-ID	P	5.0	F	
		The operator who made corresponds-				
		rection. When appended to				
		the appropriate facility code,				
		a foreign key to TB_OPERATOR_-				
		ID.				
1	CB	ANC-CORRECTION-REASON	A	2.0	F	
		The correction reason. A for-				
		eign key to TBCORREASON.				
1	CC	ANC-CORRECTION-NOTE	A	79.0	N	
		Free text entered by correc-				
		operator describing the				
		correction (optional).				
1	CD	ANC-CORRECTION-TYPE	A	1.0	F	
		Type of correction, A for add,				
		C for change, D for delete.				
1	AB	ANC-ANIMAL-NOTE-CORR-KEY-LINK-UQ	A	36.0	N	UQ
		If ANC-CORRECTION-TYPE is 'C'				
		and either TDD-EXP-TEST-ANI-				
		MAL-DATE-DE or TDD-NOTE-TIME				
		was changed, then the value is				
		built by concatenating the				
		the prior values those fields				
		with CORRECTION-DATE and COR-				
		RECTION-TIME from the trans-				
		action. A foreign key to				
		ANIMAL-NOTE-CORRECTION itself.				

ANIMAL-OBSERVATION-CORRECTION

Logical view of animal observation corrections. Master file is TD-ANIMAL-DATA. Used by EIS Report 22.

FILE...: ANIMAL-OBSERVATION-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AOC-ANIMAL-OBS-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AC	AOC-ANIMAL-OBS-CORR-KEY-UQ	A	36.0	N	UQ	
<p>A composite field formed at Update by concatenating the values of TDD-EXP-TEST-ANIMAL-DATE-DE and TDD-OBSERVATION-TIME after the record change, and CORRECION-DATE and CORREC-TION-TIME from the transac-tion. The first 21 characters are a foreign key to TD-ANI-MAL-DATA.</p>							
1	ED	AOC-OBSERVATION-OPERATOR-ID	P	5.0	N		
<p>Historic value of TDD-OBSERVA-TION-OPERATOR-ID. If AOC-COR-RECTION-TYPE is 'A', then it is the original value; other-wise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_-OPERATOR_ID.</p>							
1	DA	AOC-MANUAL-ENTRY-FLAG	A	1.0	F		
<p>Historic value of TDD-MANUAL-ENTRY-FLAG-OBS. If AOC-COR-RECTION-TYPE is 'A', then it is the original value; other-wise, it is the value prior to modification.</p>							
P 1	FA	AOC-OBSERVATION-DATA					
<p>Historic record of the TDD-OBSERVATION-DATA periodic group.</p>							
2	FB	AOC-OBSERVATION-DE	P	5.0	N		
<p>Historic values of TDD-OBSER-</p>							

FILE...: ANIMAL-OBSERVATION-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AOC-ANIMAL-OBS-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

VATION. If AOC-CORRECTION-
TYPE is 'A', then they are the
original values; otherwise,
they are the values prior to
modification. A foreign key
to TBCLOBSC.

2	FC	AOC-SITE	P	5.0	N	
---	----	----------	---	-----	---	--

Historic values of TDD-SITE.
If AOC-CORRECTION-TYPE is 'A',
then they are the original
values; otherwise, they are
the values prior to modifica-
tion. A foreign key to
TBCLOBSC.

2	FD	AOC-EIQUAL1	P	5.0	N	
---	----	-------------	---	-----	---	--

Historic values of TDD-EIQUAL1
If AOC-CORRECTION-TYPE is 'A',
then they are the original
values; otherwise, they are
the values prior to modifica-
tion. A foreign key to
TBCLOBSC.

2	FE	AOC-EIQUAL2	P	5.0	N	
---	----	-------------	---	-----	---	--

Historic values of TDD-EIQUAL2
If AOC-CORRECTION-TYPE is 'A',
then they are the original
values; otherwise, they are
the values prior to modifica-
tion. A foreign key to
TBCLOBSC.

1	CA	AOC-CORRECTION-OPERATOR-ID	P	5.0	F	
---	----	----------------------------	---	-----	---	--

The operator who made the corresponds-
rection. When appended to the
appropriate facility code, a
foreign key to TB_OPERATOR_ID.

1	CB	AOC-CORRECTION-REASON	A	2.0	F	
---	----	-----------------------	---	-----	---	--

The correction reason. A for-

FILE...: ANIMAL-OBSERVATION-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AOC-ANIMAL-OBS-CORR-KEY-UQ

T L	DB NAME	F	LENG	S	DE
	eign key to TBCORREASON.				
1	CC AOC-CORRECTION-NOTE	A	79.0	N	
	Free text entered by correc- operator describing the correction (optional).				
1	CD AOC-CORRECTION-TYPE	A	1.0	F	
	Type of correction, A for add, C for change, D for delete.				
1	AD AOC-ANIMAL-OBS-CORR-KEY-LINK-UQ	A	36.0	N	UQ
	If ANC-CORRECTION-TYPE is 'C' and either TDD-EXP-TEST-ANI- MAL-DATE-DE or TDD-OBSERVA- TION-TIME was changed, then the value is built by concate- nating the prior values of TDD-EXP-TEST-ANIMAL-DATE-DE and TDD-OBSERVATION-TIME, and CORRECTION-DATE and CORREC- TION-TIME from the transac- tion. A foreign key to ANI- MAL-OBSERVATION-CORRECTION it- self.				

ANIMAL-REMOVAL-CORRECTION

Logical view of animal removal corrections. Master file is TD-ANIMAL. Used by EIS Report 22.

FILE...: ANIMAL-REMOVAL-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ARC-ANIMAL-REM-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AE	ARC-ANIMAL-REM-CORR-KEY-UQ	A	26.0	N	UQ	
Composite field made at Update by concatenating the value of TDA-EXP-TEST-ANIMAL-UQ with CORRECTION-DATA and CORRECTION-TIME in the transaction. The first 16 characters form a foreign key to TD-ANIMAL.							
1	DE	ARC-CARCASS-IDENTIFICATION-DE	A	14.0	N		
Historic value of TDA-CARCASS-IDENTIFICATION-DE. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.							
1	ED	ARC-REASON-FOR-REMOVAL-DE	P	5.0	N		
Historic value of TDA-REASON-FOR-REMOVAL-DE. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it the value prior to modification. Foreign key to TBCLOBSC.							
1	EE	ARC-REMOVAL-DATE	P	5.0	N		
Historic value of TDA-REMOVAL-DATE. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.							
1	EF	ARC-REMOVAL-TIME	P	5.0	N		
Historic value of TDA-REMOVAL-TIME. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.							

FILE...: ANIMAL-REMOVAL-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ARC-ANIMAL-REM-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

nal value; otherwise, it is the value prior to modification.

1	EG	ARC-REMOVAL-OPERATOR-ID	P	5.0	N		
---	----	-------------------------	---	-----	---	--	--

Historic value of TDA-REMOVAL-OPERATOR-ID. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When concatenated to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	EH	ARC-DAYS-ON-EXPERIMENT	P	5.0	N		
---	----	------------------------	---	-----	---	--	--

Historic value of TDA-DAYS-ON-EXPERIMENT. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.

1	DA	ARC-MANUAL-ENTRY-FLAG-WGT	A	1.0	F		
---	----	---------------------------	---	-----	---	--	--

Historic value of TDA-MANUAL-ENTRY-FLAG-WGT. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.

1	EI	ARC-REM-WEIGHT-STATUS	P	5.0	N		
---	----	-----------------------	---	-----	---	--	--

Historic value of TDA-REM-WEIGHT-STATUS. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. Foreign key to TBWGTSTATUS.

1	DB	ARC-MANUAL-FLAG	A	1.0	F		
---	----	-----------------	---	-----	---	--	--

Historic value of TDA-MANUAL-

FILE...: ANIMAL-REMOVAL-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ARC-ANIMAL-REM-CORR-KEY-UQ

T L DB NAME	F	LENG S DE
ENTRY-FLAG-REM. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.		
1 EL ARC-REMOVAL-WEIGHT	P	11.0 N
Historic value of TDA-ANIMAL-REM-WEIGHT. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.		
P 1 FA ARC-OBSERVATION-DATA		
Historic record of the TDA-OBSERVATION-DATA periodic group.		
2 FB ARC-OBSERVATION-DE	P	5.0 N
Historic values of TDA-OBSERVATION. If ARC-CORRECTION-TYPE is 'A', then they are the original values; otherwise, they are the values prior to modification. Foreign key to TBCLOBSC.		
2 FC ARC-SITE	P	5.0 N
Historic values of TDA-SITE. If ARC-CORRECTION-TYPE is 'A', then they are the original values; otherwise, they are the values prior to modification. Foreign key to TBCLOBSC		
2 FD ARC-EIQUAL1	P	5.0 N
Historic values of TDA-EIQUAL1. If ARC-CORRECTION-TYPE is 'A', then they are the original values; otherwise, they are the values prior to modification. Foreign key to TBCLOBSC		
2 FE ARC-EIQUAL2	P	5.0 N
Historic values of TDA-EIQUAL2. If ARC-CORRECTION-TYPE is 'A', then they are the original		

FILE...: ANIMAL-REMOVAL-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ARC-ANIMAL-REM-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

		values; otherwise, they are the values prior to modification. Foreign key to TBCLOBSC				
1	DC	ARC-MANUAL-ENTRY-FLAG-OBS	A	1.0	F	
		Historic value of TDA-MANUAL-ENTRY-FLAG-OBS. If ARC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	CA	ARC-CORRECTION-OPERATOR-ID	P	5.0	F	
		Operator who made correction. When appended to appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	CB	ARC-CORRECTION-REASON	A	2.0	F	
		The correction reason. A foreign key to TBCORREASON.				
1	CC	ARC-CORRECTION-NOTE	A	79.0	N	
		Free text entered by correction operator describing the correction (optional).				
1	CD	ARC-CORRECTION-TYPE	A	1.0	F	
		Type of correction, A for add, C for change, D for delete.				
1	AF	ARC-ANIMAL-REM-CORR-KEY-LINK-UQ	A	26.0	N UQ	
		No possible value.				

ANIMAL-TRANSFER-CORRECTION

Logical view of animal transfer corrections. Master table is TD-ANIMAL. Used by EIS Report 22.

FILE...: ANIMAL-TRANSFER-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ATC-ANIMAL-TRAN-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AG	ATC-ANIMAL-TRAN-CORR-KEY-UQ	A	26.0	N	UQ	
		Composite field formed at Update by concatenating TDA-EXP-TEST-ANIMAL-UQ with CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 16 characters are a foreign key to TD-ANIMAL.					
1	ED	ATC-TRANSFER-DATE	P	5.0	N		
		Historic value of TDA-TRANSFER-DATE. If ATC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.					
1	EE	ATC-TRANSFER-TIME	P	5.0	N		
		Historic value of TDA-TRANSFER-TIME. If ATC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.					
1	EF	ATC-TRANSFER-OPERATOR-ID	P	5.0	N		
		Historic value of TDA-TRANSFER-OPERATOR-ID. If ATC-CORRECTION-TYPE is 'A', then it is the original value; otherwise it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.					
1	EG	ATC-HOSPITAL-CAGE	P	5.0	N		

FILE...: ANIMAL-TRANSFER-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ATC-ANIMAL-TRAN-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

Historic value of TDA-HOSPITAL-CAGE-DE. If ATC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.

1	CA	ATC-CORRECTION-OPERATOR-ID	P	5.0	F	
Operator who made correction. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.						
1	CB	ATC-CORRECTION-REASON	A	2.0	F	
Correction reason. Foreign key to TBCORREASON.						
1	CC	ATC-CORRECTION-NOTE	A	79.0	N	
Free text entered by correction operator describing the correction (optional).						
1	CD	ATC-CORRECTION-TYPE	A	1.0	F	
Type of correction, A for add, C for change, D for delete.						
1	AH	ATC-ANIMAL-TRAN-CORR-KEY-LINK-UQ	A	26.0	N	UQ
No value possible.						

ANIMAL-WEIGHT-CORRECTION

Logical layout of animal weight corrections. Master file is TD-ANIMAL-DATA. Used by EIS Report 22.

FILE...: ANIMAL-WEIGHT-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AWC-ANIMAL-WGT-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AI		AWC-ANIMAL-WGT-CORR-KEY-UQ	A	36.0	N	UQ
			Composite field formed at Update by concatenating the values of DD-EXP-TEST-ANIMAL-DATE-DE and TDD-WEIGHT-TIME after the record change, and CORRECTION-DATE and CORRECTION-TIME from the transaction. First 21 characters form a foreign key to TD-ANIMAL-DATA.				
1	EL		AWC-ANIMAL-WEIGHT	P	11.0	N	
			Historic value of TDD-ANIMAL-WEIGHT. If AWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	ED		AWC-WEIGHT-OPERATOR-ID	P	5.0	N	
			Historic value of TDD-WEIGHT-OPERATOR-ID. If AWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR-ID.				
1	DA		AWC-MANUAL-ENTRY-FLAG-WGT	A	1.0	F	
			Historic value of TDD-MANUAL-ENTRY-FLAG-WGT. If AWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to				

FILE...: ANIMAL-WEIGHT-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AWC-ANIMAL-WGT-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

modification.

1	EE	AWC-WEIGHT-STATUS	P	5.0	N	
Historic value of TDD-WEIGHT-STATUS. If AWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. Foreign key to TBWGTSTATUS.						
1	CA	AWC-CORRECTION-OPERATOR-ID	P	5.0	F	
ID of operator who made corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.						
1	CB	AWC-CORRECTION-REASON	A	2.0	F	
The correction reason. A foreign key to TBCORREASON.						
1	CC	AWC-CORRECTION-NOTE	A	79.0	N	
Free text entered by correction operator describing the correction (optional).						
1	CD	AWC-CORRECTION-TYPE	A	1.0	F	
Type of correction, A for add, C for change, D for delete.						
1	AJ	AWC-ANIMAL-WGT-CORR-KEY-LINK-UQ	A	36.0	N	UQ
If AWC-CORRECTION-TYPE is 'C' and either TDD-EXP-TEST-ANIMAL-DATE-DE or TDD-WEIGHT-TIME was changed, then the value is built by concatenating the prior values of TDD-EXP-TEST-ANIMAL-DATE-DE and TDD-WEIGHT-TIME from TD-ANIMAL-DATA, and CORRECTION-DATE and CORRECTION-TIME from ANIMAL-WEIGHT-						

FILE...: ANIMAL-WEIGHT-CORRECTION
TYPE...: USER VIEW
FILE-NR: 87
PRIMARY SEQUENCE FIELD: AWC-ANIMAL-WGT-CORR-KEY-UQ

T L DB NAME	F	LENG S DE
TRANSACTION.	A foreign key	
to ANIMAL-WEIGHT-CORRECTION		
itself.		

ANY-ANIMAL-CORRECTION

Logical view of animal allocation corrections. Master table is TD-ANIMAL. Used by EIS Report 22.

FILE...: ANY-ANIMAL-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AAC-ANY-ANIMAL-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AK		AAC-ANY-ANIMAL-CORR-KEY-UQ	A	26.0	N	UQ
Composite field formed at Up-date by concatenating TDA-EXP-TEST-ANIMAL-UQ, and CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 16 characters form a foreign key to TD-ANIMAL.							
1	DA		AAC-ANIMAL-IDENTIFICATION	A	1.0	F	
Historic value of TDA-ANIMAL-IDENTIFICATION. If AAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise it is the value prior to modification.							
1	ED		AAC-ANIMAL-CLASS	P	5.0	N	
Historic value of TDA-ANIMAL-CLASS. If AAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. Foreign key to TBANCLQF.							
1	DB		AAC-SEX	A	1.0	F	
Historic value of TDA-SEX. If AAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.							
1	EE		AAC-SPECIES	P	5.0	N	
Historic value of TDA-SPECIES. If AAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value							

FILE...: ANY-ANIMAL-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AAC-ANY-ANIMAL-CORR-KEY-UQ

T L DB NAME	F	LENG S DE
-------------	---	-----------

1 EF AAC-STRAIN	P	5.0 N
<p>prior to modification.</p> <p>Foreign key to TBSTRAIN.</p> <p>Historic value of TDA-STRAIN.</p> <p>If AAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.</p> <p>Foreign key to TBSTRAIN.</p>		
1 EG AAC-SUBSTRAIN	P	5.0 N
<p>Historic value of TDA-SUB-STRAIN. If AAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. Foreign key to TBSUBSTR.</p>		
1 EH AAC-ANIMAL-DATE	P	5.0 N
<p>Historic value of TDA-ANIMAL-DATE. If AAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.</p>		
1 EI AAC-ANIMAL-TIME	P	5.0 N
<p>Historic value of TDA-ANIMAL-TIME. If AAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.</p>		
1 EJ AAC-ANIMAL-OPERATOR-ID	P	5.0 N
<p>Historic value of TDA-ANIMAL-OPERATOR-ID. If AAC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.</p>		

FILE...: ANY-ANIMAL-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AAC-ANY-ANIMAL-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
		fication. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	CA	AAC-CORRECTION-OPERATOR-ID	P	5.0	F	
		ID of operator who made corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	CB	AAC-CORRECTION-REASON	A	2.0	F	
		The correction reason. A foreign key to TBCORREASON.				
1	CC	AAC-CORRECTION-NOTE	A	79.0	N	
		Free text entered by correction operator describing the correction (optional).				
1	CD	AAC-CORRECTION-TYPE	A	1.0	F	
		Type of correction, A for add, C for change, D for delete.				
1	AL	AAC-ANY-ANIMAL-CORR-KEY-LINK-UQ	A	26.0	N	UQ
		No possible value.				

ANY-MICRO-OBS-CORRECTION

Logical view of pathology microscopic observation audit changes.
Master file is TD-ORGAN-DATA. Not referenced by any TDMS
retrieval component.

FILE...: ANY-MICRO-OBS-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: MOC-MICRO-OBS-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	BA		MOC-MICRO-OBS-CORR-KEY-UQ	A	42.0	N	UQ
			Composite field formed at Up-date by concatenating the values of TDH-ORGAN-OBSERVATION-ID-DE and TDH-MICRO-OBSERVATION-SEQUENCE after the transaction has been applied, and CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 27 characters form a foreign key to TD-ORGAN-DATA.				
1	ED		MOC-MICRO-OBS-OPERATOR-ID	P	5.0	N	
			Historic value of TDH-MICRO-OBS-OPERATOR-ID. If MOC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	EE		MOC-PRIMARY-SITE	P	5.0	N	
			Historic value of TDH-PRIMARY-SITE. If MOC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. Foreign key to TBPCT.				
1	EF		MOC-OBSERVATION-MORPHOLOGY-DE	P	5.0	N	
			Historic value of TDH-OBSERVATION-MORPHOLOGY-DE. If MOC-CORRECTION-TYPE is 'A', then				

FILE...: ANY-MICRO-OBS-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: MOC-MICRO-OBS-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

it is the original value; otherwise, it is the value prior to modification. Foreign key to TBPCT.

M	1	FF	MOC-OBSERVATION-QUALIFIER	P	5.0	N
---	---	----	---------------------------	---	-----	---

Historic values of TDH-OBSERVATION-QUALIFIER. If MOC-CORRECTION-TYPE is 'A', then they are the original values; otherwise, they are the values prior to modification. Foreign key to TBPCT.

M	1	FG	MOC-OBSERVATION-SITE	P	5.0	N
---	---	----	----------------------	---	-----	---

Historic values of TDH-OBSERVATION-SITE. If MOC-CORRECTION-TYPE is 'A', then they are the original values; otherwise, they are the values prior to modification. Foreign key to TBPCT.

1	EJ		MOC-MICRO-OBS-PATHOLOGIST-ID	P	5.0	N
---	----	--	------------------------------	---	-----	---

Historic value of TDH-MICRO-OBS-PATHOLOGIST-ID. If MOC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	DA		MOC-COD-FLAG	A	1.0	F
---	----	--	--------------	---	-----	---

Historic value of TDH-COD-FLAG. If MOC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.

1	EI		MOC-TRACE-LESION-NUMBER	P	5.0	N
---	----	--	-------------------------	---	-----	---

FILE...: ANY-MICRO-OBS-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: MOC-MICRO-OBS-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
		Historic value of TDH-TRACE- LESION-NUMBER. If MOC-CORREC- TION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modi- fication.				
1	CA	MOC-CORRECTION-OPERATOR-ID	P	5.0	F	
		ID of operator who made corresponds- rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	CB	MOC-CORRECTION-REASON	A	2.0	F	
		The correction reason. A for- eign key to TBCORRREASON.				
1	CC	MOC-CORRECTION-NOTE	A	79.0	N	
		Free text entered by correc- operator describing the correction (optional).				
1	CD	MOC-CORRECTION-TYPE	A	1.0	F	
		Type of correction, A for add, C for change, D for delete.				
1	BB	MOC-MICRO-OBS-CORR-KEY-LINK-UQ	A	42.0	N	UQ
		If MOC-CORRECTION-TYPE is 'C' and either TDH-ORGAN-OBSERVA- TION-ID-DE or TDH-MICRO-OBS- SEQUENCE was changed, then the value is built by concatenat- ing the prior values of TDH- ORGAN-OBSERVATION-ID-DE and TDH-MICRO-OBS-SEQUENCE from TD-ORGAN-DATA, and CORRECTION- DATE and CORRECTION-TIME from ANY-MICRO-OBS-TRANSACTION. Foreign key to ANY-MICRO-OBS- CORRECTION itself.				

ANY-MICRO-ORGAN-CORRECTION

Logical view for tissue status observation audit changes. Master file is TD-ORGAN. Not used for retrieval by any TDMS component.

FILE...: ANY-MICRO-ORGAN-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AMC-MICRO-ORG-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AM		AMC-MICRO-ORG-CORR-KEY-UQ	A	32.0	N	UQ
			Composite field formed at Update by concatenating the value of TDO-EXP-TST-ANUM-ORG-PATH-T-UQ after the transaction has been applied, and CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 22 characters form a foreign key to TD-ORGAN.				
1	ED		AMC-ORGAN-STATUS-DATE	P	5.0	N	
			Historic value of TDO-ORGAN-STATUS-DATE. If AMC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EE		AMC-ORGAN-STATUS-TIME	P	5.0	N	
			Historic value of TDO-ORGAN-STATUS-TIME. If AMC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EF		AMC-ORGAN-STATUS-OPERATOR-ID	P	5.0	N	
			Historic value of TDO-ORGAN-STATUS-OPERATOR-ID. If AMC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				

FILE...: ANY-MICRO-ORGAN-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: AMC-MICRO-ORG-CORR-KEY-UQ

T L DB NAME	F	LENG S DE
1 EG AMC-ORGAN-STATUS Historic value of TDO-ORGAN-STATUS. If AMC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. Foreign key to TBORGSTA.	P	5.0 N
1 EH AMC-ORGAN-STAT-PATHOLOGIST-ID Historic value of TDO-ORGAN-STAT-PATHOLOGIST-ID. If AMC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.	P	5.0 N
1 CA AMC-CORRECTION-OPERATOR-ID ID of operator who made corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.	P	5.0 F
1 CB AMC-CORRECTION-REASON The correction reason. Foreign key to TBCORREASON.	A	2.0 F
1 CC AMC-CORRECTION-NOTE Free text entered by correction operator describing the correction (optional).	A	79.0 N
1 CD AMC-CORRECTION-TYPE Type of correction, A for add, C for change, D for delete.	A	1.0 F
1 AN AMC-MICRO-ORG-CORR-KEY-LINK-UQ No value is possible.	A	32.0 N UQ

BOTTLE-WEIGHT-CORRECTION

Logical view of bottle weight corrections. Master table is TD-CAGE-DATA. Used by EIS Report 22.

FILE...: BOTTLE-WEIGHT-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: BWC-BOTTLE-WGT-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AO		BWC-BOTTLE-WGT-CORR-KEY-UQ	A	31.0	N	UQ
			Composite field formed at Update by concatenating the values of TDG-EXP-TEST-CAGE-DATE and TDG-WATER-TIME after the transaction has been applied, and CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 16 characters form a foreign key to TD-CAGE-DATA.				
1	CA		BWC-BOTTLE-OPERATOR-ID	P	5.0	F	
			Historic value of TDG-WATER-OPERATOR-ID. If BWC-CORRECTION-OPERATOR-ID is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	DA		BWC-MANUAL-ENTRY-FLAG-OLD-BOTT	A	1.0	F	
			Historic value of TDG-MANUAL-ENTRY-FLAG-OLD-BOTT. If BWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EL		BWC-OLD-BOTTLE-WEIGHT	P	11.0	N	
			Historic value of TDG-OLD-BOTTLE-WEIGHT. If BWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				

FILE...: BOTTLE-WEIGHT-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: BWC-BOTTLE-WGT-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
		fication.				
1	ED	BWC-OLD-BOTTLE-STATUS	P	5.0	N	
		Historic value of TDG-OLD-BOTTLE-STATUS. If BWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. Foreign key to TBWGTSTATUS.				
1	DB	BWC-MANUAL-ENTRY-FLAG-NEW-BOTT	A	1.0	F	
		Historic value of TDG-MANUAL-ENTRY-FLAG-NEW-BOTT. If BWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EM	BWC-NEW-BOTTLE-WEIGHT	P	11.0	N	
		Historic value of TDG-NEW-BOTTLE-WEIGHT. If BWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EE	BWC-CORRECTION-OPERATOR-ID	P	5.0	N	
		ID of operator who made corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	CB	BWC-CORRECTION-REASON	A	2.0	F	
		The correction reason. A foreign key to TBCORREASON.				
1	CC	BWC-CORRECTION-NOTE	A	79.0	N	
		Free text entered by correction operator describing the correction (optional).				
1	CD	BWC-CORRECTION-TYPE	A	1.0	F	
		Type of correction, A for add,				

FILE...: BOTTLE-WEIGHT-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: BWC-BOTTLE-WGT-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
			C for change, D for delete.				
1	AP	BWC-BOTTLE-WGT-CORR-KEY-LINK-UQ	A	31.0	N	UQ	

If BWC-CORRECTION-TYPE is 'C'

and either TDG-EXP-TEST-CAGE-

DATE-DE or TDG-WATER-TIME

was changed, then the value is

built by concatenating the

prior values of TDG-EXP-TEST-

CAGE-DATE-DE and TDG-WATER-

TIME from TD-CAGE-DATA, and

CORRECTION-DATE and CORREC-

TION-TIME from BOTTLE-WEIGHT-

TRANSACTION. The first 16

characters form a foreign key

to TD-CAGE-DATA.

CAGE-ID-CORRECTION

Logical view of cage ID corrections. Master table is TD-CAGE. Used by EIS Report 22.

FILE...: CAGE-ID-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: CIC-CAGE-ID-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AQ		CIC-CAGE-ID-CORR-KEY-UQ	A	21.0	N	UQ
			Composite field formed at Up-date by concatenating the value of TDC-EXP-TEST-CAGE-UQ after the transaction has been applied, and CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 11 characters form a foreign key to TD-CAGE.				
1	ED		CIC-CAGE-OPERATOR-ID	P	5.0	N	
			Historic value of TDC-CAGE-OPERATOR-ID. If CIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	EE		CIC-CAGE-DATE	P	5.0	N	
			Historic value of TDC-CAGE-DATE. If CIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EF		CIC-CAGE-TIME	P	5.0	N	
			Historic value of TDC-CAGE-TIME. If CIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				

FILE...: CAGE-ID-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: CIC-CAGE-ID-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
1	EA	CIC-TREATMENT-NUMBER	P	3.0	N	
		Historic value of TDC-TREATMENT-NUMBER. If CIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the test number, a foreign key to PAS_TREATMNT_SEG.				
1	EG	CIC-PROCEDURE-ACTION-SET-NUM	P	5.0	N	
		Historic value of TDC-PROCEDURE-ACTION-SET-NUM. If CIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the test number, a foreign key to PAS_PROCEDUR_SEG.				
1	EH	CIC-START-DATE	P	5.0	N	
		Historic value of TDC-START-DATE. If CIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EI	CIC-AREA-NUMBER	P	5.0	N	
		Historic value of TDC-AREA-NUMBER. If CIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EJ	CIC-CORRECTION-OPERATOR-ID	P	5.0	N	
		The ID of the operator who made the correction. When appended to the appropriate facility code, a foreign key				

FILE...: CAGE-ID-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: CIC-CAGE-ID-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
		to TB_OPERATOR_ID.				
1	DD	CIC-CORRECTION-REASON	A	2.0	N	
		The correction reason. A foreign key to TBCORREASON.				
1	DH	CIC-CORRECTION-NOTE	A	79.0	N	
		Free text entered by correction operator describing the correction (optional).				
1	DA	CIC-CORRECTION-TYPE	A	1.0	F	
		Type of correction, A for add, C for change, D for delete.				
1	AR	CIC-CAGE-ID-CORR-KEY-LINK-UQ	A	21.0	N	UQ
		No possible value.				

CAGE-NOTE-CORRECTION

Logical view of cage note corrections. Master file is TD-CAGE-DATA.
Used by EIS Report 22.

FILE...: CAGE-NOTE-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: CNC-CAGE-NOTE-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AS		CNC-CAGE-NOTE-CORR-KEY-UQ	A	31.0	N	UQ
Composite field formed at Update by concatenating the value of TDG-EXP-TEST-CAGE-DATE-DE and TDG-NOTE-TIME after the transaction is applied, and CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 16 characters form a foreign key to TD-CAGE-DATA.							
1	ED		CNC-NOTE-OPERATOR-ID	P	5.0	N	
Historic value of TDC-NOTE-OPERATOR-ID. If CNC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR-ID.							
M 1	DI		CNC-CAGE-NOTE	A	80.0	N	
Historic value of TDC-CAGE-NOTE. If CNC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.							
1	DA		CNC-CAGE-NOTE-MANUAL-FLAG	A	1.0	F	
Historic value of TDG-CAGE-NOTE-MANUAL-FLAG. If CNC-CORRECTION-TYPE is 'C', then it is the original value; otherwise, it is the value prior to							

FILE...: CAGE-NOTE-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: CNC-CAGE-NOTE-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

modification.

1	CA	CNC-CORRECTION-OPERATOR-ID	P	5.0	F		
		ID of operator who made corresponds-					
		rection. When appended to the					
		appropriate facility code, a					
		foreign key to TB_OPERATOR_ID.					
1	CB	CNC-CORRECTION-REASON	A	2.0	F		
		The correction reason. A for-					
		oreign key to TBCORREASON.					
1	CC	CNC-CORRECTION-NOTE	A	79.0	N		
		Free text entered by correc-					
		operator describing the					
		correction (optional).					
1	CD	CNC-CORRECTION-TYPE	A	1.0	F		
		Type of correction, A for add,					
		C for change, D for delete.					
1	AT	CNC-CAGE-NOTE-CORR-KEY-LINK-UQ	A	31.0	N	UQ	
		If CNC-CORRECTION-TYPE is 'C'					
		and either TDG-EXP-TEST-CAGE-					
		DATE-DE or TDG-NOTE-TIME					
		was changed, then the value is					
		built by concatenating the					
		prior values of TDG-EXP-TEST-					
		CAGE-DATE-DE and TDG-NOTE-					
		TIME from TD-CAGE-DATA, and					
		CORRECTION-DATE and CORREC-					
		TION-TIME from CAGE-NOTE-					
		TRANSACTION. A foreign key					
		to CAGE-NOTE-CORRECTION it-					
		self.					

FEEDER-WEIGHT-CORRECTION

Logical view of feeder weight corrections. Master file is TD-CAGE-DATA. Used by EIS Report 22.

FILE...: FEEDER-WEIGHT-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: FWC-FEEDER-WGT-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AU		FWC-FEEDER-WGT-CORR-KEY-UQ	A	31.0	N	UQ
			Composite field formed at Update by concatenating the values of TDG-EXP-TEST-CAGE-DATE-DE and TDG-FEEDER-TIME after the transaction is applied, and CORRECTION-DATE and CORRECTION-TIME from the transaction. First 16 characters form a foreign key to TD-CAGE-DATA.				
1	ED		FWC-FEEDER-OPERATOR-ID	P	5.0	N	
			Historic value of TDG-FEEDER-OPERATOR-ID. If FWC-CORRECTION-OPERATOR-ID is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	DA		FWC-MANUAL-ENTRY-FLAG-OLD-FEED	A	1.0	F	
			Historic value of TDG-MANUAL-ENTRY-FLAG-OLD-FEED. If FWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EL		FWC-OLD-FEEDER-WEIGHT	P	11.0	N	
			Historic value of TDG-OLD-FEEDER-WEIGHT. If FWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				

FILE...: FEEDER-WEIGHT-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: FWC-FEEDER-WGT-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
		fication.				
1	EE	FWC-OLD-FEEDER-STATUS	P	5.0	N	
		Historic value of TDG-OLD-FEEDER-STATUS. If FWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. A foreign key to TBWGTSTATUS.				
1	DB	FWC-MANUAL-ENTRY-FLAG-NEW-FEED	A	1.0	F	
		Historic value of TDG-MANUAL-ENTRY-FLAG-NEW-FEED. If FWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EM	FWC-NEW-FEEDER-WEIGHT	P	11.0	N	
		Historic value of TDG-NEW-FEEDER-WEIGHT. If FWC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	CA	FWC-CORRECTION-OPERATOR-ID	P	5.0	F	
		ID of operator who made corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
1	CB	FWC-CORRECTION-REASON	A	2.0	F	
		The correction reason. A foreign key to TBCORREASON.				
1	CC	FWC-CORRECTION-NOTE	A	79.0	N	
		Free text entered by correction operator describing the correction (optional).				
1	CD	FWC-CORRECTION-TYPE	A	1.0	F	
		Type of correction, A for add,				

FILE...: FEEDER-WEIGHT-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: FWC-FEEDER-WGT-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
		C for change, D for delete.				
1	AV	FWC-FEEDER-WGT-CORR-KEY-LINK-UQ	A	31.0	N	UQ
		If FWC-CORRECTION-TYPE is 'C'				
		and either TDG-EXP-TEST-CAGE-				
		DATE-DE or TDG-FEEDER-TIME				
		was changed, then the value is				
		built by concatenating the				
		prior values of TDG-EXP-TEST-				
		CAGE-DATE-DE and TDG-FEEDER-				
		TIME from TD-CAGE-DATA, and				
		CORRECTION-DATE and CORREC-				
		TION-TIME from FEEDER-WEIGHT-				
		TRANSACTION. A foreign key				
		to FEEDER-WEIGHT-CORRECTION				
		itself.				

HISTOLOGY-NUMBER-CORRECTION

Logical view of histology number audit changes. Master file is TD-ANIMAL. Not used for retrieval by any TDMS component.

FILE...: HISTOLOGY-NUMBER-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: HIC-HISTOLOGY-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	AW		HIC-HISTOLOGY-CORR-KEY-UQ	A	26.0	N	UQ
			Composite field formed at Update by concatenating the values of TDA-EXP-TEST-ANIMAL-UQ, CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 16 characters are a foreign key to TD-ANIMAL.				
1	ED		HIC-PATHOLOGIST-IDENTIFICATION	P	5.0	N	
			Historic value of TDA-PATHOLOGIST-IDENTIFICATION. If HIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EE		HIC-HISTOLOGY-NUMBER-DATE	P	5.0	N	
			Historic value of TDA-HISTOLOGY-NUMBER-DATE. If HIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EF		HIC-HISTOLOGY-NUMBER-TIME	P	5.0	N	
			Historic value of TDA-HISTOLOGY-NUMBER-TIME. If HIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EG		HIC-HISTOLOGY-NUM-OPERATOR-ID	P	5.0	N	
			Historic value of TDA-HISTOLOGY-NUM-OPERATOR-ID. If HIC-CORRECTION-TYPE is 'A', then				

FILE...: HISTOLOGY-NUMBER-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: HIC-HISTOLOGY-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	DF	HIC-HISTOLOGY-NUMBER	A	16.0	N		
Historic value of TDA-HISTOLOGY-NUMBER. If HIC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.							
1	EH	HIC-CORRECTION-OPERATOR-ID	P	5.0	N		
ID of operator who made corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.							
1	DD	HIC-CORRECTION-REASON	A	2.0	N		
The correction reason. A foreign key to TBCORREASON.							
1	DH	HIC-CORRECTION-NOTE	A	79.0	N		
Free text entered by correction operator describing the correction (optional).							
1	DA	HIC-CORRECTION-TYPE	A	1.0	F		
Type of correction, A for add, C for change, D for delete.							
1	AX	HIC-HISTOLOGY-CORR-KEY-LINK-UQ	A	26.0	N	UQ	
No value is possible.							

MICRO-SITE-STATUS-CORRECTION

Logical view of accountable site status observation corrections. Master file is TD-ORGAN. Not used for retrieval by any TDMS component.

FILE...: MICRO-SITE-STATUS-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: SSC-MICRO-SITE-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	BC		SSC-MICRO-SITE-CORR-KEY-UQ	A	42.0	N	UQ
Composite field formed at Up-date by concatenating the value of TDO-EXP-TST-ANUM-ORG-PATHT-UQ, the values of TDO-SITE-STATUS-DATE and TDO-SITE-STATUS-TIME after the transaction is applied, and CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 22 characters form a foreign key to TD-ORGAN.							
1	ED		SSC-SITE-STATUS-PATHOLOGIST-ID	P	5.0	N	
Historic value of TDO-SITE-STATUS-PATHOLOGIST-ID. If SSC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.							
1	EE		SSC-SITE-STATUS-OPERATOR-ID	P	5.0	N	
Historic value of TDO-SITE-STATUS-OPERATOR-ID. If SSC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.							
1	EF		SSC-SITE-CODE	P	5.0	N	
Historic value of TDO-SITE-							

FILE...: MICRO-SITE-STATUS-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: SSC-MICRO-SITE-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

CODE. If SSC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. A foreign key to TBPCT.

1	EG	SSC-SITE-STATUS	P	5.0	N	
---	----	-----------------	---	-----	---	--

Historic value of TDO-SITE-STATUS. If SSC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. A foreign key to TBORGSTA.

1	CA	SSC-CORRECTION-OPERATOR-ID	P	5.0	F	
---	----	----------------------------	---	-----	---	--

ID of operator who made corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	CB	SSC-CORRECTION-REASON	A	2.0	F	
---	----	-----------------------	---	-----	---	--

The correction reason. A foreign key to TBCORREASON.

1	CC	SSC-CORRECTION-NOTE	A	79.0	N	
---	----	---------------------	---	------	---	--

Free text entered by correction operator describing the correction (optional).

1	CD	SSC-CORRECTION-TYPE	A	1.0	F	
---	----	---------------------	---	-----	---	--

Type of correction, A for add, C for change, D for delete.

1	BD	SSC-MICRO-SITE-CORR-KEY-LINK-UQ	A	42.0	N	UQ
---	----	---------------------------------	---	------	---	----

If SSC-CORRECTION-TYPE is 'C' and either TDO-SITE-STATUS-DATE or TDO-SITE-STATUS-TIME was changed, then the value is built by concatenating the prior values of TDO-EXP-TST-ANUM-ORG-PATHT-UQ, TDO-SITE-STATUS-DATE, and TDO-SITE-

FILE...: MICRO-SITE-STATUS-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: SSC-MICRO-SITE-CORR-KEY-UQ

T L DB NAME	F	LENG S DE
STATUS-TIME from TD-ORGAN, and CORRECTION-DATE and CORREC- TION-TIME from MICRO-SITE- STATUS-TRANSACTION. A for- eign key to MICRO-SITE-STATUS- CORRECTION itself.		

NOT-EXAMINED-CORRECTION

Logical view of not examined audit changes. Master file is TD-ANIMAL. Not used for retrieval by any TDMS component.

FILE...: NOT-EXAMINED-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: NEC-NO-EXAM-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	FJ	NEC-NO-EXAM-CORR-KEY-UQ	A	26.0	N	UQ	
Composite field formed at Update by concatenating the value of TDA-EXP-TEST-ANIMAL-UQ with CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 16 characters form a foreign key to TD-ANIMAL.							
1	ED	NEC-NO-EXAM-DATE	P	5.0	N		
Historic value of TDA-NO-EXAM-DATE. If NEC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.							
1	EE	NEC-NO-EXAM-TIME	P	5.0	N		
Historic value of TDA-NO-EXAM-TIME. If NEC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.							
1	EF	NEC-NO-EXAM-OPERATOR-ID	P	5.0	N		
Historic value of TDA-NO-EXAM-OPERATOR-ID. If NEC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR-ID.							
1	EG	NEC-NO-EXAM-PATH-ID	P	5.0	N		

FILE...: NOT-EXAMINED-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: NEC-NO-EXAM-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

Historic value of TDA-NO-EXAM-PATHOLOGIST. If NEC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	CA	NEC-CORRECTION-OPERATOR-ID	P	5.0	F	
ID of operator who made corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.						
1	CB	NEC-CORRECTION-REASON	A	2.0	F	
The correction reason. A foreign key to TBCORREASON.						
1	CC	NEC-CORRECTION-NOTE	A	79.0	N	
Free text entered by correction operator describing the correction (optional).						
1	CD	NEC-CORRECTION-TYPE	A	1.0	F	
Type of correction, A for add, C for change, D for delete.						
1	FK	NEC-NO-EXAM-CORR-KEY-LINK-UQ	A	26.0	N	UQ
No value is possible.						

ORGAN-NOTES-CORRECTION

Logical view of organ note audit changes. Master file is TD-ORGAN-DATA. Not used for retrieval by any TDMS component.

FILE...: ORGAN-NOTES-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ONC-ORGAN-NOTE-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	BE		ONC-ORGAN-NOTE-CORR-KEY-UQ	A	42.0	N	UQ
			Composite field formed at Up-date by concatenating the values of TDH-ORGAN-OBSERVATION-ID-DE and TDH-ORGAN-NOTE-SEQUENCE after the transaction is applied, and CORRECTION-DATE and CORRECTION-TIME from the transaction. The first 27 characters form a foreign key to TD-ORGAN-DATA.				
1	ED		ONC-ORGAN-NOTE-OPERATOR-ID	P	5.0	N	
			Historic value of TDH-ORGAN-NOTE-OPERATOR-ID. If ONC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.				
M 1	DI		ONC-ORGAN-NOTE	A	80.0	N	
			Historic value of TDH-ORGAN-NOTE. If ONC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EE		ONC-ORGAN-NOTE-PATHOLOGIST-ID	P	5.0	N	
			Historic value of TDH-ORGAN-NOTE-PATHOLOGIST-ID. If ONC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior				

FILE...: ORGAN-NOTES-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: ONC-ORGAN-NOTE-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

to modification. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	CA	ONC-CORRECTION-OPERATOR-ID	P	5.0	F	
---	----	----------------------------	---	-----	---	--

ID of operator who made corresponds-rection. When appended to the appropriate facility code, a foreign key to TB_OPERATOR_ID.

1	CB	ONC-CORRECTION-REASON	A	2.0	F	
---	----	-----------------------	---	-----	---	--

The correction reason. A foreign key to TBCORREASON.

1	CC	ONC-CORRECTION-NOTE	A	79.0	N	
---	----	---------------------	---	------	---	--

Free text entered by correction operator describing the correction (optional).

1	CD	ONC-CORRECTION-TYPE	A	1.0	F	
---	----	---------------------	---	-----	---	--

Type of correction, A for add, C for change, D for delete.

1	BF	ONC-ORGAN-NOTE-CORR-KEY-LINK-UQ	A	42.0	N	UQ
---	----	---------------------------------	---	------	---	----

If ONC-CORRECTION-TYPE is 'C' and either TDH-ORGAN-OBSERVATION-ID-DE or TDH-ORGAN-NOTE-SEQUENCE was changed, then the value is built by concatenating the prior values of TDH-ORGAN-OBSERVATION-ID-DE and TDH-ORGAN-NOTE-SEQUENCE from TD-ORGAN-DATA, and CORRECTION-DATE and CORRECTION-TIME from

ORGAN-NOTE-TRANSACTION. A

foreign key to ORGAN-NOTES-CORRECTION itself.

TEST-NOTE-CORRECTION

Logical view for test note corrections. Master file is TD-TEST-NOTE.
Used by EIS Report 22.

FILE...: TEST-NOTE-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: TNC-TEST-NOTE-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	BG	TNC-TEST-NOTE-CORR-KEY-UQ	Primary key for TEST-NOTE-CORRECTION. Effects a link with the master record in TD-TEST-NOTE. The following elements are concatenated by Update to form the field's value: the current value of TDS-EXP-TEST-DATE-TIME-UQ from TD-TEST-NOTE, and CORRECTION-DATE and CORRECTION-TIME from TEST-NOTE-TRANSACTION.	A	27.0	N	UQ
1	ED	TNC-TEST-OPERATOR-ID	Historic value of TDS-OPER-NUM. If TNC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.	P	5.0	N	
1	DJ	TNC-TEST-NOTE	Historic value of TDS-NOTE-TEXT. If TNC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.	A	160.0	N	
1	CA	TNC-CORRECTION-OPERATOR-ID	ID of operator who made corresponds-rection.	P	5.0	F	
1	CB	TNC-CORRECTION-REASON	The correction reason (TBCORR-REASON) code assigned to the correction.	A	2.0	F	
1	CC	TNC-CORRECTION-NOTE	Free text entered by correc-	A	79.0	N	

FILE...: TEST-NOTE-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: TNC-TEST-NOTE-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
		operator describing the correction (optional).				
1	CD	TNC-CORRECTION-TYPE	A	1.0	F	
		Type of correction, A for add, C for change, D for delete.				
1	BH	TNC-TEST-NOTE-CORR-KEY-LINK-UQ	A	27.0	N	UQ
		If TNC-CORRECTION-TYPE is 'C' and TDS-EXP-TEST-DATE-TIME-UQ was changed, then the value is built by concatenating the prior value of TDS-EXP-TEST- DATE-TIME-UQ from TD-TEST- NOTE, and CORRECTION-DATE and CORRECTION-TIME from TEST- NOTE-TRANSACTION. This ef- fects a link with previous corrections of the same master record.				

TREAT-NOTE-CORRECTION

Logical view of treatment note corrections. Master file is TD-TREATMENT-NOTE. Used by EIS Report 22.

FILE...: TREAT-NOTE-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87RR-KEY-UQ

PRIMARY SEQUENCE FIELD: TGC-TREAT-NOTE-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	BI	TGC-TREAT-NOTE-CORR-KEY-UQ	Primary key for TREAT-NOTE-CORRECTION. Effects a link with the master record in TD-TREAT-NOTE. The following elements are concatenated by Update to form the field's value: the current value of TDT-EXP-TRT-DATE-TIME-UQ from TD-TREAT-NOTE, and CORRECTION-DATE and CORRECTION-TIME from TREAT-NOTE-TRANSACTION.	A	30.0	N	UQ
1	ED	TGC-TREAT-OPERATOR-ID	Historic value of TDT-OPER-NUM. If TGC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.	P	5.0	N	
1	DJ	TGC-TREAT-NOTE	Historic value of TDT-NOTE-TEXT. If TGC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.	A	160.0	N	
1	CA	TGC-CORRECTION-OPERATOR-ID	ID of operator who made correspondence.	P	5.0	F	
1	CB	TGC-CORRECTION-REASON	The correction reason (TBCORREASON) code assigned to the correction.	A	2.0	F	
1	CC	TGC-CORRECTION-NOTE	Free text entered by correc-	A	79.0	N	

FILE...: TREAT-NOTE-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87RR-KEY-UQ

PRIMARY SEQUENCE FIELD: TGC-TREAT-NOTE-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

operator describing the
correction (optional).

1	CD	TGC-CORRECTION-TYPE	A	1.0	F	
---	----	---------------------	---	-----	---	--

Type of correction, A for add,
C for change, D for delete.

1	BJ	TGC-TREAT-NOTE-CORR-KEY-LINK-UQ	A	30.0	N	UQ
---	----	---------------------------------	---	------	---	----

If TGC-CORRECTION-TYPE is 'C'
and TDT-EXP-TEST-TRT-DATE-
TIME-UQ was changed, then the
value is built by concatenat-
ing the prior value of TDT-
EXP-TEST-TRT-DATE-TIME-UQ from
TD-TREAT-NOTE, and CORRECTION-
DATE and CORRECTION-TIME from
TREAT-NOTE-TRANSACTION. This
effects a link with previous
corrections of the same master
record.

UNCERTAIN-CORRECTION

Logical view of uncertain cause-of-death audit changes. Master table is TD-ANIMAL. Not used for retrieval by any TDMS component.

FILE...: UNCERTAIN-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: UCC-UNCERTAIN-CORR-KEY-UQ

T	L	DB	NAME	F	LENG	S	DE
1	FH		UCC-UNCERTAIN-CORR-KEY-UQ	A	26.0	N	UQ
			Primary key for UNCERTAIN-CORRECTION. Effects a link with the master record in TD-ANIMAL. The following elements are concatenated by Update to form the field's value: the current value of TDA-EXP-TEST-ANIMAL-DATE-UQ from TD-ANIMAL, and CORRECTION-DATE and CORRECTION-TIME from UNCERTAIN-TRANSACTION.				
1	ED		UCC-UNCERTAIN-DATE	P	5.0	N	
			Historic value of TDA-UNCERTAIN-DATE. If UCC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EE		UCC-UNCERTAIN-TIME	P	5.0	N	
			Historic value of TDA-UNCERTAIN-DATE. If UCC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EF		UCC-UNCERTAIN-OPERATOR-ID	P	5.0	N	
			Historic value of TDA-UNCERTAIN-OPERATOR-ID. If UCC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.				
1	EG		UCC-UNCERTAIN-PATH-ID	P	5.0	N	

FILE...: UNCERTAIN-CORRECTION

TYPE...: USER VIEW

FILE-NR: 87

PRIMARY SEQUENCE FIELD: UCC-UNCERTAIN-CORR-KEY-UQ

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

Historic value of TDA-UNCERTAIN-PATHOLOGIST. If UCC-CORRECTION-TYPE is 'A', then it is the original value; otherwise, it is the value prior to modification.

1	CA	UCC-CORRECTION-OPERATOR-ID	P	5.0	F	
ID of operator who made corresponds-rection.						
1	CB	UCC-CORRECTION-REASON	A	2.0	F	
The correction reason (TBCORR-REASON) code assigned to the correction.						
1	CC	UCC-CORRECTION-NOTE	A	79.0	N	
Free text entered by correction operator describing the correction (optional).						
1	CD	UCC-CORRECTION-TYPE	A	1.0	F	
Type of correction, A for add, C for change, D for delete.						
1	FI	UCC-UNCERTAIN-CORR-KEY-LINK-UQ	A	26.0	N	UQ
No value is possible.						

DESTINATION_TBL

View of valid TDMS output destinations. Used by TDMS menus to validate output destinations on TDMS screens.

FILE...: DESTINATION_TBL

TYPE...: USER VIEW

FILE-NR: 35

PRIMARY SEQUENCE FIELD: TABLE-DESTINATION

T	L	DB	NAME	F	LENG	S	DE
1	AA	TABLE-DEST-TYPE		A	1.0	F	
			Flag indicating type of description- tination. H - Lab; I - Inact- ive; V - VAX queue.				
1	AB	TABLE-DESTINATION		A	10.0	N	DE
			Valid output destination. (VAX queues longer than 10 characters must be redefined as a logical of no more than 10 characters. That logical would be entered here, and would be what users enter on TDMS screens.)				
1	AC	TABLE-FACILITY		N	5.0	N	UQ
			If a lab output destination has been assigned a TDMS fac- ility code, that would be entered here.				

4.5 TDMS Table & PCT Views

TABLES_DESCRIPTN

Logical layout of Tables Description table, which is used by NATURAL.

FILE...: TABLES_DESCRIPTN

TYPE...: USER VIEW

FILE-NR: 16

PRIMARY SEQUENCE FIELD: TABLE_NAME

T	L	DB	NAME	F	LENG	S	DE
G	1	MF	TABLES_DESCRIPTION_FILE				
	2	TN	TABLE_NAME	A	8.0	N	DE
			Primary key for table. Contains the eight-character view names for each TDMS table.				
	2	TB	TABLE_NUM	N	2.0	F	
			The file number assigned to the file when loaded into ADABAS with ADALOD or ADAFDU.				
	2	NC	NUM_COLS	N	3.0		
			The number of columns defined in the FIELD_GROUP repeating group.				
P	1	FG	FIELD_GROUP				
			Repeating group containing each field in table. This group is position-sensitive. Do not reorder elements unless you understand the TBGET retrieval mechanism.				
	2	FN	FIELD_NAME	A	2.0	N	
			The short field name assigned to the field when loaded into ADABAS.				
	2	FL	FIELD_LENGTH	N	3.0	N	
			The length of the field.				

TBACTONC

Logical layout of action table. Table lists actions that may be performed by animal caretakers as specified in protocol procedures. Information is downloaded to LDAS and is used in creating the procedure action schedule.

FILE...: TBACTONC

TYPE...: USER VIEW

FILE-NR: 29

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

-- -- -- -- --

G 1 AT ACTION-TABLE

Group embracing all fields in the table.

2	AA	ROW-NUMBER	N	3.0	UQ		
---	----	------------	---	-----	----	--	--

Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.

2	BB	DATE-ENTERED	A	8.0			
---	----	--------------	---	-----	--	--	--

The date of entry for a record. No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.

2	CC	ACTIVE-INACTIVE	A	1.0	F		
---	----	-----------------	---	-----	---	--	--

Status of record. A for active and I for inactive. Informational only; not enforced by PAS.

2	DD	CODE-VALUE	N	6.0	UQ		
---	----	------------	---	-----	----	--	--

Unique, sequential code to identify record within table. Primary key for record retrieval.

2	EE	EE-FILLER	A	1.0	F		
---	----	-----------	---	-----	---	--	--

Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.

2	FF	FF-FILLER	A	1.0	F		
---	----	-----------	---	-----	---	--	--

Field's sole purpose is to

FILE...: TBACTONC

TYPE...: USER VIEW

FILE-NR: 29

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
-------------	---	-----------

hold a column position due to
TBGET retrieval mechanism.
Contains no data.

2 GG LONG-TEXT	A	32.0 N
----------------	---	--------

Full description of action.

2 HH TYPE-ACTION	A	1.0 N
------------------	---	-------

Alpha code signifying to which
test object the action per-
tains. A for animal; C for
cage.

TBAGENCY

Logical layout of agency table. Contains a list of agencies. Table is a legacy to TDMS. Only one entry, that for NIEHS/NTP is valid. The code for that entry, which is '05', is build into the nine-digit test number entered in PAS, and is hard-coded by the menu system when EIS/PEIS reports are requested.

FILE...: TBAGENCY

TYPE...: USER VIEW

FILE-NR: 22

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
- - - - -						
G 1	AT	AGENCY-TABLE				
		Group embracing all fields in the table.				
2	AA	ROW-NUMBER	N	3.0	UQ	
		Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB	DATE-ENTERED	A	8.0		
		The date of entry for a record. No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD	CODE-VALUE	N	6.0	UQ	
		Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE	EE-FILLER	A	1.0	F	
		Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF	FF-FILLER	A	1.0	F	

FILE...: TBAGENCY

TYPE...: USER VIEW

FILE-NR: 22

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
-------------	---	-----------

Field's sole purpose is to
hold a column position due to
TBGET retrieval mechanism.
Contains no data.

2 GG CHARACTER-CODE	A	2.0 N
---------------------	---	-------

Legacy field. Appears to be
right-justified and padded
form of the code value.

2 HH LONG-TEXT	A	32.0 N
----------------	---	--------

Full description of the agency

2 II REPORT-HEADINGS	A	50.0 N
----------------------	---	--------

Legacy field. Appears to be
a long description of the
agency, centered within the
field string.

2 JJ SHORT-TEXT	A	12.0 N
-----------------	---	--------

Abbreviated description of the
agency.

TBANCLQF

Logical layout of animal class qualifier table. Animal class qualifiers are a legacy to TDMS, with no use in its current version. The PEX maintenance software requires an animal class qualifier, so the code for "Not Applicable" is entered.

FILE...: TBANCLQF

TYPE...: USER VIEW

FILE-NR: 23

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	CT ANIMAL-CLASS-QUALIFIER-TABLE				
		Group embracing all fields in the table.				
2	AA	ROW-NUMBER	N	3.0		UQ
		Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB	DATE-ENTERED	A	8.0		
		The date of entry for a record No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC	ACTIVE-INACTIVE	A	1.0		F
		Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD	CODE-VALUE	N	6.0		UQ
		Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE	EE-FILLER	A	1.0		F
		Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF	FF-FILLER	A	1.0		F
		Field's sole purpose is to				

FILE...: TBANCLQF
TYPE...: USER VIEW
FILE-NR: 23
PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
-------------	---	-----------

hold a column position due to
TBGET retrieval mechanism.
Contains no data.

2 GG LONG-TEXT	A	32.0 N
----------------	---	--------

Full description of the animal
class qualifier.

TBANMLSUPPLIER

Logical view of animal supplier table. Contains valid animal suppliers for tests. (There is a dummy record in table for tests that preceded entry of animal suppliers in test protocol. This record has a key value of 0 and text of "NOT SPECIFIED".)

FILE...: TBANMLSUPPLIER

TYPE...: USER VIEW

FILE-NR: 32

PRIMARY SEQUENCE FIELD: ANML-SUPPLR

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	FA	ANIMAL-SUPPLIER-TABLE				
			Group embracing all fields in the table.				
2	FB	ANML-SUPPLR		P	5.0	F	UQ
			Unique numeric code that identifies an animal supplier. (Code 0 exists on table.)				
2	FC	SUPPLR-TEXT		A	64.0	N	
			Full name of animal supplier. Appears on EIS/PEIS report splash pages. .				

TBCLOBSC

Logical layout of the clinical observation table, which contains clinical observations, in-life sites, in-life qualifiers, and removal reasons. Data from this table appears in nearly all EIS/PEIS reports in some manner or fashion.

FILE...: TBCLOBSC

TYPE...: USER VIEW

FILE-NR: 25

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	CT CLINICAL-OBSERVATION-TABLE				
		Group embracing all fields in the table.				
2	AA	ROW-NUMBER	N	3.0	UQ	
		Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB	DATE-ENTERED	A	8.0		
		The date of entry for a record. No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD	CODE-VALUE	N	6.0	F UQ	
		Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE	EE-FILLER	A	1.0	F	
		Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF	FF-FILLER	A	1.0	F	
		Field's sole purpose is to hold a column position due to				

FILE...: TBCLOBSC

TYPE...: USER VIEW

FILE-NR: 25

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
TBGET retrieval mechanism. Contains no data.		
2 GG SUBSET-INDICATOR	A	1.0 F
Indicator that identifies to which of four table subsets a record belongs. The possible values, and their meaning are: A - Clinical Observations B - In-life Sites C - In-life Qualifiers (Size,Count) D - Removal Reasons		
2 HH SHORT-TEXT	A	14.0 DE
Abbreviated description of the in-life term.		
2 II LONG-TEXT	A	33.0
Full description of the in-life term.		
2 JJ JJ-FILLER	A	1.0 N
Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.		
2 KK CLASS-VALUE	N	6.0 N
Observation class for clinical observation. Legacy field with no use in current version of TDMS.		
2 LL INCOMPAT-TOTAL	N	6.0 N
Total of all class values with which the current entry is incompatible. Legacy field with no meaning in current version of TDMS. (Incompatibility is built into LDAS.)		

FILE...: TBCLOBSC

TYPE...: USER VIEW

FILE-NR: 25

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
2	MM	SEX-FLAG	A	1.0	F	
Indicates gender-specific terms. "M" for male; "F" for female; blank for none.						
2	NN	MULTI-SELECT-FLAG	A	1.0	N	
Legacy field whose function was replaced by count qualifiers-fiers.						
2	OO	NON-CHG-FLAG	A	1.0	N	
Legacy field with no meaning in current version of TDMS.						
1	PP	CLNOBS-GROUP	P	5.0	N	
Clinical observation group to which a clinical observation belongs. Only valid for subset "A". Foreign key to TBCLOBSCGRP.						
M 1	QQ	EIQUAL-GROUP	P	5.0	N	
Valid in-life qualifier groups which may be associated with a clinical observation. Only valid for subset "A". Foreign key to TBQUALGRP.						
M 1	RR	EISITE-NUM	P	5.0	N	
Valid in-life sites which may be associated with a clinical observation. Only valid for subset "A". Foreign key to TBCLOBSC.						
1	TT	EIQUAL-GROUP-NUM	P	5.0	N	DE
In-life qualifier group to which an in-life qualifier belongs. Only valid for subset "C". Foreign key to TBQUALGRP.						

TBCLOBSCGRP

View of clinical observation groups, which classify entries in TBCLOBSC with SUBSET-INDICATOR equal to "A".

FILE...: TBCLOBSCGRP

TYPE...: USER VIEW

FILE-NR: 36

PRIMARY SEQUENCE FIELD: CLNOBS-GROUP

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	DA	CLINICAL-OBSERVATION-GROUP-TABLE				
			Group embracing each element in the table.				
2	DB		CLNOBS-GROUP	A	5.0	N	UQ
			Code which identifies a unique clinical observation group.				
2	DC		CLNOBS-GROUP-SHORT-TEXT	A	14.0	N	
			Abbreviated description of the clinical observation group.				
2	DD		CLNOBS-GROUP-LONG-TEXT	A	32.0	N	
			Full description of the clinical observation group.				

TBCLOBSI

Logical layout of clinical incompatibility table. This is a legacy table with no purpose in the current version of TDMS. Clinical observation incompatibility is simply enforced by LDAS software according to the following rules: unremarkable is incompatible with any other observation, hyperactive and lethargic are mutually incompatible.

FILE...: TBCLOBSI

TYPE...: USER VIEW

FILE-NR: 26

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G 1	IT	CLIN-OBS-INCOMPATIBILITY-TABLE				
		Group embracing all fields in the table.				
2	AA	ROW-NUMBER	N	3.0		UQ
		Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB	DATE-ENTERED	A	8.0		
		The date of entry for a record. No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC	ACTIVE-INACTIVE	A	1.0		F
		Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD	CODE-VALUE	N	6.0		UQ
		Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE	EE-FILLER	A	1.0		F
		Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				

FILE...: TBCLOBSI

TYPE...: USER VIEW

FILE-NR: 26

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
2 FF FF-FILLER	A	1.0 F
Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.		
2 GG LONG-TEXT	A	32.0 N
Full description of incompatibility class.		
2 HH CLASS-VALUE	N	6.0 DE
Value for incompatibility class, assigned as a unique bit switch, (.i.e, 1,2,4,8...)		
2 II TOTAL-CL-VALUE	N	6.0 N
Sum of classes with which a class is incompatible.		

TBCONDTC

Logical layout of cage condition table. This table is a legacy to TDMS and has no use in the current version of TDMS.

FILE...: TBCONDTC

TYPE...: USER VIEW

FILE-NR: 24

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	CT	CAGE-CONDITION-TABLE				
			Group embracing all field in the table.				
2	AA		ROW-NUMBER	N	3.0		UQ
			Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB		DATE-ENTERED	A	8.0		
			The date of entry for a record. No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC		ACTIVE-INACTIVE	A	1.0		F
			Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD		CODE-VALUE	N	6.0		UQ
			Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE		EE-FILLER	A	1.0		F
			Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF		FF-FILLER	A	1.0		F
			Field's sole purpose is to hold a column position due to TBGET retrieval mechanism.				

FILE...: TBCONDTC

TYPE...: USER VIEW

FILE-NR: 24

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
Contains no data.		
2 GG SHORT-TEXT	A	8.0 N
Abbreviated description of the cage condition.		
2 HH LONG-TEXT	A	32.0 N
Full description of the cage condition.		
2 II II-FILLER	A	1.0 N
Field's sole purpose is to hold a column position due to TBGET retrieval mechanism.		
Contains no data.		
2 JJ CAGE-COND-INCOMP	N	6.0 N
Incompatibility class for cage condition.		

TBCORREASON

Logical view of correction reason table. Contains the valid correction reasons that can be used to classify an entry made in the error correction system.

FILE...: TBCORREASON

TYPE...: USER VIEW

FILE-NR: 31

PRIMARY SEQUENCE FIELD: CORR-CODE

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	CA	CORRECTION-REASON-TABLE				
			Group embracing all fields in the correction reason table.				
2	CB		CORR-CODE	A	2.0	N	UQ
			Unique, two-letter, code associated with a correction reason. Displayed on E22.				
2	CC		CORR-REASON	A	30.0	N	
			Full description of correction reason.				

TBDOSRTE

Logical layout of Dose Route table. This table contains all routes by which a dose for a test may be administered. There are five categories to which current dose routes belong, which are: food, water, gavage, skin paint, and respiration.

FILE...: TBDOSRTE

TYPE...: USER VIEW

FILE-NR: 18

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE

G	1	DT DOSE-ROUTE-TABLE				
		Group embracing all fields in table.				
2	AA	ROW-NUMBER	N	3.0	UQ	
		Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB	DATE-ENTERED	A	8.0		
		The date of entry for a record No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD	CODE-VALUE	N	6.0	UQ	
		Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE	EE-FILLER	A	1.0	F	
		Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF	FF-FILLER	A	1.0	F	
		Field's sole purpose is to				

FILE...: TBDOSRTE
TYPE...: USER VIEW
FILE-NR: 18
PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
-------------	---	-----------

hold a column position due to
TBGET retrieval mechanism.
Contains no data.

2 GG LONG-TEXT	A	32.0 N
----------------	---	--------

Full description of dose
route. Used in standard
EIS/PEIS report headers.

2 HH SHORT-TEXT	A	8.0 N
-----------------	---	-------

Abbreviated description of the
dose route.

2 II GENERAL-ROUTE	N	6.0 N
--------------------	---	-------

Refers to another record in
the table, which record is
the route category to which
a specific record belongs.

TBFACLTY

Logical layout of facility table. Contains all facilities for TDMS, current and former. Used in EIS/PEIS reports splash pages, PAS validation reports, and almost every other part of the application.

FILE...: TBFACLTY

TYPE...: USER VIEW

FILE-NR: 27

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
- - - - -						
G	1	FT FACILITIES-TABLE				
		Group embracing all fields in the table.				
2	AA	ROW-NUMBER	N	3.0		UQ
		Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB	DATE-ENTERED	A	8.0		
		The date of entry for a record No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC	ACTIVE-INACTIVE	A	1.0		F
		Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD	CODE-VALUE	N	6.0		UQ
		Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE	EE-FILLER	A	1.0		F
		Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF	FF-FILLER	A	1.0		F
		Field's sole purpose is to hold a column position due to				

FILE...: TBFACLT

TYPE...: USER VIEW

FILE-NR: 27

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
		TBGET retrieval mechanism. Contains no data.				
2	GG	CHAR-ENTRY-CODE	A	5.0	N	
		Code value right-justified and padded with 0's for bar-code use. Legacy field with no use in current version of TDMS				
2	HH	LONG-TEXT	A	50.0	N	DE
		Full name of facility.				
2	II	CITY-STATE	A	25.0	N	
		City and state where facility is located. Informational only.				
2	JJ	SHORT-TEXT	A	23.0	N	
		Abbreviated title of facility.				
1	KK	LAB-CODE	A	2.0	N	DE
		CHEMTRACK code for facility.				

TBINLIFETERMS

Logical view of in-life terms descriptions table. Contains synonyms for clinical observation codes. Used as part of LDAS help system.

FILE...: TBINLIFETERMS

TYPE...: USER VIEW

FILE-NR: 33

PRIMARY SEQUENCE FIELD: TERM-NUM

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	GL	IN-LIFE-TERMS-TABLE				
			Group embracing all fields in the table.				
2	GM		TERM-NUM	P	5.0	N	UQ
			TBCLOBSC code for which there is a list of synonyms.				
2	GN		DEFINE-TEXT	A	200.0	N	
			List of synonyms or descriptions of current clinical observations.				

TBORGSTA

Logical layout of organ (tissue) status table. Contains the possible statuses by which an organ (tissue) that is being microscopically evaluated may be classified.

FILE...: TBORGSTA

TYPE...: USER VIEW

FILE-NR: 34

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	ST	ORGAN-STATUS-TABLE				
			Group embracing all fields in the organ status table.				
2	AA		ROW-NUMBER	N	3.0		UQ
			Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB		DATE-ENTERED	A	8.0		N
			The date of entry for a record. No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC		ACTIVE-INACTIVE	A	1.0		F
			Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD		CODE-VALUE	N	6.0		UQ
			Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE		EE-FILLER	A	1.0		F
			Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF		FF-FILLER	A	1.0		F
			Field's sole purpose is to hold a column position due to				

FILE...: TBORGSTA

TYPE...: USER VIEW

FILE-NR: 34

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
TBGET retrieval mechanism. Contains no data.		
2 GG SUBSET-INDICATOR	A	1.0
Legacy field with no use in current version of TDMS. For- mer use is unknown.		
2 HH LONG-TEXT	A	32.0 N
Full description of organ status. Downloaded to LDAS and appears on pathology individual animal report.		
2 II DATA-COLLECTION	A	1.0 DE
Legacy field with no use in current version of TDMS. For- mer use is unknown.		
2 JJ REPORTING	A	1.0
Symbol used in matrix-oriented pathology reports to indicate statuses. (Those reports have legends that describe these symbols. The legends are hard-coded. Any new statuses that might be added to the table would require modifica- of those reports.)		

TBPCT

Logical layout of the pathology code table. This layout is the basis on which NATURAL retrieves data for downloads and pathology reports.

FILE...: TBPCT

TYPE...: USER VIEW

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T	L	DB	NAME	F	LENG	S	DE
G	1	AA	PCT				
			Group embracing all fields in the pathology code table.				
2	AB		TOXTIME	A	8.0	N	DE
			The date upon which the record was added, or last update, in TDMS toxdate format. (Values are integers stored in an alphanumeric field.)				
2	AC		ACTIVITY-STATUS	A	1.0	F	
			Status of record. A for active; I for inactive. Not checked in PCT reports, but is checked by PAS validation reports.				
2	AD		TDMS-CODE	N	5.0	F	UQ
			Unique integer identifying a record on the table. The common practice, which is neither enforced by software nor standard operating procedures, is to assign codes sequentially, as follows: 1- 999 Morphologies 15000-15999 Qualifiers 25000-25999 Topologies				
2	AE		REORDER	N	5.0	N	
			Legacy field with no use in TDMS. Original use unknown.				
2	AF		TERMINOLOGY-PTR	N	5.0	N	
			Legacy field with no use in TDMS. Original use unknown. Contains no values.				

FILE...: TBPCT

TYPE...: USER VIEW

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T	L	DB NAME	F	LENG	S	DE
2	AG	INITIAL-ENTRY-DATE	N	5.0	N	DE
		Date on which record was first added to the table. Stored in TDMS toxdate format.				
2	AH	INITIAL-AUTHORIZATION	N	5.0	N	
		Integer code denoting the pathologist who authorized the addition of the record. The possible values are: 1001 - Dr. C. Montgomery 1002 - Dr. S. Eustis 1003 - Dr. R. Haley These codes are not operator IDs.				
2	AI	CHANGE-AUTHORIZATION	N	5.0	N	
		Code denoting the pathologist who authorizes a change to a record. (See INITIAL-AUTHORIZATION.)				
2	AJ	SUBTABLE-QUALIFIER	A	1.0	N	DE
		Letter indicating to which of three subtables a record belongs: M for morphologies, T for topographies, Q for qualifiers.				
2	AK	SYNONYM	N	5.0	N	
		Legacy field with no use in TDMS. Original use unknown. Contains no values.				
G 2	AL	ALL-HIERARCHY-LEVELS				
		Group embracing all LEVEL fields.				
3	AM	LEVEL1	N	5.0	N	DE
		Morphologies - Generic morphology type. Foreign key to TBPCT. Qualifiers - ID of the qual-				

FILE...: TBPCT

TYPE...: USER VIEW

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
<p>ifier group to which the qual- ifier belongs.</p>		
3 AN LEVEL2	N	5.0 N DE
<p>Morph. - Neoplastic flag. 4 for neoplastic; 5 for non-neoplastic.</p> <p>Qual. - Integer indicating sort sequence on pathology reports and the qualifier pull-down window on LDAS.</p>		
3 AO LEVEL3	N	5.0 N
<p>Morph. - Malignancy flag: 4 is malignant; 5 is be- nign; 6 is uncertain</p>		
3 AP LEVEL4	N	5.0 N
<p>Morph. - Metastatic flag: 4 is metastatic; 5 is non- metastatic; 6 is un- certain</p>		
3 AQ LEVEL5	N	5.0 N
<p>Legacy field. Original use unknown. Morphology terms inherited from NCTR contain the values 1 or 2.</p>		
3 AR LEVEL6	N	5.0 N
<p>Legacy field. Original use unknown. Contains no values.</p>		
3 AS LEVEL7	N	5.0 N
<p>Legacy field. Original use unknown. Contains no values.</p>		
3 AT LEVEL8	N	5.0 N
<p>Legacy field. Original use unknown. Contains no values.</p>		
G 2 AU ALL-VALIDATION-FLAGS		
<p>Group embracing all FLAG fields.</p>		
3 AV FLAG1	N	1.0 F
<p>Morph. & Top. - Gender speci- ficity. 1 is</p>		

FILE...: TBPCT

TYPE...: USER VIEW

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME

F LENG S DE

both genders,
2 is gender
dependant.

(Note: For topographies, only
organs are flagged as gender-
dependent. Sites are gender-
dependent by being located on
a gender-dependent organ.)

3 AW FLAG2 N 1.0 F

Morph. & Top. - Species speci-
ficity. 1 is
all species,
2 is gender
dependant.

(Note: For topographies, only
organs are flagged as species-
dependent. Sites are species-
dependent by being located on
a species-dependent organ.)

3 AX FLAG3 N 1.0 F

Morph. - Organ association
flag. 1 for all or-
gans; 2 for a speci-
fied list of organs;
3 for all but a
specified list of
organ. (The list of
organs for values 2
and 3 are found in
the multilple-entry
field ALL-FLAG3-
VALUES.)

3 AY FLAG4 N 1.0 F

Top. - Topography type: 1 is
a sytem, 2 is an or-
gan, 3 is a site
Qual. - Qualifier type: 1 is

FILE...: TBPCT

TYPE...: USER VIEW

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
qualifier group, 2 is a qualifier		
3 AZ FLAG5 Legacy field. Original use unknown. Contains no values.	N	1.0 F
3 BA FLAG6 Gross morphology flag. Gross morphologies are a legacy to TDMS, with no use in the cur- rent version.	N	1.0 F
3 BB FLAG7 Legacy field. Original use unknown. Contains no values.	N	1.0 F
3 BC FLAG8 Legacy field. Original use unknown. Contains no values.	N	1.0 F
2 BD REVERSE-PTR Legacy field. Original use unknown. Contains no values.	N	5.0 N
2 BE NCTR-CODE Legacy field with no use in current version of TDMS. Ap- pears to be a code that was associated with morphology terms inherited from NCTR.	A	3.0 N
2 BF SNOP-CODE Legacy field. Original use unknown.	A	4.0 N
G 2 BG SNOMED Group embracing SNOMED defini- tions. SNOMED appears to be an alternative coding system for pathology terms. Information- ational only. Only assigned to terms inherited from NCTR. Not entered for new terms.		
3 BH SNOMED-TYPE The single-letter identifying the SNOMED type. M for morph- ology; T for topography.	A	1.0 N
3 BI SNOMED-CODE Five-digit numeric code that	A	5.0 N

FILE...: TBPCT

TYPE...: USER VIEW

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T	L	DB NAME	F	LENG	S	DE
		uniquely identifies the term in the SNOMED system.				
2	BJ	ASSOCIATED-SNOMED-TERM	A	6.0	N	
		Legacy field. Original use unknown.				
2	BK	SYSTEMIC-FLAG	A	1.0	N	
		Applies to morphology terms only. Indicates whether a term is a systemic lesion.				
2	BL	GROSS-MICRO-FLAG	A	1.0	N	
		Flag indicating whether a morphology term is micro- scopically, or grossly, object- served. The possible values are: 0 - Microscopic only 1 - Gross only 2 - Both (Gross pathology is not part of the current version of TDMS.)				
2	BM	PAIRED-ORGAN-FLAG	A	1.0	N	
		Applies to organ topographies only. Indicates whether an organ is paired. The values are: 1 is not paired, 2 is paired.				
2	BN	FILLER1	A	13.0	N	
		Legacy field. Original use unknown. Contains no values.				
G 2	BO	SHORT-TEXT				
		Group which merely includes the short text field.				
3	BP	SHORT-TEXT-FIELD	A	14.0	N	DE
		Abbreviated description of the term. Used on Individual Animal Data report and LDAS				

FILE...: TBPCT

TYPE...: USER VIEW

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
pull-down screens.		
G 2 BQ LONG-TEXT		
Group which merely includes the long text field.		
3 BR LONG-TEXT-FIELD	A	60.0 N
Full description of the term.		
2 BS NUM-FLAG1-VALUES	N	5.0 N
Number of entries in ALL-FLAG1-VALUES multiple-entry field. Should have same value as MU count field for ALL-FLAG1-VALUES. Since ALL-FLAG1-VALUES can never have more than one entry, the value of this field should not be higher than one.		
2 BT NUM-FLAG2-VALUES	N	5.0 N
Number of entries in ALL-FLAG2-VALUES multiple-entry field. Should have same value as MU count field for ALL-FLAG2-VALUES. Since ALL-FLAG2-VALUES can never have more than one entry, the value of this field should not be higher than one.		
2 BU NUM-FLAG3-VALUES	N	5.0 N
Number of entries in ALL-FLAG3-VALUES multiple-entry field. Should have same value as MU count field for ALL-FLAG3-VALUES.		
2 BV NUM-FLAG4-VALUES	N	5.0 N
Number of entries in ALL-FLAG4-VALUES multiple-entry field. Should have same value as MU count field for ALL-FLAG4-VALUES.		
2 BW NUM-FLAG5-VALUES	N	5.0 N
Number of entries in ALL-FLAG5-VALUES multiple-entry field. Should have		

FILE...: TBPCT

TYPE...: USER VIEW

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
same value as MU count field for ALL-FLAG5-VALUES.		
2 BX NUM-FLAG6-VALUES Number of entries in ALL-FLAG6-VALUES multiple- entry field. Should have same value as MU count field for ALL-FLAG6-VALUES.	N	5.0 N
2 BY NUM-FLAG7-VALUES Number of entries in ALL-FLAG7-VALUES multiple- entry field. Should have same value as MU count field for ALL-FLAG7-VALUES.	N	5.0 N
2 BZ NUM-FLAG8-VALUES Number of entries in ALL-FLAG8-VALUES multiple- entry field. Should have same value as MU count field for ALL-FLAG8-VALUES.	N	5.0 N
2 CJ NUM-OF-FREE-TEXT-CHARS Number of entries in multiple- entry field FREE-TEXT-CHARS, which field is a legacy to TDMS. Number is always 0.	N	5.0 N
M 2 CA ALL-FLAG1-VALUES Gender-dependency for morph- ologies and topographies. If FLAG1 is set to 2, then the first entry in this multiple-entry field will be either M for F to indicate the gender upon which the term is dependent. There is never more than one entry in this MU field.	A	1.0 N
M 2 CB ALL-FLAG2-VALUES Indicate species dependency for morphologies and topo- graphies. If FLAG2 is set to 2, then the first entry in this MU field will be	N	5.0 N

the TBSTRAIN species code upon which the term is dependent. No more than one entry is used in this MU field.

M 2 CC ALL-FLAG3-VALUES N 5.0 N DE

Morph. - If FLAG3 is 2 or 3, the dependant organs. Foreign key to TBPCT.

Top. - If FLAG4 is 2, the system to which the organ belongs; if 3, then the organs for which the site is valid. Foreign key to TBPCT.

Qual. - If FLAG4 is 2, the qualifier groups to which the qualifier belongs. Foreign key to TBPCT.

M 2 CD ALL-FLAG4-VALUES N 5.0 N

Top. - If FLAG4 is 1, the organs included in the system; if 2, then the sites on the organ. Foreign key to TBPCT.

Qual. - If FLAG4 is 1, then the qualifiers included in the group. Foreign key to TBPCT.

M 2 CE ALL-FLAG5-VALUES N 5.0 N DE

Morph. - The IDs for the qualifier groups that are appropriate to the term.

Qual. - If FLAG4 is 2, the the IDs for the qualifier groups in which the qualifier is classified. (The LEVEL1 values for the terms in ALL-FLAG3-VALUES.)

M 2 CF ALL-FLAG6-VALUES N 5.0 N

Legacy field. Equivalent to

FILE...: TBPCT

TYPE...: USER VIEW

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
ALL-FLAG3-VALUES, but for gross terms. Gross pathology is not part of the current version of TDMS.		
M 2 CG ALL-FLAG7-VALUES	N	5.0 N
Legacy field. Equivalent to ALL-FLAG4-VALUES, but for gross terms. Gross pathology is not part of the current version of TDMS.		
M 2 CH ALL-FLAG8-VALUES	N	5.0 N
Legacy field. Equivalent to ALL-FLAG5-VALUES, but for gross terms. Gross pathology is not part of the current version of TDMS.		
M 2 CI FREE-TEXT-CHARS	A	1.0 N
Legacy field. Original use unknown. Contains no values.		

TBQUALGRP

Logical view of the in-life qualifier group table. It contains the possible groups by which current in-life qualifiers are classified.

FILE...: TBQUALGRP

TYPE...: USER VIEW

FILE-NR: 31

PRIMARY SEQUENCE FIELD: EIQUAL-GROUP

T	L	DB NAME	F	LENG	S	DE
G	1	HA INLIFE-QUALIFIER-GROUP-TABLE				
		Group embracing all fields in the table.				
2	HB	EIQUAL-GROUP	A	5.0	N	UQ
		Unique, sequentially assigned, number that identifies an in-life qualifier group.				
2	HC	EIQUAL-GROUP-SHORT-TEXT	A	14.0	N	
		Description of in-life qualifiers-fier group.				
2	HD	EIQUAL-GROUP-LONG-TEXT	A	32.0	N	
		Description of in-life qualifiers-fier group.				

TBSTAINC

Logical layout of stain table. The table contains stains that could be used for tissue section cuts. Stains are a legacy to TDMS, with no use in the current version. PEX software does not require any stain codes to be associated with a test.

FILE...: TBSTAINC

TYPE...: USER VIEW

FILE-NR: 15

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
G	1	ST	STAIN-TABLE				
			Embraces all fields in a stain record.				
2	AA	ROW-NUMBER		N	3.0	UQ	
			Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB	DATE-ENTERED		A	8.0		
			The date of entry for a record. No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC	ACTIVE-INACTIVE		A	1.0	F	
			Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD	CODE-VALUE		N	6.0	UQ	
			Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE	EE-FILLER		A	1.0	F	
			Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF	FF-FILLER		A	1.0	F	
			Field's sole purpose is to				

FILE...: TBSTAINC

TYPE...: USER VIEW

FILE-NR: 15

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
-------------	---	-----------

hold a column position due to
TBGET retrieval mechanism.
Contains no data.

2 GG LONG-TEXT	A	30.0 N
----------------	---	--------

Full description of stain.

2 HH SHORT-TEXT	A	14.0 N
-----------------	---	--------

Abbreviated description of the
stain.

TBSTRAIN

Logical layout of strain table. Table contains species and species/strains that are used in TDMS. The three current species are mice, rats, and hamsters. Strain information appears on most EIS/PEIS reports.

FILE...: TBSTRAIN

TYPE...: USER VIEW

FILE-NR: 21

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	ST	STRAIN-TABLE				
			Group embracing all fields in the table.				
2	AA	ROW-NUMBER		N	3.0	UQ	
			Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB	DATE-ENTERED		A	8.0		
			The date of entry for a record. No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC	ACTIVE-INACTIVE		A	1.0	F	
			Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD	CODE-VALUE		N	6.0	UQ	
			Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE	EE-FILLER		A	1.0	F	
			Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF	FF-FILLER		A	1.0	F	
			Field's sole purpose is to				

FILE...: TBSTRAIN

TYPE...: USER VIEW

FILE-NR: 21

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
		hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	GG	SPECIES-CODE	N	6.0	N	
		The TBSTRAIN code for the spe- cies to which a strain belongs.				
2	HH	SUBSET-INDICATOR	N	6.0	N	
		Indicates whether a record is a species or strain. 1 is a species; 2 is a strain.				
2	II	SHORT-TEXT	A	12.0	N	
		Abbreviated description of species/strain. Used in EIS/ PEIS reports.				
2	JJ	BREEDING-CODE	A	1.0	N	
		Legacy field with no meaning in TDMS. Contains no values.				
2	KK	STR-FEM-PARENT	N	6.0	N	
		Legacy field with no meaning in TDMS. Contains no values.				
2	LL	STR-MALE-PARENT	N	6.0	N	
		Legacy field with no meaning in TDMS. Contains no values.				
2	MM	LONG-TEXT	A	38.0	N	
		Full description of species/ strain.				
2	NN	NN-FILLER	A	1.0	N	
		Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				

TBSUBSTR

Logical layout of substrain table. Substrains are a legacy to TDMS, with no use in the current version. PAS software requires a substrain code for in-life test protocol, so the codes for "Not Applicable" or "Not Specified" are entered.

FILE...: TBSUBSTR

TYPE...: USER VIEW

FILE-NR: 17

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
G	1	ST	SUBSTRAIN-TABLE				
			Group embracing all fields				
			in table.				
2	AA	ROW-NUMBER		N	3.0		UQ
			Legacy field with no meaning				
			in TDMS. It was probably the				
			primary key when the file was				
			in VSAM.				
2	BB	DATE-ENTERED		A	8.0		
			The date of entry for a record				
			No enforced date format, but				
			YYYYMMDD is preferred and TDMS				
			tox date is discouraged. In-				
			formational only; not used in				
			reports.				
2	CC	ACTIVE-INACTIVE		A	1.0		F
			Status of record. A for active				
			and I for inactive. Information-				
			tional only; not enforced by				
			PAS.				
2	DD	CODE-VALUE		N	6.0		UQ
			Unique, sequential code to				
			identify record within table.				
			Primary key for record re-				
			trieval.				
2	EE	EE-FILLER		A	1.0		F
			Field's sole purpose is to				
			hold a column position due to				
			TBGET retrieval mechanism.				
			Contains no data.				
2	FF	FF-FILLER		A	1.0		F
			Field's sole purpose is to				

FILE...: TBSUBSTR

TYPE...: USER VIEW

FILE-NR: 17

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
-------------	---	-----------

hold a column position due to
TBGET retrieval mechanism.
Contains no data.

2 GG SHORT-TEXT	A	14.0 N
-----------------	---	--------

Abbreviated description of the
substrain.

2 HH LONG-TEXT	A	32.0 N
----------------	---	--------

Full description of the sub-
strain.

TBTESTYP

Logical layout of test type table. Table contains possible test types by which a test may be classified. Test type information appears on most EIS/PEIS reports, usually on the splash page.

FILE...: TBTESTYP

TYPE...: USER VIEW

FILE-NR: 28

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
G	1	TT	TEST-TYPE-TABLE				
			Group embracing all fields in the test type table.				
2	AA	ROW-NUMBER		N	3.0	UQ	
			Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.				
2	BB	DATE-ENTERED		A	8.0		
			The date of entry for a record. No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.				
2	CC	ACTIVE-INACTIVE		A	1.0	F	
			Status of record. A for active and I for inactive. Informational only; not enforced by PAS.				
2	DD	CODE-VALUE		N	6.0	UQ	
			Unique, sequential code to identify record within table. Primary key for record retrieval.				
2	EE	EE-FILLER		A	1.0	F	
			Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.				
2	FF	FF-FILLER		A	1.0	F	
			Field's sole purpose is to hold a column position due to				

FILE...: TBTESTYP

TYPE...: USER VIEW

FILE-NR: 28

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
TBGET retrieval mechanism. Contains no data.		
2 GG CHARACTER-CODE	A	2.0 N
Code value right-justified and padded with 0's for bar-code use. Legacy field with no use in current version of TDMS.		
2 HH SHORT-TEXT	A	16.0 N
Abbreviated description of the test type. Displayed in EIS/PEIS reports.		
2 II LONG-TEXT	A	32.0 N
Full description of the test type. Used in EIS/PEIS reports and downloaded to LDAS.		
2 JJ CURRENT-PERIOD-LENGTH	P	3.0 N
Defines the length, in days, of the current period for a test type. Used in some EIS reports that perform current-period reporting.		

TBTRTROLE

Logical layout of treatment role table. Contains the possible roles that a treatment group may have on a test. Some older TDMS software still has these roles hard-coded.

FILE...: TBTRTROLE

TYPE...: USER VIEW

FILE-NR: 31

PRIMARY SEQUENCE FIELD: TREAT-ROLE

T	L	DB NAME	F	LENG	S	DE

G	1	BA TREATMENT-ROLE-TABLE				
		Group embracing all fields in the treatment role table.				
2	BB	TREAT-ROLE	A	1.0	N	UQ
		Unique alphabetic code for a treatment role. Primary key for table retrieval.				
2	BC	ROLE-TEXT	A	32.0	N	
		Full description for treatment role.				

TBUNITS

Logical layout of units table. This table lists quantitative units of measure. This table is a legacy to TDMS. The only valid unit of measure in TDMS is decigrams, which unit is hard-coded into TDMS software.

FILE...: TBUNITS

TYPE...: USER VIEW

FILE-NR: 20

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

- - - - -

G 1 UT UNITS-TABLE

Group embracing all fields in table.

2	AA	ROW-NUMBER	N	3.0	UQ
---	----	------------	---	-----	----

Legacy field with no meaning in TDMS. It was probably the primary key when the file was in VSAM.

2	BB	DATE-ENTERED	A	8.0	
---	----	--------------	---	-----	--

The date of entry for a record
No enforced date format, but YYYYMMDD is preferred and TDMS tox date is discouraged. Informational only; not used in reports.

2	CC	ACTIVE-INACTIVE	A	1.0	F
---	----	-----------------	---	-----	---

Status of record. A for active and I for inactive. Informational only; not enforced by PAS.

2	DD	CODE-VALUE	N	6.0	UQ
---	----	------------	---	-----	----

Unique, sequential code to identify record within table. Primary key for record retrieval.

2	EE	EE-FILLER	A	1.0	F
---	----	-----------	---	-----	---

Field's sole purpose is to hold a column position due to TBGET retrieval mechanism. Contains no data.

2	FF	FF-FILLER	A	1.0	F
---	----	-----------	---	-----	---

Field's sole purpose is to

FILE...: TBUNITS

TYPE...: USER VIEW

FILE-NR: 20

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
-------------	---	-----------

hold a column position due to
TBGET retrieval mechanism.
Contains no data.

2 GG SHORT-TEXT	A	8.0 N
-----------------	---	-------

Abbreviated description of the
unit of measure.

2 HH LONG-TEXT	A	50.0 N
----------------	---	--------

Full description of the unit
of measure.

TBWEIGHTCURVE

Logical view of animal weight growth table. Contains weight curve upper and lower limits as a function of species, sex, and week-on-test.

FILE...: TBWEIGHTCURVE

TYPE...: USER VIEW

FILE-NR: 31

PRIMARY SEQUENCE FIELD: WEIGHT-CURVE-START-WEEK

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-

G 1 GA WEIGHT-CURVE-TABLE

Group embracing all fields in the body weight growth curve table.

2 GB WEIGHT-CURVE-START-WEEK P 9.0 N UQ

This field is a numeric string and comprises two sub-fields. The first four digits denote the species-sex; the second five, the starting test-on-week. (The ending test-on-week would be one less than the next highest starting test-on-week for the species-sex.) The possible species-sex indicators are:

0001 - Male Rats
0002 - Female Rats
0003 - Male Mice
0004 - Female Mice

2 GD MIN-WEIGHT P 5.0 N

The minimum weekly change in body weight, expressed as percentage.

2 GE MAX-WEIGHT P 5.0 N

Maximum weekly body weight change, expressed as a percentage.

TBWGTOBJ

Logical layout of Weight Object table. This table contains objects which may be weighed during the in-life portion of a test, and for which there are weight ranges. Codes from this table are entered in PAS weight items segment and included in in-life downloads.

FILE...: TBWGTOBJ

TYPE...: USER VIEW

FILE-NR: 19

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
G	1	WT	WEIGHT-OBJECT-TABLE				
			Group embracing all fields				
			in table.				
2	AA	ROW-NUMBER		N	3.0		UQ
			Legacy field with no meaning				
			in TDMS. It was probably the				
			primary key when the file was				
			in VSAM.				
2	BB	DATE-ENTERED		A	8.0		
			The date of entry for a record				
			No enforced date format, but				
			YYYYMMDD is preferred and TDMS				
			tox date is discouraged. In-				
			formational only; not used in				
			reports.				
2	CC	ACTIVE-INACTIVE		A	1.0		F
			Status of record. A for active				
			and I for inactive. Information-				
			tional only; not enforced by				
			PAS.				
2	DD	CODE-VALUE		N	6.0		UQ
			Unique, sequential code to				
			identify record within table.				
			Primary key for record re-				
			trieval.				
2	EE	EE-FILLER		A	1.0		F
			Field's sole purpose is to				
			hold a column position due to				
			TBGET retrieval mechanism.				
			Contains no data.				
2	FF	FF-FILLER		A	1.0		F
			Field's sole purpose is to				

FILE...: TBWGTOBJ

TYPE...: USER VIEW

FILE-NR: 19

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
-------------	---	-----------

hold a column position due to
TBGET retrieval mechanism.
Contains no data.

2 GG SHORT-TEXT	A	8.0 N
-----------------	---	-------

Abbreviated description of
weight object.

2 HH LONG-TEXT	A	32.0 N
----------------	---	--------

Full description of weight object-
ject. First 20 characters in-
cluded in in-life downloads.

TBWGTSTATUS

Logical view of weight status table. Contains valid statuses for objects whose scheduled weight action is being bypassed. Downloaded to LDAS as part of test protocol. Displayed in all reports that contain weight detail information.

FILE...: TBWGTSTATUS

TYPE...: USER VIEW

FILE-NR: 31

PRIMARY SEQUENCE FIELD: STATUS-NUM

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
G	1	AA	WEIGHT-STATUS-TABLE				
			Group embracing all fields in the weight status table.				
	2	AB	STATUS-NUM	P	5.0	N	UQ
			Unique, sequential, identification number for weight status used as primary index for table retrieval.				
	2	AC	STATUS-TEXT	A	30.0	N	
			Full description of weight status.				
	2	AD	STATUS-SHORT-TEXT	A	8.0	N	
			Abbreviated description of weight status. Displayed on EIS/PEIS detail reports.				
	2	AE	STATUS-FLAG	A	1.0	F	
			Flag that indicates to which type of weight action the status applies. A for animal; C for container; R for removal				

TB_OPERATOR_ID

User view for operator table. Corresponds to the LDAS "operator" table. Contains all valid operators of TDMS. Each operator is assigned a unique, sequential number within the facility with which they are associated. All users of the TDMS menu have the first eight characters of their user account entered in the STARUSER field, which the menu uses to validate users invoking the menu.

FILE...: TB_OPERATOR_ID

TYPE...: USER VIEW

FILE-NR: 14

PRIMARY SEQUENCE FIELD: OPERATOR-ID

T L DB NAME	F	LENG S DE
1 AA OPERATOR-ID Unique identifier for records in the table. The first 5 characters are the associated facility; the next four, a unique, usually sequential, number within the facility. The primary key for the file.	N	9.0 UQ
1 AC OP-ID-PASSWORD Defunct field. Used by MOD-COMP data collection system, which has been replaced by LDAS.	A	8.0 N
1 AE OP-ID-FAC-ACCESS The facility with which the user is associated. Should have the same value as the first 5 characters of the OPERATOR-ID field, and is, hence, redundant. Foreign key to TBFACILITY.	N	5.0 N
1 AF OP-ID-AGENCY The agency with which the operator is associated. This field is superfluous as NIEHS is the only valid TDMS agency. Foreign key to TBAGENCY.	N	5.0 N
1 AG OP-ID-STARUSER The first eight characters of the operator's VAX account on	A	8.0 N DE

FILE...: TB_OPERATOR_ID

TYPE...: USER VIEW

FILE-NR: 14

PRIMARY SEQUENCE FIELD: OPERATOR-ID

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

the NIEHS cluster. This value is used by the TDMS main menu to validate a user attempting access to the menu system.

P 1 AH OP-ID-USER-INFO

Repeating group for each person who's been assigned the identifier code. LDAS doesn't support reassignment of IDs so this group should only have one group entry. (MODCOMP supported reassignment.)

2	AI	OP-ID-USER-START-DATE	A	8.0	N	
---	----	-----------------------	---	-----	---	--

The date on which the operator was assigned. The date format is MMDDYY.

2	AJ	OP-ID-USER-STOP-DATE	A	8.0	N	
---	----	----------------------	---	-----	---	--

The date on which the user was removed as a TDMS operator. The date format is MMDDYY.

2	AK	OP-ID-USER-NAME	A	32.0	N	DE
---	----	-----------------	---	------	---	----

The name of the operator. The usual format is last name and first initial, separated by a comma.

4.6 ADABAS Physical Definitions

ACCT-SITES-PEIS

Physical layout of file 55. Supports the PAS_ACCT_SITE view.
PASGET uses this layout to issue direct calls for EIS/PEIS reports.

FILE...: ACCT-SITES-PEIS

TYPE...: ADABAS

FILE-NR: 55

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
G 1	TA	PAS-ACCT-SITE	Not referenced by PASGET.				
G 2	TB	PAS-KEY	Not referenced by PASGET.				
3	TC	DATA-BASE-KEY	Used by PASGET to retrieve data from the file.	A	9.0		DE
2	TE	DATA-DATE	Not referenced by PASGET.	A	8.0		N
2	TF	DATA-TIME	Not referenced by PASGET.	A	6.0		N
2	TG	DATA-FACILITY-NUMBER	Not referenced by PASGET.	B	2.0		N
2	TH	OPERATOR-ID	Not referenced by PASGET.	B	2.0		N
G 2	TK	CORRECTION	Not referenced by PASGET.				
3	TL	ORG-OR-COR	Not referenced by PASGET.	A	1.0		F
3	TM	REASON	Not referenced by PASGET.	A	78.0		N
3	TN	COR-DATE	Not referenced by PASGET.	A	8.0		N
3	TO	COR-TIME	Not referenced by PASGET.	A	6.0		N
G 1	AA	ACCT-SITE-REQ	Not referenced by PASGET.				
2	BB	ORGAN-CODE	Retrieved by PASGET for PAS validation reports and EIS/PEIS reports.	B	2.0		

FILE...: ACCT-SITES-PEIS

TYPE...: ADABAS

FILE-NR: 55

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

2	CC	SITE-GRP-CT	B	2.0	N	
---	----	-------------	---	-----	---	--

Retrieved by PASGET for PAS
validation reports and EIS/
PEIS reports.

M 2	DD	SITE-CODE	B	2.0	N	
-----	----	-----------	---	-----	---	--

Retrieved by PASGET for PAS
validation reports and EIS/
PEIS reports.

1	SD	SUPER-KEY	A	11.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-

		DATA-BASE-KEY	1	9		
--	--	---------------	---	---	--	--

		ORGAN-CODE	1	2		
--	--	------------	---	---	--	--

Used by PASGET to select data
from the file.

ACTION-TABLE

Physical layout of action table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBACTONC. TBGET uses that layout description to form direct calls to retrieve data from the table for PAS validation reports. This layout corresponds to the view TBACTONC.

FILE...: ACTION-TABLE

TYPE...: ADABAS

FILE-NR: 29

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	AT ACTION-TABLE				
2	AA	ROW-NUMBER	B	2.0		UQ
		Refer to TBACTONC. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBACTONC. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0		F
		Refer to TBACTONC. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0		UQ
		Refer to TBACTONC. Key used by TBGET to retrieve data from the table.				
2	EE	EE-FILLER	A	1.0		F
		Refer to TBACTONC. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0		F
		Refer to TBACTONC. Not referenced by TBGET.				
2	GG	LONG-TEXT	A	32.0		N
		Refer to TBACTONC. Retrived by TBGET for PAS validation reports.				
2	HH	TYPE-ACTION	A	1.0		N
		Refer to TBACTONC. Not referenced by TBGET.				

AGENCY-TABLE

Physical layout of agency table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBAGENCY. TBGET uses that layout description to form direct calls to retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBAGENCY.

FILE...: AGENCY-TABLE

TYPE...: ADABAS

FILE-NR: 22

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G 1	AT	AGENCY-TABLE				
2	AA	ROW-NUMBER	B	2.0		UQ
		Refer to TBAGENCY. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBAGENCY. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0		F
		Refer to TBAGENCY. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0		UQ
		Refer to TBAGENCY. Key used by TBGET to retrieve data from this table.				
2	EE	EE-FILLER	A	1.0		F
		Refer to TBAGENCY. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0		F
		Refer to TBAGENCY. Not referenced by TBGET.				
2	GG	CHARACTER-CODE	A	2.0		N
		Refer to TBAGENCY. Not referenced by TBGET.				
2	HH	LONG-TEXT	A	32.0		N
		Refer to TBAGENCY. TBGET retrieves this field for PAS validation reports.				
2	II	REPORT-HEADINGS	A	50.0		N
		Refer to TBAGENCY. Not referenced by TBGET.				
2	JJ	SHORT-TEXT	A	12.0		N

FILE...: AGENCY-TABLE
TYPE...: ADABAS
FILE-NR: 22
PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F L E N G S D E
-------------	-----------------

Refer to TBAGENCY. Used by
TBGET to retrieve data from
this table for EIS/PEIS re-
ports.

ANIMAL-CLASS-QUAL-TABLE

Physical layout of animal class qualifier table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBANCLQF. TBGET uses that layout description to form direct calls to retrieve data from the table for PAS validation reports. This file corresponds to the view TBANCLQF.

FILE...: ANIMAL-CLASS-QUAL-TABLE
 TYPE...: ADABAS
 FILE-NR: 23
 PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	CT ANIMAL-CLASS-QUALIFIER-TABLE				
2	AA	ROW-NUMBER	B	2.0	UQ	
		Refer to TBANCLQF. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBANCLQF. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Refer to TBANCLQF. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0	UQ	
		Refer to TBANCLQF. Key by which TBGET retrieves records from the table.				
2	EE	EE-FILLER	A	1.0	F	
		Refer to TBANCLQF. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0	F	
		Refer to TBANCLQF. Not referenced by TBGET.				
2	GG	LONG-TEXT	A	32.0	N	
		Refer to TBANCLQF. Retrieved by TBGET for PAS validation reports.				

ANIMAL-DATA-EIS

Physical layout of file 106. Supports the TD-ANIMAL-DATA view.
No TDMS software accesses this file using direct calls.

FILE...: ANIMAL-DATA-EIS

TYPE...: ADABAS

FILE-NR: 106

PRIMARY SEQUENCE FIELD: EXP-TEST-ANIMAL-DATE

T	L	DB	NAME	F	LENG	S	DE
1	AA		EXP-TEST-ANIMAL-DATE	A	21.0		DE
1	AL		BALANCE-NUMBER	U	5.0		N
1	AB		WEIGHT-TIME	U	5.0		N
1	AC		ANIMAL-WEIGHT	U	10.0		N
1	AD		WEIGHT-OPERATOR-ID	U	5.0		N
1	AE		WEIGHT-ERROR-FLAG	A	1.0		F
1	AF		ACTION-OVERRIDE-FLAG-WGT	A	1.0		F
1	AG		MANUAL-ENTRY-FLAG-WGT	A	1.0		F
1	AH		REMOVAL-ACTION-FLAG-WGT	A	1.0		F
1	AI		BACKWARD-FORWARD-FLAG-WGT	A	1.0		F
1	AJ		EXPECTED-SELECTED-FLAG-WGT	A	1.0		F
1	AK		WEIGHT-OVERRIDE	A	1.0		F
1	BB		OBSERVATION-TIME	U	5.0		N
1	AM		OBSERVATION-OPERATOR-ID	U	5.0		N
1	AN		OBSERVATION-ERROR-FLAG	A	1.0		F
1	AO		ACTION-OVERRIDE-FLAG-OBS	A	1.0		F
1	AP		MANUAL-ENTRY-FLAG-OBS	A	1.0		F
1	AQ		REMOVAL-ACTION-FLAG-OBS	A	1.0		F
1	AR		BACKWARD-FORWARD-FLAG-OBS	A	1.0		F
1	AS		EXPECTED-SELECTED-FLAG-OBS	A	1.0		F
1	AT		NOTE-OPERATOR-ID	U	5.0		N
1	BC		NOTE-TIME	U	5.0		N
1	AU		NOTE-ERROR-FLAG	A	1.0		F
1	AV		ACTION-NUMBER	U	5.0		N
1	AW		ANIMAL-NOTE	A	80.0		N
P 1	DA		OBSERVATION-DATA				
2	DB		OBSERVATION	U	5.0		N
2	DC		SITE	U	5.0		N
2	DD		EIQUAL1	U	5.0		N
2	DE		EIQUAL2	U	5.0		N
1	DF		MANUAL-ENTRY-FLAG-NOTE	A	1.0		F
1	DG		WEIGHT-STATUS	P	5.0		N

ANIMAL-EIS

Physical layout of file 103. Supports the TD-ANIMAL view. No current TDMS software accesses this file using direct calls.

FILE...: ANIMAL-EIS

TYPE...: ADABAS

FILE-NR: 103

PRIMARY SEQUENCE FIELD: EXP-TEST-ANIMAL

T	L	DB	NAME	F	LENG	S	DE
1	AA	EXP-TEST-ANIMAL	A	16.0	UQ		
1	AB	ORIGINAL-CAGE-NUMBER	U	4.0			
1	CA	CARCASS-IDENTIFICATION	A	14.0	N	DE	
1	CB	BALANCE-NUMBER	U	5.0			
1	CC	REASON-FOR-REMOVAL	U	5.0	N		
1	CD	HOSPITAL-CAGE	U	4.0	N		
1	AC	ANIMAL-IDENTIFICATION	A	1.0	F		
1	AD	ANIMAL-CLASS	U	5.0			
1	AE	SEX	A	1.0	F		
1	AF	SPECIES	U	5.0			
1	AG	STRAIN	U	5.0			
1	AH	SUBSTRAIN	U	5.0			
1	AI	ANIMAL-DATE	U	5.0			
1	AJ	ANIMAL-TIME	U	5.0			
1	AK	ANIMAL-OPERATOR-ID	U	5.0			
1	AL	ANIMAL-ERROR-FLAG	A	1.0	F		
1	AM	ACTION-OVERRIDE-FLAG-CID	A	1.0	F		
1	AN	MANUAL-ENTRY-FLAG-CID	A	1.0	F		
1	AO	REMOVAL-ACTION-FLAG-CID	A	1.0	F		
1	AP	BACKWARD-FORWARD-FLAG-CID	A	1.0	F		
1	AQ	EXPECTED-SELECTED-FLAG-CID	A	1.0	F		
1	AS	ACTION-OVERRIDE-FLAG-REM	A	1.0	F		
1	AT	MANUAL-ENTRY-FLAG-REM	A	1.0	F		
1	AU	REMOVAL-ACTION-FLAG-REM	A	1.0	F		
1	AV	BACKWARD-FORWARD-FLAG-REM	A	1.0	F		
1	AW	EXPECTED-SELECTED-FLAG-REM	A	1.0	F		
1	AX	REMOVAL-DATE	U	5.0	N		
1	CF	REMOVAL-TIME	U	5.0	N		
1	CE	REMOVAL-OPERATOR-ID	U	5.0	N		
1	CG	REMOVAL-ERROR-FLAG	A	1.0	F		
1	CM	REMOVED-DATE	U	5.0	N		
1	AY	DAYS-ON-EXPERIMENT	U	5.0			
1	AZ	ACTION-OVERRIDE-FLAG-WGT	A	1.0	F		
1	BA	MANUAL-ENTRY-FLAG-WGT	A	1.0	F		

FILE...: ANIMAL-EIS

TYPE...: ADABAS

FILE-NR: 103

PRIMARY SEQUENCE FIELD: EXP-TEST-ANIMAL

T	L	DB NAME	F	LENG	S	DE
1	BB	REMOVAL-ACTION-FLAG-WGT	A	1.0	F	
1	BC	BACKWARD-FORWARD-FLAG-WGT	A	1.0	F	
1	BD	EXPECTED-SELECTED-FLAG-WGT	A	1.0	F	
1	BE	ANIMAL-REM-WEIGHT	U	10.0	N	
P 1	CH	OLD-OBSERVATION-GROUP				
2	CI	OLD-OBSERVATION	U	5.0	N	
2	CJ	OLD-SITE-OR-QUALIFIER	U	5.0	N	
2	CK	OLD-SIZE-OR-COLOR	U	5.0	N	
1	BF	WEIGHT-OVERRIDE	A	1.0	F	
1	BI	ACTION-OVERRIDE-FLAG-OBS	A	1.0	F	
1	BJ	MANUAL-ENTRY-FLAG-OBS	A	1.0	F	
1	BK	REMOVAL-ACTION-FLAG-OBS	A	1.0	F	
1	BL	BACKWARD-FORWARD-FLAG-OBS	A	1.0	F	
1	BM	EXPECTED-SELECTED-FLAG-OBS	A	1.0	F	
M 1	BN	DISPOSITION-CODE	U	5.0	N	
1	BP	TRANSFER-DATE	U	5.0	N	
1	BQ	TRANSFER-TIME	U	5.0	N	
1	BR	TRANSFER-OPERATOR-ID	U	5.0	N	
1	BS	TRANSFER-ERROR-FLAG	A	1.0	F	
1	BT	PATHOLOGIST-IDENTIFICATION	U	5.0		
1	BU	HISTOLOGY-NUMBER-DATE	U	5.0	N	
1	BV	HISTOLOGY-NUMBER-TIME	U	5.0	N	
1	BW	HISTOLOGY-NUMBER-OPERATOR-ID	U	5.0	N	
1	BY	HISTOLOGY-NUMBER-ERROR-FLAG	A	1.0	F	
1	BX	HISTOLOGY-NUMBER	A	16.0	N	
1	BZ	PATHOLOGIST-TYPE	A	1.0	F	
1	EA	REMOVAL-WEIGHT-STATUS	P	5.0	F	
P 1	FA	OBSERVATION-DATA				
2	FB	OBSERVATION	P	3.0	N	
2	FC	SITE	P	3.0	N	
2	FD	EIQUAL1	P	3.0	N	
2	FE	EIQUAL2	P	3.0	N	
1	AR	NO-EXAM-PATHOLOGIST	P	5.0	N	
1	A0	NO-EXAM-OPERATOR-ID	P	5.0	N	
1	A1	NO-EXAM-DATE	P	5.0	N	DE
1	A2	NO-EXAM-TIME	P	5.0	N	
1	A3	UNCERTAIN-PATHOLOGIST	P	5.0	N	

FILE...: ANIMAL-EIS

TYPE...: ADABAS

FILE-NR: 103

PRIMARY SEQUENCE FIELD: EXP-TEST-ANIMAL

T	L	DB	NAME	F	LENG	S	DE
1	A4		UNCERTAIN-OPERATOR-ID	P	5.0	N	
1	A5		UNCERTAIN-DATE	P	5.0	N	
1	A6		UNCERTAIN-TIME	P	5.0	N	
1	DP		EXP-TEST-ANIMAL-SP	A	12.0		SP
			SOURCE FIELD(S) --- -START- --END-				
			EXP-TEST-ANIMAL		1		7
			EXP-TEST-ANIMAL		12		16
1	CN		EXP-TEST-HOSPITAL-CAGE	A	11.0		SP
			SOURCE FIELD(S) --- -START- --END-				
			EXP-TEST-ANIMAL		1		7
			HOSPITAL-CAGE		1		4

ANIMAL-ID-EIS

Physical layout of file 49. Supports the PAS_ANIMAL_ID view. This is the layout used by which PASGET issues direct calls for PAS validation reports.

FILE...: ANIMAL-ID-EIS

TYPE...: ADABAS

FILE-NR: 49

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA		PAS-ANIMAL-ID				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from the file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		STR-ANIMAL-ID-SEG				
			Not referenced by PASGET.				
2	BB		ID-FLAG	A	1.0		F
			Retrieved by PASGET for PAS validation reports.				
2	CC		CID-FLAG	A	1.0		F
			Retrieved by PASGET for PAS validation reports.				

FILE...: ANIMAL-ID-EIS

TYPE...: ADABAS

FILE-NR: 49

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
2 DD ID-CT Retrieved by PASGET for PAS validation reports.	B	2.0 N
P 1 EE ANY-ID-INFO Retrieved by PASGET for PAS validation reports.		
2 FF ANIMAL-ID Retrieved by PASGET for PAS validation reports.	A	1.0 F
2 GG DESCRIPTION Retrieved by PASGET for PAS validation reports.	A	16.0 N

ANIMAL-NOTE-PEIS

Physical layout of file 110. Supports TD-MICRO-ANIMAL-NOTE view. No current TDMS software accesses this file using direct calls.

FILE...: ANIMAL-NOTE-PEIS

TYPE...: ADABAS

FILE-NR: 110

PRIMARY SEQUENCE FIELD: EXP-TEST-ANIMAL-DATE

T	L	DB	NAME	F	LENG	S	DE
1	AA	EXP-TEST-ANIMAL-DATE	A	22.0	DE		
1	AB	MICRO-ANIMAL-NOTE-TIME	U	5.0			
1	AC	MICRO-ANIMAL-NOTE-OPERATOR-ID	U	5.0			
1	AD	MICRO-ANIMAL-NOTE-ERROR-FLAG	A	1.0	F		
1	AF	MICRO-ANIMAL-NOTE-PATHOLOGIST-ID	U	5.0			
M 1	AE	MICRO-ANIMAL-NOTE	A	80.0	N		
1	AG	MICRO-ANIMAL-NOTE-COD-FLAG	A	1.0	F		

ANIMAL-SUPPLIER-TABLE

Physical layout of animal supplier table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBANMLSU. TBGET uses that layout description to form direct calls to retrieve data from the table for PAS validation reports. This layout corresponds to the view TBANMLSUPPLIER.

FILE...: ANIMAL-SUPPLIER-TABLE
 TYPE...: ADABAS
 FILE-NR: 32
 PRIMARY SEQUENCE FIELD: ANML-SUPPLR

T	L	DB NAME	F	LENG	S	DE
G	1	FA ANIMAL-SUPPLIER-TABLE				
		Refer to TBANMLSUPPLIER. Not referenced by TBGET.				
2	FB	ANML-SUPPLR	P	3.0	F	UQ
		Refer to TBANMLSUPPLIER. Key used by TBGET to retrieve data from table.				
2	FC	SUPPLR-TEXT	A	64.0	N	
		Refer to TBANMLSUPPLIER. Retrieved by TBGET for PAS validation-reports.				

BALANCE

Physical layout of file 104. Supports the TD-BALANCE view. No current TDMS software accesses this file using direct calls.

FILE...: BALANCE

TYPE...: ADABAS

FILE-NR: 104

PRIMARY SEQUENCE FIELD: BALANCE-NUMBER

T	L	DB	NAME	F	LENG	S	DE
1	AA	BALANCE-NUMBER		N	5.0	F	DE
1	AB	AGENCY-NUMBER		N	2.0	F	
1	AC	WEIGHT-UNIT-CODE		N	5.0	F	

BALANCE-CALIBRATION-EIS

Physical layout of file 102. Supports the TD-BALANCE-CALIBRATION view. No current TDMS software accesses this file using direct calls.

FILE...: BALANCE-CALIBRATION-EIS

TYPE...: ADABAS

FILE-NR: 102

PRIMARY SEQUENCE FIELD: EXP-TEST-DATE

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
1	AA		EXP-TEST-DATE	U	12.0		DE
1	AB		BALANCE-NUMBER	U	5.0		
P 1	AC		CALIBRATION-WEIGHT-DATA				
2	AD		CALIBRATION-WEIGHT	U	5.0		N
2	AE		WEIGHT-OVERRIDE	A	1.0		N
1	AF		TIME	U	5.0		
1	AG		OPERATOR-ID	U	5.0		
1	AH		ERROR-FLAG	A	1.0		F
1	AI		MANUAL-ENTRY-FLAG	A	1.0		F

CAGE-CENSUS-EIS

Physical layout of file 44. Supports the PAS_CENSUS view. This is the layout by which PASGET selects data for PAS validation reports .

FILE...: CAGE-CENSUS-EIS

TYPE...: ADABAS

FILE-NR: 44

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
1	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from this file.				
1	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
1	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
1	TG		DATA-FACILITY-NUMBER	B	2.0		F
			Not referenced by PASGET.				
1	TH		OPERATOR-ID	B	2.0		F
			Not referenced by PASGET.				
G 1	TK		CORRECTION				
			Not referenced by PASGET.				
2	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
2	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
2	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
1	BB		STARTING-CAGE	B	2.0		F
			Retrieved by PASGET for PAS validation reports.				
1	CC		ENDING-CAGE	B	2.0		F
			Retrieved by PASGET for PAS validation reports.				
1	DD		TREATMENT-ID	A	3.0		N
			Retrieved by PASGET for PAS validation reports.				
1	EE		PROCEDURE-ID	A	3.0		N
			Retrieved by PASGET for PAS validation reports.				

FILE...: CAGE-CENSUS-EIS

TYPE...: ADABAS

FILE-NR: 44

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
1	FF	WGT-ITEM	U	3.0	N	
		Retrieved by PASGET for PAS validation reports.				
1	GG	RELATIVE-START	B	2.0	F	
		Retrieved by PASGET for PAS validation reports.				
P 1	HH	ANY-ANIMAL-IDS				
		Retrieved by PASGET for PAS validation reports.				
2	II	SSSC-ID	B	2.0	F	
		Retrieved by PASGET for PAS validation reports.				
2	JJ	NO-ANIMALS	U	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
1	SD	SUPER-KEY	A	11.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY		1	9	
		STARTING-CAGE		1	2	
		Used by PASGET to select data from the file.				
1	SL	AGENCY-EXP-TEST-ENDING-CAGE	A	11.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY		1	9	
		ENDING-CAGE		1	2	
		Not used by PASGET.				

CAGE-CONDITION-TABLE

Physical layout of cage condition table. Layout is stored in the Tables Description table with TABLE_NAME set to TBCONDTC. TBGET uses that information to form direct calls which retrieve data from this table for PAS validation reports. This layout correspond to the view TBCONDTC.

FILE...: CAGE-CONDITION-TABLE

TYPE...: ADABAS

FILE-NR: 24

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	CT CAGE-CONDITION-TABLE				
2	AA	ROW-NUMBER	B	2.0	UQ	
		Refer to TBCONDTC. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBCONDTC. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Refer to TBCONDTC. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0	UQ	
		Refer to TBCONDTC. Key used by TBGET to retrieve data for PAS validation reports.				
2	EE	EE-FILLER	A	1.0	F	
		Refer to TBCONDTC. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0	F	
		Refer to TBCONDTC. Not referenced by TBGET.				
2	GG	SHORT-TEXT	A	8.0	N	
		Refer to TBCONDTC. Not referenced by TBGET.				
2	HH	LONG-TEXT	A	32.0	N	
		Refer to TBCONDTC. Retrieved by TBGET for PAS validation reports.				
2	II	II-FILLER	A	1.0	N	
		Refer to TBCONDTC. Not referenced by TBGET.				

FILE...: CAGE-CONDITION-TABLE

TYPE...: ADABAS

FILE-NR: 24

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
2 JJ CAGE-COND-INCOMP	B	2.0 N
Refer to TBCONDTC. Not referenced by TBGET.		

CAGE-DATA-EIS

Physical layout for file 107. Supports the TD-CAGE-DATA view. No current TDMS software accesses this file using direct calls.

FILE...: CAGE-DATA-EIS

TYPE...: ADABAS

FILE-NR: 107

PRIMARY SEQUENCE FIELD: EXP-TEST-CAGE-DATE

T	L	DB	NAME	F	LENG	S	DE
1	AA		EXP-TEST-CAGE-DATE	U	16.0		DE
1	AY		FEEDER-BALANCE-NUMBER	U	5.0		N
1	BV		BOTTLE-BALANCE-NUMBER	U	5.0		N
1	AB		FEEDER-TIME	U	5.0		N
1	AC		FEEDER-OPERATOR-ID	U	5.0		N
1	AD		FEEDER-ERROR-FLAG	A	1.0		F
1	AE		FEED-CONSUMED-BY-CAGE	U	10.0		N
1	AF		HOURS-FEEDER-IN-CAGE	U	5.0		N
1	AG		AVG-FEED-PER-ANIMAL-HOUR	U	10.0		N
1	AH		FEEDER-ANIMAL-HOURS-USED	U	5.0		N
1	AI		ACTION-OVERRIDE-FLAG-OLD-FEEDER	A	1.0		F
1	AJ		MANUAL-ENTRY-FLAG-OLD-FEEDER	A	1.0		F
1	AK		REMOVAL-ACTION-FLAG-OLD-FEEDER	A	1.0		F
1	AL		BACKWARD-FORWARD-FLAG-OLD-FEEDER	A	1.0		F
1	AM		EXPECTED-SELECTED-FLAG-OLD-FEED	A	1.0		F
1	AN		OLD-FEEDER-MIX-NUMBER	U	10.0		N
1	AO		OLD-FEEDER-WEIGHT	U	10.0		N
1	AP		OLD-FEEDER-WEIGHT-OVERRIDE	A	1.0		F
1	AQ		ACTION-OVERRIDE-FLAG-NEW-FEEDER	A	1.0		F
1	AR		MANUAL-ENTRY-FLAG-NEW-FEEDER	A	1.0		F
1	AS		REMOVAL-ACTION-FLAG-NEW-FEEDER	A	1.0		F
1	AT		BACKWARD-FORWARD-FLAG-NEW-FEEDER	A	1.0		F
1	AU		EXPECTED-SELECTED-FLAG-NEW-FEED	A	1.0		F
1	AV		NEW-FEEDER-MIX-NUMBER	U	10.0		N
1	AW		NEW-FEEDER-WEIGHT	U	10.0		N
1	AX		NEW-FEEDER-WEIGHT-OVERRIDE	A	1.0		F
1	CL		WATER-TIME	U	5.0		N
1	AZ		WATER-OPERATOR-ID	U	5.0		N
1	BA		WATER-ERROR-FLAG	A	1.0		F
1	BB		WATER-CONSUMED-BY-CAGE	U	10.0		N
1	BC		AVG-WATER-PER-ANIMAL-HOUR	U	10.0		N
1	BD		HOURS-BOTTLE-IN-CAGE	U	5.0		N
1	BE		BOTTLE-ANIMAL-HOURS-USED	U	5.0		N
1	BF		ACTION-OVERRIDE-FLAG-OLD-BOTTLE	A	1.0		F

FILE...: CAGE-DATA-EIS

TYPE...: ADABAS

FILE-NR: 107

PRIMARY SEQUENCE FIELD: EXP-TEST-CAGE-DATE

T	L	DB	NAME	F	LENG	S	DE
1	BG		MANUAL-ENTRY-FLAG-OLD-BOTTLE	A	1.0	F	
1	BH		REMOVAL-ACTION-FLAG-OLD-BOTTLE	A	1.0	F	
1	BI		BACKWARD-FORWARD-FLAG-OLD-BOTTLE	A	1.0	F	
1	BJ		EXPECTED-SELECTED-FLAG-OLD-BOTT	A	1.0	F	
1	BK		OLD-BOTTLE-MIX-NUMBER	U	10.0	N	
1	BL		OLD-BOTTLE-WEIGHT	U	10.0	N	
1	BM		OLD-BOTTLE-WEIGHT-OVERRIDE	A	1.0	F	
1	BN		ACTION-OVERRIDE-FLAG-NEW-BOTTLE	A	1.0	F	
1	BO		MANUAL-ENTRY-FLAG-NEW-BOTTLE	A	1.0	F	
1	BP		REMOVAL-ACTION-FLAG-NEW-BOTTLE	A	1.0	F	
1	BQ		BACKWARD-FORWARD-FLAG-NEW-BOTTLE	A	1.0	F	
1	BR		EXPECTED-SELECTED-FLAG-NEW-BOTT	A	1.0	F	
1	BS		NEW-BOTTLE-MIX-NUMBER	U	10.0	N	
1	BT		NEW-BOTTLE-WEIGHT	U	10.0	N	
1	BU		NEW-BOTTLE-WEIGHT-OVERRIDE	A	1.0	F	
1	CM		OBSERVATION-TIME	U	5.0	N	
1	BW		OBSERVATION-OPERATOR-ID	U	5.0	N	
1	BX		OBSERVATION-ERROR-FLAG	A	1.0	F	
1	BY		ACTION-OVERRIDE-FLAG-OBS	A	1.0	F	
1	BZ		MANUAL-ENTRY-FLAG-OBS	A	1.0	F	
1	CA		REMOVAL-ACTION-FLAG-OBS	A	1.0	F	
1	CB		BACKWARD-FORWARD-FLAG-OBS	A	1.0	F	
1	CC		EXPECTED-SELECTED-FLAG-OBS	A	1.0	F	
1	CE		MIS-ANIMAL-ID-TYPE	A	1.0	F	
1	CN		NOTE-TIME	U	5.0	N	
1	CF		NOTE-OPERATOR-ID	U	5.0	N	
1	CG		NOTE-ERROR-FLAG	A	1.0	F	
1	CH		NOTE-ACTION-NUMBER	U	5.0	N	
1	CI		CAGE-NOTE	A	80.0	N	
M 1	CJ		CAGE-OBSERVATION	U	5.0	N	
M 1	CK		MIS-ANIMAL-IDENTIFICATION	U	5.0	N	
1	CO		MISID-ANIMAL-ERROR-FLAG	A	1.0	F	
1	CP		MISID-ANIMAL-TIME	U	5.0	N	
1	CQ		MISID-ANIMAL-OPERATOR-ID	U	5.0	N	
1	EA		OLD-BOTTLE-STATUS	P	5.0	N	
1	EB		NEW-BOTTLE-STATUS	P	5.0	N	
1	EC		CAGE-NOTE-MANUAL-FLAG	A	1.0	N	

FILE...: CAGE-DATA-EIS

TYPE...: ADABAS

FILE-NR: 107

PRIMARY SEQUENCE FIELD: EXP-TEST-CAGE-DATE

T	L	DB NAME	F	LENG	S	DE
1	DA	FEEDER-SUPERKEY	N	21.0		SP
		SOURCE FIELD(S) --- -START- --END-				
		EXP-TEST-CAGE-DATE		1		16
		FEEDER-TIME		1		5
1	DB	WATER-SUPERKEY	N	21.0		SP
		SOURCE FIELD(S) --- -START- --END-				
		EXP-TEST-CAGE-DATE		1		16
		WATER-TIME		1		5
1	DC	NOTES-SUPERKEY	N	21.0		SP
		SOURCE FIELD(S) --- -START- --END-				
		EXP-TEST-CAGE-DATE		1		16
		NOTE-TIME		1		5

CAGE-EIS

Physical layout of file 105. Supports the TD-CAGE view. No current TDMS software accesses this file using direct calls.

FILE...: CAGE-EIS

TYPE...: ADABAS

FILE-NR: 105

PRIMARY SEQUENCE FIELD: EXP-TEST-CAGE

T	L	DB	NAME	F	LENG	S	DE
1	AA	EXP-TEST-CAGE		U	11.0	F	DE
1	AB	CAGE-OPERATOR-ID		U	5.0	N	
1	AC	CAGE-ERROR-FLAG		A	1.0	F	
1	AD	CAGE-DATE		U	5.0	N	
1	AE	CAGE-TIME		U	5.0	N	
1	AF	TREATMENT-NUMBER		U	3.0	F	
1	AG	PROCEDURE-ACTION-SET-NUMBER		U	5.0	N	
1	AH	START-DATE		U	5.0	N	
1	AI	AREA-NUMBER		U	5.0	N	
1	AJ	ABORT-DATE		U	5.0	N	
1	AK	ABORT-TIME		U	5.0	N	
1	AL	ABORT-OPERATOR-ID		U	5.0	N	
1	AM	ABORT-ERROR-FLAG		A	1.0	F	
1	AN	ABORTED-DATE		U	5.0	N	
1	AO	EXP-TEST-TREATMENT		U	10.0		SP
		SOURCE FIELD(S) --- -START- --END-					
		EXP-TEST-CAGE			1		7
		TREATMENT-NUMBER			1		3
1	AP	EXP-TEST-TREATMENT-CAGE		U	14.0		SP
		SOURCE FIELD(S) --- -START- --END-					
		EXP-TEST-CAGE			1		7
		TREATMENT-NUMBER			1		3
		EXP-TEST-CAGE			8		11
1	AQ	EXP-TEST		N	7.0		SP
		SOURCE FIELD(S) --- -START- --END-					
		EXP-TEST-CAGE			1		7

CALIB-WGTS-EIS

Physical layout of file 60. Supports the PAS_CALI_WGT view.

PASGET uses this layout to issue direct calls for PAS validation reports.

FILE...: CALIB-WGTS-EIS

TYPE...: ADABAS

FILE-NR: 60

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA	PAS-CALI-WGT	Not referenced by PASGET.				
G 2	TB	PAS-KEY	Not referenced by PASGET.				
3	TC	DATA-BASE-KEY	Used by PASGET to retrieve data from the file.	A	9.0		DE
2	TE	DATA-DATE	Not referenced by PASGET.	A	8.0		N
2	TF	DATA-TIME	Not referenced by PASGET.	A	6.0		N
2	TG	DATA-FACILITY-NUMBER	Not referenced by PASGET.	B	2.0		N
2	TH	OPERATOR-ID	Not referenced by PASGET.	B	2.0		N
G 2	TK	CORRECTION	Not referenced by PASGET.				
3	TL	ORG-OR-COR	Not referenced by PASGET.	A	1.0		F
3	TM	REASON	Not referenced by PASGET.	A	78.0		N
3	TN	COR-DATE	Not referenced by PASGET.	A	8.0		N
3	TO	COR-TIME	Not referenced by PASGET.	A	6.0		N
G 1	AA	STR-CALIBRATION-SEG	Not referenced by PASGET.				
2	BB	NO-OF-WEIGHINGS	Retrieved by PASGET for PAS validation reports.	B	2.0		N
2	CC	CALIBRATION-CT	Retrieved by PASGET for PAS validation reports.	B	2.0		N

FILE...: CALIB-WGTS-EIS

TYPE...: ADABAS

FILE-NR: 60

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

P 1 DD CALIBRATION-WT		
-----------------------	--	--

Retrieved by PASGET for PAS
validation reports.

2 EE CAL-WT-DES	A	32.0 N
-----------------	---	--------

Retrieved by PASGET for PAS
validation reports.

2 FF MAX-WT	B	4.0 N
-------------	---	-------

Retrieved by PASGET for PAS
validation reports.

2 GG MIN-WT	B	4.0 N
-------------	---	-------

Retrieved by PASGET for PAS
validation reports.

CLIN-OBS-GROUP-TABLE

Physical layout of file 36. Supports the TBCLOBSCGRP view. No TDMS software uses this layout to issue direct calls.

FILE...: CLIN-OBS-GROUP-TABLE

TYPE...: ADABAS

FILE-NR: 36

PRIMARY SEQUENCE FIELD: CLNOBS-GROUP

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	DA	CLINICAL-OBSERVATION-GROUP-TABLE				
	2	DB	CLNOBS-GROUP	A	5.0	N	UQ
	2	DC	CLNOBS-GROUP-SHORT-TEXT	A	14.0	N	
	2	DD	CLNOBS-GROUP-LONG-TEXT	A	32.0	N	

CLIN-OBS-INCOMP-TABLE

Physical layout of clinical observation incompatibility table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBCLOBSI. TBGET uses that layout description to form direct calls to retrieve data from the table for PAS validation reports. This layout corresponds to the view TBCLOBSI.

FILE...: CLIN-OBS-INCOMP-TABLE

TYPE...: ADABAS

FILE-NR: 26

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	IT CLIN-OBS-INCOMPATIBILITY-TABLE				
2	AA	ROW-NUMBER	B	2.0	UQ	
		Refer to TBCLOBSI. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBCLOBSI. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Refer to TBCLOBSI. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0	UQ	
		Refer to TBCLOBSI. Not referenced by TBGET.				
2	EE	EE-FILLER	A	1.0	F	
		Refer to TBCLOBSI. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0	F	
		Refer to TBCLOBSI. Not referenced by TBGET.				
2	GG	LONG-TEXT	A	32.0	N	
		Refer to TBCLOBSI. Retrieved by TBGET for PAS validation reports.				
2	HH	CLASS-VALUE	B	2.0	DE	
		Refer to TBCLOBSI. Key by which TBGET retrieves data from this table.				
2	II	TOTAL-CL-VALUE	B	2.0	N	
		Refer to TBCLOBSI. Not referenced by TBGET.				

CLIN-OBSERV-EIS

Physical layout of file 48. Supports the PAS_CLINICAL_OBS view.
This is the layout by which PASGET selects data for PAS validation reports.

FILE...: CLIN-OBSERV-EIS

TYPE...: ADABAS

FILE-NR: 48

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA		PAS-CLINICAL-OBS				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from the file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		STR-CLINICAL-OBS-SEG				
			Not referenced by PASGET.				
2	BB		PAS-OBS-CODE-CT	B	2.0		
			Retrieved by PASGET for PAS validation reports.				
2	CC		PAS-SELECT-SITE-CT	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				

FILE...: CLIN-OBSERV-EIS

TYPE...: ADABAS

FILE-NR: 48

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
2	DD	PAS-SELECT-SIZE-CT	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	EE	PAS-REM-REASON-CT	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	FF	PAS-CAGE-COND-CT	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
P 1	GG	ANY-OBSERVATION				
		Retrieved by PASGET for PAS validation reports.				
2	HH	OBSERVATION-CODE	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	II	SITE-REQ	A	1.0		
		Retrieved by PASGET for PAS validation reports.				
2	JJ	SIZE-REQ	A	1.0		
		Retrieved by PASGET for PAS validation reports.				
M 1	LL	SITE-CODE	B	2.0		
		Retrieved by PASGET for PAS validation reports.				
M 1	MM	SIZE-CODE	B	2.0		
		Retrieved by PASGET for PAS validation reports.				
M 1	NN	REMOVAL-CODE	B	2.0		
		Retrieved by PASGET for PAS validation reports.				
M 1	OO	CAGE-COND-CODE	B	2.0		
		Retrieved by PASGET for PAS validation reports.				

CLINICAL-OBS-TABLE

Physical layout of clinical observation table. The layout is also stored in the Tables Description table with the field TABLE_NAME set to TBCLOBSC. That information is used by TBGET to form direct calls against this table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBCLOBSC.

FILE...: CLINICAL-OBS-TABLE

TYPE...: ADABAS

FILE-NR: 25

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	CT CLINICAL-OBSERVATION-TABLE				
2	AA	ROW-NUMBER	B	2.0		UQ
		Refer to TBCLOBSC. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBCLOBSC. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0		F
		Refer to TBCLOBSC. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0	F	UQ
		Refer to TBCLOBSC. Key used by TBGET to retrieve data from table.				
2	EE	EE-FILLER	A	1.0		F
		Refer to TBCLOBSC. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0		F
		Refer to TBCLOBSC. Not referenced by TBGET.				
2	GG	SUBSET-INDICATOR	A	1.0		F
		Refer to TBCLOBSC. Not referenced by TBGET.				
2	HH	SHORT-TEXT	A	14.0		DE
		Refer to TBCLOBSC. Not referenced by TBGET.				
2	II	LONG-TEXT	A	33.0		
		Refer to TBCLOBSC. Retrieved by TBGET for EIS/PEIS and PAS validation reports.				

FILE...: CLINICAL-OBS-TABLE

TYPE...: ADABAS

FILE-NR: 25

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
2 JJ JJ-FILLER Refer to TBCLOBSC. Not referenced by TBGET.	A	1.0 N
2 KK CLASS-VALUE Refer to TBCLOBSC. Not referenced by TBGET.	B	2.0 N
2 LL INCOMPAT-TOTAL Refer to TBCLOBSC. Not referenced by TBGET.	B	2.0 N
2 MM SEX-FLAG Refer to TBCLOBSC. Not referenced by TBGET.	A	1.0 F
2 NN MULTI-SELECT-FLAG Refer to TBCLOBSC. Not referenced by TBGET.	A	1.0 N
2 OO NON-CHG-FLAG Refer to TBCLOBSC. Not referenced by TBGET.	A	1.0 N
1 PP CLNOBS-GROUP Refer to TBCLOBSC. Not referenced by TBGET.	P	3.0 N
M 1 QQ EIQUAL-GROUP Refer to TBCLOBSC. Not referenced by TBGET.	P	3.0 N
M 1 RR EISITE-NUM Refer to TBCLOBSC. Not referenced by TBGET.	P	3.0 N
1 TT EIQUAL-GROUP-NUM Refer to TBCLOBSC. Not referenced by TBGET.	P	3.0 N DE

CORRECTION-DATA

This is the physical layout of correction data. The physical file supports the views for each of the 20 correction segments. In-life correction segments are reported in EIS report 22; pathology correction segments are not reported at all. File contains only correction data since the introduction of LDAS. Prior correction data is found in TD-CORRECTION-REASON.

FILE...: CORRECTION-DATA

TYPE...: ADABAS

FILE-NR: 87

PRIMARY SEQUENCE FIELD:

T	L	DB	NAME	F	LENG	S	DE
1	AA	ANC-ANIMAL-NOTE-CORR-KEY-UQ	Primary key for ANIMAL-NOTE-CORRECTION.	A	36.0	N	UQ
1	AB	ANC-ANIMAL-NOTE-CORR-KEY-LINK-UQ	Link key for ANIMAL-NOTE-CORRECTION.	A	36.0	N	UQ
1	AC	AOC-ANIMAL-OBS-CORR-KEY-UQ	Primary key for ANIMAL-OBSERVATION-CORRECTION.	A	36.0	N	UQ
1	AD	AOC-ANIMAL-OBS-CORR-KEY-LINK-UQ	Link key for ANIMAL-OBSERVATION-CORRECTION.	A	36.0	N	UQ
1	AE	ARC-ANIMAL-REM-CORR-KEY-UQ	Primary key for ANIMAL-REMOVAL-CORRECTION.	A	26.0	N	UQ
1	AF	ARC-ANIMAL-REM-CORR-KEY-LINK-UQ	Link key for ANIMAL-REMOVAL-CORRECTION.	A	26.0	N	UQ
1	AG	ATC-ANIMAL-TRAN-CORR-KEY-UQ	Primary key for ANIMAL-TRANSFER-CORRECTION.	A	26.0	N	UQ
1	AH	ATC-ANIMAL-TRAN-CORR-KEY-LINK-UQ	Link key for ANIMAL-TRANSFER-CORRECTION.	A	26.0	N	UQ
1	AI	AWC-ANIMAL-WGT-CORR-KEY-UQ	Primary key for ANIMAL-WEIGHT-CORRECTION.	A	36.0	N	UQ
1	AJ	AWC-ANIMAL-WGT-CORR-KEY-LINK-UQ	Link key for ANIMAL-WEIGHT-	A	36.0	N	UQ

FILE...: CORRECTION-DATA

TYPE...: ADABAS

FILE-NR: 87

PRIMARY SEQUENCE FIELD:

T	L	DB NAME	F	LENG	S	DE
		CORRECTION.				
1	AK	AAC-ANY-ANIMAL-CORR-KEY-UQ	A	26.0	N	UQ
		Primary key for ANY-ANIMAL-CORRECTION.				
1	AL	AAC-ANY-ANIMAL-CORR-KEY-LINK-UQ	A	26.0	N	UQ
		Link key for ANY-ANIMAL-CORRECTION.				
1	AM	AMC-MICRO-ORG-CORR-KEY-UQ	A	32.0	N	UQ
		Primary key for ANY-MICRO-ORGAN-CORRECTION.				
1	AN	AMC-MICRO-ORG-CORR-KEY-LINK-UQ	A	32.0	N	UQ
		Link key for ANY-MICRO-ORGAN-CORRECTION.				
1	AO	BWC-BOTTLE-WGT-CORR-KEY-UQ	A	31.0	N	UQ
		Primary key for BOTTLE-WEIGHT-CORRECTION.				
1	AP	BWC-BOTTLE-WGT-CORR-KEY-LINK-UQ	A	31.0	N	UQ
		Link key for BOTTLE-WEIGHT-CORRECTION.				
1	AQ	CIC-CAGE-ID-CORR-KEY-UQ	A	21.0	N	UQ
		Primary key for CAGE-ID-CORRECTION.				
1	AR	CIC-CAGE-ID-CORR-KEY-LINK-UQ	A	21.0	N	UQ
		Link key for CAGE-ID-CORRECTION.				
1	AS	CNC-CAGE-NOTE-CORR-KEY-UQ	A	31.0	N	UQ
		Primary key for CAGE-NOTE-CORRECTION.				
1	AT	CNC-CAGE-NOTE-CORR-KEY-LINK-UQ	A	31.0	N	UQ
		Link key for CAGE-NOTE-CORRECTION.				
1	AU	FWC-FEEDER-WGT-CORR-KEY-UQ	A	31.0	N	UQ
		Primary key for FEEDER-WEIGHT-CORRECTION.				
1	AV	FWC-FEEDER-WGT-CORR-KEY-LINK-UQ	A	31.0	N	UQ
		Link key for FEEDER-WEIGHT-CORRECTION.				

FILE...: CORRECTION-DATA

TYPE...: ADABAS

FILE-NR: 87

PRIMARY SEQUENCE FIELD:

T	L	DB	NAME	F	LENG	S	DE
1	AW	HIC-HISTOLOGY-CORR-KEY-UQ	Primary key for HISTOLOGY-NUMBER-CORRECTION.	A	26.0	N	UQ
1	AX	HIC-HISTOLOGY-CORR-KEY-LINK-UQ	Link key for HISTOLOGY-NUMBER-CORRECTION.	A	26.0	N	UQ
1	AY	MAC-MICRO-ANI-NOTE-CORR-KEY-UQ	Primary key for ANIMAL-MICRO-NOTES-CORRECTION.	A	37.0	N	UQ
1	AZ	MAC-MICRO-ANI-NTE-COR-KEY-LNK-UQ	Link key for ANIMAL-MICRO-NOTES-CORRECTION.	A	37.0	N	UQ
1	BA	MOC-MICRO-OBS-CORR-KEY-UQ	Primary key for ANY-MICRO-OBS-CORRECTION.	A	42.0	N	UQ
1	BB	MOC-MICRO-OBS-CORR-KEY-LINK-UQ	Link key for ANY-MICRO-OBS-CORRECTION.	A	42.0	N	UQ
1	BC	SSC-MICRO-SITE-CORR-KEY-UQ	Primary key for MICRO-SITE-STATUS-CORRECTION.	A	42.0	N	UQ
1	BD	SSC-MICRO-SITE-CORR-KEY-LINK-UQ	Link key for MICRO-SITE-STATUS-CORRECTION.	A	42.0	N	UQ
1	BE	ONC-ORGAN-NOTE-CORR-KEY-UQ	Primary key for ORGAN-NOTES-CORRECTION.	A	42.0	N	UQ
1	BF	ONC-ORGAN-NOTE-CORR-KEY-LINK-UQ	Link key for ORGAN-NOTES-CORRECTION.	A	42.0	N	UQ
1	BG	TNC-TEST-NOTE-CORR-KEY-UQ	Primary key for TEST-NOTE-CORRECTION.	A	27.0	N	UQ
1	BH	TNC-TEST-NOTE-CORR-KEY-LINK-UQ	Link key for TEST-NOTE-CORRECTION.	A	27.0	N	UQ
1	BI	TGC-TREAT-NOTE-CORR-KEY-UQ	Primary key for TREAT-NOTE-CORRECTION.	A	30.0	N	UQ
1	BJ	TGC-TREAT-NOTE-CORR-KEY-LINK-UQ	Link key for TREAT-NOTE-CORRECTION.	A	30.0	N	UQ

FILE...: CORRECTION-DATA

TYPE...: ADABAS

FILE-NR: 87

PRIMARY SEQUENCE FIELD:

T L DB NAME	F	LENG S DE
1 CA CORRECTION-OPERATOR-ID ID of operator who made correction. Common to all correction segments.	P	3.0 F
1 CB CORRECTION-REASON Correction reason code assigned to correction. Common to all segments.	A	2.0 F
1 CC CORRECTION-NOTE Description of correction. Common to all segments.	A	79.0 N
1 CD CORRECTION-TYPE Type of correction, A for add, C for change, D for delete. Common to all segments.	A	1.0 F
1 DA FLAG1 Field used in ANY-ANIMAL-, CAGE-ID-, ANIMAL-OBSERVATION-, ANIMAL-WEIGHT-, ANIMAL-NOTE-, FEEDER-WEIGHT-, CAGE-NOTE-, BOTTLE-WEIGHT-, ANIMAL-REMOVAL-, HISTOLOGY-NUMBER-, ANY-MICRO-OBS-, and MICRO-ANIMAL-NOTES-CORRECTION segments.	A	1.0 F
1 DB FLAG2 Field used in ANY-ANIMAL-, ANIMAL-REMOVAL-, BOTTLE-WEIGHT-, and FEEDER-WEIGHT-CORRECTION segments.	A	1.0 F
1 DC FLAG3 Field used in ANIMAL-REMOVAL-CORRECTION.	A	1.0 F
1 DD CORR-REASON Correction reason code for CAGE-ID-, and HISTOLOGY-NUMBER-CORRECTION.	A	2.0 N

FILE...: CORRECTION-DATA

TYPE...: ADABAS

FILE-NR: 87

PRIMARY SEQUENCE FIELD:

T L DB NAME	F	LENG S DE
1 DE CARCASS-IDENTIFICATION Used in ANIMAL-REMOVAL-CORRECTION.	A	14.0 N
1 DF HISTOLOGY-NUMBER Used in HISTOLOGY-NUMBER-CORRECTION.	A	16.0 N
1 DH CORR-NOTE Used in CAGE-ID-, and HISTOLOGY-NUMBER-CORRECTION.	A	79.0 N
M 1 DI SHORT-NOTE Used in ANIMAL-NOTE-, CAGE-NOTE-, MICRO-ANIMAL-NOTES-, and ORGAN-NOTE-CORRECTION.	A	80.0 N
1 DJ LONG-NOTE Used in TEST-NOTE-, and TREAT-NOTE-CORRECTION.	A	160.0 N
1 EA TREATMENT-NUMBER Used in CAGE-ID-CORRECTION.	P	3.0 N
1 ED NUMBER1 Common to all correction segments.	P	3.0 N
1 EE NUMBER2 Used in ANY-ANIMAL-, ANY-MICRO-ORGAN-, ANIMAL-REMOVAL-, ANIMAL-TRANSFER-, ANIMAL-WEIGHT-, BOTTLE-WEIGHT-, CAGE-ID-, FEEDER-WEIGHT-, HISTOLOGY-NUMBER-, MICRO-ANIMAL-NOTES-, ANY-MICRO-OBS-, NOT-EXAMINED-, ORGAN-NOTE-, MICRO-SITE-STATUS-, and UNCERTAIN-CORRECTION.	P	3.0 N
1 EF NUMBER3 Used in ANY-ANIMAL-, ANY-MICRO-ORGAN-, ANIMAL-REMOVAL-, ANIMAL-TRANSFER-, CAGE-ID-, HISTOLOGY-NUMBER-, ANY-MICRO-OBS-, NOT-EXAMINED-, MICRO-SITE-STATUS-, and UNCERTAIN-CORRECTION.	P	3.0 N
1 EG NUMBER4 Used in ANY-ANIMAL-, ANY-MIC-	P	3.0 N

FILE...: CORRECTION-DATA

TYPE...: ADABAS

FILE-NR: 87

PRIMARY SEQUENCE FIELD:

T L DB NAME	F	LENG S DE
RO-ORGAN-, ANIMAL-REMOVAL-, ANIMAL-TRANSFER-, CAGE-ID-, HISTOLOGY-NUMBER-, NOT-EXAM- INED-, MICRO-SITE-STATUS-, and UNCERTAIN-CORRECTION.		
1 EH NUMBER5 Used by ANY-ANIMAL-, ANY-MIC- RO-ORGAN-, ANIMAL-REMOVAL-, CAGE-ID-, and HISTOLOGY-NUM- BER-CORRECTION.	P	3.0 N
1 EI NUMBER6 Used by ANY-ANIMAL-, ANIMAL- REMOVAL-, CAGE-ID-, and ANY- MICRO-OBS-CORRECTION.	P	3.0 N
1 EJ NUMBER7 Used by ANY-ANIMAL-, CAGE-ID-, and ANY-MICRO-OBS-CORRECTION.	P	3.0 N
1 EL WEIGHT1 Used by ANIMAL-REMOVAL-, ANI- MAL-WEIGHT, BOTTLE-WEIGHT-, and FEEDER-WEIGHT-CORRECTION.	P	7.0 N
1 EM WEIGHT2 Used by BOTTLE-WEIGHT- and FEEDER-WEIGHT-CORRECTION.	P	7.0 N
P 1 FA OBSERVATION-DATA Used by ANIMAL-OBSERVATION- and ANIMAL-REMOVAL-CORRECTION.		
2 FB OBSERVATION Used by ANIMAL-OBSERVATION- and ANIMAL-REMOVAL-CORRECTION.	P	3.0 N
2 FC SITE Used by ANIMAL-OBSERVATION- and ANIMAL-REMOVAL-CORRECTION.	P	3.0 N
2 FD EIQUAL1 Used by ANIMAL-OBSERVATION- and ANIMAL-REMOVAL-CORRECTION.	P	3.0 N
2 FE EIQUAL2 Used by ANIMAL-OBSERVATION-	P	3.0 N

FILE...: CORRECTION-DATA

TYPE...: ADABAS

FILE-NR: 87

PRIMARY SEQUENCE FIELD:

T	L	DB NAME	F	LENG	S	DE
		and ANIMAL-REMOVAL-CORRECTION.				
M	1	FF PATHOLOGY-QUALIFIER	P	3.0	N	
		Used in ANY-MICRO-OBS-CORRECTION.				
M	1	FG PATHOLOGY-SITE	P	3.0	N	
		Used in ANY-MICRO-OBS-CORRECTION.				
1	FH	UCC-UNCERTAIN-CORR-KEY-UQ	A	26.0	N	UQ
		Primary key for UNCERTAIN-CORRECTION.				
1	FI	UCC-UNCERTAIN-CORR-KEY-LINK-UQ	A	26.0	N	UQ
		Link key for UNCERTAIN-CORRECTION.				
1	FJ	NEC-NO-EXAM-CORR-KEY-UQ	A	26.0	N	UQ
		Primary key for NOT-EXAMINED-CORRECTION.				
1	FK	NEC-NO-EXAM-CORR-KEY-LINK-UQ	A	26.0	N	UQ
		Link key for NOT-EXAMINED-CORRECTION.				

DATABASE-TRANSACTION

Physical/logical layout of transaction super-table, which embraces 21 user views. Each view corresponds to one of the transaction types that UPDATE will process. This file is updated by both RECEIVE and ECS. The file design is not programmer-friendly as it maintains the older VSAM file structure.

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
- - - - -						
G	1	AC HEADER				
		Group embracing fields that are common to most segments.				
2	AA	SEGMENT-ID	A	2.0	F	
		Two-letter segment identifier. Not referenced by Update; informational only.				
2	AB	MICRO-OR-ECS-FLAG	A	1.0	F	
		Identifier for source of transaction. A - LDAS, E - Mainframe ECS. Used by Update to generate Update statistics report.				
2	AE	TOX-MICRO-DATE	B	2.0	F	
		The date, in TDMS TOX format, on which the source set for the transaction was created. Not referenced by Update.				
2	AF	TOX-MICRO-TIME	B	2.0	F	
		The time of day, in TDMS TOX time format, on which the source set for the transaction was created. Not referenced by Update.				
2	AG	MICRO-TRANS-SET	B	2.0	F	
		The number of the transaction's source set. Not referenced by Update.				
2	AI	TOX-MAINFRAME-DATE	B	2.0	F	
		The date, in TDMS TOX date				

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

format, on which the source set was received. Not referenced by Update.

2	AJ	TOX-MAINFRAME-TIME	B	2.0	F	
---	----	--------------------	---	-----	---	--

The time of day, in TDMS TOX time format, on which the source set was received. Not referenced by Update.

2	AL	AGENCY-NUMBER	B	2.0	F	
---	----	---------------	---	-----	---	--

Always 5, for NIH/NIEHS/NTP Not referenced by Update.

2	AP	ENTRY-OPERATOR-NUMBER	B	2.0	F	
---	----	-----------------------	---	-----	---	--

When the transaction is a correspondence, this indicates the operator ID for the person making the correction.

2	AS	SOFTWARE-VERSION	B	2.0	F	
---	----	------------------	---	-----	---	--

The version of LDAS under which the source set was generated. Not referenced by Update.

2	AT	TOX-DATA-DATE	B	2.0	F	
---	----	---------------	---	-----	---	--

The date, in TDMS TOX date format, on which the transaction was created.

2	AU	TOX-DATA-TIME	B	2.0	F	
---	----	---------------	---	-----	---	--

The time of day, in TDMS TOX time format, at which the transaction was created.

2	AV	DATA-FACILITY-NUMBER	B	2.0	F	
---	----	----------------------	---	-----	---	--

The facility associated with the machine on which the source set was created. Not referenced by Update.

2	AW	OPERATOR-OR-PATHOLOGIST	B	2.0	F	
---	----	-------------------------	---	-----	---	--

For an original in-life transaction, the ID of the operator who created the transaction. For an original pathology transaction, this is the operator ID of the pathologist

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

who is assigned to the test.
For corrections, this is the
ID of the operator who origi-
nally made the record being
corrected.

2	AX	PROCESS-FLAG	A	2.0	F		
---	----	--------------	---	-----	---	--	--

Two-letter identifier for the
type of transaction. The first
digit can be A, C, or D, which
signifies add, change, or del-
ete, respectively. The second
digit is either blank, indica-
ting an original transaction,
or 'A', indicating an audit
transaction.

2	AZ	DATA-BASE-KEY-LENGTH	B	2.0	F		
---	----	----------------------	---	-----	---	--	--

The number of digits being
used in the DATA-BASE-KEY
field. Not referenced by Up-
date.

2	A0	NUMBER-OF-TREE-KEY-LEVELS	B	2.0	F		
---	----	---------------------------	---	-----	---	--	--

The number of two-byte binary
segments of the SEGMENT-TREE-
KEY field being used.

2	A2	DATA-BASE-INDICATOR	A	1.0	F		
---	----	---------------------	---	-----	---	--	--

Indicates sub-system. 2 for
in-life, 3 for pathology, and
1 for balance calibration.

2	A3	DATA-BASE-KEY	A	14.0	F		
---	----	---------------	---	------	---	--	--

Indicates the master record
to which the transaction ap-
plies. For example, for an
animal observation transac-
, this field would con-
tain the test and animal num-
ber.

2	A4	SEGMENT-TREE-KEY	B	20.0	F		
---	----	------------------	---	------	---	--	--

This field is parsed into 10
2-byte binary sub-fields.
Different segments use these
sub-fields differently. The
number being used for a trans-

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

action is defined in the
NUMBER-OF-TREE-KEY-LEVELS
field.

G 1 A5 CAGE-ID-DATA

Group embracing fields speci-
fic to the cage ID segment.
Identified by SEGMENT-ID = CI.

2 A6 TREATMENT-NUMBER B 2.0 N

The number of the treatment
group to which the cage is
assigned. Loaded by Update
into the TDC-TREATMENT-NUMBER
field in the TD-CAGE file.

2 A7 CONTROL-TREATMENT-FLAG A 1.0 N

Alphabetic identifier of the
role of the cage's treatment
group. Not processed by Up-
date.

2 A8 PROCEDURE-ACTION-SET-NUMBER B 2.0 N

The procedure action set that
applies to the cage, as de-
fined in test protocol. Load-
ed by Update into the field
TDC-PROCEDURE-ACTION-SET-NUM
of file TD-CAGE.

2 A9 TEST-TYPE B 2.0 N

The test type for the test.
Loaded by Update into the
field TDT-TEST-TYPE of the
TD-TEST file.

2 BA DOSE-DATE B 2.0 N

The first dosing date for the
cage, in TDMS TOX date format.
Loaded by Update into the
field TDC-START-DATE of file
TD-CAGE.

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
2 BB AREA-NUMBER The lab room where the cage is located, as defined in test protocol. Loaded by Update into the field TDC-AREA-NUMBER of file TD-CAGE.	B	2.0 N
G 1 BC ANY-ANIMAL-DATA Group embracing all fields specific to the animal allocation segment. Identified by SEGMENT-ID = AA.		
2 BE ANIMAL-NUMBER The number of the animal being allocated. Redundant field. Update uses same value, which is stored in SEGMENT-TREE-KEY.	B	2.0 N
2 BF ANIMAL-ID Not loaded by RECEIVE. Used by mainframe ECS to correct old TDA-ANIMAL-IDENTIFICATION values in TD-ANIMAL. However, that field was never assigned values by any version of TDMS.	A	1.0 N
2 BG SEX Sex of animal as defined in test protocol. Loaded by Update into the TDA-SEX field of the TD-ANIMAL file. M for male; F for female.	A	1.0 N
2 BH SPECIES The species code of the animal as defined in test protocol. Loaded by Update into the TDA-SPECIES field of the TD-ANIMAL file. 4 for mice; 10000 for rats; 20000 for hamsters.	B	2.0 N
2 BI STRAIN The strain code of the animal as defined in test protocol. Loaded by Update into the TDA-STRAIN field of the TD-ANIMAL	B	2.0 N

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME F LENG S DE

file. 5-99, mice strains;
10001-10099, rat strains;
20001-20099, hamster strains.

2 BJ SUBSTRAIN B 2.0 N

The substrain code of the animal-
mal as defined in test proto-
col. Loaded by Update into
the TDA-SUBSTRAIN field of the
TD-ANIMAL file.

2 BK ANIMAL-CLASS B 2.0 N

The animal class code as de-
fined in test protocol. Load-
ed by Update into the TDA-
ANIMAL-CLASS fields of the
TD-ANIMAL file.

G 1 BL BALANCE-CALIBRATION-DATA

Group embracing segment-speci-
fic fields of balance calibra-
transactions. Identified
by SEGMENT-ID = BC.

2 BN BALANCE-NUMBER B 2.0 N

Number of the balance whose
calibration weights are re-
corded in the transaction.
Loaded by Update into the
TDF-BALANCE-NUMBER-DE field of
the TD-BALANCE-CALIBRATION
file and, if it is a new
balance, into the TDB-BALANCE-
NUMBER-UQ field of the TD-
BALANCE file.

2 BP CALIBRATION-WEIGHT-MANUAL A 1.0 N

Flag indicating whether the
calibration weights were en-
tered manually. Loaded by
Update into the TDF-CALI-
MANUAL-FLAG field of the TD-
BALANCE-CALIBRATION file.

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

P	1	BQ	BALANCE-CALIBRATION-WEIGHT-DATA				
---	---	----	---------------------------------	--	--	--	--

Repeating group, which is not a group at all, containing only calibration weights.

	2	BR	CALIBRATION-WEIGHT	B	2.0	N	
--	---	----	--------------------	---	-----	---	--

The calibration weights recorded automatically or manually on LDAS. Loaded into the TDF-CALIBRATION-WEIGHT field of the TD-BALANCE-CALIBRATION file.

G	1	BT	OBSERVATION-DATA				
---	---	----	------------------	--	--	--	--

Group embracing segment-specific, non-periodic group fields of animal observation transactions. Identified by SEGMENT-ID = AO. (There is only one field in this group.)

	2	BV	MANUAL-ANIMAL-OBSERVATION-FLAG	A	1.0	N	
--	---	----	--------------------------------	---	-----	---	--

Flag indicating whether the observation was entered during a manual entry session. 0 - No; 1 - Yes. Loaded by Update into the TDD-MANUAL-ENTRY-FLAG-OBS field of the TD-ANIMAL-DATA file.

P	1	BY	ALL-ANIMAL-OBSERVATION-SETS				
---	---	----	-----------------------------	--	--	--	--

Repeating group embracing segment-specific transactions of both the animal observation and animal removal segments., identified by SEGMENT-ID = AO and AR, respectively.

	2	BZ	ANIMAL-OBSERVATION-NUM	B	2.0	N	
--	---	----	------------------------	---	-----	---	--

Code for clinical sign observed on animal. Loaded by Update into the TDD-OBSERVATION and TDA-OBSERVATION fields of the TD-ANIMAL-DATA

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

and TD-ANIMAL files, respectively.

2	B0	ANIMAL-OBSERVATION-SITE	B	2.0	N		
---	----	-------------------------	---	-----	---	--	--

Code for site of clinical sign when that sign requires that a site be specified. Loaded by Update into the TDD-SITE in the TD-ANIMAL-DATA file and the TDA-SITE field of the TD-ANIMAL file.

2	B1	ANIMAL-OBSERVATION-QUALIFIER-1	B	2.0	N		
---	----	--------------------------------	---	-----	---	--	--

Code for the size of the clinical sign when that sign requires that a size be specified. Loaded by Update into the TDD-EIQUAL1 field of the TD-ANIMAL-DATA file and the TDA-EIQUAL1 of the TD-ANIMAL file.

2	B2	ANIMAL-OBSERVATION-QUALIFIER-2	B	2.0	N		
---	----	--------------------------------	---	-----	---	--	--

Code for the count of the clinical sign when that sign requires that count be specified. Loaded by Update into the TDD-EIQUAL2 field of the TD-ANIMAL-DATA file and the TDA-EIQUAL2 of the TD-ANIMAL file.

G	1	B3	ANIMAL-WEIGHT-DATA				
---	---	----	--------------------	--	--	--	--

Group embracing segment-specific fields of animal weight transactions. Identified by SEGMENT-ID = AW.

2	B5	ANIMAL-MANUAL-WEIGHT-FLAG	A	1.0	N		
---	----	---------------------------	---	-----	---	--	--

Flag indicating whether animal weight was manually entered. 0 - Automatically recorded from attached scale 1 - Manually entered. Loaded by Update into the TDD-MANUAL-ENTRY-FLAG-WGT field of the TD-ANIMAL-DATA file.

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB NAME	F	LENG	S	DE
2	B7	ANIMAL-WEIGHT-STATUS	P	3.0	N	
<p>When no weight is recorded, status code that explains why weight was not recorded. Loaded by Update into the field TDD-WEIGHT-STATUS of the TD-ANIMAL-DATA file.</p>						
2	B8	ANIMAL-WEIGHT	B	4.0	N	
<p>Weight of animal, expressed in decigrams. Loaded by Update into the TDD-ANIMAL-WEIGHT of the TD-ANIMAL-DATA file.</p>						
G 1	B9	CONTAINER-WEIGHT-DATA				
<p>Group embracing segment-specific fields of feeder and bottle weight transactions. Identified by SEGMENT-ID = FW and BW.</p>						
2	CB	OLD-CONTAINER-MANUAL-WEIGHT-FLAG	A	1.0	N	
<p>Flag indicating whether used container weight was recorded manually. 0 - Automatically recorded by attached scale; 1 - Entered manually. Loaded by Update into the TDG-MANUAL-ENTRY-FLAG-OLD-BOTT and TDG-MANUAL-ENTRY-FLAG-OLD-FEED fields in the TD-CAGE-DATA file.</p>						
2	CD	OLD-CONTAINER-WEIGHT	B	4.0	N	
<p>Weight, in decigrams, of used container. Loaded by Update into the TDG-OLD-FEEDER-WEIGHT and TDG-OLD-BOTTLE-WEIGHT fields in the TD-CAGE-DATA file.</p>						

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

2	CF		NEW-CONTAINER-MANUAL-WEIGHT-FLAG	A	1.0		N
---	----	--	----------------------------------	---	-----	--	---

Flag indicating whether full container weight was recorded manually. 0 - Automatically recorded by attached scale; 1 - Recorded manually. Loaded by Update into the TDG-MANUAL-ENTRY-FLAG-NEW-BOTT and TDG-MANUAL-ENTRY-FLAG-NEW-FEED in the TD-CAGE-DATA files.

2	CH		NEW-CONTAINER-WEIGHT	B	4.0		N
---	----	--	----------------------	---	-----	--	---

Weight, in decigrams, of the full container. Loaded by Update into the TDG-NEW-FEEDER-WEIGHT and TDG-NEW-BOTTLE-WEIGHT in the TD-CAGE-DATA file.

2	CJ		OLD-CONTAINER-STATUS	P	3.0		N
---	----	--	----------------------	---	-----	--	---

Weight status code indicating reason why used weight was not recorded. Loaded by Update into the TDG-OLD-FEEDER-STATUS and TDG-OLD-FEEDER-STATUS in the TD-CAGE-DATA file.

G	1	CT	INLIFE-NOTE-DATA				
---	---	----	------------------	--	--	--	--

Group embracing segment-specific fields of in-life animal and cage note transactions. Identified by SEGMENT-ID = AN and CN.

2	CV		INLIFE-NOTE	A	80.0		N
---	----	--	-------------	---	------	--	---

Free text entered by operator modifying or explaining animal or cage data. Loaded by Update into the TDD-ANIMAL-NOTE of the TD-ANIMAL-DATA file and the TDG-CAGE-NOTE field of the TD-CAGE-DATA file.

2	CW		INLIFE-NOTE-MANUAL-FLAG	A	1.0		N
---	----	--	-------------------------	---	-----	--	---

Flag indicating whether note was entered as part of a manual entry session. 0 - En-

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

tered as part of an entry session of schedule actions;
1 - Entered as part of a manual entry session. Loaded by Update into the TDD-MANUAL-ENTRY-FLAG-NOTE and TDG-CAGE-NOTE-MANUAL-FLAG of the TD-ANIMAL-DATA and TD-CAGE-DATA files, respectively.

G 1 CX ANIMAL-TRANSFER-DATA

Group embracing segment-specific fields in animal transfer transactions, identified by SEGMENT-ID = AT.

2 CY NEW-CAGE-NUMBER P 3.0 N

Number of cage to which the animal is being transferred. Loaded by Update into the TDA-HOSPITAL-CAGE-DE field of the TD-ANIMAL file.

2 CZ TOX-ALLOCATION-DATE B 2.0 N

The date, in TDMS TOX date format, on which the animal was transferred to the new cage. Loaded by Update into the TDA-TRANSFER-DATE field of the TD-ANIMAL file.

2 C0 TOX-ALLOCATION-TIME B 2.0 N

The time of day, in TDMS TOX time format, at which the animal was transferred. Loaded by Update into the TDA-TRANSFER-TIME of the TD-ANIMAL file.

G 1 C1 ANIMAL-REMOVAL-DATA

Group embracing segment-specific, non-repeating group, fields in animal removal

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

transactions. Identified by
SEGMENT-ID = AR.

2	C3	CARCASS-IDENTIFICATION	A	14.0	N		
---	----	------------------------	---	------	---	--	--

Alternative identifier for
animal. The first 2 digits are
05 for agency; next 7 digits
are the test number; the last
5 are animal number. (For
very old tests, the last 5
digits may be constructed as
the 4-digit cage number plus
a sequential cage identifier
within the cage, up to 5.)
Loaded by Update into the
TDA-CARCASS-IDENTIFICATION-DE
field of the TD-ANIMAL file.

2	C5	MANUAL-REMOVAL-FLAG	A	1.0	N		
---	----	---------------------	---	-----	---	--	--

Flag indicating whether remov-
al was recorded as part of a
manual entry session. 0 - No;
1 - Yes. Loaded by Update
into TDA-MANUAL-ENTRY-FLAG-
REM field of the TD-ANIMAL
file.

2	C8	DAYS-ON-EXPERIMENT	B	2.0	N		
---	----	--------------------	---	-----	---	--	--

Dynamically calculated value
equal to the removal date min-
us the start date of the animal-
mals cage. For LDAS tests,
this value is the same as
days-on-dose. Loaded by Up-
date into the TDA-DAYS-ON-
EXPERIMENT field of the TD-
ANIMAL file.

2	DA	MANUAL-REMOVAL-WEIGHT-FLAG	A	1.0	N		
---	----	----------------------------	---	-----	---	--	--

Flag indicating whether remov-
al weight was recorded as part
of a manual entry session. 0-

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

No; 1 - Yes. Loaded by Update into TDA-MANUAL-ENTRY-FLAG-WGT field of the TD-ANIMAL file.

2	DC	REMOVAL-WEIGHT-STATUS	P	3.0	N		
---	----	-----------------------	---	-----	---	--	--

Weight status code indicating the reason for which a removal weight was not taken. Loaded by Update into the TDA-REM-WEIGHT-STATUS field of the TD-ANIMAL file.

2	DD	REMOVAL-WEIGHT	B	4.0	N		
---	----	----------------	---	-----	---	--	--

Weight, in decigrams, of the animal at removal. Loaded by Update into the TDA-ANIMAL-REM-WEIGHT field of the TD-ANIMAL file.

2	DF	REMOVAL-REASON	B	2.0	N		
---	----	----------------	---	-----	---	--	--

Code indicating reason for removal. Loaded by Update into the TDA-REASON-FOR-REMOVAL-DE field of the TD-ANIMAL file.

2	DH	MANUAL-REMOVAL-OBSERVATION-FLAG	A	1.0	N		
---	----	---------------------------------	---	-----	---	--	--

Flag indicating whether removal observations were recorded as part of a manual entry session. 0 - No; 1 - Yes. Loaded by Update into the TDA-MANUAL-ENTRY-FLAG-OBS field of the TD-ANIMAL file.

G	1	DJ	HISTOLOGY-NUMBER-DATA				
---	---	----	-----------------------	--	--	--	--

Group embracing segment-specific fields of histology number transactions. Identified by SEGMENT-ID = HN.

2	DK	HISTOLOGY-OPERATOR	B	2.0	N		
---	----	--------------------	---	-----	---	--	--

ID for operator assigning histology number. Loaded by Update into the TDA-HISTOLOGY-NUM-OPERATOR-ID field of the TD-ANIMAL file.

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
2	DL		HISTOLOGY-NUMBER	A	16.0	N	
			Lab-specific identification number for carcasses. Loaded by Update into TDA-HISTOLOGY-NUMBER field of TD-ANIMAL file				
G 1	DR		ANY-MICRO-ORGAN-DATA				
			Group embracing segment-specific fields in organ status transactions. Identified by SEGMENT-ID = AM.				
2	DS		MICRO-ORGAN-OPERATOR	B	2.0	N	
			ID of operator assigning organ status. Loaded by Update into the TDO-ORGAN-STATUS-OPERATOR-ID of the TD-ORGAN file.				
2	DU		STATUS-NUMBER	B	2.0	N	
			Status of organ being microscopically evaluated. Loaded by Update into the TDO-ORGAN-STATUS field of the TD-ORGAN file.				
G 1	DV		MICRO-OBSERVATION-DATA				
			Group embracing segment-specific fields of microscopic observation transactions. Identified by SEGMENT-ID = MO.				
2	DW		MICRO-OBSERVATION-OPERATOR	B	2.0	N	
			ID of operator who recorded the microscopic observation. Loaded by Update into the TDH-MICRO-OBS-OPERATOR-ID field of the TD-ORGAN-DATA file.				
2	DY		PRIMARY-SITE	B	2.0	N	
			Organ from which the microscopic observation metastasized.				

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

sized. (Blank indicates that organ is original site.)
Loaded by Update into the TDH-PRIMARY-SITE of the TD-ORGAN-DATA file.

2	FR		TRACE-LESION-NUMBER	B	2.0	N	
---	----	--	---------------------	---	-----	---	--

No meaning in TDMS. However, it is referenced by ECS and Update.

2	D1		MORPHOLOGY-NUMBER	B	2.0	N	
---	----	--	-------------------	---	-----	---	--

Morphology being microscopically observed. Loaded by Update into the TDH-OBSERVATION-MORPHOLOGY-DE field of the TD-ORGAN-DATA file.

M	2	D2	QUALIFIER-NUMBER	B	2.0	N	
---	---	----	------------------	---	-----	---	--

Multiple-entry field containing qualifiers that describe morphology being observed. Loaded by Update into the multiple-entry field TDH-OBSERVATION-QUALIFIER of the TD-ORGAN-DATA file. Up to 4 entries per field.

M	2	D3	SITE-NUMBER	B	2.0	N	
---	---	----	-------------	---	-----	---	--

Multiple-entry field containing sites on the organ being observed where the morphology was detected. Loaded by Update into the TDH-OBSERVATION-SITE field of the TD-ORGAN-DATA file. Up to 3 entries per field.

2	D4		CAUSE-OF-DEATH-FLAG	A	1.0	N	
---	----	--	---------------------	---	-----	---	--

Flag indicating whether the observation caused the death of the animal. C for contributory; P for primary. Loaded by Update into the TDH-COD-FLAG field of the TD-ORGAN-DATA file.

G	1	D5	PATHOLOGY-NOTE-DATA				
---	---	----	---------------------	--	--	--	--

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME F LENG S DE

Group embracing segment-specific fields of organ and pathology animal note transactions. Identified by SEGMENT-ID = ON or MN, respectively.

2 D6 PATHOLOGY-NOTE-OPERATOR B 2.0 N

ID of operator recording the note. Loaded by Update into the TDH-ORGAN-NOTE-OPERATOR-ID field of the TD-ORGAN-DATA file or the TDM-OPERATOR-IDENTIFICATION field of the TD-MICRO-ANIMAL-NOTE field.

M 2 D8 PATHOLOGY-NOTE A 80.0 N

Free text describing a microscopically evaluated animal or organ. Loaded by Update into the TDH-ORGAN-NOTE field of the TDH-ORGAN-DATA field or the TDM-MICRO-ANIMAL-NOTE field of the TD-MICRO-ANIMAL-NOTE file.

G 1 EL MICRO-SITE-STATUS-DATA

Legacy group embracing all segment-specific fields in site status transactions. Site status transactions are not generated by LDAS. Any MODCOMP site status information is required to be converted to organ status before modification. Identified by SEGMENT-ID = SS.

2 EM MICRO-SITE-STATUS-OPERATOR-ID B 2.0 N

Id for operator assigning site status. Loaded by Update into the TDO-SITE-STAT-

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
			US-OPERATOR-ID field of the TD-ORGAN file.				
2	EN		MICRO-SITE	B	2.0	N	
			Organ-Site for which a status is being assigned. Loaded by Update into the TDO-SITE-CODE field of the TD-ORGAN file.				
2	EO		MICRO-SITE-STATUS	B	2.0	N	
			Status of the organ-site being evaluated. Loaded by Update into the TDO-SITE-STATUS field of the TD-ORGAN file.				
G 1	FF		NO-EXAM-UCOD-DATA				
			Group embracing segment-specific fields in "not examined" or "uncertain cause-of-death" transactions. Identified by SEGMENT-ID = NE or UC, respectively.				
2	FG		NO-EXAM-UCOD-STUDY	A	7.0	N	
			TDMS test number to which the observation applies.				
2	FH		NO-EXAM-UCOD-ANIMAL-NUM	P	3.0	N	
			Number of animal being classified as either not examined or having uncertain cause-of-death.				
2	FI		NO-EXAM-UCOD-PATHOLOGIST	P	3.0	N	
			Operator ID of the pathologist assigned to the test. Loaded by Update into the TDA-NO-EXAM-PATHOLOGIST or TDA-UNCERTAIN-PATHOLOGIST in the TD-ANIMAL file.				
G 1	FK		TEST-NOTE-DATA				
			Group embracing segment-specific fields in test note transactions. Identified by SEGMENT-ID = TN.				
2	FL		TEST-STUDY	A	7.0	N	
			The TDMS test number to which the note applies.				

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

2	FM		TEST-NOTE	A	160.0	N	
---	----	--	-----------	---	-------	---	--

Free text commenting on some aspect of a test. Loaded by Update into the TDS-NOTE-TEXT field of the TD-TEST-NOTE file

G	1	FN	TREAT-NOTE-DATA				
---	---	----	-----------------	--	--	--	--

Group embracing segment-specific fields of treatment group note transactions. Identified by SEGMENT-ID = TG

2	FO		TREAT-STUDY	A	7.0	N	
---	----	--	-------------	---	-----	---	--

The study to which the note applies.

2	FP		TREAT-GROUP	A	3.0	N	
---	----	--	-------------	---	-----	---	--

The ID of the treatment group to which the note applies.

2	FQ		TREAT-NOTE	A	160.0	N	
---	----	--	------------	---	-------	---	--

Free text commenting on some aspect of a treatment group. Loaded by Update into the TDT-NOTE-TEXT field of the TD-TREATMENT-NOTE file.

G	1	EP	CORRECTION-PORION				
---	---	----	-------------------	--	--	--	--

Group embracing segment-common correction fields.

2	EQ		CORRECTION-DATE	P	3.0	N	
---	----	--	-----------------	---	-----	---	--

Date, in TDMS TOX date format, on which correction was made. Used by Update, along with CORRECTION-TIME, in forming primary key of correction record.

2	ER		CORRECTION-TIME	P	3.0	N	
---	----	--	-----------------	---	-----	---	--

Time of day, in TDMS TOX time format, at which correction was made. Used by Update, along with CORRECTION-DATE, in forming primary key of correction record.

2	ES		CORRECTION-REASON	A	2.0	N	
---	----	--	-------------------	---	-----	---	--

Code for correction reason.

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME F LENG S DE

Loaded by Update into corresponds-
ponding field in correction
record.

2 ET CORRECTION-NOTE A 79.0 N

Free text describing the corresponds-
rection. Loaded by Update in-
to corresponding fields in
correction segments.

2 EU CORRECTION-TYPE A 1.0 N

Single letter indicating type
of correction. A-Add, C-Change
D-Delete

2 EV CAGE-CORRECTION-KEY A 21.0 N DE

Key for Cage Correction segment-
ment. Elements of which the
key is comprised are: the
value that will be the seven-
digit test number, the four-
digit cage number, the
CORRECTION-DATE field, and the
CORRECTION-TIME field. If the
cage number is being changed,
then this value will be load-
ed by Update into the CIC-
CAGE-ID-CORR-LINK-UQ field of
the CAGE-ID-CORRECTION file.

2 EX ANIMAL-CORRECTION-KEY A 26.0 N DE

Key for animal allocation,
animal transfer, animal remov-
al, histology number, not ex-
amined, and uncertain cause-
of-death correction records.
The elements of the key are:
test number (7), cage (4), and
animal number (5), CORRECTION-
DATE, and CORRECTION-TIME. If
the non-CORRECTION elements

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

are being changed, Update will load the field into the CORR-LINK field appropriate for the segment.

2	EY		TEST-NOTE-CORRECTION-KEY	A	27.0	N	DE
---	----	--	--------------------------	---	------	---	----

Key for test note correction records. The elements of the key are: test number (7), note date (5), note time (5), CORRECTION-DATE, and CORRECTION TIME. If the non-CORRECTION elements are being changed, Update will load the load the TNC-TEST-NOTE-CORR-KEY-LINK-UQ field of the TEST-NOTE-CORRECTION file. segment.

2	EZ		TREAT-NOTE-CORRECTION-KEY	A	30.0	N	DE
---	----	--	---------------------------	---	------	---	----

Key for treatment group note corrections. Elements of key are: test number (7), treatment group ID (3), note date (5), note time (5), CORRECTION-DATE, and CORRECTION TIME. If the non-CORRECTION elements are being changed, Update will load the key into the TGC-TREAT-NOTE-CORR-LINK-KEY-UQ field of the TREAT-NOTE-CORRECTION file.

2	FA		CAGE-DATA-CORRECTION-KEY	A	31.0	N	DE
---	----	--	--------------------------	---	------	---	----

Key for feeder weight, bottle weight, and cage note corrections. The elements of the key are: test number (7), cage number (4), observation date (5), observation time (5), CORRECTION-DATE, and CORRECTION TIME. If the non-CORRECTION fields are being changed, then Update will load the key into the CORR-KEY-LINK field

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	L	E	N	G	S	D	E
---	---	----	------	---	---	---	---	---	---	---	---

appropriate for the segment.

2	FB	ORGAN-CORRECTION-KEY	A	32.0	N	D	E
---	----	----------------------	---	------	---	---	---

Key for organ status and site status correction segments.

The elements of the key are:

test number (7), cage number (4), animal number (5), organ code (5), pathologist code (1), CORRECTION-DATE, and CORRECTION-TIME. If any of the non-CORRECTION elements are being changed, then Update loads the key into the appropriate CORR-KEY-LINK field for the segment.

2	FC	ANIMAL-DATA-CORRECTION-KEY	A	36.0	N	D	E
---	----	----------------------------	---	------	---	---	---

The key for the animal observation, animal weight, and animal note correction segments. The elements of the key are: test number (7), cage number (4), animal number (5), observation date (5), observation time (5), CORRECTION-DATE, and CORRECTION-TIME.

If any of the non-CORRECTION elements are changed, then Update loads the key into the CORR-KEY-LINK field of the appropriate correction segment.

2	FD	PATH-NOTE-CORRECTION-KEY	A	37.0	N	D	E
---	----	--------------------------	---	------	---	---	---

Key for pathology animal note corrections. Elements of the key are: test number (7), cage number (4), animal number (5), pathologist code (1), note

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T	L	DB	NAME	F	LENG	S	DE
---	---	----	------	---	------	---	----

date (5), note time (5), CORRECTION-DATE, and CORRECTION-TIME. If the non-CORRECTION elements are changed, then Update loads the key into the MAC-MICRO-ANI-NTE-COR-KEY-LNK-UQ field of the ANIMAL-MICRO-NOTE-CORRECTION file.

2	FE	ORGAN-DATA-CORRECTION-KEY	A	42.0	N	DE
---	----	---------------------------	---	------	---	----

Key used for micro observation and organ note corrections. Elements of the field are: test number (7), cage number (4), animal number (5), organ code (5), pathologist code (1), observation date (5), observation time (5), CORRECTION-DATE (5), and CORRECTION-TIME. If any of the non-CORRECTION elements are being changed, then Update will load the key into the CORRECTION-KEY-LINK field of the appropriate correction segment.

1	EW	RECEIVE-UPDATE-KEY	A	54.0	N	DE
---	----	--------------------	---	------	---	----

Unique key for database transactions. Elements of key are: SEGMENT-ID, MICRO-OR-ECS-FLAG, applicable segment key (27), observation date (5), observation time (5), receive date (5), receive time (5), and sequence number (4).

1	GA	ANY-ANIMAL-KEY	A	9.0	N	SP
---	----	----------------	---	-----	---	----

SOURCE FIELD(S) --- -START- --END-						
DATA-BASE-KEY	3	9				
ANIMAL-NUMBER	1	2				

Super-key whose elements are the test number (7) and the animal (B2).

FILE...: DATABASE-TRANSACTION

TYPE...: ADABAS

FILE-NR: 86

PRIMARY SEQUENCE FIELD: RECEIVE-UPDATE-KEY

T L DB NAME	F	LENG S DE
1 GB BAL-CAL-SUPERKEY	A	37.0 N SP
SOURCE FIELD(S) --- -START- --END-		
DATA-BASE-KEY	3	9
BALANCE-NUMBER	1	2
TOX-DATA-DATE	1	2
TOX-DATA-TIME	1	2
RECEIVE-UPDATE-KEY	31	54

Super-key used to identify
whether there is an existing
record to which the calibra-
should be appended.

DESTINATION-TABLE

Physical layout of file 35. Supports the DESTINATION_TBLE view.
No TDMS software uses physical layout to issue direct calls.

FILE...: DESTINATION-TABLE

TYPE...: ADABAS

FILE-NR: 35

PRIMARY SEQUENCE FIELD: TABLE-DESTINATION

T	L	DB	NAME	F	LENG	S	DE
1	AA	TABLE-DEST-TYPE		A	1.0	F	
1	AB	TABLE-DESTINATION		A	10.0	N	DE
1	AC	TABLE-FACILITY		U	5.0	N	UQ

DOSE-ROUTE-TABLE

Physical layout of Dose Route table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBDOSRTE. TBGET uses that layout description to form direct calls to retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBDOSRTE.

FILE...: DOSE-ROUTE-TABLE

TYPE...: ADABAS

FILE-NR: 18

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
G	1	DT	DOSE-ROUTE-TABLE				
2	AA	ROW-NUMBER		B	2.0	UQ	
			Refer to TBDOSRTE. Not referenced by TBGET.				
2	BB	DATE-ENTERED		A	8.0		
			Refer to TBDOSRTE. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE		A	1.0	F	
			Refer to TBDOSRTE. Not referenced by TBGET.				
2	DD	CODE-VALUE		B	2.0	UQ	
			Refer to TBDOSRTE. Key by which TBGET retrieves data from this table.				
2	EE	EE-FILLER		A	1.0	F	
			Refer to TBDOSRTE. Not referenced by TBGET.				
2	FF	FF-FILLER		A	1.0	F	
			Refer to TBDOSRTE. Not referenced by TBGET.				
2	GG	LONG-TEXT		A	32.0	N	
			Refer to TBDOSRTE. Used by TBGET for EIS/PEIS and PAS validation reports.				
2	HH	SHORT-TEXT		A	8.0	N	
			Refer to TBDOSRTE. Not referenced by TBGET.				
2	II	GENERAL-ROUTE		B	2.0	N	
			Refer to TBDOSRTE. Used by TBGET for PAS validation reports.				

FACILITY-CHEMICAL-EIS

Physical layout of file 59. Supports the PAS_FAC_CHEM view.
 PASGET uses this layout to issue direct calls for PAS validation and
 EIS/PEIS reports.

FILE...: FACILITY-CHEMICAL-EIS

TYPE...: ADABAS

FILE-NR: 59

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
G 1	TA		PAS-FAC-CHEM				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Not referenced by PASGET.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		STR-FAC-CHEM-SEG				
			Not referenced by PASGET.				
2	BB		TEST-TYPE	B	2.0		N DE
			Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	CC		DOSE-RT	B	2.0		N
			Retrieved by PASGET for PAS validation reports and EIS/				

FILE...: FACILITY-CHEMICAL-EIS

TYPE...: ADABAS

FILE-NR: 59

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
		PEIS reports.				
2	DD	FACILITY-ID	B	2.0	N	DE
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	EE	FACILITY-NAME	A	54.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	FF	PRIN-INV	A	32.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	GG	CONTRACT-1	A	12.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	HH	CONTRACT-2	A	12.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
M 2	II	HAZARD-TXT	A	64.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	JJ	TEST-COMPOUND-CT	B	2.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
P 1	KK	ANY-TEST-COMPOUND				
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	LL	NTP	A	6.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				

FILE...: FACILITY-CHEMICAL-EIS

TYPE...: ADABAS

FILE-NR: 59

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
2 MM CAS Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	A	12.0 N
2 NN COMPOUND Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	A	54.0 N DE
1 OO LOCK-DATE Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	P	3.0 N
1 PP LOCK-ENTRY-DATE Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	P	3.0 N
1 QQ LOCK-OPER-ID Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	P	3.0 N
1 RR ANML-SUPPLR Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	P	3.0 N
1 SS PROTO-REQ-TIS-FLAG Not referenced by PASGET.	A	1.0 N

FACILITY-TABLE

Physical layout of facility table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBFACLT. TBGET uses that layout description to form direct calls to retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBFACLT.

FILE...: FACILITY-TABLE

TYPE...: ADABAS

FILE-NR: 27

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	FT FACILITIES-TABLE				
2	AA	ROW-NUMBER	B	2.0	UQ	
		Refer to TBFACLT. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBFACLT. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Refer to TBFACLT. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0	UQ	
		Refer to TBFACLT. Key used by TBGET to retrieve data from table.				
2	EE	EE-FILLER	A	1.0	F	
		Refer to TBFACLT. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0	F	
		Refer to TBFACLT. Not referenced by TBGET.				
2	GG	CHAR-ENTRY-CODE	A	5.0	N	
		Refer to TBFACLT. Not referenced by TBGET.				
2	HH	LONG-TEXT	A	50.0	N	DE
		Refer to TBFACLT. Retrieved by TBGET for PAS validation reports.				
2	II	CITY-STATE	A	25.0	N	
		Refer to TBFACLT. Not referenced by TBGET.				

FILE...: FACILITY-TABLE

TYPE...: ADABAS

FILE-NR: 27

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
2 JJ SHORT-TEXT	A	23.0 N
Refer to TBFACTY. Retrieved by TBGET for EIS/PEIS reports.		
1 KK LAB-CODE	A	2.0 N DE
Refer to TBFACTY. Not refer- enced by TBGET.		

IN-LIFE-TERMS-TABLE

Physical layout of in-life terms table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBINLIFE. TBGET can use that layout description to form direct calls to retrieve data from the table, but does not. This layout corresponds to the view TBINLIFETERMS.

FILE...: IN-LIFE-TERMS-TABLE

TYPE...: ADABAS

FILE-NR: 33

PRIMARY SEQUENCE FIELD: TERM-NUM

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	GL	IN-LIFE-TERMS-TABLE				
			Refer to TBINLIFETERMS. Not referenced by TBGET.				
2	GM	TERM-NUM		P	3.0	N	UQ
			Refer to TBINLIFETERMS. Not referenced by TBGET.				
2	GN	DEFINE-TEXT		A	200.0	N	
			Refer to TBINLIFETERMS.. Not referenced by TBGET.				

LOAD-WEIGHT-EIS

Physical layout of file 47. Supports the PAS_LOAD_WT view. This is the layout by which PASGET selects data for PAS validation reports.

FILE...: LOAD-WEIGHT-EIS

TYPE...: ADABAS

FILE-NR: 47

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA		PAS-LOAD-WT				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from this file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		STR-LOAD-WT-SEG				
			Not referenced by PASGET.				
2	BB		SSSC-ID	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				
2	CC		SPECIES-CODE	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				

FILE...: LOAD-WEIGHT-EIS

TYPE...: ADABAS

FILE-NR: 47

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
2	DD	STRAIN-CODE	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	EE	SUBSTRAIN-CODE	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	FF	CLASS-CODE	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	GG	SEX-OF-ANIMAL	A	1.0	F	
		Retrieved by PASGET for PAS validation reports.				
2	HH	NO-ANIMALS	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	II	DESIRED-AGE	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	JJ	HI-ALLOCATION-WEIGHT	B	4.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	KK	LO-ALLOCATION-WEIGHT	B	4.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	LL	WT-DEVIATION	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	MM	MAX-TEST-WT	B	4.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	NN	MIN-TEST-WT	B	4.0	N	
		Retrieved by PASGET for PAS validation reports.				
1	SD	SUPER-KEY	A	11.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY	1	9		
		SSSC-ID	1	2		

FILE...: LOAD-WEIGHT-EIS
TYPE...: ADABAS
FILE-NR: 47
PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F LENG S DE
Used by PASGET to select data from the file.	

MASTER

Physical layout of file 40. Supports the PAS_MASTER view. This is the layout by which PASGET selects data for PAS validation reports and EIS/PEIS reports.

FILE...: MASTER

TYPE...: ADABAS

FILE-NR: 40

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	TA	PAS-MASTER				
Not referenced by PASGET.							
G	2	TB	PAS-KEY				
Not referenced by PASGET.							
	3	TC	DATA-BASE-KEY	A	9.0		DE
Selected and retrieved by PAS-GET for PAS validation and EIS/PEIS reports.							
P	1	TD	DOWNLOAD-INFO				
Not referenced by PASGET.							
	2	TE	DOWNLOAD-DATE	A	8.0		N
Not referenced by PASGET.							
	2	TF	DOWNLOAD-TIME	A	6.0		N
Not referenced by PASGET.							

MICRO-PROC-PEIS

Physical layout of file 56. Supports the PAS_MICRO_PROC view. No TDMS software uses this layout to issue direct calls.

FILE...: MICRO-PROC-PEIS

TYPE...: ADABAS

FILE-NR: 56

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
G	1	TA	PAS-MICRO-PROC				
			Not referenced by PASGET.				
G	2	TB	PAS-KEY				
			Not referenced by PASGET.				
	3	TC	DATA-BASE-KEY	A	9.0		DE
			Not referenced by PASGET.				
	2	TE	DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
	2	TF	DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
	2	TG	DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
	2	TH	OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G	2	TK	CORRECTION				
			Not referenced by PASGET.				
	3	TL	ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
	3	TM	REASON	A	78.0		N
			Not referenced by PASGET.				
	3	TN	COR-DATE	A	8.0		N
			Not referenced by PASGET.				
	3	TO	COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G	1	AA	STR-MICRO-PROC-SEG				
			Not referenced by PASGET.				
	2	BB	MICRO-PROC-CT	B	2.0		N
			Not referenced by PASGET.				
M	2	CC	MICRO-PROC-CHAR	A	65.0		N
			Not referenced by PASGET.				

MICRO-TABLE-PEIS

Physical layout of file 58. Supports the PAS_MCRO_PTH_TBL view.
PASGET uses this layout to issue direct calls for PAS validation reports.

FILE...: MICRO-TABLE-PEIS

TYPE...: ADABAS

FILE-NR: 58

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA		PAS-MICRO-PATH-TABLES				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from the file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		STR-MICRO-PATH-TABLES-SEQ				
			Not referenced by PASGET.				
2	BB		P-STAIN-CT	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				
2	CC		P-COD-COM-CT	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				

FILE...: MICRO-TABLE-PEIS

TYPE...: ADABAS

FILE-NR: 58

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
2 DD P-QUALIFIER-CT Retrieved by PASGET for PAS validation report.	B	2.0 N
M 2 EE STAIN-CODE Retrieved by PASGET for PAS validation reports.	B	2.0 N
M 2 FF COD-COM-CODE Retrieved by PASGET for PAS validation reports.	B	2.0 N
M 2 GG QUALIFIER-CODE Retrieved by PASGET for PAS validation reports.	B	2.0 N

MODCOMP-CORRECTION-DATA

Physical layout of file 109. Supports the TD-CORRECTION-REASON view. No TDMS software accesses this file using direct calls.

FILE...: MODCOMP-CORRECTION-DATA

TYPE...: ADABAS

FILE-NR: 109

PRIMARY SEQUENCE FIELD: CORRECTION-KEY

T	L	DB	NAME	F	LENG	S	DE
1	AA		CORRECTION-KEY	A	29.0		DE
1	AC		CORRECTION-SEQUENCE-NUMBER	U	2.0		
M 1	AD		CORRECTION-REASON	A	80.0	N	
1	AE		CORRECTION-DATE	U	5.0		
1	AF		CORRECTION-TIME	U	5.0		
1	AG		CORRECTION-OPERATOR-ID	U	5.0		
1	AH		EXP-TEST-CAGE	A	11.0		SB
SOURCE FIELD(S) --- -START- --END-							
			CORRECTION-KEY		3		13

MORPHOLOGIES-PEIS

Physical layout of file 53. Supports the PAS_MORPHOLOGIES view.
 PASGET uses this layout to issue direct calls for PAS validation reports.

FILE...: MORPHOLOGIES-PEIS

TYPE...: ADABAS

FILE-NR: 53

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA		PAS-MORPHOLOGIES				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from the file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		STR-MORPHOLOGIES-SEG				
			Not referenced by PASGET.				
2	BB		ORGAN-CODE	B	2.0		
			Retrieved by PASGET for PAS validation reports.				
2	CC		NON-NEO-MOPRH-CT	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				

FILE...: MORPHOLOGIES-PEIS

TYPE...: ADABAS

FILE-NR: 53

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
2	DD	NEO-MORPH-CT	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
P 1	EE	ANY-NON-NEOPLASTIC				
		Retrieved by PASGET for PAS validation reports.				
2	FF	MORPHOLOGY-CODE-NON	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	GG	QUALIFIER-FLAG-NON	A	1.0	F	
		Retrieved by PASGET for PAS validation reports.				
P 1	HH	ANY-NEOPLASTIC-GRP				
		Retrieved by PASGET for PAS validation reports.				
2	II	MORPHOLOGY-CODE	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
2	JJ	QUALIFIER-FLAG	A	1.0	F	
		Retrieved by PASGET for PAS validation reports.				
1	SD	SUPER-KEY	A	11.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY		1	9	
		ORGAN-CODE		1	2	

Used by PASGET to select data from the file.

OPERATOR-ID-TABLE

Physical layout for operator table. Supports TB_OPERATOR_ID user view.

FILE...: OPERATOR-ID-TABLE

TYPE...: ADABAS

FILE-NR: 14

PRIMARY SEQUENCE FIELD: OPERATOR-ID

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
1	AA		OPERATOR-ID	U	9.0		UQ
1	AC		OP-ID-PASSWORD	A	8.0	N	
1	AE		OP-ID-FAC-ACCESS	B	2.0	N	
1	AF		OP-ID-AGENCY	B	2.0	N	
1	AG		OP-ID-STARUSER	A	8.0	N	DE
P 1	AH		OP-ID-USER-INFO				
2	AI		OP-ID-USER-START-DATE	A	8.0	N	
2	AJ		OP-ID-USER-STOP-DATE	A	8.0	N	
2	AK		OP-ID-USER-NAME	A	32.0	N	DE

ORGAN-DATA-PEIS

Physical layout of file 108. Supports the TD-ORGAN-DATA view. No current TDMS software accesses this file using direct calls.

FILE...: ORGAN-DATA-PEIS

TYPE...: ADABAS

FILE-NR: 108

PRIMARY SEQUENCE FIELD: EXP-TEST-ANIMAL-ORGAN-DATE

T	L	DB	NAME	F	LENG	S	DE
1	AA		EXP-TEST-ANIMAL-ORGAN-DATE	A	27.0		DE
1	AB		MICRO-OBS-TIME	U	5.0		N
1	AC		MICRO-OBS-OPERATOR-ID	U	5.0		N
1	AD		MICRO-OBS-ERROR-FLAG	A	1.0		F
1	AE		MICRO-OBS-PRIMARY-SITE	U	5.0		N
1	AF		MICRO-OBS-MORPHOLOGY	U	5.0		N
1	AV		MICRO-OBS-MORPH-GENERIC-TYPE	U	5.0		N
1	AW		MICRO-OBS-NEOPLASTIC-FLAG	A	1.0		F
M 1	AG		MICRO-OBS-QUALIFIER	U	5.0		N
1	AX		MICRO-OBS-TRACE-LESION-NUMBER	U	5.0		N
M 1	AH		MICRO-OBS-SITE	U	5.0		N
1	AI		ORGAN-NOTE-OPERATOR-ID	U	5.0		N
1	AR		ORGAN-NOTE-TIME	U	5.0		N
1	AJ		ORGAN-NOTE-ERROR-FLAG	A	1.0		F
M 1	AK		ORGAN-NOTE	A	80.0		N
1	AL		ORGAN-STAIN-OPERATOR-ID	U	5.0		N
1	AS		ORGAN-STAIN-TIME	U	5.0		N
1	AM		ORGAN-STAIN-ERROR-FLAG	A	1.0		F
P 1	AT		ORGAN-STAIN-INFORMATION				
2	AN		ORGAN-STAIN	U	5.0		N
2	AU		ORGAN-STAIN-CONTRIBUTORY-FLAG	A	1.0		N
1	AO		MICRO-OBS-PATHOLOGIST-ID	U	5.0		N
1	AP		ORGAN-NOTE-PATHOLOGIST-ID	U	5.0		N
1	AQ		ORGAN-STAIN-PATHOLOGIST-ID	U	5.0		N
1	BA		MICRO-OBS-COD-FLAG	A	1.0		F
1	UB		EXP-TEST-ORGAN	A	12.0		SP
SOURCE FIELD(S) --- -START- --END-							
			EXP-TEST-ANIMAL-ORGAN-DATE		1		7
			EXP-TEST-ANIMAL-ORGAN-DATE		17		21

ORGAN-PEIS

Physical layout of file 111. Supports the TD-ORGAN view. No TDMS software accesses this file using direct calls.

FILE...: ORGAN-PEIS

TYPE...: ADABAS

FILE-NR: 111

PRIMARY SEQUENCE FIELD: EXP-TEST-ANIMAL-ORGAN-PATHT

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
1	AA	EXP-TEST-ANIMAL-ORGAN-PATHT	A	22.0	UQ		
1	AB	ORGAN-STATUS-DATE	U	5.0	N		
1	AC	ORGAN-STATUS-TIME	U	5.0	N		
1	AD	ORGAN-STATUS-OPERATOR-ID	U	5.0	N		
1	AE	ORGAN-STATUS-ERROR-FLAG	A	1.0	F		
1	AF	ORGAN-STATUS	U	5.0	N		
1	AG	ORGAN-STATUS-PATHOLOGIST-ID	U	5.0	N		
P 1	AH	SITE-STATUS-INFORMATION					
2	AK	SITE-STATUS-DATE	U	5.0	N		
2	AL	SITE-STATUS-TIME	U	5.0	N		
2	AM	SITE-STATUS-PATHOLOGIST-ID	U	5.0	N		
2	AN	SITE-STATUS-OPERATOR-ID	U	5.0	N		
2	AO	SITE-STATUS-ERROR-FLAG	A	1.0	N		
2	AI	SITE-CODE	U	5.0	N		
2	AJ	SITE-STATUS	U	5.0	N		

ORGAN-STATUS-TABLE

Physical layout of organ (tissue) status table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBORGSTA. TBGET uses that layout description to form direct calls to retrieve data from the table for pathology reports. This layout corresponds to the view TBORGSTA.

FILE...: ORGAN-STATUS-TABLE

TYPE...: ADABAS

FILE-NR: 34

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	ST ORGAN-STATUS-TABLE				
		Refer to TBORGSTA. Not referenced by TBGET.				
2	AA	ROW-NUMBER	B	2.0	UQ	
		Refer to TBORGSTA. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0	N	
		Refer to TBORGSTA. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Refer to TBORGSTA. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0	UQ	
		Refer to TBORGSTA. Key used by TBGET to retrieve data from the table.				
2	EE	EE-FILLER	A	1.0	F	
		Refer to TBORGSTA. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0	F	
		Refer to TBORGSTA. Not referenced by TBGET.				
2	GG	SUBSET-INDICATOR	A	1.0		
		Refer to TBORGSTA. Not referenced by TBGET.				
2	HH	LONG-TEXT	A	32.0	N	
		Refer to TBORGSTA. Not referenced by TBGET.				
2	II	DATA-COLLECTION	A	1.0	DE	
		Refer to TBORGSTA. Not referenced by TBGET.				

FILE...: ORGAN-STATUS-TABLE

TYPE...: ADABAS

FILE-NR: 34

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
-------------	---	-----------

2 JJ REPORTING	A	1.0
----------------	---	-----

Refer to TBORGSTA. Retrieved
by TBGET for use in pathology
reports.

ORGANS-PEIS

Physical layout of file 52. Supports the PAS_ORGANS view. PASGET uses this layout to issue direct calls for PAS validation reports.

FILE...: ORGANS-PEIS

TYPE...: ADABAS

FILE-NR: 52

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA		PAS-ORGANS				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from the file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		STR-ORGANS-SEG				
			Not referenced by PASGET.				
2	BB		SSSC-ID-CT	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				
2	CC		ORGAN-GROUP-CT	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				

FILE...: ORGANS-PEIS

TYPE...: ADABAS

FILE-NR: 52

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
-------------	---	-----------

M 2 DD SSSC-ID	B	2.0
----------------	---	-----

Retrieved by PASGET for PAS
validation reports.

M 2 EE ORGAN-CODE	B	2.0
-------------------	---	-----

Retrieved by PASGET for PAS
validation reports.

1 FF STUDY-ORGAN	A	11.0	SP
------------------	---	------	----

SOURCE FIELD(S) --- -START- --END-

DATA-BASE-KEY	1	9
---------------	---	---

ORGAN-CODE	1	2
------------	---	---

Not used by PASGET.

PATHOLOGY-CODE-TABLE

Physical layout of the pathology code table. This is layout is the basis on which PCTGET retrieves data for PCT and pathology reports. (PCTGET retrieves all fields in the table for pathology code table reports.) This layout corresponds to the view TBPCT.

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T	L	DB	NAME	F	L	E	N	G	S	D	E
-	-	-	-	-	-	-	-	-	-	-	-
G	1	AA	PCT								
			Group embracing all fields								
			in the pathology code table.								
2	AB	TOXTIME		A	8.0	N	D	E			
			The date upon which the record								
			was added, or last update, in								
			TDMS toxdate format. (Values								
			are integers stored in an al-								
			phanumeric field.)								
2	AC	ACTIVITY-STATUS		A	1.0	F					
			Status of record. A for ac-								
			tive; I for inactive. Not								
			checked in PCT reports, but is								
			checked by PAS validation re-								
			ports.								
2	AD	TDMS-CODE		B	2.0	F	U	Q			
			Unique integer identifying								
			a record on the table. The								
			common practice, which is								
			neither enforced by software								
			nor standard operating proce-								
			dures, is to assign codes								
			sequentially, as follows:								
			1- 999 Morphologies								
			15000-15999 Qualifiers								
			25000-25999 Topologies								
2	AE	REORDER		B	2.0	N					
			Legacy field with no use in								
			TDMS. Original use unknown.								
2	AF	TERMINOLOGY-PTR		B	2.0	N					
			Legacy field with no use in								
			TDMS. Original use unknown.								

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
Contains no values.		
2 AG INITIAL-ENTRY-DATE	B	2.0 N DE
Date on which record was first added to the table. Stored in TDMS toxdate format.		
2 AH INITIAL-AUTHORIZATION	B	2.0 N
Integer code denoting the pathologist who authorized the addition of the record. The possible values are:		
1001 - Dr. C. Montgomery		
1002 - Dr. S. Eustis		
1003 - Dr. R. Haley		
These codes are not operator IDs.		
2 AI CHANGE-AUTHORIZATION	B	2.0 N
Code denoting the pathologist who authorizes a change to a record. (See INITIAL-AUTHORIZATION.)		
2 AJ SUBTABLE-QUALIFIER	A	1.0 N DE
Letter indicating to which of three subtables a record belongs: M for morphologies, T for topographies, Q for qualifiers.		
2 AK SYNONYM	B	2.0 N
Legacy field with no use in TDMS. Original use unknown. Contains no values.		
G 2 AL ALL-HIERARCHY-LEVELS		
Group embracing all LEVEL fields.		
3 AM LEVEL1	B	2.0 N DE
Morphologies - TDMS code for the generic morphology		

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
<p>type.</p> <p>Qualifiers - ID of the qualifier group to which the qualifier belongs.</p>		
3 AN LEVEL2	B	2.0 N DE
<p>Morph. - Neoplastic flag. 4 for neoplastic; 5 for non-neoplastic.</p> <p>Qual. - Integer indicating sort sequence on pathology reports and the qualifier pull-down window on LDAS.</p>		
3 AO LEVEL3	B	2.0 N
<p>Morph. - Malignancy flag: 4 is malignant; 5 is benign; 6 is uncertain</p>		
3 AP LEVEL4	B	2.0 N
<p>Morph. - Metastatic flag: 4 is metastatic; 5 is non-metastatic; 6 is uncertain</p>		
3 AQ LEVEL5	B	2.0 N
<p>Legacy field. Original use unknown. Morphology terms inherited from NCTR contain the values 1 or 2.</p>		
3 AR LEVEL6	B	2.0 N
<p>Legacy field. Original use unknown. Contains no values.</p>		
3 AS LEVEL7	B	2.0 N
<p>Legacy field. Original use unknown. Contains no values.</p>		
3 AT LEVEL8	B	2.0 N
<p>Legacy field. Original use unknown. Contains no values.</p>		
G 2 AU ALL-VALIDATION-FLAGS		
<p>Group embracing all FLAG fields.</p>		
3 AV FLAG1	N	1.0 F

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME

F LENG S DE

Morph. & Top. - Gender specificity. 1 is both genders, 2 is gender dependant.

(Note: For topographies, only organs are flagged as gender-dependent. Sites are gender-dependent by being located on a gender-dependent organ.)

3 AW FLAG2

N 1.0 F

Morph. & Top. - Species specificity. 1 is all species, 2 is gender dependant.

(Note: For topographies, only organs are flagged as species-dependent. Sites are species-dependent by being located on a species-dependent organ.)

3 AX FLAG3

N 1.0 F

Morph. - Organ association flag. 1 for all organs; 2 for a specified list of organs; 3 for all but a specified list of organ. (The list of organs for values 2 and 3 are found in the multiple-entry field ALL-FLAG3-VALUES.)

3 AY FLAG4

N 1.0 F

Top. - Topography type: 1 is

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
a sytem, 2 is an or- gan, 3 is a site Qual. - Qualifier type: 1 is qualifier group, 2 is a qualifier		
3 AZ FLAG5 Legacy field. Original use unknown. Contains no values.	N	1.0 F
3 BA FLAG6 Gross morphology flag. Gross morphologies are a legacy to TDMS, with no use in the cur- rent version.	N	1.0 F
3 BB FLAG7 Legacy field. Original use unknown. Contains no values.	N	1.0 F
3 BC FLAG8 Legacy field. Original use unknown. Contains no values.	N	1.0 F
2 BD REVERSE-PTR Legacy field. Original use unknown. Contains no values.	B	2.0 N
2 BE NCTR-CODE Legacy field with no use in current version of TDMS. Ap- pears to be a code that was associated with morphology terms inherited from NCTR.	A	3.0 N
2 BF SNOP-CODE Legacy field. Original use unknown.	A	4.0 N
G 2 BG SNOMED Group embracing SNOMED defini- tions. SNOMED appears to be an alternative coding system for pathology terms. Infor- mational only. Only assigned to terms inherited from NCTR. Not entered for new terms.		
3 BH SNOMED-TYPE The single-letter identifying the SNOMED type. M for morph-	A	1.0 N

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T	L	DB NAME	F	LENG	S	DE
		ology; T for topography.				
3	BI	SNOMED-CODE	A	5.0	N	
		Five-digit numeric code that uniquely identifies the term in the SNOMED system.				
2	BJ	ASSOCIATED-SNOMED-TERM	A	6.0	N	
		Legacy field. Original use unknown.				
2	BK	SYSTEMIC-FLAG	A	1.0	N	
		Applies to morphology terms only. Indicates whether a term is a systemic lesion.				
2	BL	GROSS-MICRO-FLAG	A	1.0	N	
		Flag indicating whether a morphology term is microscopically, or grossly, observed. The possible values are:				
		0 - Microscopic only				
		1 - Gross only				
		2 - Both				
		(Gross pathology is not part of the current version of TDMS.)				
2	BM	PAIRED-ORGAN-FLAG	A	1.0	N	
		Applies to organ topographies only. Indicates whether an organ is paired. The values are: 1 is not paired, 2 is paired.				
2	BN	UNKNOWN-FIELD	A	13.0	N	
		Legacy field. Original use unknown. Contains no values.				
G 2	BO	SHORT-TEXT				
		Group which merely includes the short text field.				

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
3 BP SHORT-TEXT-FIELD	A	14.0 N DE
Abbreviated description of the term. Used on Individual Animal Data report and LDAS pull-down screens.		
G 2 BQ LONG-TEXT		
Group which merely includes the long text field.		
3 BR LONG-TEXT-FIELD	A	60.0 N
Full description of the term.		
2 BS NUM-FLAG1-VALUES	B	2.0 N
Number of entries in ALL-FLAG1-VALUES multiple-entry field. Should have same value as MU count field for ALL-FLAG1-VALUES. Since ALL-FLAG1-VALUES can never have more than one entry, the value of this field should not be higher than one.		
2 BT NUM-FLAG2-VALUES	B	2.0 N
Number of entries in ALL-FLAG2-VALUES multiple-entry field. Should have same value as MU count field for ALL-FLAG2-VALUES. Since ALL-FLAG2-VALUES can never have more than one entry, the value of this field should not be higher than one.		
2 BU NUM-FLAG3-VALUES	B	2.0 N
Number of entries in ALL-FLAG3-VALUES multiple-entry field. Should have same value as MU count field for ALL-FLAG3-VALUES.		
2 BV NUM-FLAG4-VALUES	B	2.0 N
Number of entries in ALL-FLAG4-VALUES multiple-entry field. Should have same value as MU count field for ALL-FLAG4-VALUES.		

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
2 BW NUM-FLAG5-VALUES Number of entries in ALL-FLAG5-VALUES multiple- entry field. Should have same value as MU count field for ALL-FLAG5-VALUES.	B	2.0 N
2 BX NUM-FLAG6-VALUES Number of entries in ALL-FLAG6-VALUES multiple- entry field. Should have same value as MU count field for ALL-FLAG6-VALUES.	B	2.0 N
2 BY NUM-FLAG7-VALUES Number of entries in ALL-FLAG7-VALUES multiple- entry field. Should have same value as MU count field for ALL-FLAG7-VALUES.	B	2.0 N
2 BZ NUM-FLAG8-VALUES Number of entries in ALL-FLAG8-VALUES multiple- entry field. Should have same value as MU count field for ALL-FLAG8-VALUES.	B	2.0 N
2 CJ NUM-OF-FREE-TEXT-CHARS Number of entries in multiple- entry field FREE-TEXT-CHARS, which field is a legacy to TDMS. Number is always 0.	B	2.0 N
M 2 CA ALL-FLAG1-VALUES Gender-dependency for morph- ologies and topographies. If FLAG1 is set to 2, then the first entry in this multiple-entry field will be either M for F to indicate the gender upon which the term is dependent. There is never more than one entry in this MU field.	A	1.0 N
M 2 CB ALL-FLAG2-VALUES Indicate species dependency	B	2.0 N

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME

F LENG S DE

for morphologies and topographies. If FLAG2 is set to 2, then the first entry in this MU field will be the TBSTRAIN species code upon which the term is dependent. No more than one entry is used in this MU field.

M 2 CC ALL-FLAG3-VALUES

B 2.0 N DE

Morph. - If FLAG3 is set to 2 or 3, the organ codes on which the term is dependent.

Top. - If FLAG4 is set to 2, the code for the system to which the organ belongs. If FLAG4 is set to 3, then the codes for the organs for which the term is a site.

Qual. - If FLAG4 is 2, then the codes for the groups to which the qualifier belongs.

M 2 CD ALL-FLAG4-VALUES

B 2.0 N

Top. - If FLAG4 is 1, then the values are the codes for the organs included in the system. If FLAG4 is 2, then the values are the codes for the sites on the organ.

Qual. - If FLAG4 is 1, then the values are the codes for the qualifiers included in

FILE...: PATHOLOGY-CODE-TABLE

TYPE...: ADABAS

FILE-NR: 30

PRIMARY SEQUENCE FIELD: TDMS-CODE

T L DB NAME	F	LENG S DE
the group.		
M 2 CE ALL-FLAG5-VALUES	B	2.0 N DE
Morph. - The IDs for the qual- ifier groups that are appropriate to the term.		
Qual. - If FLAG4 is 2, the the IDs for the qualifier groups in which the qualifier is classified. (The LEVEL1 values for the terms in ALL-FLAG3- VALUES.)		
M 2 CF ALL-FLAG6-VALUES	B	2.0 N
Legacy field. Equivalent to ALL-FLAG3-VALUES, but for gross terms. Gross pathology is not part of the current version of TDMS.		
M 2 CG ALL-FLAG7-VALUES	B	2.0 N
Legacy field. Equivalent to ALL-FLAG4-VALUES, but for gross terms. Gross pathology is not part of the current version of TDMS.		
M 2 CH ALL-FLAG8-VALUES	B	2.0 N
Legacy field. Equivalent to ALL-FLAG5-VALUES, but for gross terms. Gross pathology is not part of the current version of TDMS.		
M 2 CI FREE-TEXT-CHARS	A	1.0 N
Legacy field. Original use unknown. Contains no values.		
2 CL UNKNOWN-FIELD-1	A	1.0 N DE
Legacy field. Original use unknown. Contains no values.		

PREPARATION-EIS

Physical layout of file 43. Supports the PAS_PREPARTN_SEG view.
This is the layout by which PASGET selects data for PAS validation reports and EIS/PEIS reports.

FILE...: PREPARATION-EIS

TYPE...: ADABAS

FILE-NR: 43

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
G 1	TA	PAS-PREP	Not referenced by PASGET.				
G 2	TB	PAS-KEY	Not referenced by PASGET.				
3	TC	DATA-BASE-KEY	Used by PASGET to select data from the file.	A	9.0		DE
2	TE	DATA-DATE	Not referenced by PASGET.	A	8.0		N
2	TF	DATA-TIME	Not referenced by PASGET.	A	6.0		N
2	TG	DATA-FACILITY-NUMBER	Not referenced by PASGET.	B	2.0		N
2	TH	OPERATOR-ID	Not referenced by PASGET.	B	2.0		N
G 2	TK	CORRECTION	Not referenced by PASGET.				
3	TL	ORG-OR-COR	Not referenced by PASGET.	A	1.0		F
3	TM	REASON	Not referenced by PASGET.	A	78.0		N
3	TN	COR-DATE	Not referenced by PASGET.	A	8.0		N
3	TO	COR-TIME	Not referenced by PASGET.	A	6.0		N
G 1	AA	STR-PREP-SEG	Not referenced by PASGET.				
2	BB	ID-NO	Retrieved by PASGET for PAS validation reports and EIS/PEIS reports.	P	3.0		N
2	CC	COMPOUND-CT	Retrieved by PASGET for PAS	B	2.0		N

FILE...: PREPARATION-EIS

TYPE...: ADABAS

FILE-NR: 43

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

		validation reports and EIS/ PEIS reports.				
P	1	DD ANY-COMPOUND-INFO Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
	2	EE NTP-NO Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	A	6.0	N	
	2	FF COMPOUND-NAME Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	A	54.0	N	
	2	GG COMP-TYPE Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	A	1.0	F	
	2	HH COMP-PERCENT Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	P	5.0	N	
	2	II MEAS-TYPE Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.	A	2.0	N	
	1	SD SUPER-KEY SOURCE FIELD(S) --- -START- --END- DATA-BASE-KEY ID-NO	A	11.0	SP	
				1		9
				1		2

Used by PASGET to select data
from the file.

PROCEDURE-EIS

Physical layout of file 42. Supports the PAS_PROCEDUR_SEG view.
This is the layout by which PASGET selects data for PAS validation reports and EIS/PEIS reports.

FILE...: PROCEDURE-EIS

TYPE...: ADABAS

FILE-NR: 42

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
G 1	TA	PAS-PROCEDURE	Not referenced by PASGET.				
G 2	TB	PAS-KEY	Not referenced by PASGET.				
3	TC	DATA-BASE-KEY	Used by PASGET to select data from this file.	A	9.0		DE
2	TE	DATA-DATE	Not referenced by PASGET.	A	8.0		N
2	TF	DATA-TIME	Not referenced by PASGET.	A	6.0		N
2	TG	DATA-FACILITY-NUMBER	Not referenced by PASGET.	B	2.0		N
2	TH	OPERATOR-ID	Not referenced by PASGET.	B	2.0		N
G 2	TK	CORRECTION	Not referenced by PASGET.				
3	TL	ORG-OR-COR	Not referenced by PASGET.	A	1.0		F
3	TM	REASON	Not referenced by PASGET.	A	78.0		N
3	TN	COR-DATE	Not referenced by PASGET.	A	8.0		N
3	TO	COR-TIME	Not referenced by PASGET.	A	6.0		N
G 1	AA	STR-PROCEDURE-SEG	Not referenced by PASGET.				
2	BB	ID-NO	Retrieved by PASGET for PAS validation reports.	B	2.0		N
2	CC	PROCEDURE-GRP-CNT	Retrieved by PASGET for PAS validation reports.	B	2.0		N

FILE...: PROCEDURE-EIS

TYPE...: ADABAS

FILE-NR: 42

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

P 1 DD ANY-PROCEDURE-GROUP

Retrieved by PASGET for PAS
validation reports.

2 EE ACTION-CODE B 2.0 N DE

Retrieved by PASGET for PAS
validation reports.

2 FF REL-START I 2.0 F

Retrieved by PASGET for PAS
validation reports.

2 GG REL-STOP I 2.0 F

Retrieved by PASGET for PAS
validation reports.

2 HH WEEKLY-INTERVALS B 2.0 N

Retrieved by PASGET for PAS
validation reports.

2 II REL-OR-SPEC A 1.0

Retrieved by PASGET for PAS
validation reports.

2 JJ DAY-OF-WEEK A 7.0 N

Retrieved by PASGET for PAS
validation reports.

1 SD SUPER-KEY A 11.0 SP

SOURCE FIELD(S) --- -START- --END-

DATA-BASE-KEY 1 9

ID-NO 1 2

Used by PASGET to select data
from the file.

1 SE PROCEDURE-ACTION A 11.0 SP

SOURCE FIELD(S) --- -START- --END-

DATA-BASE-KEY 1 9

ACTION-CODE 1 2

PROTO-REQ-TISS-PEIS

Physical layout of file 84. Supports the PAS_PROTO_REQ_TIS view.
No direct calls are issued based on this layout.

FILE...: PROTO-REQ-TISS-PEIS

TYPE...: ADABAS

FILE-NR: 84

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
G	1	TA	PAS-PROTO-REQ-TIS				
			Not referenced by PASGET.				
G	2	TB	PAS-KEY				
			Not referenced by PASGET.				
	3	TC	DATA-BASE-KEY	A	9.0		DE
			Not referenced by PASGET.				
	2	TE	DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
	2	TF	DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
	2	TG	DATA-FACILITY-NUMBER	U	5.0		N
			Not referenced by PASGET.				
	2	TH	OPERATOR-ID	U	5.0		N
			Not referenced by PASGET.				
G	2	TK	CORRECTION				
			Not referenced by PASGET.				
	3	TL	ORG-OR-COR	A	1.0		N
			Not referenced by PASGET.				
	3	TM	REASON	A	78.0		N
			Not referenced by PASGET.				
	3	TN	COR-DATE	A	8.0		N
			Not referenced by PASGET.				
	3	TO	COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G	1	AA	REQ-TIS-SEQ				
			Not referenced by PASGET.				
	2	BB	TREAT-NUM	U	3.0		
			Not referenced by PASGET.				
	2	CC	BEG-REMOVAL-DATE	U	5.0		
			Not referenced by PASGET.				
	2	DD	END-REMOVAL-DATE	U	5.0		
			Not referenced by PASGET.				
	2	EE	PROTO-REQ-TIS-PATHOLOGIST	U	5.0		N
			Not referenced by PASGET.				

FILE...: PROTO-REQ-TISS-PEIS

TYPE...: ADABAS

FILE-NR: 84

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
M	2	FF REASON-NUM	U	5.0	N	
		Not referenced by PASGET.				
M	2	GG TISSUE-NUM	P	3.0	N	
		Not referenced by PASGET.				
	2	HH ALL-TREAT-FLAG	A	1.0	F	
		Not referenced by PASGET.				
	2	II ALL-REM-REASON-FLAG	A	1.0	F	
		Not referenced by PASGET.				
	1	AB SYSTEM-DATE-TIME	A	13.0	DE	
		Not referenced by PASGET.				
	1	SD PROTO-REQ-TIS-SUPERKEY	A	22.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY		1	9	
		BEG-REMOVAL-DATE		1	5	
		END-REMOVAL-DATE		1	5	
		TREAT-NUM		1	3	
		Not used by PASGET.				
	1	SE AGENCY-EXP-TEST-TREAT-REASON	A	17.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY		1	9	
		TREAT-NUM		1	3	
		REASON-NUM		1	5	
		Not used by PASGET.				

PTH-AST-OPR-PEIS

Physical layout of file 57. Supports the PAS_PH_P_AST_OPR view.
PASGET uses this layout to issue direct calls for PAS validation reports.

FILE...: PTH-AST-OPR-PEIS

TYPE...: ADABAS

FILE-NR: 57

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
-	-	-	-	-	-	-	-
G	1	TA	PAS-PATH-P-ASSIT-OPER				
			Not referenced by PASGET.				
G	2	TB	PAS-KEY				
			Not referenced by PASGET.				
	3	TC	DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from the table.				
	2	TE	DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
	2	TF	DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
	2	TG	DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
	2	TH	OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G	2	TK	CORRECTION				
			Not referenced by PASGET.				
	3	TL	ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
	3	TM	REASON	A	78.0		N
			Not referenced by PASGET.				
	3	TN	COR-DATE	A	8.0		N
			Not referenced by PASGET.				
	3	TO	COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G	1	AA	STR-PATH-P-ASSIT-OPER-SEG				
			Not referenced by PASGET.				
	2	BB	PERSONNEL-ID-CT	B	2.0		N
			Not referenced by PASGET.				
P	1	CC	ANY-ID-GRP				
			Not referenced by PASGET.				
	2	DD	NAME	A	32.0		N
			Retrieved by PASGET for PAS				

FILE...: PTH-AST-OPR-PEIS

TYPE...: ADABAS

FILE-NR: 57

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

validation reports.

2	EE	PERSONNEL-CODE	A	1.0	F	
---	----	----------------	---	-----	---	--

Retrieved by PASGET for PAS

validation reports.

2	FF	PATH-CODE	A	12.0	N	
---	----	-----------	---	------	---	--

Retrieved by PASGET for PAS

validation reports.

1	GG	SUPER-KEY	A	22.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-

		DATA-BASE-KEY	1	9		
--	--	---------------	---	---	--	--

		PERSONNEL-CODE	1	1		
--	--	----------------	---	---	--	--

		PATH-CODE	1	12		
--	--	-----------	---	----	--	--

Not used by PASGET.

RACK-CAGE-EIS

Physical layout of file 51. Supports the PAS_RACK_CAGE view.

PASGET uses this layout to issue direct calls for PAS validation reports.

FILE...: RACK-CAGE-EIS

TYPE...: ADABAS

FILE-NR: 51

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA		PAS-RACK-CAGE				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from the file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		RACK-CAGE-SEQ				
			Not referenced by PASGET.				
2	BB		C-AREA	U	5.0		
			Retrieved by PASGET for PAS validation reports.				
2	CC		RACK-NO	U	2.0		
			Retrieved by PASGET for PAS validation reports.				

FILE...: RACK-CAGE-EIS

TYPE...: ADABAS

FILE-NR: 51

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
2 DD C-SIDE Retrieved by PASGET for PAS validation reports.	A	1.0 F
2 EE RACK-ROW Retrieved by PASGET for PAS validation reports.	U	2.0 N
2 FF RACK-COL Retrieved by PASGET for PAS validation reports.	U	2.0 N
P 1 GG CAGES Retrieved by PASGET for PAS validation reports.		
2 HA COLUMN-1 Retrieved by PASGET for PAS validation reports.	B	2.0 N
2 HB COLUMN-2 Retrieved by PASGET for PAS validation reports.	B	2.0 N
2 HC COLUMN-3 Retrieved by PASGET for PAS validation reports.	B	2.0 N
2 HD COLUMN-4 Retrieved by PASGET for PAS validation reports.	B	2.0 N
2 HE COLUMN-5 Retrieved by PASGET for PAS validation reports.	B	2.0 N
2 HF COLUMN-6 Retrieved by PASGET for PAS validation reports.	B	2.0 N
2 HG COLUMN-7 Retrieved by PASGET for PAS validation reports.	B	2.0 N
2 HH COLUMN-8 Retrieved by PASGET for PAS validation reports.	B	2.0 N
2 HI COLUMN-9	B	2.0 N

FILE...: RACK-CAGE-EIS

TYPE...: ADABAS

FILE-NR: 51

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
		Retrieved by PASGET for PAS validation reports.				
2	HJ	COLUMN-10	B	2.0	N	
		Retrieved by PASGET for PAS validation reports.				
1	SD	SUPER-KEY	A	17.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY		1	9	
		C-AREA		1	5	
		RACK-NO		1	2	
		C-SIDE		1	1	

Used by PASGET to select data
from the file.

RECEIVE-STATUS

Physical layout of file 38. Supports the STATUS_FILE view. No TDMS software accesses the file using direct calls.

FILE...: RECEIVE-STATUS

TYPE...: ADABAS

FILE-NR: 38

PRIMARY SEQUENCE FIELD: STATUS_EXP_AND_TEST_NUMBER

T	L	DB	NAME	F	LENG	S	DE
G	1	AA	TRANSMISSION_STATUS_REC				
G	2	BB	STATUS_TRANSACTION_SET_ID				
G	3	CC	STATUS_DISKETTE_ID				
	4	CD	STATUS_AGENCY_FACILT_DISK_SUBSYS	A	11	0	
G	4	HH	STATUS_EXPERIMENT_ID				
	5	IJ	STATUS_EXP_AND_TEST_NUMBER	U	7	0	DE
	3	KK	STATUS_TRANSACTION_SET_NUM	B	2	0	
	3	LL	STATUS_HEADER_TOX_DATE	A	8	0	
	3	MM	STATUS_HEADER_TOX_TIME	A	6	0	
G	2	NN	STATUS_INFO				
	3	OO	STATUS_RECEIVE_DATE	A	8	0	N
	3	PP	STATUS_RECEIVE_TIME	A	6	0	N
	3	QQ	STATUS_RECEIVE_INDICATOR	A	1	0	F
	3	RR	STATUS_VERIFY_INDICATOR	A	1	0	F
	3	SS	STATUS_CONVERT_INDICATOR	A	1	0	F
	3	TT	STATUS_NOTIFY_INDICATOR	A	1	0	N DE
	3	UU	STATUS_TRANSACTION_COUNT	B	2	0	N
	1	KE	STATUS_KEY	A	34	0	SP
			SOURCE FIELD(S) --- -START- --END-				
			STATUS_AGENCY_FACILT_DISK_SUBSYS	1	11		
			STATUS_EXP_AND_TEST_NUMBER	1	7		
			STATUS_TRANSACTION_SET_NUM	1	2		
			STATUS_HEADER_TOX_DATE	1	8		
			STATUS_HEADER_TOX_TIME	1	6		
	1	EE	STATUS_FACILITY_CODE	A	5	0	SB
			SOURCE FIELD(S) --- -START- --END-				
			STATUS_AGENCY_FACILT_DISK_SUBSYS	3	7		
	1	VV	STATUS_EXP_TEST_TRANS_SET_NUM	B	9	0	SP
			SOURCE FIELD(S) --- -START- --END-				
			STATUS_EXP_AND_TEST_NUMBER	1	7		
			STATUS_TRANSACTION_SET_NUM	1	2		

REMOVAL-ACTION-EIS

Physical layout of file 50. Supports the PAS_REMVL_ACTION view.
PASGET uses this layout to issue direct calls for PAS validation reports.

FILE...: REMOVAL-ACTION-EIS

TYPE...: ADABAS

FILE-NR: 50

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
G	1	TA	PAS-REMOVAL-ACTIONS				
			Not referenced by PASGET.				
G	2	TB	PAS-KEY				
			Not referenced by PASGET.				
	3	TC	DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from the file.				
	2	TE	DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
	2	TF	DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
	2	TG	DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
	2	TH	OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G	2	TK	CORRECTION				
			Not referenced by PASGET.				
	3	TL	ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
	3	TM	REASON	A	78.0		N
			Not referenced by PASGET.				
	3	TN	COR-DATE	A	8.0		N
			Not referenced by PASGET.				
	3	TO	COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G	1	AA	STR-REMOVAL-ACTIONS-SEQ				
			Not referenced by PASGET.				
	2	BB	REMOVAL-ACT-CT	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				
M	2	CC	REMOVAL-ACT-CODE	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				

RESP-PARTY-EIS

Physical layout of file 45. Supports the PAS_RESP_PARTIES view.
This is the layout by which PASGET selects data for PAS validation reports.

FILE...: RESP-PARTY-EIS

TYPE...: ADABAS

FILE-NR: 45

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
G 1	TA		PAS-RESP-PARTIES				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to select data from this file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		STR-RESP-PARTIES-SEG				
			Not referenced by PASGET.				
2	BB		RESP-PARTY	A	32.0		N
			Retrieved by PASGET for PAS validation reports.				
2	CC		JOB-POSITION	A	2.0		N
			Retrieved by PASGET for PAS validation reports.				

FILE...: RESP-PARTY-EIS

TYPE...: ADABAS

FILE-NR: 45

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T L DB NAME	F	LENG S DE
2 DD FACILITY-CODE Retrieved by PASGET for PAS validation reports.	B	2.0 N
2 EE EFFECTIVE-DATE Retrieved by PASGET for PAS validation reports.	U	6.0 N
2 FF PHONE-NUMBER Retrieved by PASGET for PAS validation reports.	U	10.0 N
2 GG EXTENSION Retrieved by PASGET for PAS validation reports.	U	4.0 N
2 HH FACILITY-ADDRESS Retrieved by PASGET for PAS validation reports.	A	128.0 N

SITES-PEIS

Physical layout of file 54. Supports the PAS_SITES view. PASGET uses this layout to issue direct calls for PAS validation reports.

FILE...: SITES-PEIS

TYPE...: ADABAS

FILE-NR: 54

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA		PAS-SITES				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from the file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		STR-SITES-SEG				
			Not referenced by PASGET.				
2	BB		ORGAN-CODE	B	2.0		
			Retrieved by PASGET for PAS validation reports.				
2	CC		SITE-GRP-CT	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				

FILE...: SITES-PEIS

TYPE...: ADABAS

FILE-NR: 54

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

M	2	DD SITE-CODE	B	2.0	N	
---	---	--------------	---	-----	---	--

Retrieved by PASGET for PAS
validation reports.

1	SD	SUPER-KEY	A	11.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-

DATA-BASE-KEY	1	9
---------------	---	---

ORGAN-CODE	1	2
------------	---	---

Used by PASGET to select data
from the file.

STAIN-TABLE

Physical layout of stain table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBSTAINC. TBGET uses that layout description to form direct calls to retrieve data from the table for PCT and PAS validation reports. This layout corresponds to the view TBSTAINC. (Stain data is not collected or reported by the current version of TDMS.)

FILE...: STAIN-TABLE

TYPE...: ADABAS

FILE-NR: 15

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG	S	DE
- - - - -				
G 1 ST STAIN-TABLE				
Not referenced by TBGET.				
2 AA ROW-NUMBER	B	2.0	UQ	
Refer to TBSTAINC. Not referenced by TBGET.				
2 BB DATE-ENTERED	A	8.0		
Refer to TBSTAINC. Not referenced by TBGET.				
2 CC ACTIVE-INACTIVE	A	1.0	F	
Refer to TBSTAINC. Used by TBGET for PCT reports.				
2 DD CODE-VALUE	B	2.0	UQ	
Refer to TBSTAINC. Key by which TBGET retrieves data from this table.				
2 EE EE-FILLER	A	1.0	F	
Refer to TBSTAINC. Not referenced by TBGET.				
2 FF FF-FILLER	A	1.0	F	
Refer to TBSTAINC. Not referenced by TBGET.				
2 GG LONG-TEXT	A	30.0	N	
Refer to TBSTAINC. Used by TBGET for PAS validation reports.				
2 HH SHORT-TEXT	A	14.0	N	
Refer to TBSTAINC. Used by TBGET for PCT reports.				

STRAIN-TABLE

Physical layout of strain table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBSTRAIN. TBGET uses that layout description to form direct calls to retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBSTRAIN.

FILE...: STRAIN-TABLE

TYPE...: ADABAS

FILE-NR: 21

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	ST STRAIN-TABLE				
2	AA	ROW-NUMBER	B	2.0	UQ	
		Refer to TBSTRAIN. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBSTRAIN. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Refer to TBSTRAIN. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0	UQ	
		Refer to TBSTRAIN. Key by which TBGET retrieves data from this table.				
2	EE	EE-FILLER	A	1.0	F	
		Refer to TBSTRAIN. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0	F	
		Refer to TBSTRAIN. Not referenced by TBGET.				
2	GG	SPECIES-CODE	B	2.0	N	
		Refer to TBSTRAIN. Not referenced by TBGET.				
2	HH	SUBSET-INDICATOR	B	2.0	N	
		Refer to TBSTRAIN. Used by TBGET for PAS validation reports.				
2	II	SHORT-TEXT	A	12.0	N	
		Refer to TBSTRAIN. Used by TBGET for EIS/PEIS and PAS validation reports.				

FILE...: STRAIN-TABLE

TYPE...: ADABAS

FILE-NR: 21

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
2 JJ BREEDING-CODE Refer to TBSTRAIN. Not referenced by TBGET.	A	1.0 N
2 KK STR-FEM-PARENT Refer to TBSTRAIN. Not referenced by TBGET.	B	2.0 N
2 LL STR-MALE-PARENT Refer to TBSTRAIN. Not referenced by TBGET.	B	2.0 N
2 MM LONG-TEXT Refer to TBSTRAIN. Used by TBGET for EIS/PEIS reports.	A	38.0 N
2 NN NN-FILLER Refer to TBSTRAIN. Not referenced by TBGET.	A	1.0 N

SUBSTRAIN-TABLE

Physical layout of substrain table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBSUBSTR. TBGET uses that layout description to form direct calls to retrieve data from the table for PAS validation reports. This layout corresponds to the view TBSUBSTR. Substrain data is not used in TDMS. Dummy values are entered in PAS.

FILE...: SUBSTRAIN-TABLE

TYPE...: ADABAS

FILE-NR: 17

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
G	1	ST	SUBSTRAIN-TABLE				
	2	AA	ROW-NUMBER	B	2.0		UQ
			Refer to TBSUBSTR. Not referenced by TBGET.				
	2	BB	DATE-ENTERED	A	8.0		
			Refer to TBSUBSTR. Not referenced by TBGET.				
	2	CC	ACTIVE-INACTIVE	A	1.0		F
			Refer to TBSUBSTR. Not referenced by TBGET.				
	2	DD	CODE-VALUE	B	2.0		UQ
			Refer to TBSUBSTR. Key by which TBGET retrieves data from this table				
	2	EE	EE-FILLER	A	1.0		F
			Refer to TBSUBSTR. Not used by TBGET.				
	2	FF	FF-FILLER	A	1.0		F
			Refer to TBSUBSTR. Not used by TBGET.				
	2	GG	SHORT-TEXT	A	14.0		N
			Refer to TBSUBSTR. Used by TBGET for PAS validation reports.				
	2	HH	LONG-TEXT	A	32.0		N
			Refer to TBSUBSTR. Used by TBGET for PAS validation reports.				

TABLES_DESCRIPTION

Physical layout of Tables Description table. This layout is used by TBGET when it retrieves data from file. File contains physical layouts of TDMS tables used by TBGET.

FILE...: TABLES_DESCRIPTION

TYPE...: ADABAS

FILE-NR: 16

PRIMARY SEQUENCE FIELD: TABLE_NAME

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G	1	MF	TABLES_DESCRIPTION_FILE				
	2	TN	TABLE_NAME	A	8.0	N	DE
			Primary key for table. Contains the eight-character view names for each TDMS table.				
	2	TB	TABLE_NUM	B	1.0	F	
			The file number assigned to the file when loaded into ADABAS with ADALOD or ADAFDU.				
	2	NC	NUM_COLS	B	2.0		
			The number of columns defined in the FIELD_GROUP repeating group.				
P	1	FG	FIELD_GROUP				
			Repeating group containing each field in table. This group is position-sensitive. Do not reorder elements unless you understand the TBGET retrieval mechanism.				
	2	FN	FIELD_NAME	A	2.0	N	
			The short field name assigned to the field when loaded into ADABAS.				
	2	FL	FIELD_LENGTH	B	2.0	N	
			The length of the field.				

TBLDAS

Physical layout of super-table that embraces five tables, those tables being the weight status, treatment role, correction reason, in-life qualifier group, and weight curve tables. Each of these tables is defined in the Tables Description table with the field TABLE_NAME set to TBWGTSTATUS, TBTRTROLE, TBCORRREASON, TBQUALGRP, and TBWGTCURVE, respectively. TBGET uses the descriptions to retrieve data via direct calls. Supports views named as TABLE_NAME.

FILE...: TBLDAS

TYPE...: ADABAS

FILE-NR: 31

PRIMARY SEQUENCE FIELD:

T	L	DB	NAME	F	LENG	S	DE
G	1	AA	WEIGHT-STATUS-TABLE				
			Refer to TBWGTSTATUS. Not referenced by TBGET.				
	2	AB	STATUS-NUM	P	5.0	N	UQ
			Refer to TBWGTSTATUS. Not referenced by TBGET.				
	2	AC	STATUS-TEXT	A	30.0	N	
			Refer to TBWGTSTATUS. Not referenced by TBGET.				
	2	AD	STATUS-SHORT-TEXT	A	8.0	N	
			Refer to TBWGTSTATUS. Not referenced by TBGET.				
	2	AE	STATUS-FLAG	A	1.0	F	
			Refer to TBWGTSTATUS. Not referenced by TBGET.				
G	1	BA	TREATMENT-ROLE-TABLE				
			Refer to TBTRTROLE. Not referenced by TBGET.				
	2	BB	TREAT-ROLE	A	1.0	N	UQ
			Refer to TBTRTROLE. Key used by TBGET to retrieve data for PAS validation reports.				
	2	BC	ROLE-TEXT	A	32.0	N	
			Refer to TBTRTROLE. Retrieved by TBGET for PAS validation reports.				
G	1	CA	CORRECTION-REASON-TABLE				
			Refer to TBCORRREASON. Not				

FILE...: TBLDAS

TYPE...: ADABAS

FILE-NR: 31

PRIMARY SEQUENCE FIELD:

T L DB NAME	F	LENG S DE
referenced by TBGET.		
2 CB CORR-CODE	A	2.0 N UQ
Refer to TBCORREASON. Not referenced by TBGET.		
2 CC CORR-REASON	A	30.0 N
Refer to TBCORREASON. Not referenced by TBGET.		
G 1 GA WEIGHT-CURVE-TABLE		
Refer to TBWEIGHTCURVE. Not referenced by TBGET.		
2 GB WEIGHT-CURVE-START-WEEK	P	5.0 N UQ
Refer to TBWEIGHTCURVE. Not referenced by TBGET.		
2 GD MIN-WEIGHT	P	5.0 N
Refer to TBWEIGHTCURVE. Not referenced by TBGET.		
2 GE MAX-WEIGHT	P	5.0 N
Refer to TBWEIGHTCURVE. Not referenced by TBGET.		
G 1 HA INLIFE-QUALIFIER-GROUP-TABLE		
Refer to TBQUALGRP. Not referenced by TBGET.		
2 HB EIQUAL-GROUP	A	5.0 N UQ
Refer to TBQUALGRP. Not referenced by TBGET.		
2 HC EIQUAL-GROUP-SHORT-TEXT	A	14.0 N
Refer to TBQUALGRP. Not referenced by TBGET.		
2 HD EIQUAL-GROUP-LONG-TEXT	A	32.0 N
Refer to TBQUALGRP. Not referenced by TBGET.		

TD-NOTES

Physical layout of file 114. Supports TD-TEST-NOTE and TD-TREATMENT-NOTE views. No TDMS software accesses this file using direct calls.

FILE...: TD-NOTES

TYPE...: ADABAS

FILE-NR: 114

PRIMARY SEQUENCE FIELD:

T	L	DB	NAME	F	LENG	S	DE
1	AA	EXP-TEST-TRT-DATE-TIME-UQ	A	20.0	N	UQ	
1	AB	EXP-TEST-DATE-TIME-UQ	A	17.0	N	UQ	
1	AC	OPER-NUM	P	3.0	F		
1	AD	NOTE-TEXT	A	160.0	N		
1	BB	TEST	A	7.0	SP		
SOURCE FIELD(S) --- -START- --END-							
		EXP-TEST-TRT-DATE-TIME-UQ		1		7	

TEST

Physical layout for file 112. Supports the TD-TEST view. No TDMS software accesses this file using direct calls.

FILE...: TEST

TYPE...: ADABAS

FILE-NR: 112

PRIMARY SEQUENCE FIELD: EXPERIMENT-TEST

T	L	DB	NAME	F	LENG	S	DE
1	AA	EXPERIMENT-TEST		U	7.0	UQ	
1	AC	AGENCY-NUMBER		U	2.0		
1	AD	TEST-TYPE		U	5.0		

TEST-TYPE-TABLE

Physical layout of test type table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBTESTYP.

TBGET uses that layout description to form direct calls to retrieve data from the table for EIS/PEIS and PAS validation reports. This layout corresponds to the view TBTESTYP.

FILE...: TEST-TYPE-TABLE

TYPE...: ADABAS

FILE-NR: 28

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G 1	TT	TEST-TYPE-TABLE				
2	AA	ROW-NUMBER	B	2.0		UQ
		Refer to TBTESTYP. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBTESTYP. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0		F
		Refer to TBTESTYP. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0		UQ
		Refer to TBTESTYP. Key used by TBGET to retrieve data from table.				
2	EE	EE-FILLER	A	1.0		F
		Refer to TBTESTYP. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0		F
		Refer to TBTESTYP. Not referenced by TBGET.				
2	GG	CHARACTER-CODE	A	2.0		N
		Refer to TBTESTYP. Not referenced by TBGET.				
2	HH	SHORT-TEXT	A	16.0		N
		Refer to TBTESTYP. Retrieved by TBGET for EIS/PEIS reports.				
2	II	LONG-TEXT	A	32.0		N
		Refer to TBTESTYP. Retrieved by TBGET for PAS validation reports.				

FILE...: TEST-TYPE-TABLE

TYPE...: ADABAS

FILE-NR: 28

PRIMARY SEQUENCE FIELD: CODE-VALUE

T L DB NAME	F	LENG S DE
2 JJ CURRENT-PERIOD-LENGTH	P	3.0 N
Refer to TBTESTYP. Not referenced by TBGET.		

TREATMENT-EIS

Physical layout of file 41. Supports the PAS_TREATMNT_SEG view.
This is the layout by which PASGET selects data for PAS validation reports and EIS/PEIS reports.

FILE...: TREATMENT-EIS

TYPE...: ADABAS

FILE-NR: 41

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
G 1	TA	PAS-TREATMENT	Not referenced by PASGET.				
G 2	TB	PAS-KEY	Not referenced by PASGET.				
3	TC	DATA-BASE-KEY	Used by PASGET to select data from the file.	A	9.0		DE
2	TE	DATA-DATE	Not referenced by PASGET.	A	8.0		N
2	TF	DATA-TIME	Not referenced by PASGET.	A	6.0		N
2	TG	DATA-FACILITY-NUMBER	Not referenced by PASGET.	B	2.0		N
2	TH	OPERATOR-ID	Not referenced by PASGET.	B	2.0		N
G 2	TK	CORRECTION	Not referenced by PASGET.				
3	TL	ORG-OR-COR	Not referenced by PASGET.	A	1.0		F
3	TM	REASON	Not referenced by PASGET.	A	78.0		N
3	TN	COR-DATE	Not referenced by PASGET.	A	8.0		N
3	TO	COR-TIME	Not referenced by PASGET.	A	6.0		N
G 1	AA	STR-TREATMENT-SEG	Not referenced by PASGET.				
2	BB	ID-NO	Retrieved by PASGET for PAS validation reports and EIS/PEIS reports.	B	2.0		N
2	CC	CONTROL-DES	Retrieved by PASGET for PAS	A	1.0		F

FILE...: TREATMENT-EIS

TYPE...: ADABAS

FILE-NR: 41

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
		validation reports and EIS/ PEIS reports.				
2	DD	TEXT-ID	A	16.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	EE	REGIMEN-CT	B	2.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
P 1	FF	ANY-REGIMEN				
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	GG	PREP-ID	B	2.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	HH	ROUTE	B	2.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	II	PREP-VOL	P	3.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	JJ	DOSE-CALC	A	1.0		
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	KK	REL-START	B	2.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	LL	REL-STOP	B	2.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				

FILE...: TREATMENT-EIS

TYPE...: ADABAS

FILE-NR: 41

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
		PEIS reports.				
2	MM	INTERVALS	B	2.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	NN	REL-OR-SPEC	A	1.0		
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
2	OO	DAY-OF-WEEK	A	7.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
1	PP	MIN-WATER	P	5.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
1	QQ	MAX-WATER	P	5.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
1	RR	MIN-FOOD	P	5.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
1	SS	MAX-FOOD	P	5.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
1	TT	CURVE-NO	P	5.0	N	
		Retrieved by PASGET for PAS validation reports and EIS/ PEIS reports.				
1	SD	SUPER-KEY	A	11.0	SP	
		SOURCE FIELD(S) --- -START- --END-				
		DATA-BASE-KEY		1	9	
		ID-NO		1	2	
		Used by PASGET to select data from the file.				

TREATMENT-PEIS

Physical layout of file 100. Supports TD-TREATMENT view. No current TDMS software accesses this file using direct calls.

FILE...: TREATMENT-PEIS

TYPE...: ADABAS

FILE-NR: 100

PRIMARY SEQUENCE FIELD: EXP-TEST-TREATMENT

T	L	DB	NAME	F	LENG	S	DE
1	AA	EXP-TEST-TREATMENT		U	10.0	UQ	
1	AB	TREATMENT-CODE		A	1.0	F	
M 1	AD	PROTOCOL-REQUIRED-TISSUE-CODE		U	5.0	N	
1	AE	PATHOLOGIST-IDENTIFICATION		U	5.0		
1	AF	OPERATOR-IDENTIFICATION		U	5.0		
1	AG	DATE		U	5.0		
1	AH	SEQUENCE		U	5.0		
1	AI	ERROR-FLAG		A	1.0	F	

UNITS-TABLE

Physical layout of units table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBUNITS. TBGET can use that layout description to form direct calls to retrieve data from the table, but does not since the table has no meaning in TDMS. This file corresponds to the view TBUNITS.

FILE...: UNITS-TABLE

TYPE...: ADABAS

FILE-NR: 20

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB NAME	F	LENG	S	DE
G	1	UT UNITS-TABLE				
2	AA	ROW-NUMBER	B	2.0	UQ	
		Refer to TBUNITS. Not referenced by TBGET.				
2	BB	DATE-ENTERED	A	8.0		
		Refer to TBUNITS. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE	A	1.0	F	
		Refer to TBUNITS. Not referenced by TBGET.				
2	DD	CODE-VALUE	B	2.0	UQ	
		Refer to TBUNITS. Not used by TBGET.				
2	EE	EE-FILLER	A	1.0	F	
		Refer to TBUNITS. Not referenced by TBGET.				
2	FF	FF-FILLER	A	1.0	F	
		Refer to TBUNITS. Not referenced by TBGET.				
2	GG	SHORT-TEXT	A	8.0	N	
		Refer to TBUNITS. Not used by TBGET.				
2	HH	LONG-TEXT	A	50.0	N	
		Refer to TBUNITS. Not referenced by TBGET.				

WEIGHT-ITEM-EIS

Physical layout of file 46. Supports the PAS_WGT_ITEMS view. This is the layout by which PASGET selects data for PAS validation reports.

FILE...: WEIGHT-ITEM-EIS

TYPE...: ADABAS

FILE-NR: 46

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB	NAME	F	LENG	S	DE
- - - - -							
G 1	TA		PAS-WGT-ITEMS				
			Not referenced by PASGET.				
G 2	TB		PAS-KEY				
			Not referenced by PASGET.				
3	TC		DATA-BASE-KEY	A	9.0		DE
			Used by PASGET to retrieve data from this file.				
2	TE		DATA-DATE	A	8.0		N
			Not referenced by PASGET.				
2	TF		DATA-TIME	A	6.0		N
			Not referenced by PASGET.				
2	TG		DATA-FACILITY-NUMBER	B	2.0		N
			Not referenced by PASGET.				
2	TH		OPERATOR-ID	B	2.0		N
			Not referenced by PASGET.				
G 2	TK		CORRECTION				
			Not referenced by PASGET.				
3	TL		ORG-OR-COR	A	1.0		F
			Not referenced by PASGET.				
3	TM		REASON	A	78.0		N
			Not referenced by PASGET.				
3	TN		COR-DATE	A	8.0		N
			Not referenced by PASGET.				
3	TO		COR-TIME	A	6.0		N
			Not referenced by PASGET.				
G 1	AA		WGT-ITEMS-SEG				
			Not referenced by PASGET.				
2	BB		WGT-GRP-ID	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				
2	CC		WGT-CT	B	2.0		N
			Retrieved by PASGET for PAS validation reports.				

FILE...: WEIGHT-ITEM-EIS

TYPE...: ADABAS

FILE-NR: 46

PRIMARY SEQUENCE FIELD: DATA-BASE-KEY

T	L	DB NAME	F	LENG	S	DE
---	---	---------	---	------	---	----

P	1	DD ANY-WGT-ITEM				
---	---	-----------------	--	--	--	--

Retrieved by PASGET for PAS
validation reports.

2	EE	ITEM-ID	B	2.0	N	
---	----	---------	---	-----	---	--

Retrieved by PASGET for PAS
validation reports.

2	FF	LOW-WEIGHT	B	4.0	N	
---	----	------------	---	-----	---	--

Retrieved by PASGET for PAS
validation reports.

2	GG	HI-WEIGHT	B	4.0	N	
---	----	-----------	---	-----	---	--

Retrieved by PASGET for PAS
validation reports.

1	SD	SUPER-KEY	A	11.0	SP	
---	----	-----------	---	------	----	--

SOURCE FIELD(S) --- -START- --END-

		DATA-BASE-KEY	1	9		
--	--	---------------	---	---	--	--

		WGT-GRP-ID	1	2		
--	--	------------	---	---	--	--

Used by PASGET to select data
from the file.

WEIGHT-OBJECT-TABLE

Physical layout of Weight Object table. This layout is also stored in the Tables Description table with the field TABLE_NAME set to TBWGTOBJ. TBGET uses that layout description to form direct calls to retrieve data from the table for PAS validation reports. This layout corresponds to the view TBWGTOBJ.

FILE...: WEIGHT-OBJECT-TABLE

TYPE...: ADABAS

FILE-NR: 19

PRIMARY SEQUENCE FIELD: CODE-VALUE

T	L	DB	NAME	F	LENG	S	DE
G	1	WT	WEIGHT-OBJECT-TABLE				
2	AA	ROW-NUMBER		B	2.0	UQ	
			Refer to TBWGTOBJ. Not referenced by TBGET.				
2	BB	DATE-ENTERED		A	8.0		
			Refer to TBWGTOBJ. Not referenced by TBGET.				
2	CC	ACTIVE-INACTIVE		A	1.0	F	
			Refer to TBWGTOBJ. Not referenced by TBGET.				
2	DD	CODE-VALUE		B	2.0	UQ	
			Refer to TBWGTOBJ. Key by which TBGET retrieves data from this table.				
2	EE	EE-FILLER		A	1.0	F	
			Refer to TBWGTOBJ. Not referenced by TBGET.				
2	FF	FF-FILLER		A	1.0	F	
			Refer to TBWGTOBJ. Not referenced by TBGET.				
2	GG	SHORT-TEXT		A	8.0	N	
			Refer to TBWGTOBJ. Used by TBGET for PAS validation reports.				
2	HH	LONG-TEXT		A	32.0	N	
			Refer to TBWGTOBJ. Not referenced by TBGET.				

Appendix A

A-1: Legend for File Listing

Heading	Description	Contents
Name	File Name	32-character alphanumeric string
Typ	File Type	U - User View A - ADABAS File (Physical)
NR	File Number	ADABAS File Number
Description	File Description	Free text describing file

A-2: Legend for Listing of Files with Fields

Heading	Description	Values
T	Type	(blank) - Elementary field G - Group field P - Periodic Group M - Multiple Field
L	Level	1-9
DB	ADABAS Field Name	A 2-character alphanumeric string
Name	Field Name	32-character alphanumeric string
F	Field Format	B - Binary A - Alphanumeric I - Integer N or U - Unpacked Numeric P - Packed Numeric F - Floating Point
Length	Field Length	Number; For numeric strings, length is the sum of values on each side of decimal point, with each value representing the places on the respective sides of the decimal point.
S	Suppression Option	(blank) - Default N - Null Suppression F - Fixed Length
DE	Descriptor	DE - Non-unique Descriptor UQ - Unique Descriptor SB - Sub-Descriptor SP - Super-Descriptor

****End of Document****